axiom^{x3}[®]

Technical and scientific review

Bone preservation

BONE IS GOLD

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KEY CHAPTERS

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REFERENCES

- [1] Anthogyr. Axiom X3® international product launch survey - Report T22025, 2022.
- [2] Misch CE, Judy KW. Classification of partially edentulous arches for implant dentistry. *The International Journal of oral implantology* 1987;4(2):7–13.
- [3] Olivier C, Maire E, Richard H, et al. Bone behavior around a self-tapping dental implant during its insertion: European Synchrotron Radiation Facility 2025.
- [4] Lahens B, Neiva R, Tovar N, et al. Biomechanical and histologic basis of osseodensification drilling for endosteal implant placement in low density bone. An experimental study in sheep. *Journal of the Mechanical Behavior of Biomedical Materials* 2016;63:56–65.
- [5] Gaspar J, Proença L, Botelho J, et al. Implant Stability of Osseodensification Drilling Versus Conventional Surgical Technique: A Systematic Review. *Int J Oral Maxillofac Implants* 2021;36(6):1104–10.
- [6] Mercier F, Bartala M, Ella B. Evaluation of the Osseodensification Technique in Implant Primary Stability: Study on Cadavers. *Int J Oral Maxillofac Implants* 2022;37(3):593–600.
- [7] Shah FA, Palmquist A. Evidence that Osteocytes in Autogenous Bone Fragments can Repair Disrupted Canalicular Networks and Connect with Osteocytes in de novo Formed Bone on the Fragment Surface. *Calcified Tissue International* 2017;101(3):321–27.
- [8] Schmitt CM, Nogueira-Filho G, Tenenbaum HC, et al. Performance of conical abutment (Morse Taper) connection implants: A systematic review. *Journal of Biomedical Materials Research Part A* 2014;102(2):552–74.
- [9] Zipprich H, Weigl P, Lange B, et al. Micromovements at the implant-abutment interface: Measurement, causes, and consequences (in German). *Implantologie* 2007;1:31–46.
- [10] Assenza B, Tripodi D, Scarano A, et al. Bacterial leakage in implants with different implant-abutment connections: an in vitro study. *J Periodontol*, 2012, 83(4):491-7.
- [11] Adrien, J. et al (2016). Comparative in-situ loading analysis of narrow diameter implants by high resolution X-Ray tomography, Poster presented at the 25th EAO's Annual Scientific Congress, Paris, France.
- [12] ISO 5832-3. Medical devices — Implants for surgery — Metallic materials — Part 3: Wrought titanium 6-aluminium 4-vanadium alloy.
- [13] ASTM F136. Standard Specification for Wrought Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications (UNS R56401).
- [14] Liens A, Fabregue D, Courtois N, et al. Multi-scale microstructural and mechanical characterization of titanium and titanium alloys used in dentistry. Poster presented at the 24th EAO's Annual Scientific Congress, Stockholm, Sweden. 2015.
- [15] Imam AY, Moshaverinia A, McGlumphy EA. Implant-abutment interface: A comparison of the ultimate force to failure among narrow-diameter implant systems. *Journal of Prosthetic Dentistry* 2014;112(2):136–42
- [16] Hirata R, Bonfante EA, Machado LS, et al. Mechanical Evaluation of Two Grades of Titanium Used in Implant Dentistry. *INT J ORAL MAXILLOFAC IMPLANTS* 2015;30(4):800–05.

- [17] Velasco-Ortega E, Jos A, Camean AM, et al. In vitro evaluation of cytotoxicity and genotoxicity of a commercial titanium alloy for dental implantology. *Mutation research* 2010;702(1):17–23.
- [18] Chacun D, Lafon A, Courtois N, et al. Histologic and histomorphometric evaluation of new zirconia-based ceramic dental implants: A preclinical study in dogs. *Dental materials official publication of the Academy of Dental Materials* 2021;37(9):1377–89.
- [19] Hoornaert A, Vidal L, Besnier R, et al. Biocompatibility and osseointegration of nanostructured titanium dental implants in minipigs. *Clinical oral implants research* 2020;31(6):526–35.
- [20] D.Duddeck, Department for oral surgery and implantology, University of Köln, Germany. *Implant Study 2014 /2015 Quantitative and qualitative element analysis of implant surfaces by SEM/EDX.*
- [21] Wennerberg A, Albrektsson T. Suggested guidelines for the topographic evaluation of implant surfaces. *The International journal of oral & maxillofacial implants* 2000;15(3):331–44.
- [22] Albrektsson T, Wennerberg A. Oral implant surfaces: Part 1--review focusing on topographic and chemical properties of different surfaces and in vivo responses to them. *International Journal of Prosthodontics* 2004;17(5):536–43.
- [23] Albrektsson T, Wennerberg A. Oral Implant Surfaces: Part2-Review Focusing on Clinical Knowledge of different surfaces. *International Journal of Prosthodontics* 2004;17:544–64.
- [24] Le Guehennec L, Lopez-Heredia M-A, Enkel B, et al. Osteoblastic cell behaviour on different titanium implant surfaces. *Acta biomaterialia* 2008;4(3):535–43.
- [25] Le Guehennec L, Goyenvalle E, Lopez-Heredia MA, et al. Histomorphometric analysis of the osseointegration of four different implant surfaces in the femoral epiphyses of rabbits. *Clinical oral implants research* 2008;19(11):1103–10.
- [26] Habibovic P, Li J, van der Valk CM, et al. Biological performance of uncoated and octacalcium phosphate-coated Ti6Al4V. *Biomaterials* 2005;26(1):23–36.
- [27] Shams N, Ghasemi M, Sadatmansouri S, et al. Morphology and differentiation of MG63 osteoblast cells on saliva contaminated implant surfaces. *Journal of Dentistry (Tehran, Iran)* 2015;12(6):424.
- [28] Anthogyr. Internal data.
- [29] Albrektsson T, Eriksson AR, Friberg B, et al. Histologic investigations on 33 retrieved Nobelpharma implants. *Clinical materials* 1993;12(1):1–9.
- [30] Chacun D, Gritsch K, Droguet G, et al. Retrospective cohort study of bone level implants up to 7years of follow-up: Poster presented at Europerio 2022. *Journal of Clinical Periodontology* 2022;49(S23):143–288.

C225_EN - 2023-06

Photo credits: Anthogyr. All rights reserved

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