M series

by  amanngirrbach

 **straumann**
simply doing more

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User Manual

Mode d'emploi

Istruzioni d'uso

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Deutsch 3 - 42

English 43 - 78

Français 79 - 117

Italiano 118 - 156

Español 157 - 196

Português 197 - 237

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Tab. 6

- Translation of the original Instruction Manual -

Table of Contents

1	Explanation of Symbols	44	7	Handling the Coolant/Lubricant (C/L) ..	67
2	General Safety Instructions	45	7.1	C/L analysis	67
3	Suitable Personnel	45	7.2	Removing the C/L container	67
4	Machine Specifications	46	7.3	Emptying the C/L container	69
4.1	Milling machine	46	7.4	Cleaning the C/L container	69
4.2	Supply unit	50	7.5	Filling in C/L	69
5	Installation	53	7.6	Installing the C/L container	70
5.1	Milling machine	53	7.7	Cleaning the C/L system	70
5.2	Maintenance unit	57	8	Cleaning and Maintenance	71
5.3	Installation of the software	58	8.1	Milling machine	71
5.4	Establishing the Network Connection .	59	8.2	Supply unit	75
5.5	Starting-up the milling machine ...	60	9	Malfunctions, Repairs and Warranty ..	76
6	Application and Operation	61	9.1	Malfunctions	76
6.1	Switching between wet and dry processing	62	9.2	Repairs	76
6.2	Changing the blank holder	64	9.3	Warranty	76
6.3	Checking the Model Area	64	10	Environmental Protection	76
6.4	Function of the buttons	65	11	Technical Data and Accessories/Spare Parts	77
6.5	Reference run	65			
6.6	Program execution	65			
6.7	Checking the workpiece during a running program	66			
6.8	Moving the the service position ...	66			
6.9	Inspection-window cleaning	66			
6.10	Moving to the change position for the blank holder	66			

1 Explanation of Symbols

Warning indications

Warning indications in the text are marked with a colour-backed triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- _ **NOTE** means that property damage can occur.
- _ **CAUTION** means that light to fairly serious personal injury can occur.
- _ **WARNING** means that serious personal injury can occur.
- _ **DANGER** means that serious personal injury can occur.

Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
—	Item of a list
▪	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

Tab. 1

Other symbols on the machine

Symbol	Meaning
	USB port (socket type B)
	Network connection (Ethernet)
	Slot for SD card
	Control cable connection, extraction
	Fuse
	Compressed-air connection
	Control cable connection, supply unit
	Supply connection, coolant/lubricant

Tab. 2

2 General Safety Instructions

When installing, starting-up and operating the machine, always observe the following safety instructions:

CAUTION:

A defective machine can lead to malfunctions!
On detection of damage or a functional defect of the machine:

- ▷ Label the machine as defective.
- ▷ Prevent further operation until the machine has been repaired.

NOTE:

Possible machine damage through insufficient extraction during dry processing!

- ▷ Operate the machine only with the extraction system recommended by the manufacturer or another compatible extraction system (see page 58).

NOTE:

Possible machine damage through insufficient cooling during grinding operations!

- ▷ For grinding operations, operate the machine only with the supply unit specially intended for this purpose (see page 50ff).

NOTE:

Damage through escaping coolant/lubricant!

- ▷ Regularly check the coolant/lubricant connection for tightness against leaks.

NOTE:

- ▷ Switch the machine off when not in use or unsupervised for longer periods, e.g., overnight. This method of saving electrical energy also benefits the environment.

3 Suitable Personnel

NOTE:

Starting-up and operation of the machine may only be carried out by trained specialised personnel.

4 Machine Specifications

4.1 Milling machine

This User Manual describes different models of the milling machine. Therefore, the representations can deviate from the actual appearance of the machine.

4.1.1 Delivery Scope

- _ Milling machine
- _ Mains cable
- _ USB cable (5 m)
- _ Network cable (5 m)
- _ Control cable, extraction
- _ Compressed-air hoses (250 mm and 2000 mm)
- _ Cleaning hose with brush head
- _ Maintenance unit
- _ Service set for Jäger spindle
- _ Blank template, zirconium oxide opaque
- _ Test blank
- _ CNC cutter, 2.5 (milling tool mounted in collet)
- _ Screw set
- _ Torque wrench, 4 Nm
- _ Blank holder 71
- _ Extractor adapter
- _ Extraction hood
- _ Drain cover
- _ Slide caliper

The milling machine is delivered with a test blank inserted in the blank holder. Prior to delivery, a test object was factory-milled in order to check the zero point settings. Possibly remaining milled shavings/chips are not a reason for rejection. The test blank is provided for test milling and calibration of the milling machine.

- ▷ After unpacking, check the machine for completeness and possible transport damages. Please claim any transport damages immediately with your supplier.

4.1.2 Intended Use

The milling machine is a PC-controlled milling and grinding machine for the fabrication of dentures by means of:

- _ Dry processing of blanks made of
 - pre-sintered zirconium oxide
 - unsintered non-precious metals (NPMs)
 - wax
 - plastics
- _ Wet processing of blanks made of
 - sintered glass ceramics
 - plastics

Using blanks and tools not approved by the manufacturer can damage the machine and make the product unusable. For such cases, the manufacturer shall assume no liability whatsoever.

Unauthorised modifications/alterations of the machine shall also void the warranty.

A compressed-air connection is required for operation of milling machine (see page 57).

Dry processing is permitted only in conjunction with an appropriate extraction system (for technical requirements, see page 58).

Wet milling and grinding is permitted only in conjunction with the supply unit (accessory).

4.1.3 CE Declaration of Conformity

In terms of design and performance, this product complies with the European Directives and the supplementary national requirements. Conformity has been confirmed with the CE marking.

The product's Declaration of Conformity can be requested from the manufacturer.

4.1.4 Components and Interfaces


Fig. 1 Machine overview, front

- 1 Measuring device, tool length
- 2 C/L nozzles
- 3 Milling spindle
- 4 Extraction hood
- 5 Service button
- 6 Operation button
- 7 C/L drain
- 8 Tool holder
- 9 Blank holder


Fig. 2 Machine overview, rear

- 1 USB data port 
- 2 Ethernet interface 
- 3 Control interface for supply unit 
- 4 Slot for SD memory card 
- 5 Control interface for extraction 
- 6 Compressed-air inlet
- 7 Compressed-air hose, 250 mm
- 8 Maintenance unit
- 9 Compressed-air connection, milling spindle
- 10 Connection, C/L supply line 
- 11 Connection, C/L drain 
- 12 Fuse compartment
- 13 Power supply connection
- 14 Main switch
- 15 Extraction-hose connection

4.1.5 Accessories

Supply unit

The supply unit (see section 4.2) is mandatory for wet processing. It provides for the necessary coolant/lubricant circuit and includes the coolant/lubricant filtration. Additionally, it offers space for the extraction.

Extraction

The extraction recommended by the manufacturer is required for dry processing with the milling machine. It is equipped with a special filter bag and a HEPA micro filter, which filters up to 99.97 % of fine dust particles (i.e. zirconium oxide dust particles) and is classified in the fine dust particle category M.

The separate socket for the control cable allows for automatic operation.

The extraction can be inserted into the supply unit.

Milling and grinding tools

Various cutting and grinding tools are available for the milling machine (see catalogue).

4.2 Supply unit

4.2.1 Delivery Scope

- _ Supply unit with
 - Container for coolant/lubricant (C/L)
 - Filter for (C/L)
- _ Control cable
- _ C/L supply line
- _ Piping kit
- _ Mains cable
- _ Swarf sieve
- _ Splash guard
- _ 2 litres of coolant/lubricant
- _ Analytic test strips
 - Multifunction additive test
 - pH-value
 - Nitrite test
 - Test record

4.2.2 Intended Use

The supply unit is for wet processing with the milling machine. It ensures the application-efficient supply, drainage and filtering of the required coolant/lubricant.

Operation with other machines than the milling machine is not permitted.

4.2.3 CE Declaration of Conformity

In terms of design and performance, this product complies with the European Directives and the supplementary national requirements. Conformity has been confirmed with the CE marking.

The product's Declaration of Conformity can be requested from the manufacturer.

4.2.4 Components and Interfaces


Fig. 3 Machine overview, front

- 1** C/L drain for set-up under a table
- 2** C/L drain for set-up directly below the milling machine
- 3** Retainer for adjustable legs of milling machine


Fig. 4 Machine overview, rear

- 1 Control interface for supply unit
- 2 Power supply connection
- 3 Control interface for extraction
- 4 C/L supply line
- 5 C/L drain for set-up directly below the milling machine
- 6 C/L drain for set-up under a table
- 7 C/L supply line
- 8 C/L container
- 9 C/L drain cock

5 Installation

5.1 Milling machine

Countries with 100-115 V mains voltage

The machine is preset for a mains voltage of 230 V. When the mains voltage at the set-up location is 100-115 V:

- ▷ Pull out the fuse compartment ([11] in Fig. 2 on page 49) below the mains plug.
- ▷ Replace both T3.15 A fuses with the provided T6.3 A fuses.

5.1.1 Set-up of the milling machine without supply unit

Please refer to the separate milling machine set-up instructions.

- _ The machine is intended exclusively for use within dry, closed rooms.
 - _ A clearance of at least 200 mm is to be observed sideways, to the rear and upwards. Minimum space requirements for the milling machine (incl. connections):
 - Width: 970 mm
 - Depth: 850 mm
 - Height: 940 mm
 - _ The empty weight of the machine is between 75 kg and 78 kg. The set-up surface must withstand an appropriate load.
 - _ During operation, the room temperature must be between 18 °C and 30 °C; avoid large variations in temperature.
- ▷ Set up the machine on a sturdy work table or work bench (off the ground and clear of walls).
 - ▷ Adjust the height of the machine by screwing the feet in or out, ensuring that the machine is positioned firmly on all four feet.

- ▷ Insert the extraction hose of the extraction into connection [7] on the rear side of the milling machine.
- ▷ Connect the compressed-air hose on the rear side of the milling machine.

5.1.2 Set-up of the milling machine with supply unit

Please refer to the separate set-up instructions for milling machine and supply unit.

In conjunction with the supply unit, the milling machine is set up either directly on the supply unit (tower version) or on a table above the supply unit (table version).

Installation of the extraction in the supply unit

When dry processing with air extraction via the extraction, it can be installed in the supply unit (see extraction set-up instructions). This simplifies the cabling.


Fig. 5

- 1 Supply unit
- 2 Extraction

- ▷ Set the extraction to the operating mode AUTO.
- ▷ Insert the extraction with the extraction opening facing toward the rear into the supply unit.

Malfunction!

- ▷ Do **not** plug in the control cable of the extraction to the milling machine!

- ▷ Plug in the control cable of the extraction at the rear side of the supply unit.
- ▷ Insert the extraction hose of the extraction into the milling machine.

Set-up of the milling machine directly on the supply unit (tower version)

A clearance of at least 200 mm is to be observed sideways, to the rear and upwards. Minimum space requirements for the milling machine on the supply unit (incl. connections):

- _ Height: 1591 mm
- _ Width: 970 mm
- _ Depth: 860 mm

The set-up surface must be level.


Fig. 6

- 1 Milling machine
- 2 Supply unit

NOTE:

Possible damage of the electric lines and the C/L supply line:

- ▷ When moving the supply unit, pay attention that no lines are driven over.

- Place the milling machine onto the supply unit, ensuring that the legs engage into the retainers.
- Mount the service unit (see page 57).
- Remove the C/L container of the supply unit (see section 7.2 on page 67).
- Fill coolant/lubricant into the C/L container (see section 7.5 on page 69).
- Reinstall the C/L container (see section 7.6 on page 70).
- Connect the C/L supply line to the milling machine and the supply unit.


Fig. 7

- 1 Milling machine
- 2 Supply unit
- 3 C/L drain pipe

- Plug the C/L drain pipe [3] together and connect it to the milling machine and the supply unit (see set-up instructions).
- Move the supply unit with the milling machine to the set-up location.
- Lock the casters of the supply unit with the brakes.
- Insert the extraction hose of the extraction to the extraction-hose connection ([15] in Fig. 2 on page 49) on the rear side of the milling machine.

Milling machine set-up on a table above the supply unit (table version)

Under the table, a clearance of at least 200 mm is to be observed toward the sides and the rear, and at least 10 mm upwards. Minimum space requirements for the supply unit below the table:

- _ Height below the table: 675 mm
- _ Width below the table: 970 mm
- _ Depth below the table: 860 mm

The set-up surface must be level.

For the set-up dimensions of the milling machine, see page 53.


Fig. 8

- 1 Milling machine
- 2 Supply unit

In order for the C/L drain pip to be mounted, the table must be set up with a certain clearance to the wall or have an appropriate cut-out for the piping of the CL/ drain (see set-up instructions).

- Position the milling machine on the table.

- ▷ Adjust the height of the machine by screwing the feet in or out, ensuring that the machine is positioned firmly on all four feet.
- ▷ Mount the service unit (see page 57).
- ▷ Remove the C/L container of the supply unit (see section 7.2 on page 67).
- ▷ Fill coolant/lubricant into the C/L container (see section 7.5 on page 69).
- ▷ Reinstall the C/L container (see section 7.6 on page 70).

NOTE:

Possible damage of the electric lines and the C/L supply line:

- ▷ When moving the supply unit, pay attention that no lines are driven over.

- ▷ Move the supply unit below the table in such a manner that both units are aligned with respect to each other.


Fig. 9

- 1** Milling machine
- 2** Supply unit
- 3** Connection for C/L return pipe
- 4** C/L drain pipe

- ▷ Connect the C/L supply line to the milling machine and the supply unit.
- ▷ Plug the C/L drain pipe [4] together and connect it to the milling machine and the supply unit (see set-up instructions).
- ▷ Insert the extraction hose of the extraction to the extraction-hose connection ([15] in Fig. 2 on page 49) on the rear side of the milling machine.

5.2 Maintenance unit

The Jäger spindle is equipped with a sealing-air feature. This sealing air prevents shavings/chips and dust from entering the spindle. The service unit filters possible contamination in the sealing air and thus prevents costly machine damage.

Installing the service unit


Fig. 10 Maintenance unit

- 1 Knurled screws
- 2 Compressed-air hose, 250 mm
- 3 Compressed-air inlet
- 4 Water separator
- 5 Compressed-air connection, milling spindle
- 6 Connection, C/L drain pipe
- 7 Connection, C/L supply line

- Mount the service unit on the rear side of the machine using the knurled screws [1].
- Connect the 250 mm compressed-air hose [2] to the right-hand connection of the service unit and to the compressed-air connection [5] of the milling machine.

- Connect the 2,000 mm compressed-air hose to the compressed-air inlet [3] of the service unit and the other end to the compressed-air supply of the laboratory.

Requirements for the compressed-air supply

For operation, the compressed air must fulfil the following minimum requirements:

- _ Flow rate: 50 l/min (1.77 CFM)
- _ Pressure: 6 bar / 87 psi

Requirements on air purity:

- _ Solid contaminants: Class 3;
better than 5 µm for solids
- _ Water content: Class 4;
Max. pressure dew point +3°C
- _ Total oil content: Class 2;
Max. oil content 0.1 mg/m³

In continuous operation, the pressure and the air-flow rate are monitored by the milling machine. When the required values are fallen below, the milling program stops. Once the correct air values are available again, the milling program can be continued after fault acknowledgement in the software.

Requirements on air purity for the Jäger spindle

When another service unit than the one supplied is to be used, the following connection conditions according to

DIN-ISO 8573-1 must be ensured:

- _ Solid contaminants: Class 3
(particle size max. 5 µm, particle density max. 5 mg/m³)
- _ Water content: Class 4
(max. pressure dew point +3 °C, water content max. 6000 mg/m³)
- _ Total oil content: Class 2;
(max. oil content 0.1 mg/m³)

Extraction requirements

When another extraction system is to be used, it must meet the following requirements:

- _ Suction capacity: 56.6 l/min
- _ Suitability for zirconium-oxide dusts
- _ HEPA microfilter (97.97 %), filter class H12, dust category M

The connection of the extraction hose to the machine is made with an adapter (in delivery scope). The adapter fits extraction hoses with a 38 mm interior diameter.

5.3 Installation of the software

NOTE:

Malfunction!

► Operate the milling machine only with the PCs recommended by the manufacturer.

- Insert the data carrier into the drive of the PC.
The data carrier contains an installation guide.
- Start the installation guide and follow the instructions.

Connecting the interfaces


Fig. 11 Interfaces on the rear side of the milling machine

- 1 USB data port
 - 2 Ethernet interface
 - 3 Control interface for supply unit
 - 4 Slot for SD memory card
 - 5 Control interface for extraction
 - 6 Main switch
 - 7 Power supply connection
- Make sure that the main switch [6] of the milling machine is in the 0 position.
 - Check if the SD memory card is properly seated in the slot [4] on the rear side of the milling machine.
 - Connect the mains cable to power supply connection [7] of the milling machine as well as to a mains socket outlet.

For installation without supply unit:

- Connect the control cable to the control interface for extraction [5] on the rear side of the milling machine and to the interface on the rear side of the extraction system.

For installation with supply unit:


Fig. 12 Interfaces on the rear side of the supply unit

- 1 Power supply connection
- 2 Connection for C/L supply line
- 3 Control interface for supply unit
- 4 Control interface for extraction

- ▷ Connect the control cable to the control interface for the supply unit on the rear side of the milling machine ([3] in Fig. 11) and to the rear side of the supply unit ([3] in Fig. 12).
- ▷ Connect the control cable of the extraction system to the control interface for extraction on the rear side of the supply unit ([4] in Fig. 12).
- ▷ Connect the mains cable to power supply connection ([1] in Fig. 12) of the supply unit as well as to a mains socket outlet.

Connecting to the PC

- ▷ Connect a USB cable to the USB interface ([1] in fig. 11) of the milling machine and to a free USB port of the PC.
- The maximum permitted length of the USB cable is 5 m.

5.4 Establishing the Network Connection

To establish the network connection, the milling machine and the PC must first be connected via USB.

- ▷ Connect the Ethernet interface [2] of the milling machine with the PC or with the network using a network cable.
- The maximum permitted length of the network cable to the PC or to the network router is 100 m.
- ▷ Switch the milling machine on.
 - ▷ Switch the supply unit on.
 - ▷ Start the milling machine software.
 - ▷ After initializing the milling machine in the software under *Settings > Communication*, set the connection type from USB to Ethernet and adapt the IP address, if required (see description of the milling machine software).
 - ▷ Switch the milling machine off and then on again.
 - ▷ Initialize the milling machine in the software.
 - ▷ After the software has detected the Ethernet connection, remove the USB cable.

When operating via a network, up to eight machines can be connected and controlled via the CAM software.
The PC and the milling machine must both be in the same sub-network.

5.5 Starting-up the milling machine

The three standard cutters in tool stations 1-3 are mandatory for operation.

In addition to standard tools, sister tools can optionally be inserted in the tool holder. When a standard tool has reached its wear limit, the milling machine automatically accesses the respective sister tool.


Fig. 13 Tool holder

- ▷ Insert the CNC cutter standard tools into the tool holder of the milling machine. Position the tool shank facing toward the outside.
 - CNC cutter 2.5 is in the spindle.
 - Set CNC cutter 1.0 to position 2.
 - Set CNC cutter 0.6 to position 3.
- ▷ Optionally: Insert the sister tools into the tool holder of the milling machine. Position the tool shank facing toward the outside.
 - Set CNC cutter 2.5 to position 4.
 - Set CNC cutter 1.0 to position 5.
 - Set CNC cutter 0.6 to position 6.

- ▷ Insert the selected blank into the holder in the milling machine. The high/thick side of the blank must face the spindle.


Fig. 14 Inserting a blank

- ▷ Fasten the blank using four screws and washers.

NOTE:

Do not damage the blank!

- ▷ Do not overtighten the screws! Avoid mechanical stress!

- ▷ Tighten the screws with the Allen key provided.


Fig. 15 Tightening a blank

- ▷ Mount the respective hood.
 - Extraction hood for dry processing
 - Splash guard for wet processing

For wet processing:

- Insert the swarf sieve with the wide side facing toward the rear.

Without the swarf sieve, chips/swarf will be flushed into the container. The filter bag will fill up more quickly.

The swarf sieve is included in the delivery scope of the supply unit. It can remain inserted in the machine during dry processing.

- Shut the door of the milling machine.
- Switch the milling machine on via the main switch.
- Switch the PC on.
- Start the milling machine software.

When the software is started, a machine initialization takes place.

If the milling machine is not detected:

- Restart the PC.

After re-positioning the milling machine, it is necessary to check the system axes.

- For this, observe the instructions for milling the test bodies.

6 Application and Operation**WARNING:**

Danger from shavings/chips or tool fragments being thrown from the machine!

- Make sure that the door of the milling machine is always closed during processing!

**NOTE:**

Insufficient milling results and increased wear when dry processing!

- Do not operate the milling machine without extraction!

**NOTE:**

Insufficient milling/grinding results and increased wear when wet processing!

- Do not operate the milling machine without supply unit!

6.1 Switching between wet and dry processing

The milling machine can be used both for wet and dry processing.

In delivery condition, it is set up for dry processing.

6.1.1 Switching from dry to wet processing


Fig. 16

- 1 Extraction hood
- 2 Lever
- 3 Swarf sieve
- 4 Splash guard
- 5 Extraction hood
- B Service button

- ▷ Disengage lever [2] and pull off the white extraction hood [1] toward the front.
- ▷ Connect the extraction hose to the extraction opening in the milling chamber.
- ▷ With the door open, briefly press service button [B] on the front side of the milling machine. The extraction starts.
- ▷ Vacuum out the milling chamber.

When the milling chamber is clean:

- ▷ Press service button [B] again.

The extraction stops.

- ▷ Pull off the extraction hose.

- ▷ Insert the cleaned swarf sieve [3].

The wide side faces toward the rear.

For wet processing, it is mandatory to use the swarf sieve [3]. The sieve serves as the first coolant/lubricant stage for coarse particles.

- ▷ Check if the outlet holes of the C/L nozzles are unobstructed.
- ▷ Check if the door gasket is clean and free of swarf/chips.

▷ Check if the C/L drain pipe is free of obstructions and if the piping on the outside is securely connected.

▷ Mount the transparent splash guard [4] to the extraction opening and lock it in place with lever [5].

▷ Transfer the milling file for wet processing.

6.1.2 Switching from wet to dry processing


Fig. 17

- 1** Extraction hood
- 2** Lever of the extraction hood
- 3** Swarf sieve
- 4** Splash guard
- 5** Lever of the splash guard
- B** Service button

- ▷ Disengage lever [5] and pull off the transparent splash guard [4] toward the front.
- ▷ Remove the swarf sieve [3].
- ▷ Clean and dry the milling chamber and the swarf sieve (see page 71).

- ▷ Insert the blank and insert or change the tooling, as required.
- ▷ Mount the white extraction hood [1] to the extraction opening and lock it in place with lever [2].
- ▷ Check the air extraction.
- ▷ Transfer the milling file for dry processing.

The swarf sieve [3] is not required for dry processing.

6.2 Changing the blank holder

In delivery condition, the milling machine is equipped with a 71 blank holder. For processing varying materials, it may be required to use a different blank holder.

To change the blank holder:

- ▷ With the door closed, press operation button [A] for at least three seconds.

The current tool is placed down. The blank holder [2] moves to the change position.


Fig. 18 Changing-position of the blank holder

1 Fastening screw for blank holder

2 Blank holder

- ▷ Loosen fastening screw [1].
- ▷ Depending on the model, pull off the blank holder toward the right or toward the bottom.
- ▷ Clean the locating bore and surfaces.
- ▷ Check the pin of the blank holder and clean if required.
- ▷ Insert the blank holder with the pin to the stop in the locating bore.

NOTE:

Inaccurate working results or damage to the blank-holder fixture possible!

- ▷ Tighten fastening screw [1] with a torque of 4 Nm.

The torque wrench included in the delivery scope is appropriately set.

- ▷ Tighten fastening screw [1] with the supplied torque wrench.

- ▷ Briefly press operation button [A].

The machine moves to the starting position.

6.3 Checking the Model Area

The model area, on which the work is to be machined, must fit within the contour of the enclosed template. Otherwise, the work cannot be machined in the blank.

The template takes into account the shrinkage of the zirconium oxide in the final sintering.


Fig. 19 Template for zirconium oxide opaque

- ▷ For work in wax or plastic, measure out the model area.

6.4 Function of the buttons


Fig. 20 Buttons

Key	Door shut	Door open
Operation button [A]	Starts or stops the reference run; starts or stops the current program; activates or deactivates the control position; activates the change position for the blank holder	Activates or deactivates the extrac-tion system
Service button [B]	Interrupts the current program; activates the service position; activates the tool-changing position; starts a cleaning run of the inspection window	Activates or deactivates the sealing air

Tab. 3

Instead of with the buttons, these functions can also be actuated via the software (see instructions of the milling machine software).

6.5 Reference run

After switching on, a reference run must be carried out to determine the mechanical zero points. Operation button [A] flashes.

To start the reference run:

- ▷ Briefly press operation button [A].

The reference run is started. Operation button [A] lights up.

The reference run can be stopped:

- ▷ Briefly press operation button [A].

The reference run is stopped. Operation button [A] flashes.

To continue the reference run:

- ▷ Briefly press operation button [A].

The reference run is continued at the position where it was stopped. Operation button [A] lights up.

On completion of the reference run, the machine returns to the starting position. Operation button [A] goes out.

6.6 Program execution

Before executing a program, it must first be loaded (see CAM software). The respective program (milling file) is generated via the CAM software.

- ▷ Briefly press operation button [A].

The program is started. Operation button [A] lights up.

The program can be stopped:

- ▷ Briefly press operation button [A] or service button [B].

-or-

- ▷ Open the door.

The program is interrupted. Operation button [A] flashes.

To continue the program:

- ▷ With the door closed, briefly press operation button [A].

The program is continued at the position where it was stopped. Operation button [A] lights up.

On completion of the program, the machine returns to the starting position. Operation button [A] goes out.

6.7 Checking the workpiece during a running program

- ▷ Push and hold operation button [A] for at least three seconds.

The program is interrupted. The spindle returns to the rear, the current tool is placed down and the workpiece is turned around. While travelling, operation button [A] is lit; after reaching the position, it flashes.

- ▷ Visually check the workpiece.
- ▷ Briefly press operation button [A].
The tool is picked up again. The spindle travels to the position, where the program was interrupted. While travelling, operation button [A] is lit; after reaching the position, it flashes.
- ▷ Briefly press operation button [A].
The program is continued. Operation button [A] lights up.

6.8 Moving the the service position

The service position can be moved to from the starting position or from an interrupted program.

- ▷ Briefly press service button [B].
The current tool is placed down and the spindle returns to the service position. The spindle is released; the motor is blocked. While travelling, operation button [A] is lit; after reaching the position, it flashes.

To exit the service position:

- ▷ Briefly press operation button [A].

The tool is picked up again. The spindle travels back to the starting position or to the position, where the program was interrupted. While travelling, operation button [A] is lit; after reaching the position, it flashes.

6.9 Inspection-window cleaning

- ▷ With the door closed, press and hold operation button [B] longer than three seconds.

As long as the button is pressed, an air nozzle positioned behind the inspection window blows off the window and cleans it from dust and chips.

The inspection-window cleaning is also carried out automatically when a tool change takes place and after completion of a program. In the course of each automatic cleaning run, three short air-blasts take place.

6.10 Moving to the change position for the blank holder

The change position for the blank holder can only be moved to from the starting position.

- ▷ Push and hold operation button [A] for at least three seconds.

The current tool is placed down. The blank holder is turned in such a manner that the fastening screw becomes accessible. While travelling, operation button [A] is lit; after reaching the position, it flashes.

To exit the change position for the blank holder:

- ▷ Briefly press operation button [A].

The machine moves to the starting position.

7 Handling the Coolant/Lubricant (C/L)

CAUTION:

Health hazards are possible through improper handling of C/L!

- ▷ Use only C/L recommended by the manufacturer!
- ▷ Observe the safety instructions and handling regulations for C/L!

7.1 C/L analysis

For legal and technical reasons, the coolant/lubricant must be analysed weekly.

- ▷ Observe country-specific laws and guidelines on coolant/lubricant!

Regular analysis of the C/L is required for proper function of the milling machine. The results of the analyses must be documented and archived.

- ▷ Check the following coolant/lubricant properties weekly:
 - Concentration
 - pH-value
 - Nitrite content

The analysis is to be carried out with test strips. These test strips are simple to handle and enable the pH-value, the nitrite content and the C/L concentration to be determined.

- ▷ Remove the blank from the holder.
- ▷ Mount the splash guard.
- ▷ Place a clean, flat receptacle in the centre of the interior below the splash guard.
- ▷ Shut the door.
- ▷ Push and hold operation button [A] for approx. eight seconds.

After three seconds, the C/L flow is started.

- ▷ Open the door and remove the receptacle.
- ▷ Insert the reaction zone of the test strip into the C/L.
- ▷ Shake off excess fluid.
- ▷ After one minute, compare the reaction zone with the colour scale and read off the measured value.

The reference values and the information on measures in case of exceeding or falling below the reference values can be found on an information sheet on the inside of the container-access cover in the supply unit.

The analysis protocol can also be placed in the self-adhesive transparent folder located there.

7.2 Removing the C/L container

Only for table version:

- ▷ Remove the C/L drain pipe [1] on the side of the supply unit.


Fig. 21

1 C/L drain pipe

NOTE:

Possible damage of the electric lines and the C/L supply line:

- ▷ When moving the supply unit, pay attention that no lines are driven over.

- ▷ Move the supply unit out from under the table.

For both versions:

- ▷ Remove the left side wall of the supply unit.
- ▷ Pull the C/L drain pipe [1] and the C/L supply line [3] off of the C/L container top.
- ▷ Disconnect the control cable [2] at the partition wall (between the chamber for the extraction and the container).


Fig. 22 Disconnecting the interfaces at the C/L container (tower version)


Fig. 23 Disconnecting the interfaces at the C/L container (table version)

Legend to figure 22 and figure 23:

- 1 C/L drain pipe
- 2 Control cable
- 3 C/L supply line

NOTE:

Possible damage to the drain cock!

- ▷ Move the C/L container only by its handles!

- ▷ Carefully lift the C/L container [2] out of the supply unit.


Fig. 24

- 1 Supply unit
- 2 C/L container

7.3 Emptying the C/L container

NOTE:

Possible damage to the drain cock!

▷ Move the C/L container only by its handles!

- ▷ Place the C/L container [1] with its rear side facing frontward on a table.


Fig. 25

- 1** C/L container
2 Drain cock
3 Disposal receptacle

- ▷ Position a disposal receptacle [3] below the drain cock [2].
- ▷ If required, mount a hose to the drain cock and conduct it into the disposal receptacle.
- ▷ Open the drain cock.
- The C/L flows out of the C/L container into the disposal receptacle.

When the C/L container is empty:

- ▷ Shut the drain cock again.

NOTE:

Dispose of the C/L in accordance with the applicable regulations and laws. Proper disposal of the C/L must be documented.

7.4 Cleaning the C/L container

NOTE:

The filter and the C/L container may only be cleaned with water.

▷ Do not use cleaning agents.

- ▷ Unscrew and remove the four star-knob nuts on the top side of the C/L container.
- ▷ Remove the cover with the drain pipe.
- ▷ Remove the filter bag and empty it.
- ▷ Clean the filter bag.
- ▷ Rinse and clean the C/L container.

7.5 Filling in C/L

NOTE:

The water for the coolant/lubricant must have drinking-water quality. The optimal hardness range is between 10 - 20 °dH (corresponds to 1.8 - 3.6 mmol total hardness per litre). If the hardness is far beyond this range, we recommend the addition of distilled water.

- ▷ Fill 20 litres of water and 1 litre of coolant/lubricant concentrate into a separate container and mix well.
- ▷ Insert the filter bag into the C/L container.
- ▷ Fill the mixed C/L into the C/L container.
- ▷ Reattach the cover with the drain pipe and secure it with the four star-knob nuts.

7.6 Installing the C/L container

- ▷ Carefully lift the C/L container into the supply unit.
- ▷ Connect the control cable at the partition wall (between the chamber for the extraction and the container).
- ▷ connect the return pipe and the C/L supply line at the top of the C/L container.
- ▷ Insert the left side wall of the supply unit.

NOTE:

Danger of damaging cables and hoses!

- ▷ When moving the supply unit, do not pinch or drive over any cables or hoses.

For the tower version:

- ▷ Connect the C/L drain pipe and the C/L supply line of the supply unit to the milling machine.

For the table version:

- ▷ Move the supply unit back under the table and align it with the milling machine.
- ▷ Insert the drain pipe on the side to the stop into the supply unit.

Through this, the control system recognizes the drain pipe and enables the C/L pump.

- ▷ Connect the C/L supply line of the supply unit to the milling machine.

7.7 Cleaning the C/L system

- ▷ Remove and clean the swarf sieve.
- ▷ Clean the milling chamber.
- ▷ Reinsert the swarf sieve.
- ▷ Mount the splash guard.
- ▷ Empty the C/L container (see page 69).
- ▷ Clean the C/L container (see page 69).
- ▷ Fill 5 litres of water and 250 ml of coolant/lubricant concentrate into a separate container and mix well.
- ▷ Insert the filter bag.
- ▷ Fill the mixed C/L into the C/L container.
- ▷ Fully connect the supply unit again.
- ▷ Activate "Clean coolant/lubricant circuit" in the control service functions.
- ▷ Start the cleaning in the software

-or-

- ▷ briefly press operation button [A] on the milling machine
- C/L is circulated for 15 minutes to clean the system.

The cleaning can be interrupted and restarted by briefly pressing operation button [A].

When the cleaning is completed:

- ▷ Empty the C/L container (see page 69).

After the cleaning, the C/L used for the cleaning must be disposed of.

- ▷ Refill the C/L container (see page 69).

8 Cleaning and Maintenance

8.1 Milling machine

8.1.1 Cleaning

The milling machine must be cleaned after each working day.

- ▷ Empty the swarf sieve frequently, especially when wet grinding plastic.

When the swarf sieve is too full, it is possible that the coolant/lubricant will not drain off.

Cleaning after dry processing

NOTE:

Zirconia dust is abrasive!

- ▷ Clean the window in the door carefully, using only a soft cloth; do not scrub or scour.

- ▷ Open the door of the milling machine.
- ▷ Replace the extraction hood with the cleaning hose with brush head.
- ▷ Press operation button [A] on the milling machine.
The extraction system is activated.. The operation button lights up.
- ▷ Press service button [B] on the milling machine.
The sealing air at the spindle is activated. The service button lights up.
- ▷ If required, remove the swarf sieve and dispose of accumulated milling dust/shavings away from the machine.
In this manner, the filter bag in the extraction system will not fill up as quickly and can be used longer.
- ▷ Thoroughly vacuum out the milling chamber with the cleaning hose.

Residual zirconia dust from dry processing can clog the filter bag in case of subsequent wet processing.

- ▷ Briefly press service button [B].
The sealing air is deactivated.
- ▷ Briefly press operation button [A].
The extraction system is deactivated.

Cleaning after wet processing

NOTE:

The extraction is not designed for wet vacuuming. The extraction may not be used for cleaning the interior after wet processing!

- ▷ Remove shavings/chips manually.

NOTE:

The milling chamber may only be cleaned with water.

- ▷ Do not use cleaning agents.

Dried-in coolant/lubricant is difficult to remove.

After wet processing:

- ▷ Clean and dry the milling chamber prior to downtimes lasting several hours.

- ▷ Remove the swarf sieve and dispose of accumulated milling dust/shavings.
- ▷ Clean the milling chamber and wipe it dry using a soft cloth.

8.1.2 Regular inspections

- _ Condition of the bellow
- _ Condition of hoses and drain pipes
- ▷ Keep the door gasket and the stopping surface on the door clean. Ensure not to damage the sealing lips.

8.1.3 Weekly maintenance

NOTE:

- Danger of damaging the machine!
- ▷ Do not use compressed air, ultrasonic or a steam jet for the cleaning!

Prior to maintenance:

Before maintaining the spindle and the C/L nozzles, the spindle sealing air must be activated:

- ▷ Briefly press service button [B].
The current tool is placed down and the spindle returns to the service position. The spindle is released; the motor is blocked.
- ▷ Open the door of the milling machine.
- ▷ Press service button [B] on the milling machine.
The sealing air is activated. The service button lights up.

After the maintenance

When the maintenance is completed:

- ▷ Shut the door of the milling machine.
- ▷ Briefly press operation button [A].
The tool is picked up again. The spindle travels back to the starting position or to the position, where the program was interrupted.

Spindle

The provided Jäger spindle service set and a milling tool of your choice are required for maintenance of the spindle.


Fig. 26 Service set for Jäger spindle

- 1 Collet grease
- 2 Collet brush
- 3 Felt cone
- 4 Ejector pin
- 5 Spindle wrench

In order to clean the spindle:

- ▷ Remove any dust from the cap, the C/L pipes and the spindle using a dry, soft cloth. Do not use cleaning agents!

To remove the collet:

- ▷ Apply the black spindle wrench [5] over the collet.


Fig. 27

- ▷ Insert a milling tool into the open collet.

- ▷ Using the spindle wrench, unscrew the collet out of the spindle.


Fig. 28

- ▷ Clean the collet holder using the felt cone [3].
The collet holder must be free of shavings/chips and debris.


Fig. 29

- ▷ Clean the inside of the collet with the fine collet brush [2].


Fig. 30

- ▷ Apply a light coat of the collet grease ([1] in Fig. 26 on page 72) from the service set to the outer cone of the collet.

This improves the sliding properties and increases the clamping force of the collet.

NOTE:

Only use the provided collet grease for lubrication of the collet cone. Do not use other types of grease.

- ▷ Screw the cleaned and lubricated collet in reverse order to the stop back into the spindle.

CAUTION:

Danger of tool breakage!

- ▷ Make sure to remove the milling tool from the collet chuck!

When the milling tool was from the milling machine's tool holder:

- ▷ Place the tool properly positioned back into the tool holder.

C/L nozzles

Do not clean the C/L nozzles with the col-
let brush!

- ▷ Inspect the C/L nozzles weekly.
- ▷ Clean the outlet holes as required.

When the outlet holes are clogged:

- ▷ Loosen the Allen screw on the side of the spin-
dle cap.


Fig. 31

- ▷ Pull off the cap.


Fig. 32

- ▷ Blow out the holes from the inside of the cap
using a steam cleaner.

CAUTION:

The cap heats up!

- ▷ Apply a coat of grease to the C/L pipes.
- ▷ Clean the seat of the cap on the spindle.
- ▷ Mount the spindle cap again with the screw at
the right-hand side. Pay attention to the correct
connection of the C/L pipes.

Service unit

- ▷ Visually check the service unit for:

- Visible particles, debris or particulate matter
within the acrylic glass (signs for contami-
nated compressed air)
- Amber discolouration of the filter element
(sign of oil in the compressed air)
- Water in the filter (sign of water in the com-
pressed air)

Once any of the listed impairments is given:

- ▷ Correct the problem with the compressed air.
- ▷ Replace the complete service unit.

8.1.4 Monthly Maintenance

Checking the zero point

It is recommended to check the zero point once every month.

▷ Mill a test cube from the test blank.

The test blank is to be used exclusively for test milling and calibration of the milling machine.

Jäger spindle

The spindle must be in the starting position.

▷ Manually turn the shaft of the spindle through (at least ten times).


Fig. 33 Turning the Jäger spindle through

8.1.5 External maintenance

External maintenance is necessary after 1000 operating hours. The milling machine software will inform you appropriately.

8.2 Supply unit

8.2.1 Inspections and maintenance measures

Regular inspections:

- _ Condition of hoses and drain pipes
- _ C/L filling level in the supply unit.

When the filling level in the C/L container has fallen below the minimum level, the milling program is interrupted and an error message is indicated.

- ▷ Fill the C/L container (see page 69).
- ▷ Acknowledge the error in the software.

8.2.2 Weekly maintenance

▷ Analyse the C/L (see separate analysis instructions with test protocol).

8.2.3 Half-yearly maintenance

▷ Clean the C/L system (see page 70).

The half-yearly cleaning is indicated by the program. Upon activation of the system cleaning, the indication is blanked out again.

9 Malfunctions, Repairs and Warranty

9.1 Malfunctions

In case of malfunctions:

- ▷ Restart the software.
- ▷ Restart the milling machine.
- ▷ Restart the PC.

9.2 Repairs

Repairs may be carried out only by trained, specialised personnel.

9.3 Warranty

The warranty complies with the statutory provisions. For more information, please refer to our General Terms of Business.

10 Environmental Protection

Packaging

In terms of packaging, the manufacturer participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally friendly and recyclable.

Waste equipment

Waste equipment contain materials that should be sorted for environmental-friendly recycling.

- ▷ After the service life has elapsed, dispose of the machine via the public disposal systems in an environmentally-friendly manner.

The plastic parts are labelled. Thus, they can be sorted for recycling or disposal.

Coolant/Lubricant (C/L)

- ▷ Dispose of the C/L in accordance with the applicable regulations and laws.

Proper disposal of the C/L must be documented.

11 Technical Data and Accessories/Spare Parts

 Subject to changes.

Technical data of the milling machine

	Unit	Value
Art. No.	–	179250-STM
Dimensions (D × W × H)	mm	650 × 570 × 740
Weight:		75 - 78
Electrical connections	V/Hz	100-230/50-60
Power output	W	250
Fuse (time-delayed)	A	3.15 / 6.3
Axes		4 -5
Accuracy	μm	< 10
Allowable temperature range (operation)	°C	18 - 30
Spindle torque	Ncm	4
Spindle speed	rpm	60,000
Collet diameter	mm	3
Sound pressure level	dB(A)	60
Interfaces	–	USB / Ethernet / slot for SD card / controlling of supply unit / controlling of extraction

Tab. 4

	Unit	Value
Compressed-air connection	–	Dry, clean compressed air
– Pressure	bar	6
– Flow rate	l/min	50
Air purity		
– Solid contaminants	–	Class 3; better than 5 μm for solids
– Water content	–	Class 4; Max. pressure dew point +3 °C
– Total oil content:	–	Class 2; max. oil content 0.1 mg/m ³
Extraction	–	Suitable for zirconium oxide dust particles, HEPA microfilter (97.97 %), filter class H12, dust category M
– Diameter of extraction connection	mm	38
– Suction capacity	l/min	56.6

Tab. 4

Technical data of the supply unit

	Unit	Value
Art. No.	–	178630-STM
Dimensions (D × W × H)	mm	656 × 572 × 651
Weight	kg	36
Electrical connections	V/Hz	100-230/50-60
Power output	W	150
Fuse (time-delayed)	A	1.0
Max. pump flow rate	l/min	3.8
Capacity of C/L container	l	50
Mesh size of filter	µm	50
Interfaces	–	Controlling of supply unit/ Controlling of extraction

Tab. 5

Accessories/Spare parts

Art. No.	Designation
179210	Maintenance unit
179211	Screw set for blank holder
179216	Threaded inserts for blank holder
179217	Test specimen for test cube
179218	Service set, Jäger spindle
179251	Cover, tool holder
179230	Collet, 3 mm
179252	Blank holder 71
179253	Swarf sieve
179254	Extraction hood
179255	Splash guard
179256	Filter bag
179257	Extraction hose with brush
178650	Coolant/lubricant concentrate
178651	Multi-function additive test, 100 pce. test strips
178652	pH-value test strips, 100 pce.
178653	Nitrite test, 100 pce. test strips
1156901	Fuse, milling machine (T 3.15 A / 250 V)
172334	Fuse, supply unit (T 1.0 A / 250 V)
179186	Blank holder 71 Milling machine (4X)
179184	Extraction hood, milling machine (4X)
179188	Splash guard, Milling machine (4X)
179189	Control system, complete Milling machine (4X)

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Amann Girrbach AG
Institut Straumann AG