




- 1 CD2818001 Position Indicator
- 2 CD2818001 Position Indicator
- 3 CD2818001 Position Indicator
- 4 CD2818001 Position Indicator
- 5 CD0420008 Implant Driver IF NP for Ratchet Long
- 6 CD0420010 Implant Driver IF SP for Ratchet Long
- 7 --- ---
- 8 --- ---
- 9 CD0326004 Countersink Drill 3.5
- 10 CD0326005 Countersink Drill 3.75
- 11 CD0326006 Countersink Drill 4.3
- 12 CD0326007 Countersink Drill 5.0
- 13 CD0420007 Implant Driver IF NP for Ratchet Short
- 14 CD0420009 Implant Driver IF SP for Ratchet Short
- 15 --- ---
- 16 --- ---

- 17 CD0312001 Drill Extender
- 18 CD0324007 Tapered Drill 3.5 +
- 19 CD0324011 Tapered Drill 3.75 +
- 20 CD0324015 Tapered Drill 4.3 +
- 21 CD0324019 Tapered Drill 5.0 +
- 22 CD0422004 Implant Driver IF NP for Hand Piece
- 23 CD0422005 Implant Driver IF SP for Hand Piece
- 24 --- ---
- 25 CD0521003 Screw Driver Hex Manual 1.2
- 26 CD0324020 Initial Drill
- 27 CD0324001 Tapered Drill 2.0
- 28 CD0324004 Tapered Drill 3.5
- 29 CD0324008 Tapered Drill 3.75
- 30 CD0324012 Tapered Drill 4.3
- 31 CD0324016 Tapered Drill 5.0
- 32 CD0521004 Adapter - Hand Piece to Manual
- 33 --- ---
- 34 CD2920001 Torque Wrench / Ratchet
- 35 CD2911002 Depth Probe



Cleaning

Clean this product according to the following instructions after each use:

1. Remove the instruments from the Cassette and clean separately (see specific instrument instructions).
2. Disassemble the three parts of the Cassette. Soak the lid, tray and base for at least 1 min in the cleaning solution (Enzol® enzymatic detergent, prepared per manufacturer directions) so that the components are sufficiently covered. Pay attention that there is no contact between the components. Assist cleaning by careful brushing with a soft brush, paying particular attention to the interface between the silicone block inserts and the tray or base. Do not use metal brushes or steel wool.
3. Soak the lid, tray and base for 15 min in the cleaning solution (Enzol® enzymatic detergent, prepared per manufacturer directions) with ultrasonic treatment so that the components are sufficiently covered. Pay attention that there is no contact between the components.

4. Remove the Cassette components from the cleaning solution and rinse them at least 3 times intensively (at least 1 min) under running water (local potable water is sufficient).
5. Repeat the cleaning steps of components that remain visibly contaminated.
6. Check all components after cleaning for signs of corrosion or damage. Do not continue to use corroded or damaged components.
7. Dry the components thoroughly before storage and sterilization; the accumulation of moisture on the components and instruments may cause corrosion that can detrimentally affect instrument performance.
8. Insert the cleaned, disinfected and dried instruments in the corresponding locations of the Cassette tray and base. Assemble the Cassette by inserting the tray into the base and attaching the lid.
9. Protect the assembled Cassette from moisture and dust if it will not be used for an extended time.

Presentation and Sterilization

For traceability, each Cassette presents on its tray the laser engraving of the UDI code (Unique Device Identification).

This product is reusable and supplied non-sterile, being unitarily packaged. This product must be correctly cleaned and sterilized before each use. Sterilize the products on the previous day or on the day of the procedure. ATTENTION: This product cannot be autoclaved in its original packaging. Please use for sterilization only the steam sterilization according to the following parameters:


	Fractionated Vacuum / Dynamic Air Removal ¹	Gravity Displacement
Sterilization Time	4 minutes	15 minutes
Sterilization Temperature ²	132°C	132°C
Drying Time	At least 20 minutes ²	At least 20 minutes ²

¹ At least three vacuum steps.

² The effectiveness required in drying time depends directly on the parameters of sole responsibility of the user (density and load configuration, sterilizing conditions, which must be determined by the user). Nevertheless, a drying time shorter than 20 minutes cannot be applied.

Note:

1. Please store the Cassette after sterilization in the sterilization packagings at a dry and dust-free place.
2. The flash/immediate use sterilization procedure cannot be used.
3. Do not use dry heat sterilization, radiation sterilization, formaldehyde and ethylene oxide sterilization, as well as plasma sterilization



 **straumann®**

© Nuvo® 2020. All rights reserved. Nuvo® and/or other trademarks and logos from Nuvo® that are mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

© All other trademarks are the property of their respective owners. All rights reserved.

ifu.nuvoimplants.com/en www.nuvoimplants.com

Straumann North American Headquarters
Straumann USA, LLC
60 Minuteman Road | Andover, MA 01810
Phone 833-BUYNUVO (833-289-6886)
Fax 978/747 2490
www.straumann.us

USLIT.2074 7/20 V1 PMR