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
**Schincaglia GP, Hebert E, Farina R, Simonelli A, Trombelli L. (2015).**

Single versus double flap approach in periodontal regenerative treatment.  
*J Clin Periodontol. Jun;42(6):557-66.*

**Farina R, Simonelli A, Minenna L, Rasperini G, Trombelli L. (2014).**

Single-flap approach in combination with enamel matrix derivative in the treatment of periodontal intraosseous defects.  
*Int J Periodontics Restorative Dent. Jul-Aug; 34(4):497-506.*

**Farina R, Simonelli A, Rizzi A, Pramstraller M, Cucchi A, Trombelli L. (2013).**

Early postoperative healing following buccal single flap approach to access intraosseous periodontal defects.  
*Clin Oral Investig. Jul; 17(6):1573-83.*

**Trombelli L, Simonelli A, Schincaglia GP, Cucchi A, Farina R. (2012).**

Single-flap approach for surgical debridement of deep intraosseous defects: a randomized controlled trial.  
*J Periodontol. Jan; 83(1):27-35.*

**Trombelli L, Farina R (2011).**

Flap designs for periodontal healing.  
*Endodontic Topics; 25:4-15.*

**Trombelli L. (2010).**

Flap design and suturing technique to optimize reconstructive outcomes.  
*In Sculean A., Periodontal regenerative therapy, Berlin: Quintessenz.*

**Trombelli L, Simonelli A, Pramstraller M, Wikesjö UM, Farina R. (2010).**

Single flap approach with and without guided tissue regeneration and a hydroxyapatite biomaterial in the management of intraosseous periodontal defects.  
*J Periodontol. Sep;81(9):1256-63.*

**Trombelli L, Farina R, Franceschetti G, Calura G. (2009).**

Single-flap approach with buccal access in periodontal reconstructive procedures.  
*J Periodontol. Feb;80(2):353-60.*

**Trombelli L & Farina,R.(2008).**

Clinical outcomes with bioactive agents alone or in combination with grafting or guided tissue regeneration.  
*J Clin Periodontol 35 Suppl 8, 117–135.*

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**Guida L, Annunziata M, Belardo S, Farina R, Scabbia A, Trombelli L. (2007).**

Effect of autogenous cortical bone particulate in conjunction with enamel matrix derivative in the treatment of periodontal intraosseous defects.  
*J Periodontol. Feb;78(2):231-8.*

**Trombelli L, Annunziata M, Belardo S, Farina R, Scabbia A, Guida L. (2006).**

Autogenous bone graft in conjunction with enamel matrix derivative in the treatment of deep periodontal intra-osseous defects: a report of 13 consecutively treated patients.

*J Clin Periodontol. Jan;33(1):69–75.*

**Trombelli L. (2005).**

Which reconstructive procedures are effective for treating the periodontal intraosseous defect?

*Periodontol 2000. 37:88-105.*

**Trombelli L, Heitz-Mayfield LJ, Needleman I, Moles D, Scabbia A. (2002).**

A systematic review of graft materials and biological agents for periodontal intraosseous defects.

*J Clin Periodontol.;29 Suppl 3:117-35; discussion 160-2.*

**Trombelli L, Bottega S, Zucchelli G. (2002).**

Supracrestal soft tissue preservation with enamel matrix proteins in treatment of deep intrabony defects.

*J Clin Periodontol. May;29(5):433-9.*

#### LITERATUR

1 Tonetti et al. Enamel matrix proteins in the regenerative therapy of deep intrabony defects – A multicentre randomized controlled trial *J Clin Periodontology* 2002;29;317-325. 2 Miron RJ, Dard M, Weinreb M. Enamel matrix derivative, inflammation and soft tissue wound healing. *J Periodontol Res.* 2014 Nov 23. 3 Ozcelik O, Haytac MC, Seydaoglu G. Immediate post-operative effects of different periodontal treatment modalities on oral health-related quality of life: a randomized clinical trial. *J Clin Periodontol.* 2007 Sep; 34(9):788-96. 4 Trombelli L, Farina F, Franceschetti G & Minenna L. (2007) Parametri fondamentali della diagnosi parodontale. *Dental Cosmos* 15–25.

**International Headquarters**

Institut Straumann AG  
Peter Merian-Weg 12  
CH-4002 Basel, Switzerland  
Phone +41 (0)61 965 11 11  
[www.straumann.com](http://www.straumann.com)

**National Distributor**

Institut Straumann AG  
Straumann Schweiz  
Peter Merian-Weg 12  
Postfach  
CH-4002 Basel  
[www.straumann.ch](http://www.straumann.ch)

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Tel. Bestellungen: 0800 810 812  
Tel. Hotline: 0800 810 814  
Tel. CARES®: 0800 810 816  
E-Mail: [ch.sales@straumann.com](mailto:ch.sales@straumann.com)