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
# Early loading (2 or 6 weeks) of sandblasted and acid-etched (SLA®) ITI implants in the posterior mandible. A 1-year randomized controlled clinical trial

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# A 3-arm study of early loading of rough-surface implants in the completely edentulous maxilla and in the edentulous posterior maxilla and mandible: results after 1 year of loading

Nordin T, Nilsson R, Frykholm A, Hallman M.  
Int J Oral Maxillofac Implants 2004;19(6):880–886.

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
# Immediate loading of implants with 3-unit fixed partial dentures: a 12-month clinical study

Cornelini R, Cangini F, Covani U, Barone A, Buser D.  
Int J Oral Maxillofac Implants 2006;21(6):914–918.

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# A five-year life table analysis on wide neck ITI implants with prosthetic evaluation and radiographic analysis: results from a private practice

Bischof M, Nedir R, Abi Najm S, Szmukler-Moncler S, Samson J.  
Clin Oral Implants Res 2006;17(5):512–520.

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# Clinical field trial examining an implant with a sand-blasted, acid-etched surface

Cochran D, Oates T, Morton D, Jones A, Buser D, Peters F.  
J Periodontol 2007;78(6):974–982.

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# Clinical performance of wide-body implants with a sandblasted and acid-etched (SLA®) surface: results of a 3-year follow-up study in a referral clinic

Bornstein MM, Harnisch H, Lussi A, Buser D.  
Int J Oral Maxillofac Implants 2007;22(4):631–638.

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

# Five-year results from a randomized, controlled trial on early and delayed loading of implant supporting full-arch prosthesis in the edentulous maxilla

Fischer K, Stenberg T, Hedin M, Sennerby L.  
Clin Oral Implants Res 2008;19(5):433–441.

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
# Performance of dental implants after staged sinus floor elevation procedures: 5-year results of a prospective study in partially edentulous patients

Bornstein MM, Chappuis V, von Arx T, Buser D.  
Clin Oral Implants Res 2008;19(10):1034–1043.

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# Clinical and radiographic study of implant treatment outcome in periodontally susceptible and non-susceptible patients: a prospective long-term study

De Boever AL, Quirynen M, Coucke W, Theuniers G, De Boever JA.  
*Clin Oral Implants Res* 2009;20(12):1341–1350.

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# Marginal bone level changes and prosthetic maintenance of mandibular overdentures supported by 2 implants: a 5-year randomized clinical trial

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
# 10-year outcome of SLA® implants in the edentulous maxilla

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
Arlin ML. Short dental implants as a treatment option: results from an observational study in a single private practice. *Int J Oral Maxillofac Implants* 2006;21(5):769-76.

Arlin ML. Survival and success of sandblasted, large-grit, acid-etched and titanium plasma-sprayed implants: a retrospective study. *J Can Dent Assoc* 2007;73(9):821.

Arvidson K, Esselin O, Felle-Persson E, Jonsson G, Smedberg JJ, Soderstrom U. Early loading of mandibular full-arch bridges screw retained after 1 week to four to five Monotype implants: 3-year results from a prospective multicentre study. *Clin Oral Implants Res* 2008;19(7):693-703.

Aydin M, Yilmaz A, Kâtibo lu B, Tunç EP. ITI implants and Dolder bars in the treatment of large traumatic defect of mandible: a clinical report. *Dent Traumatol* 2004;20(6):348-52.

Barewal RM, Oates TW, Meredith N, Cochran DL. Resonance frequency measurement of implant stability in vivo on implants with a sandblasted and acid-etched surface. *Int J Oral Maxillofac Implants* 2003;18(5):641-51.

Bergkvist G, Sahlholm S, Karlsson U, Nilner K, Lindh C. Immediately loaded implants supporting fixed prostheses in the edentulous maxilla: a preliminary clinical and radiologic report. *Int J Oral Maxillofac Implants* 2005;20(3):399-405.

Bergkvist G, Nilner K, Sahlholm S, Karlsson U, Lindh C. Immediate loading of implants in the edentulous maxilla: use of an interim fixed prosthesis followed by a permanent fixed prosthesis: a 32-month prospective radiological and clinical study. *Clin Implant Dent Relat Res* 2009;11(1):1-10.

Bernard JP, Schatz JP, Christou P, Belser U, Kiliaridis S. Long-term vertical changes of the anterior maxillary teeth adjacent to single implants in young and mature adults. A retrospective study. *J Clin Periodontol* 2004;31(11):1024-8.

Bischof M, Nedir R, Szmukler-Moncler S, Bernard JP, Samson J. Implant stability measurement of delayed and immediately loaded implants during healing. *Clin Oral Implants Res* 2004;15(5):529-39.

Bischof M, Nedir R, Abi Najm S, Szmukler-Moncler S, Samson J. A five-year life-table analysis on wide neck ITI implants with prosthetic evaluation and radiographic analysis: results from a private practice. *Clin Oral Implants Res* 2006;17(5):512-20.

Bornstein MM, Lussi A, Schmid B, Belser UC, Buser D. Early loading of nonsubmerged titanium implants with a sandblasted and acid-etched (SLA) surface: 3-year results of a prospective study in partially edentulous patients. *Int J Oral Maxillofac Implants* 2003;18(5):659-666.

Bornstein MM, Schmid B, Belser UC, Lussi A, Buser D. Early loading of non-submerged titanium implants with a sandblasted and acid-etched surface. 5-year results of a prospective study in partially edentulous patients. *Clin Oral Implants Res* 2005;16(6):631-638.

Bornstein MM, Harnisch H, Lussi A, Buser D. Clinical performance of wide-body implants with a sandblasted and acid-etched (SLA) surface: results of a 3-year follow-up study in a referral clinic. *Int J Oral Maxillofac Implants* 2007;22(4):631-8.

Bornstein MM, Chappuis V, von Arx T, Buser D. Performance of dental implants after staged sinus floor elevation procedures: 5-year results of a prospective study in partially edentulous patients. *Clin Oral Implants Res* 2008;19(10):1034-43.

Carinci F, Guidi R, Franco M, Viscioni A, Rigo L, De Santis B, Tropina E. Implants inserted in fresh-frozen bone: a retrospective analysis of 88 implants loaded 4 months after insertion. *Quintessence Int* 2009;40(5):413-9.

Carlino P, Pepe V, Pollice G, Grassi FR. Immediate transmucosal implant placement in fresh maxillary and mandibular molar extraction sockets: description of technique and preliminary results. *Minerva Stomatol* 2008;57(10):471-83.

Caruso G, Cattaneo A. Removable prosthesis supported by implants according to the Cagliari modified conometry technique: case report. *Int J Periodontics Restorative Dent* 2007;27(3):259-65.

Cehreli MC, Uysal S, Akca K. Marginal bone level changes and prosthetic maintenance of mandibular overdentures supported by 2 implants: a 5-year randomized clinical trial. *Clin Implant Dent Relat Res* 2009;12(1):114-121.

Chaddad K, Ferreira AF, Geurs N, Reddy MS. Influence of surface characteristics on survival rates of mini-implants. *Angle Orthod* 2008;78(1):107-13.

Chiapasco M, Lang NP, Bosshardt DD. Quality and quantity of bone following alveolar distraction osteogenesis in the human mandible. *Clin Oral Implants Res* 2006;17(4):394-402.

Cochran DL, Buser D, ten Bruggenkate CM, Weingart D, Taylor TM, Bernard JP, Peters F, Simpson JP. The use of reduced healing times on ITI implants with a sandblasted and acid-etched (SLA) surface: early results from clinical trials on ITI SLA implants. *Clin Oral Implants Res* 2002;13(2):144-53.

Cochran D, Oates T, Morton D, Jones A, Buser D, Peters F. Clinical field trial examining an implant with a sand-blasted, acid-etched surface. *J Periodontol* 2007;78(6):974-82.

Cordaro L, Torsello F, Mirisola Di Torresanto V, Rossini C. Retrospective evaluation of mandibular incisor replacement with narrow neck implants. *Clin Oral Implants Res* 2006;17(6):730-5.

Cornellini R, Cangini F, Covani U, Barone A, Buser D. Immediate restoration of single-tooth implants in mandibular molar sites: a 12-month preliminary report. *Int J Oral Maxillofac Implants* 2004;19(6):855-60.

- Cornellini R, Cangini F, Covani U, Barone A, Buser D. Immediate loading of implants with 3-unit fixed partial dentures: a 12-month clinical study. *Int J Oral Maxillofac Implants* 2006;21(6):914-8.
- Crismani AG, Bernhart T, Schwarz K, Celar AG, Bantleon HP, Watzek G. Ninety percent success in palatal implants loaded 1 week after placement: a clinical evaluation by resonance frequency analysis. *Clin Oral Implants Res* 2006;17(4):445-50.
- De Boever AL, De Boever JA. Guided bone regeneration around non-submerged implants in narrow alveolar ridges: a prospective long-term clinical study. *Clin Oral Implants Res* 2005;16(5):549-56.
- De Boever AL, Quirynen M, Coucke W, Theuniers G, De Boever JA. Clinical and radiographic study of implant treatment outcome in periodontally susceptible and non-susceptible patients: a prospective long-term study. *Clin Oral Implants Res* 2009;20(12):1341-50.
- Eliasson A, Blomqvist F, Wennerberg A, Johansson A. A retrospective analysis of early and delayed loading of full-arch mandibular prostheses using three different implant systems: clinical results with up to 5 years of loading. *Clin Implant Dent Relat Res* 2009;11(2):134-48.
- Feldmann I, Bondemark L. Anchorage capacity of osseointegrated and conventional anchorage systems: a randomized controlled trial. *Am J Orthod Dentofacial Orthop* 2008;133(3):339.e19-28.
- Feloutzis A, Lang NP, Tonetti MS, Bürgin W, Brägger U, Buser D, Duff GW, Kornman KS. IL-1 gene polymorphism and smoking as risk factors for peri-implant bone loss in a well-maintained population. *Clin Oral Implants Res* 2003;14(1):10-7.
- Ferrigno N, Laureti M, Fanali S. Inferior alveolar nerve transposition in conjunction with implant placement. *Int J Oral Maxillofac Implants* 2005;20(4):610-20.
- Ferrigno N, Laureti M, Fanali S. Dental implants placement in conjunction with osteotome sinus floor elevation: a 12-year life-table analysis from a prospective study on 588 ITI implants. *Clin Oral Implants Res* 2006;17(2):194-205. Erratum in: *Clin Oral Implants Res* 2006;17(4):479.
- Fischer K, Stenberg T. Early loading of ITI implants supporting a maxillary full-arch prosthesis: 1-year data of a prospective, randomized study. *Int J Oral Maxillofac Implants* 2004;19(3):374-81.
- Fischer K, Stenberg T. Three-year data from a randomized, controlled study of early loading of single-stage dental implants supporting maxillary full-arch prostheses. *Int J Oral Maxillofac Implants* 2006;21(2):245-52.
- Fischer K, Stenberg T, Hedin M, Sennerby L. Five-year results from a randomized, controlled trial on early and delayed loading of implants supporting full-arch prosthesis in the edentulous maxilla. *Clin Oral Implants Res* 2008;19(5):433-41.
- Frenken JW, Bouwman WF, Bravenboer N, Zijdeveld SA, Schulten EA, ten Bruggenkate CM. The use of Straumann Bone Ceramic in a maxillary sinus floor elevation procedure: a clinical, radiological, histological and histomorphometric evaluation with a 6-month healing period. *Clin Oral Implants Res* 2010;21(2):201-8.
- Friedmann A, Kaner D, Leonhardt J, Bernimoulin JP. Immediate substitution of central incisors with an unusual enamel parapsia by a newly developed titanium implant: a case report. *Int J Periodontics Restorative Dent* 2005;25(4):393-9.
- Fugazzotto PA, Vlassis J, Butler B. ITI implant use in private practice: clinical results with 5,526 implants followed up to 72+ months in function. *Int J Oral Maxillofac Implants* 2004;19(3):408-12.
- Fugazzotto PA. Guided bone regeneration at immediate implant insertion and loading: a case report. *Implant Dent* 2004;13(3):223-7.
- Fugazzotto PA. Implant placement at the time of maxillary molar extraction: technique and report of preliminary results of 83 sites. *J Periodontol* 2006;77(2):302-9.
- Gallucci GO, Mavropoulos A, Bernard JP, Belser UC. Influence of immediate implant loading on peri-implant soft tissue morphology in the edentulous maxilla. *Int J Oral Maxillofac Implants* 2007;22(4):595-602.
- Ganeles J. Early loading with the ITI SLA surface as a predictable, routine procedure. *J Indiana Dent Assoc* 2002;81(3):15-6.
- Ganeles J, Rosenberg MM, Holt RL, Reichman LH. Immediate loading of implants with fixed restorations in the completely edentulous mandible: report of 27 patients from a private practice. *Int J Oral Maxillofac Implants* 2001;16(3):418-26.
- Giancotti A, Greco M, Docimo R, Arcuri C. Extraction treatment using a palatal implant for anchorage. *Aust Orthod J* 2003;19(2):87-90.
- Giancotti A, Greco M, Mampieri G, Arcuri C. Clinical management in extraction cases using palatal implant for anchorage. *J Orthod* 2004;31(4):288-94.
- Gokcen-Rohlig B, Yaltirik M, Ozer S, Tuncer ED, Evlioglu G. Survival and success of ITI implants and prostheses: retrospective study of cases with 5-year follow-up. *Eur J Dent* 2009;3(1):42-9.
- Hänggi MP, Hänggi DC, Schoolfield JD, Meyer J, Cochran DL, Hermann JS. Crestal bone changes around titanium implants. Part I: A retrospective radiographic evaluation in humans comparing two non-submerged implant designs with different machined collar lengths. *J Periodontol* 2005;76(5):791-802.
- Hayakawa T, Kiba H, Yasuda S, Yamamoto H, Nemoto K. A histologic and histomorphometric evaluation of two types of retrieved human titanium implants. *Int J Periodontics Restorative Dent* 2002;22(2):164-71.
- Hoffmann O, Beaumont C, Tatakis DN, Zafiropoulos GG. Telescopic crowns as attachments for implant supported restorations: a case series. *J Oral Implantol* 2006;32(6):291-9.
- Hoshino K, Miura H, Morikawa O, Kato H, Okada D, Shinki T. Influence of occlusal height for an implant prosthesis on the periodontal tissues of the antagonist. *J Med Dent Sci* 2004;51(4):187-96.
- Huwiler MA, Pjetursson BE, Bosshardt DD, Salvi GE, Lang NP. Resonance frequency analysis in relation to jawbone characteristics and during early healing of implant installation. *Clin Oral Implants Res* 2007;18(3):275-80.
- Isono H, Kaida K, Hamada Y, Kokubo Y, Ishihara M, Hirashita A, Kuwahara Y. The reconstruction of bilateral clefts using endosseous implants after bone grafting. *Am J Orthod Dentofacial Orthop* 2002;121(4):403-10.
- Jackson A, Lemke R, Hatch J, Salome N, Gakunga P, Cochran D. A comparison of stability between delayed versus immediately loaded orthodontic palatal implants. *J Esthet Restor Dent* 2008;20(3):174-84.

- Jaquière C, Rohner D, Kunz C, Bucher P, Peters F, Schenk RK, Hammer B. Reconstruction of maxillary and mandibular defects using prefabricated microvascular fibular grafts and osseointegrated dental implants – a prospective study. *Clin Oral Implants Res* 2004;15(5):598-606.
- Jebreen SE, Khraisat A. Multicenter retrospective study of ITI implant-supported posterior partial prosthesis in Jordan. *Clin Implant Dent Relat Res* 2007;9(2):89-93.
- Jung BA, Wehrbein H, Hopfenmüller W, Harzer W, Gedrange T, Diedrich P, Kunkel M. Early loading of palatal implants (ortho-type II) a prospective multicenter randomized controlled clinical trial. *Trials* 2007;8:24.
- Jung BA, Kunkel M, Gollner P, Liechti T, Wehrbein H. Success rate of second-generation palatal implants. *Angle Orthod* 2009;79(1):85-90.
- Jung UW, Choi JY, Kim CS, Cho KS, Chai JK, Kim CK, Choi SH. Evaluation of mandibular posterior single implants with two different surfaces: a 5-year comparative study. *J Periodontol* 2008;79(10):1857-63.
- Kato S, Kato M. 'Moderately roughened- and roughened-surface implants used as rigid orthodontic anchorage: a case series. *Clin Implant Dent Relat Res* 2006;8(2):87-94.
- Kermalli JY, Deporter DA, Lai JY, Lam E, Atenafu E. Performance of threaded versus sintered porous-surfaced dental implants using open window or indirect osteotomy-mediated sinus elevation: a retrospective report. *J Periodontol* 2008;79(4):728-36.
- Kessler-Liechti G, Zix J, Mericske-Stern R. Stability measurements of 1-stage implants in the edentulous mandible by means of resonance frequency analysis. *Int J Oral Maxillofac Implants* 2008;23(2):353-8.
- Kinsel RP, Liss M. Retrospective analysis of 56 edentulous dental arches restored with 344 single-stage implants using an immediate loading fixed provisional protocol: statistical predictors of implant failure. *Int J Oral Maxillofac Implants* 2007;22(5):823-30.
- Lai HC, Zhang ZY, Wang F, Zhuang LF, Liu X, Pu YP. Evaluation of soft-tissue alteration around implant-supported single-tooth restoration in the anterior maxilla: the pink esthetic score. *Clin Oral Implants Res* 2008;19(6):560-4.
- Lai HC, Zhang ZY, Zhuang LF, Wang F, Liu X, Pu YP. Early loading of ITI implants supporting maxillary fixed full-arch prostheses. *Clin Oral Implants Res* 2008;19(11):1129-34.
- Landes CA, Kovács AF. Comparison of early telescope loading of non-submerged ITI implants in irradiated and non-irradiated oral cancer patients. *Clin Oral Implants Res* 2006;17(4):367-74.
- Lang NP, Tonetti MS, Suvan JE, Pierre Bernard J, Botticelli D, Fourmousis I, Hallund M, Jung R, Laurell L, Salvi GE, Shafer D, Weber HP; European Research Group on Periodontology. Immediate implant placement with transmucosal healing in areas of aesthetic priority. A multicentre randomized-controlled clinical trial I. Surgical outcomes. *Clin Oral Implants Res* 2007;18(2):188-96.
- Levine RA, Ganeles J, Jaffin RA, Clem DS 3rd, Beagle JR, Keller GW. Multicenter retrospective analysis of wide-neck dental implants for single molar replacement. *Int J Oral Maxillofac Implants* 2007;22(5):736-42.
- Longoni S, Sartori M, Apruzzese D, Baldoni M. Preliminary clinical and histologic evaluation of a bilateral 3-dimensional reconstruction in an atrophic mandible: a case report. *Int J Oral Maxillofac Implants* 2007;22(3):478-83.
- Luongo G, Di Raimondo R, Filippini P, Gualini F, Paoleschi C. Early loading of sandblasted, acid-etched implants in the posterior maxilla and mandible: a 1-year follow-up report from a multicenter 3-year prospective study. *Int J Oral Maxillofac Implants* 2005;20(1):84-91.
- Mannaï C. Early implant loading in severely resorbed maxilla using xenograft, autograft, and platelet-rich plasma in 97 patients. *J Oral Maxillofac Surg* 2006;64(9):1420-6.
- Männchen R, Schätzle M. Success rate of palatal orthodontic implants: a prospective longitudinal study. *Clin Oral Implants Res* 2008;19(7):665-9.
- Nedir R, Bischof M, Szmukler-Moncler S, Bernard JP, Samson J. Predicting osseointegration by means of implant primary stability. *Clin Oral Implants Res* 2004;15(5):520-528.
- Nedir R, Bischof M, Briaux JM, Beyer S, Szmukler-Moncler S, Bernard JP. A 7-year life table analysis from a prospective study on ITI implants with special emphasis on the use of short implants. Results from a private practice. *Clin Oral Implants Res* 2004;15(2):150-7.
- Nedir R, Bischof M, Vazquez L, Szmukler-Moncler S, Bernard JP. Osteotomy sinus floor elevation without grafting material: a 1-year prospective pilot study with ITI implants. *Clin Oral Implants Res* 2006;17(6):679-86.
- Nedir R, Bischof M, Szmukler-Moncler S, Belser UC, Samson J. Prosthetic complications with dental implants: from an up-to-8-year experience in private practice. *Int J Oral Maxillofac Implants* 2006;21(6):919-28.
- Nedir R, Bischof M, Vazquez L, Nurdin N, Szmukler-Moncler S, Bernard JP. Osteotomy sinus floor elevation technique without grafting material: 3-year results of a prospective pilot study. *Clin Oral Implants Res* 2009;20(7):701-7.
- Nixon KC, Chen ST, Ivanovski S. A retrospective analysis of 1,000 consecutively placed implants in private practice. *Aust Dent J* 2009;54(2):123-9.
- Nordin T, Nilsson R, Frykholm A, Hallman M. A 3-arm study of early loading of rough-surfaced implants in the completely edentulous maxilla and in the edentulous posterior maxilla and mandible: results after 1 year of loading. *Int J Oral Maxillofac Implants* 2004;19(6):880-6.
- Nordin T, Graf J, Frykholm A, Helldén L. Early functional loading of sandblasted and acid-etched (SLA) Straumann implants following immediate placement in maxillary extraction sockets. Clinical and radiographic result. *Clin Oral Implants Res* 2007;18(4):441-51.
- Oates TW, Valderrama P, Bischof M, Nedir R, Jones A, Simpson J, Toutenburg H, Cochran DL. Enhanced implant stability with a chemically modified SLA surface: a randomized pilot study. *Int J Oral Maxillofac Implants* 2007;22(5):755-60.
- Palattella P, Torsello F, Cordaro L. Two-year prospective clinical comparison of immediate replacement vs. immediate restoration of single tooth in the esthetic zone. *Clin Oral Implants Res* 2008;19(11):1148-53.
- Pan YH, Ramp LC, Liu PR. Patient responses to dental implant-retained mandibular overdenture therapy: a 6-year clinical study. *Chang Gung Med J* 2007;30(4):363-9.
- Park JC, Ha SR, Kim SM, Kim MJ, Lee JB, Lee JH. A randomized clinical 1-year trial comparing two types of non-submerged dental implants. *Clin Oral Implants Res* 2010;21(2):228-36.

- Payne AG, Tawse-Smith A, Duncan WD, Kumara R. Conventional and early loading of unsplinted ITI implants supporting mandibular overdentures. *Clin Oral Implants Res* 2002;13(6):603-9.
- Payne AG, Tawse-Smith A, Thompson WM, Kumara R. Early functional loading of unsplinted roughened surface implants with mandibular overdentures 2 weeks after surgery. *Clin Implant Dent Relat Res* 2003;5(3):143-53.
- Peñarrocha M, Palomar M, Sanchis JM, Guarinos J, Balaguer J. Radiologic study of marginal bone loss around 108 dental implants and its relationship to smoking, implant location, and morphology. *Int J Oral Maxillofac Implants* 2004;19(6):861-7.
- Pinholt EM. Brånemark and ITI dental implants in the human bone-grafted maxilla: a comparative evaluation. *Clin Oral Implants Res* 2003;14(5):584-92.
- Quinlan P, Nummikoski P, Schenk R, Cagna D, Mellonig J, Higginbottom F, Lang K, Buser D, Cochran D. Immediate and early loading of SLA ITI single-tooth implants: an in vivo study. *Int J Oral Maxillofac Implants* 2005;20(3):360-70.
- Roccuzzo M, Bunino M, Prioglio F, Bianchi SD. Early loading of sandblasted and acid-etched (SLA) implants: a prospective split-mouth comparative study. *Clin Oral Implants Res* 2001;12(6):572-8.
- Roccuzzo M, Wilson T. A prospective study evaluating a protocol for 6 weeks' loading of SLA implants in the posterior maxilla: one year results. *Clin Oral Implants Res* 2002;13(5):502-7.
- Roccuzzo M, Aglietta M, Bunino M, Bonino L. Early loading of sandblasted and acid-etched implants: a randomized-controlled double-blind split-mouth study. Five-year results. *Clin Oral Implants Res* 2008;19(2):148-52.
- Rodrigo D, Aracil L, Martin C, Sanz M. Diagnosis of implant stability and its impact on implant survival: a prospective case series study. *Clin Oral Implants Res* 2010;21(3):255-261.
- Romeo E, Ghisolfi M, Rozza R, Chiapasco M, Iops D. Short (8-mm) dental implants in the rehabilitation of partial and complete edentulism: a 3- to 14-year longitudinal study. *Int J Prosthodont* 2006;19(6):586-92.
- Salvi GE, Gallini G, Lang NP. Early loading (2 or 6 weeks) of sandblasted and acid-etched (SLA) ITI implants in the posterior mandible. A 1-year randomized controlled clinical trial. *Clin Oral Implants Res* 2004;15(2):142-9.
- Schätzle M, Männchen R, Balbach U, Hämmerle CH, Toutenburg H, Jung RE. Stability change of chemically modified sandblasted/acid-etched titanium palatal implants. A randomized-controlled clinical trial. *Clin Oral Implants Res* 2009;20(5):489-95.
- Siciliano VI, Salvi GE, Matarasso S, Cafiero C, Blasi A, Lang NP. Soft tissues healing at immediate transmucosal implants placed into molar extraction sites with buccal self-contained dehiscences. A 12-month controlled clinical trial. *Clin Oral Implants Res* 2009;20(5):482-8.
- Stricker A, Voss PJ, Gutwald R, Schramm A, Schmelzeisen R. Maxillary sinus floor augmentation with autogenous bone grafts to enable placement of SLA-surfaced implants: preliminary results after 15-40 months. *Clin Oral Implants Res* 2003;14(2):207-12.
- Stricker A, Gutwald R, Schmelzeisen R, Gellrich NG. Immediate loading of 2 interforaminal dental implants supporting an overdenture: clinical and radiographic results after 24 months. *Int J Oral Maxillofac Implants* 2004;19(6):868-72.
- Tortamano Neto P, Camargo LO. Prospective clinical evaluation of dental implants with sand-blasted, large-grit, acid-etched surfaces loaded 6 weeks after surgery. *Quintessence Int* 2004;35(9):717-22.
- Tortamano P, Orii TC, Yamanochi J, Nakamae AE, Guarnieri Tde C. Outcomes of fixed prostheses supported by immediately loaded endosseous implants. *Int J Oral Maxillofac Implants* 2006;21(1):63-70.
- Van de Velde T, Sennerby L, De Bruyn H. The clinical and radiographic outcome of implants placed in the posterior maxilla with a guided flapless approach and immediately restored with a provisional rehabilitation: a randomized clinical trial. *Clin Oral Implants Res* 2010. [Epub ahead of print].
- Vitkov L, Gellrich NC, Hannig M. Sinus floor elevation via hydraulic detachment and elevation of the Schneiderian membrane. *Clin Oral Implants Res* 2005;16(5):615-21.
- Vogel RC. Enhancing implant esthetics with ideal provisionalization. *J Indiana Dent Assoc* 2002;81(3):11-4.
- Wehrbein H, Merz BR, Hämmerle CH, Lang NP. Bone-to-implant contact of orthodontic implants in humans subjected to horizontal loading. *Clin Oral Implants Res* 1998;9(5):348-53.
- West JD, Oates TW. Identification of stability changes for immediately placed dental implants. *Int J Oral Maxillofac Implants* 2007;22(4):623-30.
- Wilson TG Jr, Carnio J, Schenk R, Cochran D. Immediate implants covered with connective tissue membranes: human biopsies. *J Periodontol* 2003;74(3):402-9.
- Zarone F, Sorrentino R, Vaccaro F, Russo S. Prosthetic treatment of maxillary lateral incisor agenesis with osseointegrated implants: a 24-39-month prospective clinical study. *Clin Oral Implants Res* 2006;17(1):94-101.
- Zhou W, Han C, Li D, Li Y, Song Y, Zhao Y. Endodontic treatment of teeth induces retrograde peri-implantitis. *Clin Oral Implants Res* 2009;20(12):1326-32.
- Zhou W, Han C, Yunming L, Li D, Song Y, Zhao Y. Is the osseointegration of immediately and delayed loaded implants the same?—comparison of the implant stability during a 3-month healing period in a prospective study. *Clin Oral Implants Res* 2009;20(12):1360-6.
- Zix J, Kessler-Liechti G, Mericske-Stern R. Stability measurements of 1-stage implants in the maxilla by means of resonance frequency analysis: a pilot study. *Int J Oral Maxillofac Implants* 2005;20(5):747-52.

# LITERATUR

- <sup>1</sup> Cochran DL. A comparison of endosseous dental implant surfaces. *J Periodontol* 1999;70(12):1523-1539.
- <sup>2</sup> Cochran DL, Nummikoski PV, Higginbottom FL, Hermann JS, Makins SR, Buser D. Evaluation of an endosseous titanium implant with a sandblasted and acid-etched surface in the canine mandible: radiographic results. *Clin Oral Implants Res* 1996;7(3):240-252.
- <sup>3</sup> Cochran DL, Schenk RK, Lussi A, Higginbottom FL, Buser D. Bone response to loaded and unloaded titanium implants with a sandblasted and acid-etched surface: a histometric study in the canine mandible. *J Biomed Mater Res* 1998;40(1):1-11.
- <sup>4</sup> Buser D, Nydegger T, Oxland T, Cochran DL, Schenk RK, Hirt HP, Snétivy D, Nolte LP. Interface shear strength of titanium implants with a sandblasted and acid-etched surface: a biomechanical study in the maxilla of miniature pigs. *J Biomed Mater Res* 1999;45(2):75-83.
- <sup>5</sup> Buser D, Nydegger T, Hirt HP, Cochran DL, Nolte LP. Removal torque values of titanium implants in the maxilla of miniature pigs. *Int J Oral Maxillofac Implants* 1998;13(5):611-619.
- <sup>6</sup> Li D, Ferguson SJ, Beutler T, Cochran DL, Sittig C, Hirt HP, Buser D. Biomechanical comparison of the sandblasted and acid-etched and the machined and acid-etched titanium surface for dental implants. *J Biomed Mater Res* 2002;60(2):325-332.
- <sup>7</sup> Buser D, Schenk RK, Steinemann S, Fiorellini JP, Fox CH, Stich H. Influence of surface characteristics on bone integration of titanium implants. A histomorphometric study in miniature pigs. *J Biomed Mater Res* 1991;25(7):889-902.
- <sup>8</sup> Shalabi MM, Gortemaker A, Van't Hof MA, Jansen JA, Creugers NH. Implant surface roughness and bone healing. A systematic review. *J Dent Res* 2006;85(6):496-500.
- <sup>9</sup> Albrektsson T, Wennerberg A. Oral implant surfaces: Part 1 – review focusing on topographic and chemical properties of different surfaces and in vivo responses to them. *Int J Prosthodont* 2004;17(5):536-543.
- <sup>10</sup> Boyan BD, Bonewald LF, Paschalis EP, Lohmann CH, Rosser J, Cochran DL, Dean DD, Schwartz Z, Boskey AL. Osteoblast-mediated mineral deposition in culture is dependent on surface microtopography. *Calcif Tissue Int* 2002;71(16):519-529.
- <sup>11</sup> Lössdörfer S, Schwartz Z, Wang L, Lohmann CH, Turner JD, Wieland M, Cochran DL, Boyan BD. Microrough implant surface topographies increase osteogenesis by reducing osteoclast formation and activity. *J Biomed Mater Res A* 2004;70(3):361-369.
- <sup>12</sup> Cochran DL, Buser D, ten Bruggenkate CM, Weingart D, Taylor TM, Bernard JP, Peters F, Simpson JP. The use of reduced healing times on ITI implants with a sandblasted and acid-etched (SLA) surface: early results from clinical trials on SLA implants. *Clin Oral Implants Res* 2002;13(2):144-153.
- <sup>13</sup> Rocuzzo M, Wilson T. A prospective study evaluating a protocol for 6 weeks' loading of SLA implants in the posterior maxilla: one year results. *Clin Oral Implants Res* 2002;13(5):502-507.
- <sup>14</sup> Salvi GE, Gallini G, Lang NP. Early loading (2 or 6 weeks) of sandblasted and acid-etched (SLA) ITI implants in the posterior mandible. A 1-year randomized controlled clinical trial. *Clin Oral Implants Res* 2004;15(2):142-149.
- <sup>15</sup> Nordin T, Nilsson R, Frykholm A, Hallman M. A 3-arm study of early loading of rough-surfaced implants in the completely edentulous maxilla and in the edentulous posterior maxilla and mandible: results after 1 year of loading. *Int J Oral Maxillofac Implants* 2004;19(6):880-886.
- <sup>16</sup> Luongo G, Di Raimondo R, Filippini P, Gualini F, Paoleschi C. Early loading of sandblasted, acid-etched implants in the posterior maxilla and mandible: a 1-year follow-up report from a multicenter 3-year prospective study. *Int J Oral Maxillofac Implants* 2005;20(1):84-91.
- <sup>17</sup> Fischer K, Stenberg T, Hedin M, Sennerby L. Five-year results from a randomized, controlled trial on early and delayed loading of implants supporting full-arch prosthesis in the edentulous maxilla. *Clin Oral Implants Res* 2008;19(5):433-441.
- <sup>18</sup> Cornellini R, Cangini F, Covani U, Barone A, Buser D. Immediate restoration of single-tooth implants in mandibular molar sites: a 12-month preliminary report. *Int J Oral Maxillofac Implants* 2004;19(6):855-860.
- <sup>19</sup> Cornellini R, Cangini F, Covani U, Barone A, Buser D. Immediate loading of implants with 3-unit fixed partial dentures. A 12-month clinical study. *Int J Oral Maxillofac Implants* 2006;21(6):914-918.
- <sup>20</sup> Cochran D, Oates T, Morton D, Jones A, Buser D, Peters F. Clinical field trial examining an implant with a sand-blasted, acid-etched surface. *J Periodontol* 2007;78(6):974-982.
- <sup>21</sup> Bornstein MM, Chappuis V, von Arx T, Buser D. Performance of dental implants after staged sinus elevation procedures: 5-year results of a prospective study in partially edentulous patients. *Clin Oral Implants Res* 2008;19(10):1034-1043.
- <sup>22</sup> Cehreli MC, Uysal S, Akca K. Marginal bone level changes and prosthetic maintenance of mandibular overdentures supported by 2 implants: a 5-year randomized clinical trial. *Clin Implant Dent Relat Res* 2010;12(2):114-121.
- <sup>23</sup> Jung JH, Choi SH, Cho KS, Kim CS. Bone-added osteotome sinus floor elevation with simultaneous placement of non-submerged sand blasted with large grit and acid etched implants: a 5-year radiographic evaluation. *J Periodontal Implant Sci* 2010;40(2):69-75.
- <sup>24</sup> Fischer K. 10-year outcome of SLA implants in the edentulous maxilla. ITI World Symposium, Geneva, Switzerland, 15-17 Apr 2010.
- <sup>25</sup> Zarone F, Sorrentino R, Vaccaro F, Russo S. Prosthetic treatment of maxillary lateral incisor agenesis with osseointegrated implants: a 24-39-month prospective clinical study. *Clin Oral Implants Res* 2006;17(1):94-101.
- <sup>26</sup> Bischof M, Nedir R, Abi Najm S, Szmukler-Moncler S, Samson J. A five-year life-table analysis on wide neck ITI implants with prosthetic evaluation and radiographic analysis: results from a private practice. *Clin Oral Implants Res* 2006;17(5):512-520.
- <sup>27</sup> Lang NP, Tonetti MS, Suvan JE, Pierre Bernard J, Botticelli D, Fourmouls I, Hallund M, Jung R, Laurell L, Salvi GE, Shafer D, Weber HP, European Research Group on Periodontology. Immediate implant placement with transmucosal healing in areas of aesthetic priority. A multicentre randomized-controlled clinical trial I. Surgical outcomes. *Clin Oral Implants Res* 2007;18(2):188-196.
- <sup>28</sup> Arlin ML. Short dental implants as a treatment option: results from an observational study in a single private practice. *Int J Oral Maxillofac Implants* 2006;21(5):769-776.
- <sup>29</sup> Bornstein MM, Lussi A, Schmid B, Belser UC, Buser D. Early loading on nonsubmerged titanium implants with a sandblasted and acid-etched surface: 3-year results of a prospective study in partially edentulous patients. *Int J Oral Maxillofac Implants* 2003;18(5):659-666.
- <sup>30</sup> Bornstein MM, Harnisch H, Lussi A, Buser D. Clinical performance of wide-body implants with a sandblasted and acid-etched (SLA) surface: results of a 3-year follow-up study in a referral clinic. *Int J Oral Maxillofac Implants* 2007;22(4):631-638.
- <sup>31</sup> Fischer K, Stenberg T. Early loading of ITI implants supporting a maxillary full-arch prosthesis: 1-year data from a prospective, randomized study. *Int J Oral Maxillofac Implants* 2004;19(3):374-381.
- <sup>32</sup> Fischer K, Stenberg T. Three-year data from a randomized, controlled study of early loading of single-stage dental implants supporting full-arch maxillary prostheses. *Int J Oral Maxillofac Implants* 2006;21(2):245-252.

## NOTIZEN

**International Headquarters**

Institut Straumann AG  
Peter Merian-Weg 12  
CH-4002 Basel, Switzerland  
Phone +41 (0)61 965 11 11  
[www.straumann.com](http://www.straumann.com)

**National Distributor**

Institut Straumann AG  
Straumann Schweiz  
Peter Merian-Weg 12  
Postfach  
CH-4002 Basel  
[www.straumann.ch](http://www.straumann.ch)

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Tel. Bestellungen: 0800 810 812  
Tel. Hotline: 0800 810 814  
Tel. CARES®: 0800 810 816  
E-Mail: [ch.sales@straumann.com](mailto:ch.sales@straumann.com)