Periodontal Treatment Guide

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Lang, Lindhe, Clinical Periodontology and Implant Dentistry (5. Ausg.) 2008(2) 655-673. **Lang** et al., Qualitätssicherung in der Parodontologie. Interdisziplinäre Fortbildungswoche IWF zur Qualitätssicherung in der Zahnmedizin der Schweizerischen Zahnärzte-Gesellschaft SSO, 1999 **Saxer UP, Muhlemann HR**, Motivation and education. SSO Schweiz Monatsschr Zahnheilkd 1975, 85, 905-919 **Lang N, Tan WC, Krähenmann MA, Zwahlen M**, A systematic review of the effects of full-mouth debridement with and without antiseptics in patients with chronic periodontitis. 6. Europäischer Workshop über Parodontologie, Feb. 2008, Ittingen, Thurgau, Schweiz **Heitz Mayfield LJA, Trombelli L, Heitz F, Needleman I, Moles D**, A systematic review of the effect of surgical debridement vs. non-surgical debridement for the treatment of chronic periodontitis. J Clin Periodontol 2002, 29 (3) 92-102 **Herrera D, Alonso B, Leon R, Roldan S, Sanz M**, Antimicrobial therapy in periodontitis: the use of systemic antimicrobials against the subgingival biofilm. 6. Europäischer Workshop über Parodontologie, Feb. 2008, Ittingen, Thurgau, Schweiz **Gaunt F, Devine M, Pennington M, Vernazza C, Gwynett E, Steen N, Heasman P**, The cost-effectiveness of supportive periodontal care for patients with chronic periodontitis. 6. Europäischer Workshop über Parodontologie, Feb. 2008, Ittingen, Thurgau, Schweiz **Schwarz F, Aoki A, Becker J, Sculean A**, Laser application in non-surgical periodontal therapy: a systematic review. 6. Europäischer Workshop über Parodontologie, Feb. 2008, Ittingen, Thurgau, Schweiz **Claffey N, Nylund K, Kiger R, Garrett S, Egelberg J**, Diagnostic predictability of scores of plaque, bleeding, suppuration and probing depth for probing attachment loss. 3.5 years of observation following initial periodontal therapy. J Clin Periodontol 1990, 17 (2) 108-114 **Eberhard J, Jervoe-Storm PM, Needleman I, Worthington H, Jepsen S**, Full-mouth treatment concepts for chronic periodontitis: a systematic review, J Clin Periodontol 2008, 35 (7) 591-604 **Lang N** et al., A systematic review of the effects of full-mouth debridement with and without antiseptics in patients with chronic periodontitis, J Clin Periodontol 2008, 35 (8) 8-21 **Quirinen M, Mongardini C, de Soete M**, The role of chlorhexidine in the one-stage full-mouth disinfection treatment of patients with advanced adult periodontitis. J Clin Periodontol 2000, 27 (8) 578-589 **Swierkot K, Flores de Jacoby L, Mengel R** et al., One-stage full-mouth disinfection versus quadrant and full-mouth root planing, J Clin Periodontol 2009, 36, (3) 240-249 **Quirinen M** et al., Benefit of „one-stage full-mouth disinfection“ is explained by disinfection and root planing within 24 hours: a randomized controlled trial. J Periodontol, 33, (9) 639-647 **Wang D** et al., Antibody response after single-visit full-mouth ultrasonic debridement versus quadrant-wise therapy. J Clin Periodontol 2006, 9, 632-638 **Lindhe J** et al., Special issue of the Journal of Clinical Periodontology on enamel matrix Proteins. J Clin Periodontol 1997 Sept. 24(9)

Bosshardt D et al., Effects of enamel matrix proteins on tissue formation along the roots of human teeth. J Periodontol. res 2004, 40, 158 **Heijl, Heden** et al., Enamel matrix derivative (Straumann® Emdogain) in the treatment of intrabony periodontal defects. J Clin Periodontology 1997; 24; 705-714 **Pontoriero** et al., The use of barrier membranes and enamel matrix proteins in the treatment of angular bone defects. J Clin Periodontol. 1999; 26(12): 833-40 **Heden, Wennström** et al., Five-Year Follow-up of Regenerative Periodontal Therapy with Enamel Matrix Derivative at Sites with Angular Bone Defects. J Periodontol 2006; 295-301 **Sculean** et al., Treatment of Intrabony Defects With an Enamel Matrix Protein Derivative or Bioabsorbable Membrane: A 8-Year Follow-up Split-Mouth Study. J Periodontol 2006; 77(11), 1879-1886 **McGuire MK, Nunn M**, Evaluation of Human recession defect treated with coronally advanced flaps and either Enamel Matrix Derivative or Connective Tissue. J Periodontol 2003; 74: 1110-1125 **McGuire MK, Cochran DL**, Evaluation of Human recession defect treated with coronally advanced flaps and either Enamel Matrix Derivative or Connective Tissue. J Periodontol 2003; 74; 1126-1135 **Cueva MA, Boltchi FE, Nunn ME, Rivera-Hidalgo F, Rees T**, A comparative study of coronally advanced flaps with and without the addition of enamel matrix derivative in the treatment of marginal tissue recession. J Periodontol 2004 Juli. 75(7):949-56 **Castellanos A, de la Rosa M, de la Garza M, Caffesse RG**, Enamel matrix derivative and coronal flaps to cover marginal tissue recessions. J Periodontol. 2006 Jan. 77(1):7-14 **Jepsen, Meyle** et al., A randomized clinical trial comparing enamel matrix derivative and membrane treatment of buccal Class II furcation involvement in mandibular molars. Part I: Study design and results for primary outcomes. Part II: secondary outcomes. J Periodontol. 2004, Aug. 75(8): 1150-60 **Meyle, Jepsen** et al., A randomized clinical trial comparing enamel matrix derivative and membrane treatment of buccal Class II furcation involvement in mandibular molars. Part i: Study design and results for primary outcomes. Part ii: secondary outcomes. J Periodontol. 2004, Sept. 75(9): 1188-95 **Lindhe J** et al., Special issue of the Journal of Clinical Periodontology on Enamel Matrix Proteins. J Clin Periodontol 1997 Sept. 24(9) **Bosshardt D** et al., Effects of enamel matrix proteins on tissue formation along the roots of human teeth. J Periodontol. res 2004, 40, 158 **Cortellini P, Tonetti MS**, microsurgical approach to periodontal regeneration. initial evaluation in a case cohort. J Periodontol 2001, 72, 559-569 **De Sanctis M, Zucchelli G, Clauser C**, Bacterial colonization of barrier material and periodontal regeneration. J Clin Periodontol 1996, 23, 1039-1046 **Esposito M, Grusovin MG, Coulthard P, Worthington HV**, Enamel matrix derivative (Emdogain) for periodontal tissue regeneration in intrabony defects. Cochrane Datenbank Syst. Rev. 2005, CD003875

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© Institut Straumann AG, 2018. Alle Rechte vorbehalten.
Institut Straumann AG
Straumann GmbH