

Clinical Performance of PURE Ceramic Implants in Single-Tooth Gaps: 5-Year Results¹

Aim

- To evaluate the mid-term clinical performance of a Straumann® PURE Ceramic Implant Monotype in single-tooth gaps in the maxilla and mandible
- To confirm the previously published results² on hard and soft tissue parameters after 5 years of follow-up

Study duration



5 Years

Patients



15

21

Mean age

53 Years

Study design

Indication

Single-tooth gaps

Protocol

Transmucosal implant placement in healed extraction sites

Provisionalization

After 12–14 weeks

Final restoration


After 28 weeks

Study implant



Straumann® PURE Ceramic Implant Monotype

STABLE bone levels



After an initial remodeling phase (0–6 months), crestal bone level remained stable, 0.04 mm change between 6 months and 5 years.

EXCELLENT performance

Implant success and survival after 5 years follow-up vs. meta-analysis⁴

Survival rate PURE



Survival and success rates of single-tooth PURE Ceramic Monotype implants after 5 years

Success rate PURE



Meta-analysis titanium implants



General survival rate of single tooth titanium implants after 5 years⁴

Conclusions

- Straumann® PURE Ceramic Implant Monotype shows excellent survival and success rates of 97.2% and stable bone levels after 5 years of follow-up.
- The PURE Ceramic Implant Monotype offers a reliable and successful treatment alternative to Titanium implants.

References

¹ A Prospective Clinical Study to Evaluate the Performance of Zirconium Dioxide Dental Implants in Single Tooth Gaps in the Maxilla and Mandible: 5-Year Results. Manuscript in preparation ² Bormann KH, Gellrich NC, Kniha H, Schild S, Weingart D, Gahlert A. A Prospective Clinical Study to Evaluate the Performance of Zirconium Dioxide Dental Implants in Single Tooth Gaps in the Maxilla and Mandible: 3-Year Results. BMC Oral Health. 2018 Nov 1;18(1):181 ³ Gahlert M, Kniha H, Weingart D, Schild S, Gellrich NC, Bormann KH. A prospective clinical study to evaluate the performance of zirconium dioxide dental implants in single-tooth gaps. Clin Oral Implants Res. 2016; 27(12): e176-e184. ⁴ Salinas T, Eckert S. Implant-supported single crowns predictably survive to five years with limited complications. J Evid Based Dent Pract. 2010 Mar;10(1):56-7.

Ihr Widerspruchsrecht: Wenn Sie der Verarbeitung Ihrer Daten für Werbezwecke widersprechen oder eine erteilte Einwilligung widerrufen möchten, genügt jederzeit eine Nachricht an unseren Datenschutzbeauftragten per E-Mail an datenschutz.de@straumann.com oder per Post an Straumann GmbH, Datenschutzbeauftragter, Heinrich-von-Stephan-Straße 21, 79100 Freiburg. Dies gilt ebenso, wenn Sie aus Gründen, die sich aus Ihrer besonderen Situation ergeben, der Verarbeitung Ihrer Daten widersprechen wollen.