







Straumann® Variobase® for Crown AS (Angled Solution)

Quick Guide

<p>1</p> 	<p>2</p> 	<p>3</p> 	<p>4</p> 	<p>5</p> 
<p>Step 1 – Attach the Burn-out Coping Base to the Variobase® AS (snap-on retention):</p> <ul style="list-style-type: none"> • Check for the correct alignment between the Variobase® AS and the Burn-out Coping Base. • Check for the proper fit and the absence of any rotational or vertical movement. 	<p>Step 2 – Place the Variobase® AS and the Burn-out Coping Base onto the implant analog:</p> <ul style="list-style-type: none"> • Screw the assembly of the Variobase® AS and the Burn-out Coping Base into the implant analog in the master cast (hand-tight). • Use only the Screw AS and the Screwdriver AS which are both color coded in green. 	<p>Step 3 – Assemble the Burn-out Coping Top onto the Burn-out Coping Base (friction retention):</p> <ul style="list-style-type: none"> • Check that the screw channel is centered with the cut-out of the Variobase® AS. • Check for proper fit between the components. • Rotate the Burn-out Coping Top in the optimal position for the final restoration (within the $\pm 45^\circ$ rotational range). • Wax-up together the Burn-out Coping Base and Top to avoid any rotation. 	<p>Step 4 – Modify the Burn-out Coping Top according to the individual anatomical circumstances:</p> <ul style="list-style-type: none"> • If necessary, shorten the upper part (blue area) of the Burn-out Coping Top according to the individual anatomical circumstances. • Shortening the lower part (red area) of the Burn-out Coping Top may prevent the possibility to remove the screw. 	<p>Step 5 – Fabricate the prosthetic restoration:</p> <ul style="list-style-type: none"> • Contour a wax-up shape according to the individual anatomical situation. • Fabricate the prosthetic restoration following the standard procedures to either press or cast the prosthetic restoration.

6



Step 6 – Finalize the prosthetic restoration before bonding:

- If the framework is veneered, ensure that the veneering material's thermal expansion coefficient matches the coping material's thermal expansion coefficient.

7



Step 7 – Final insertion:

- Position the sterilized restoration into the implant and tighten the screw to 35 Ncm using the AS Screwdriver.
- Close the SCS screw channel with cotton and sealing compound.

Please Note:

- Screws AS & Screwdrivers AS are not compatible with the standard SCS and Createch® screws and screwdrivers.
- The Burn-out Coping of the Variobase® for Crown AS allows for a fixed screw channel angulation of 25°.
- Both the Burn-out Coping Base and Top parts have rotation-indexing elements to limit the rotation of the Burn-out Coping Top to a maximum of 90° around the abutment axis (±45°).
- An incorrect alignment of the Burn-out Coping Top may prevent removal of the screw after the crown is finalized.
- Only the exit hole of the Burn-out Coping Top can be shortened according to the anatomical situation.
- Shortening the lower part of the Burn-out Coping Top may prevent the possibility to remove the screw once the restoration is finalized.
- Make sure that the wax layer on the abutment is sufficiently thick (at least 0.15 mm) to provide space for the Burn-out Coping to expand during heating.
- Respect the minimal wall thickness of the dental material used according to the manufacturer's instructions.
- The framework can be produced with reduced anatomical design or with full-contour crown.
- **The Screw AS should be tightened at 35 Ncm.** Applying a torque >35 Ncm could damage the Screw AS which may prevent the possibility to unscrew it.

Instruments for the Straumann® AS Angled Solution



Straumann® Variobase® for Crown AS (incl. Screw AS)



Straumann® Variobase® for Crown AS Burn-out-Coping



Straumann North American Headquarters

Straumann USA, LLC
60 Minuteman Road
Andover, MA 01810
Phone 800/448 8168 (US) • 800/363 4024 (CA)
Fax 978/747 2490
www.straumann.us • www.straumann.ca

ifu.straumann.com

© Straumann USA, LLC 2018. All rights reserved.
Straumann® and/or other trademarks and logos from Straumann® that are mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates. All rights reserved.