

How to Send a Case to a Lab and a Patient.

After you have scanned and validated a case, you can now send it to the lab by following the steps below:


Step 1. First select the **Send** icon from the workflow at the top of the screen. A 3D overview of the case will be shown in the center of the screen.

Step 2. You can examine the order in more detail by clicking **View** in the lower left corner to see the scans and patient details

Step 3. To send the case to the patient, press on the **Send to Patient** option.

Note: To be able to use this option, remember to save the patient's email address in the patient profile.

Step 3. Before being sent to the designated lab the scans will be postprocessed.



Roger Connery
Send to: Smith Dental Practice
Delivery date: Wednesday, September 6, 2017, 1:00 PM
Case no.: 84593_20170904_1141_50
Crown 3, Emax, B2

View

Order from settings: H40101	Default settings
Order from settings: H40101	Order from settings: H40101
Order ID:	84593
Scan file:	File name: 84593_20170904_1141_50
Scan file:	Scan location: 84593_20170904_1141_50
Scan file:	Scan method: 84593_20170904_1141_50
Scan file:	Product: 84593_20170904_1141_50

View **Send to Patient**

Send

45% Processing scans

Step 4. When the case has been successfully sent a green tick will appear.

Step 5. You can monitor the status of the case by using 3shape Communicate.

