



Centralized  
milling



# zerion<sup>®</sup> ML

Multi-layered high translucent zirconia manufactured from KATANA<sup>™</sup> ML discs by Kuraray Noritake Dental Inc.

Zirconia with natural color gradient for up to full-arch full-contour restorations



## EFFICIENCY

Post-processing time is reduced with full-contour restorations



## STRENGTH

zerion<sup>®</sup> ML is a durable material indicated for up to full-arch restorations




## FLEXIBILITY

A broad range of indications and shades provide high flexibility to dental professionals

# Available shades

- The zerion® ML restorations (available through our CARES® centralized milling solution) are available in five multi-layered high translucent shades: A Light, A Dark, B Light, C Light and D Light



# Available indications











- For final restorations
- For up to full-arch

Recommended indication









Possible indication

Not possible indication (n/a)

## TOOTH-BASED

				
Inlay	Onlay	Veneer	Partial crown	Coping
				
Full-contour crown	Telescopic crown	Bridge/bar framework	Full-contour bridge	Vestibular veneering

## ABUTMENT-BASED (AVAILABLE WITHIN THE CARES® X-STREAM™ SOLUTION)

			
Coping on Ti-Base	Full-contour crown on Ti-Base	Bridge/bar framework on Ti-Base	Full-contour bridge on Ti-Base
			
Coping on CAD/CAM abutment	Full-contour crown on CAD/CAM abutment	Bridge/bar framework on CAD/CAM abutments	Full-contour bridge on CAD/CAM abutments

# Workflow

STEPS	OPTIONS
Nesting	Position the restoration in the desired disc/block shade transition
Finishing	Polishing Staining & glazing Layering
Cementing	Adhesive Self-adhesive Conventional <i>Recommended self-adhesive cement:</i> <i>Panavia™ SA Cement Plus</i> <i>Recommended pre-treatment cement:</i> <i>Panavia™ V5</i>

## Chemical composition

ELEMENTS	WEIGHT
ZrO <sub>2</sub> + HfO <sub>2</sub>	90-95 %
Y <sub>2</sub> O <sub>3</sub>	5-8 %
Other oxides	0-2 %

## Physical properties

PROPERTIES	VALUE
Material	Y-TZP
Density [g/cm <sup>3</sup> ]	≥ 6.0
Biaxial flexural strength [MPa]*	≥ 900
CTE (25-500°C) [10 <sup>-6</sup> K <sup>-1</sup> ]	9.9 ± 0.2

\*According to ISO 6872


## Preparation guidelines

- ▶ The preparation shall not have angles or sharp edges
- ▶ The shoulder preparation shall be with rounded inner edges and/or chamfer
- ▶ The preparation shall have retentive surfaces for conventional cementation

# Wall thickness

- The cross-sections shall be adjusted in relation to the total size of the restoration
- The transition from a connector to a crown or coping shall be rounded

## MINIMAL RESTORATION DESIGN



# Sterilization

- For restorations bonded to abutments before being placed in the patient's mouth, the product must be sterilized before use. Sterilization according to the following parameters or per instructions for use of the abutment if different.

METHOD	CONDITIONS	DRYING TIME
Moist Heat (Autoclave) Pre-Vacuum Displacement	132°C / 134°C (270°F / 273°F) for 3 min	Local practice
Moist Heat (Autoclave) Gravity Displacement	132°C / 134°C (270°F / 273°F) for 5 min	Local practice

## International Headquarters

Institut Straumann AG  
 Peter Merian-Weg 12  
 CH-4002 Basel, Switzerland  
 Phone +41 (0)61 965 11 11  
 Fax +41 (0)61 965 11 01  
[www.straumann.com](http://www.straumann.com)

KATANA™ is a trademark of Noritake Co., Ltd.  
 Panavia™ is a trademark of Kuraray Co., Ltd..

© Institut Straumann AG, 2018. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.