



M/D series

In-house
milling



Sintron®

Cobalt-Chromium sinter metal

High-quality Cobalt-Chromium sinter metal for dry milling
and veneered restorations up to full-arch



HANDLING

Sintron® can be dry milled due to its wax like texture ensuring optimum edge stability



FLEXIBILITY

The Cobalt-Chromium material is veneerable with any standard non-precious porcelain



STRENGTH

The high-quality ensures homogeneous, distortion-free restorations for up to full-arch frameworks

Available discs

SINTRON® DISCS (EACH ARTICLE CONTAINS 1 DISC)

Disc height	10 mm	12 mm	14 mm	16 mm	18 mm	20 mm	25 mm
Round Disc (98 mm)	761132-STM	761130-STM	761128-STM	761126-STM	761124-STM	761122-STM	761120-STM

Available indications











- For final restorations
- For up to full-arch

Recommended indication









Possible indication

Not possible indication (n/a)

TOOTH-BASED

 [n/a] Inlay	 [n/a] Onlay	 [n/a] Veneer	 Partial crown	 Coping
 Full-contour crown	 Telescopic crown	 Bridge/bar framework	 Full-contour bridge	 Vestibular veneering

ABUTMENT-BASED

 Coping on Ti-Base	 Full-contour crown on Ti-Base	 Bridge/bar framework on Ti-Base	 Full-contour bridge on Ti-Base
 Coping on CAD/CAM abutment	 Full-contour crown on CAD/CAM abutment	 Bridge/bar framework on CAD/CAM abutments	 Full-contour bridge on CAD/CAM abutments

Workflow

STEPS	OPTIONS
Nesting	Position the restoration in the desired disc position
Milling	Dry milling with D and M series
Sintering	Recommendation: Straumann® CARES® Argotherm
Finishing	Layering Polishing
Cementing	Adhesive Self-adhesive Conventional

Chemical composition

ELEMENTS	WEIGHT
Cobalt (Co)	66 %
Chromium (Cr)	28 %
Molybdenum (Mo)	5 %
Further Elements (Mn, Si, Fe)	< 1 %
Further Elements (C)	< 0.1 %
Organic binder (for blanks in blank condition)	1–2 %

Physical properties

PROPERTIES	VALUE
Tensile strength (R_m) [MPa]	900
E-Module [GPa]	200
CTE (25 – 500 °C) [$10^{-6}K^{-1}$]	14.5


Preparation guidelines

- ▶ The preparation shall not have angles or sharp edges
- ▶ The shoulder preparation shall be with rounded inner edges and/or chamfer
- ▶ The preparation shall have retentive surfaces for conventional cementation


Wall thickness

- ▶ The cross-sections shall be adjusted in relation to the total size of the restoration
- ▶ The transition from a connector to a crown or coping shall be rounded

MINIMAL RESTORATION DESIGN



Restoration wall thickness



Connector cross-section

International Headquarters

Institut Straumann AG

Peter Merian-Weg 12

CH-4002 Basel, Switzerland

Phone +41 (0)61 965 11 11

Fax +41 (0)61 965 11 01

www.straumann.com

Sintron® is a registered trademark from Amann Girrbach AG, Austria

© Institut Straumann AG, 2021. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.