Straumann[®] n!ce[®] Glaze paste and 3D transpa paste

Straumann[®] n!ce[®] Glaze paste und 3D transpa paste

Straumann[®] n!ce[®] Glaze paste et 3D transpa paste

Straumann[®] n!ce[®] Glaze paste y 3D transpa paste

Glaze paste e 3D transpa paste **Straumann**[®] n!ce[®]

Gebrauchsanweisung | Instruction for use | Mode d'emploi |

Modo de empleo | Modalità d'uso

- *EN* • **Straumann**[®] n!ce[®] Glaze paste and 3D transpa paste

1. Description

The Straumann® n!ce® Glaze paste and 3D paste set contains 2 pastes (Glaze paste and 3D transpa paste) as well as 2 dedicated liquids (Refreshing liquid and Diluting liquid) and an accessory (1 brush). The pastes have been specially developed for the finishing of n!ce® fully anatomical restorations and small corrections. They are compatible with Straumann® n!ce® high strength glass-ceramic.

The Straumann® n!ce® Glaze and 3D paste set (REF 256-001-00) contains all of the individual components.

Description	Contents	REF
n!ce® Glaze paste	3g	256-010-03
n!ce® 3D transpa paste	3g	256-030-03
n!ce® Diluting liquid	10ml	256-012-10
n!ce® Refreshing liquid	10ml	256-014-10
n!ce® Brush	1 piece	256-906-20

The n!ce® pastes are manufactured by Dentaaurum using a special high-tech process. The particle size of the pastes have been adapted to suit the requirements for the staining technique and to offer the possibility of creating an individual morphology in the incisal area and in the occlusal surface area (3D paste). The pastes' consistency allows a homogenous wetting of the restoration surface. The thixotropic effect produces perfect surface integrity and enables small corrections. All n!ce® pastes are translucent, thereby retaining the light dynamic effect of the framework material. They give the fully anatomical restorations a natural fluorescence.

1.1. n!ce® Glaze paste

n!ce® Glaze paste is a ready to use fluorescent glaze to seal the surface of the restorations. Only one glaze firing is required (a homogeneous glaze level is already obtained after the first glaze firing). The paste's fluorescence was adapted to the fluorescence of natural teeth.

1.2. n!ce® 3D transpa paste

- With the 3D transpa paste, it is possible to make small corrections and to add contact points precisely.
- The 3D transpa paste has a firing stability which guarantees that the desired morphology is maintained.
- The paste's fluorescence was adapted to the fluorescence of natural teeth.

1.3. n!ce® liquids

The Diluting liquid and the Refreshing liquid are compatible with all n!ce® pastes. However their behavior during processing varies.

1.3.1. Diluting liquid

The Diluting liquid can be used to reduce the paste's viscosity. The Glaze paste can then be applied in a thinner layer giving a glazed effect.

The Diluting liquid is preferably used for applying very thin layers of Glaze paste.

! Please, stir the Glaze paste in the pot before using the liquids with a glass instrument.

! Please use the Diluting liquid on the mixing plate only and not in the pot.

! Using too much Diluting liquid may impair the stability of the build-up.

! Please make sure you close the pots and the bottle of Diluting liquid after use.

1.3.2. Refreshing liquid

The Refreshing liquid is the preferable solution for the 3D transpa paste since it does not affect the paste's modeling qualities.

If the 3D transpa paste is too dry, add Refreshing liquid.

The Refreshing liquid can be used to adapt the viscosity of the 3D pastes without changing their thixotropic behavior.

Usually used to refresh 3D pastes, Refreshing liquid can be an alternative to the diluting liquid to dilute the Glaze paste keeping stability.

! Do not use the Refreshing liquid directly in the pot: Add Refreshing liquid to the paste on the working plate not in the pot. Mix homogeneously on the plate.

! Please make sure you close the pots and the bottle of Refreshing liquid after use.

2. Use and handling

2.1. Milling / Preparation of the n!ce® restoration

Please follow the Straumann® n!ce® instructions for use.

Do not apply with the pastes until the surface of the restoration has been steam-cleaned.

2.2. Application of the Glaze paste

- Use a clean and dry brush.
- Use Diluting or Refreshing liquid to adapt the viscosity of the Glaze paste if needed.
- Apply a thin layer of the Glaze paste on the whole restoration.

2.3. Application of the 3D transpa paste

- Use a clean and dry brush.
- Before applying the 3D transpa paste, wet the surface of the crown with a very thin layer of Glaze paste.
- Apply the 3D transpa paste to create precise contact points or small corrections.

! The thickness of 3D transpa has to be thinner in any case than the framework

2.5. Firing Chart

Please use the following firing chart:

Table 1	Base Temperature (°C)	Drying time (min)		Heating rate (°C/min)	Vacuum start (°C)	Vacuum end (°C)	Final temperature (°C)	Holding time (min)	Cooling temperature (°C)	Cooling rate (°C/min)
		Preheating	Drying							
Glaze Paste	450	5	5	55	450	790	790	1 (without vacuum)	400	25
3D transpa Paste	450	5/6*	5/6*	55	450	790	790	1 (without vacuum)	400	25

* Increase the drying time depending on the thickness

3. Technical data

3.1. Classification

n!ce[®] pastes are dental ceramics Type I Class 1 (according to DIN EN ISO 6872:2015) for the finishing of restorations made of n!ce[®] glass-ceramic.

3.2. Indications

Staining technique, characterization, glazing and small corrections of n!ce[®] glass-ceramic.

3.3. Contraindications

n!ce[®] pastes are contraindicated for all applications that are not listed under the heading "Indications".

n!ce[®] pastes may not be used if there is a known intolerance to any of the ingredients.

3.4. Physical data

Physico-chemical properties of Glaze paste / 3D transpa paste

	Test standard	Test requirement	Values
Coefficient of thermal expansion	ISO 6872:2015	2 bakes, 4 bakes (25-500°C) [10 ⁻⁶ K ⁻¹]	Glaze paste 2 bakes = 7.0 4 bakes = 7.0 3D transpa paste 2 bakes = 7.2 4 bakes = 7.1
Glass temperature (T _g)	ISO 6872:2015	-	≈ 535°C

4. Symbols and general warnings

Consult instructions for use. Please follow the link to the e-IFU
www.ifu.straumann.com

Batch number

Storable until (Date : Year – Month)

CE mark

Rx only

Sale and use only to/ by a qualified technician

Manufacturer

Warning notices

Only for dental use by a qualified technician. The user must ensure that the intended use corresponds to the manufacturer's specifications.

Measures for protection of equipment and products

Firing ovens should be maintained, cleaned, and calibrated regularly. When grinding the objects or blasting the objects, wear safety glasses and breathing protection and use a suction unit. There is a risk of combustion when firing in the furnace. Therefore, wear gloves and use tongs.

General information

To protect the product against contamination, choose a suitable work area and use clean instruments and tools. In the processing of the medical products, products in question should be stored between 10 °C and 35 °C (50 °F and 95 °F). These non-toxic products can be disposed of without an issue if the statutory provisions are upheld.

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Dentaurum GmbH & Co KG | Turnstr. 31 | 75228 Ispringen