STRAUMANN® IMMEDIAXY™

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REFERÊNCIAS

1 van Velzen FJJ, Ofec R, Schulten EAJM, Ten Bruggenkate CM. 10-year survival rate and the incidence of peri-implant disease of 374 titanium dental implants with a SLA surface: a prospective cohort study in 177 fully and partially edentulous patients. *Clinical oral implants research*. 2015;26(10):1121-1128. 2 Kim S, Jung UW, Cho KS, Lee JS. Retrospective radiographic observational study of 1692 Straumann tissue-level dental implants over 10 years: I. Implant survival and loss pattern. *Clin Implant Dent Relat Res*. Oct 2018;20(5):860-866. 3 Nicolau P, Guerra F, Reis R, Kraft T, Benz K, Jackowski J. 10-year outcomes with immediate and early loaded implants with a chemically modified SLA surface. *Quintessence Int*. Jan 25 2019;50(2):114-124. 4 Grandin HM, Berner S, Dard M. A Review of Titanium Zirconium (TiZr) Alloys for Use in Endosseous Dental Implants. *Materials*. 2012;5(8):1348-1360. 5 Altuna P, Lucas-Taulé E, Gargallo-Albiol J, Figueras-Álvarez O, Hernández-Alfaro F, Nart J. Clinical evidence on titanium-zirconium dental implants: a systematic review and meta-analysis. *Int J Oral Maxillofac Surg*. Jul 2016;45(7):842-850. 6 Herrmann J, Hentschel A, Glauche I, Vollmer A, Schlegel KA, Lutz R. Implant survival and patient satisfaction of reduced diameter implants made from a titanium-zirconium alloy: A retrospective cohort study with 550 implants in 311 patients. *J Craniomaxillofac Surg*. Dec 2016;44(12):1940-1944. 7 Calvo-Guirado JL, López Torres JA, Dard M, Javed F, Pérez-Albacete Martínez C, Maté Sánchez de Val JE. Evaluation of extrashort 4-mm implants in mandibular edentulous patients with reduced bone height in comparison with standard implants: a 12-month results. *Clin Oral Implants Res*. Jul 2016;27(7):867-874. 8 Chen ST, Wilson TG, Jr., Hämmerle CH. Immediate or early placement of implants following tooth extraction: review of biologic basis, clinical procedures, and outcomes. *Int J Oral Maxillofac Implants*. 2004;19 Suppl:12-25. 9 Ganeles J, Zollner A, Jackowski J, ten Bruggenkate C, Beagle J, Guerra F. Immediate and early loading of Straumann implants with a chemically modified surface (SLActive) in the posterior mandible and maxilla: 1-year results from a prospective multicenter study. *Clin Oral Implants Res*. Nov 2008;19(11):1119-1128. 10 De Macedo D, Madeira H, Jarry C. Single molar replacement immediately loaded with final crown through Straumann BLX and chairside digital workflow—One tooth one time technique. *Clinical Oral Implants Research*. 2019;30:242-242. 11 Eckert SE, Huebel G, Sandler N, Elkattah R, McNeil DC. Immediately Loaded Fixed Full-Arch Implant-Retained Prosthesis: Clinical Analysis When Using a Moderate Insertion Torque. *Int J Oral Maxillofac Implants*. May/June 2019;34(3):737-744. 12 Eskin MA, Uzel G, Yilmaz S. A fixed reconstruction of fully edentulous patients with immediate function using an apically tapered implant design: a retrospective clinical study. *Int J Implant Dent*. Nov 23 2020;6(1):77. 13 Coskunses FM, Tak Ö. Clinical performance of narrow-diameter titanium-zirconium implants in immediately loaded fixed full-arch prostheses: a 2-year clinical study. *Int J Implant Dent*. Apr 16 2021;7(1):30. 14 Chmielewski K, Ryncarz W, Yüksel O, et al. Image analysis of immediate full-arch prosthetic rehabilitations guided by a digital workflow: assessment of the discrepancy between planning and execution. *International journal of implant dentistry*. 2019;5(1):26-26. 15 Arcuri C, Barlattani A, Mazzetti V, Carosi P, Ferrigno N. Immediate restoration of a single upper central integrating digital workflow and a novel dental implant: a case report. *Journal of biological regulators and homeostatic agents*. 2019;33(6):171-177. 16 Al-Nawas B, Domagala P, Fragola G, et al. A Prospective Noninterventional Study to Evaluate Survival and Success of Reduced Diameter Implants Made From Titanium-Zirconium Alloy. *J Oral Implantol*. Aug 2015;41(4):e118-125. 17 Chiapasco M, Casentini P, Zaniboni M, Corsi E, Anello T. Titanium-zirconium alloy narrow-diameter implants (Straumann Roxolid®) for the rehabilitation of horizontally deficient edentulous ridges: prospective study on 18 consecutive patients. *Clin Oral Implants Res*. Oct 2012;23(10):1136-1141. 18 Corcuera-Flores JR, Pérez-Fierro M, Blanco-Carrión A, Torres-Lagares D, Castellanos-Cosano L, Machuca-Portillo G. Bone loss around narrow implants versus standard diameter implants: Retrospective 2-years case-control study. *J Clin Exp Dent*. Jan 2020;12(1):e79-e84. 19 Caramés J, Pinto AC, Caramés G, Francisco H, Fialho J, Marques D. Survival Rate of 1008 Short Dental Implants with 21 Months of Average Follow-Up: A Retrospective Study. *J Clin Med*. Dec 5 2020;9(12). 20 Gašperšič R, Dard M, Linder S, Oblak Č. One-Year Results Assessing the Performance of Prosthetic Rehabilitations in the Posterior Maxilla Supported by 4-mm Extrashort Implants Splinted to 10-mm Implants: A Prospective Case Series. *Int J Oral Maxillofac Implants*. Mar-Apr 2021;36(2):371-378. 21 Schimmel M, Janner SFM, Joda T, Wittneben JG, McKenna G, Brägger U. Mandibular implant-supported fixed complete dental prostheses on implants with ultrashort and standard length: A pilot treatment. *J Prosthet Dent*. Jul 28 2020. 22 Slotte C, Grønningssæter A, Halmøy AM, et al. Four-millimeter implants supporting fixed partial dental prostheses in the severely resorbed posterior mandible: two-year results. *Clin Implant Dent Relat Res*. May 2012;14 Suppl 1:e46-58. 23 Brizuela-Velasco A, Perez-Pevida E, Jimenez-Garrido A, et al. Mechanical Characterisation and Biomechanical and Biological Behaviours of Ti-Zr Binary-Alloy Dental Implants. *Biomed Res Int*. 2017;2017:2785863. 24 Data on file. 25 Medvedev AE, Molotnikov A, Lapovok R, et al. Microstructure and mechanical properties of Ti-15Zr alloy used as dental implant material. *J Mech Behav Biomed Mater*. Sep 2016;62:384-398. 26 Altinci P, Can G, Gunes O, Ozturk C, Eren H. Stability and Marginal Bone Level Changes of SLActive Titanium-Zirconium Implants Placed with Flapless Surgery: A Prospective Pilot Study. *Clin Implant Dent Relat Res*. Dec 2016;18(6):1193-1199. 27 Buser D, Broggini N, Wieland M, et al. Enhanced bone apposition to a chemically modified SLA titanium surface. *J Dent Res*. Jul 2004;83(7):529-533. 28 Lang NP, Salvi GE, Huynh-Ba G, Ivanovski S, Donos N, Bosshardt DD. Early osseointegration to hydrophilic and hydrophobic implant surfaces in humans. *Clin Oral Implants Res*. Apr 2011;22(4):349-356. 29 El Chaar E, Zhang L, Zhou Y, et al. Osseointegration of Superhydrophilic Implants Placed in Defect Grafted Bones. *Int J Oral Maxillofac Implants*. March/April 2019;34(2):443-450. 30 Heberer S, Kilic S, Hossamo J, Raguse JD, Nelson K. Rehabilitation of irradiated patients with modified and conventional sandblasted acid-etched implants: preliminary results of a split-mouth study. *Clin Oral Implants Res*. May 2011;22(5):546-551. 31 Nack C, Raguse JD, Stricker A, Nelson K, Nahles S. Rehabilitation of irradiated patients with chemically modified and conventional SLA implants: five-year follow-up. *J Oral Rehabil*. Jan 2015;42(1):57-64. 32 Kopf BS, Ruch S, Berner S, Spencer ND, Maniura-Weber K. The role of nanostructures and hydrophilicity in osseointegration: In-vitro protein-adsorption and blood-interaction studies. *J Biomed Mater Res A*. Aug 2015;103(8):2661-2672. 33 Cabrera-Domínguez JJ, Castellanos-Cosano L, Torres-Lagares D, Pérez-Fierro M, Machuca-Portillo G. Clinical performance of titanium-zirconium implants with a hydrophilic surface in patients with controlled type 2 diabetes mellitus: 2-year results from a prospective case-control clinical study. *Clin Oral Investig*. Jul 2020;24(7):2477-2486.

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© Institut Straumann AG, 2023. Todos os direitos reservados.
Institut Straumann AG
3Shape® e TRIOS® são marcas comerciais registradas de 3Shape A/S, Copenhagen, Dinamarca.
Medit® é uma marca registrada da Medit Corp., Seul, Coreia do Sul.
Jason®, mucoderm® e cerabone® são marcas registradas da botiss biomaterials GmbH, Zossen, Alemanha.