
Pilar Straumann® CARES®  
de dióxido de zircônio

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Broggini N, McManus LM, Hermann JS, Medina RU, Oates TW, Schenk RK et al. (2003). Persistent acute inflammation at the implant-abutment interface. *J Dent Res* 82(3):232-237. — Cosgarea R, Gasparik C, Dudea D, Culic B, Dannewitz B, Sculean A (2015). Peri-implant soft tissue colour around titanium and zirconia abutments: a prospective randomized controlled clinical study. *Clin Oral Implants Res* 26(5):537-544. — de Medeiros RA, Vecchiato-Filho AJ, Pellizzer EP, Mazarro JV, dos Santos DM, Goiato MC (2013). Analysis of the peri-implant soft tissues in contact with zirconia abutments: an evidence-based literature review. *J Contemp Dent Pract* 14(3):567-572. — Degidi M, Artese L, Scarano A, Perrotti V, Gehrke P, Piattelli A (2006). Inflammatory infiltrate, microvessel density, nitric oxide synthase expression, vascular endothelial growth factor expression, and proliferative activity in peri-implant soft tissues around titanium and zirconium oxide healing caps. *J Periodontol* 77(1):73-80. — Dhir S (2013). Biofilm and dental implant: The microbial link. — Gigandet M, Gianni B, Francisco F, Bürgin W, Brägger U (2012). Implants with Original and Non-Original Abutment Connections. *Clinical Implant Dentistry and Related Research*: n/a. — Glauser R, Sailer I, Wohlwend A, Studer S, Schibli M, Scharer P (2004). Experimental zirconia abutments for implant-supported single-tooth restorations in esthetically demanding regions: 4-year results of a prospective clinical study. *Int J Prosthodont* 17(3):285-290. — Joda T, Burki A, Bethge A, Bragger U, Zysset P (2015). Stiffness, strength and failure modes of implant-supported monolithic lithium-disilicate (LS2) crowns: influence of titanium and zirconia abutments. *The International Journal of Oral & Maxillofacial Implants* submitted. — Jung RE, Holderegger C, Sailer I, Khraisat A, Suter A, Hammerle CH (2008). The effect of all-ceramic and porcelain-fused-to-metal restorations on marginal peri-implant soft tissue color: a randomized controlled clinical trial. *Int J Periodontics Restorative Dent* 28(4):357-365. — Kajiwara N, Masaki C, Mukaibo T, Kondo Y, Nakamoto T, Hosokawa R (2015). Soft tissue biological response to zirconia and metal implant abutments compared with natural tooth: microcirculation monitoring as a novel bioindicator. *Implant Dent* 24(1):37-41. — Keilig L, Berg J, Söhnchen P, Kocherovskaya BC (2013). Micro-mobility of the implant/abutment interface for original and third-party abutments—a combined experimental and numerical study (abstract). Poster EAO Ref. no. 346. — Kim SK, Koak JY, Heo SJ, Taylor TD, Ryoo S, Lee SY (2012). Screw loosening with interchangeable abutments in internally connected implants after cyclic loading. *Int J Oral Maxillofac Implants* 27(1):42-47. — Klotz MW, Taylor TD, Goldberg AJ (2011). Wear at the titanium-zirconia implant-abutment interface: a pilot study. *Int J Oral Maxillofac Implants* 26(5):970-975. — Lindquist LW, Carlsson GE, Jemt T (1996). A prospective 15-year follow-up study of mandibular fixed prostheses supported by osseointegrated implants. Clinical results and marginal bone loss. *Clin Oral Implants Res* 7(4):329-336. — Manicone PF, Rossi IP, Raffaelli L (2007). An overview of zirconia ceramics: basic properties and clinical applications. *J Dent* 35(11):819-826. — Nakamura K, Kanno T, Milleding P, Ortengren U (2010). Zirconia as a dental implant abutment material: a systematic review. *Int J Prosthodont* 23(4):299-309. — Nascimento CC, Albuquerque RF (2011). Bacterial Leakage Along the Implant-Abutment Interface. — Passos SP, Torrealba Y, Major P, Linke B, Flores-Mir C, Nychka JA (2014). In Vitro Wear Behavior of Zirconia Opposing Enamel: A Systematic Review. *J Prosthodont*. — Quiryren M, van Steenberghe D (1993). Bacterial colonization of the internal part of two-stage implants. An in vivo study. *Clin Oral Implants Res* 4(3):158-161. — Rimondini L, Cerroni L, Carrassi A, Torricelli P (2002). Bacterial colonization of zirconia ceramic surfaces: an in vitro and in vivo study. *Int J Oral Maxillofac Implants* 17(6):793-798. — Saidin S, Abdul Kadir MR, Sulaiman E, Abu Kasim NH (2012). Effects of different implant-abutment connections on micromotion and stress distribution: prediction of microgap formation. *J Dent* 40(6):467-474. — Sailer I, Philipp A, Zembic A, Pjetursson BE, Hammerle CH, Zwahlen M (2009). A systematic review of the performance of ceramic and metal implant abutments supporting fixed implant reconstructions. *Clin Oral Implants Res* 20 Suppl 4:4-31. — Scarano A, Piattelli M, Caputi S, Favero GA, Piattelli A (2004). Bacterial adhesion on commercially pure titanium and zirconium oxide disks: an in vivo human study. *J Periodontol* 75(2):292-296. — Stimmelmayer M, Edelhoff D, Guth JF, Erdelt K, Happe A, Beuer F (2012). Wear at the titanium-titanium and the titanium-zirconia implant-abutment interface: a comparative in vitro study. *Dent Mater* 28(12):1215-1220. — Vagkopoulou T, Koutayas SO, Koidis P, Strub JR (2009). Zirconia in dentistry: Part 1. Discovering the nature of an upcoming bioceramic. *Eur J Esthet Dent* 4(2):130-151. — Welander M, Abrahamsson I, Berglundh T (2008). The mucosal barrier at implant abutments of different materials. *Clin Oral Implants Res* 19(7):635-641. — Zembic A, Kim S, Zwahlen M, Kelly JR (2014). Systematic review of the survival rate and incidence of biologic, technical, and esthetic complications of single implant abutments supporting fixed prostheses. *Int J Oral Maxillofac Implants* 29 Suppl:99-116.

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© Institut Straumann AG, 2015. Todos os direitos reservados.
Institut Straumann AG