

Stage 3 | Prosthetic procedures

# Step 2

## Fabrication of the final prosthesis

# Prosthetic procedures

Step 2 | Fabrication of the final prosthesis

Overview



## Assessment and treatment planning



Step 1 | Patient's expectations,  
history and examination



Step 2 | Treatment planning



Step 3 | Consultation and consent



Step 4 | Fabrication of the surgical drill template

## Surgical procedures



Step 1 | Implant surgery



Step 2 | Post-operative review and suture removal

7–10 days

6–8 weeks

## Prosthetic procedures



Step 1 | Impression-taking



Step 2 | Fabrication of the final prosthesis



Step 3 | Insertion of the final prosthesis

2 weeks

## Aftercare and maintenance



Step 1 | Review visit



Step 2 | Maintenance visit

3–6 months  
(or as necessary)



*In clinic with patient*



*Office / Lab work*


## Contents

<u>Introduction</u>	4
<u>Learning objectives</u>	4
<u>Materials for the dental laboratory</u>	5
<u>Checklist for the fabrication of the final prosthesis</u>	9





## Introduction



Understand what information to provide to your lab technician.

This step explains the necessary information which should be communicated to your lab to help your dental technician fabricate the final fixed dental prosthesis (FDP) required.

We highly recommend that you work with [lab technicians](#) who are experienced in implant prosthetic restorations and familiar with the Straumann® Dental Implant System (SDIS).

Work with lab technicians who are familiar with the Straumann® Dental Implant System.

## Learning objectives

- Understand which information you have to provide to the dental technician together with the impressions, in order to obtain the desired final FDP.
- Understand that the technician should first create the framework of the 3-unit FDP for try-in before final processing.



## Materials for the dental laboratory

The following items must be prepared before communicating or sending to the dental laboratory:

1. Impression(s)
2. Bite registration (if taken)
3. Impression Posts (if the [NC or RC Impression Posts for closed tray](#) were used)

**⚠ Caution:** Items 1-3 should be disinfected and packed appropriately to protect against any damage during transportation to the laboratory.

Disinfect and securely pack the impression(s), bite registration and Impression Posts before sending them to your lab technician.

4. Lab prescription form with clear instructions on the following:

4.1 Implant platform (e.g., RN, WN, NC or RC)

4.2 Tooth shade for the final prosthesis

Inform your lab technician about the implant platforms and desired tooth shade of the final restoration.

Tips for the color assessment.



It is helpful to indicate which shade guide along with the chosen shade was used during the color assessment. The determination of the tooth shade may also be done by the dental technician if there is a dental laboratory on site, or the patient may also visit the dental technician's office for the color assessment.





# Prosthetic procedures

Step 2 | Fabrication of the final prosthesis

Materials for the dental laboratory




## 4.3 Type of restoration

### Screw-retained



[Straumann® Variobase® for Bridge/Bar](#) (for [SP](#) and [BLT Implants](#))


### Cement-retained



[synOcta® Cementable Abutments](#) (for SP Implants)  
[BL Cementable Abutments](#) (for BLT Implants)

Inform your lab technician if you desire a screw- or cement-retained restoration.

**⚠ Caution:** For cement-retained FDPs, the cement margin must lie no deeper than 2.0 mm below the gingiva.



Cement-retained restorations:

Consider that the cement margin must not be deeper than 2.0 mm below the gingiva.



Please click [here](#) for more information about cement-retained and screw-retained restorations.

Learn more about cement-retained and screw-retained restorations.



# Prosthetic procedures


Step 2 | Fabrication of the final prosthesis

Materials for the dental laboratory




## 4.4 Material for framework and final prosthesis

### Zirconia (all-ceramic)



### Porcelain fused to metal



Inform your lab technician about the desired material for the FDP.

## 4.5 Framework design and try-in

Communicate with your dental technician to provide the uncompleted prosthesis or framework for a try-in before final processing.

The try-in of the framework and uncompleted prosthesis should be carried out at these following time-points:

- After framework construction, before veneering
- At “biscuit-bake”, before final glazing
- Before final insertion of the prosthesis

**Try-in sessions should be planned with your dental technician at the following time-points:**

- After framework construction, before veneering
- At “biscuit-bake”, before final glazing
- Before final insertion of the prosthesis

**⚠ Caution:** A **PASSIVE FIT** of the cement-retained or screw-retained framework on the dental implants is necessary to avoid any physical stress or tension in the framework, or on the osseointegrated implants. This helps to prevent any future biomechanical complications.

A passive fit of the framework is mandatory to avoid any physical stress or tension.



## Prosthetic procedures

Step 2 | Fabrication of the final prosthesis


Materials for the dental laboratory



### 4.6 Pontic design

The design of the pontic for a 3-unit FDP should help to facilitate cleanability and oral hygiene procedures by the patient.

The pontic should have an ovate or a modified ridge lap design to facilitate oral hygiene procedures.



- Avoid a total ridge lap (a) design as this can cause food impaction and subsequent soft tissue inflammation.
- To allow proper cleaning, the pontic should be seated in flush contact with the underlying soft tissues, with an ovate (b) or modified ridge lap (c) design.

Please read the next step about [insertion of the final prosthesis](#) for a more detailed description of the procedure.

### 5. Intra-oral photographs


If you have taken intra-oral photographs, these can be of great help to the dental technician in designing the framework and final prosthesis, especially when the dental technician is not able to see the patient in person.

Intra-oral photographs might be helpful for the dental technician to design the FDP.



**CHECKLIST FOR THE FABRICATION OF THE FINAL PROSTHESIS**

Disinfect and prepare the following materials to be sent to the dental lab:

- ☐ Impression(s)
- ☐ Bite registration (if taken)
- ☐ Impression Posts (if the  [NC or RC Impression Posts for closed tray](#) were used)
- ☐ Lab prescription form with clear instructions on the following:
  - ☐ Implant platform (e.g., RN, WN, NC or RC)
  - ☐ Tooth shade for the final prosthesis
  - ☐ Type of restoration: cement-retained or screw-retained
  - ☐ Material for framework and final prosthesis
  - ☐ Framework design and try-in
  - ☐ Pontic design
- ☐ Intra-oral photographs (if taken)



# Prosthetic procedures

## Step 2 | Fabrication of the final prosthesis

### DISCLAIMER

Straumann® Smart is a blended training and education program focused on the education of general dentists who want to become surgically active in the field of dental implantology. The program is limited to information pertaining to straightforward implant cases and focuses on a reduced portfolio of products that are suitable for the treatment of such cases.

All clinical Straumann® Smart content – such as texts, medical record forms, pictures and videos – was created in collaboration with Prof. Dr. Christoph Hämmerle, Prof. Dr. Ronald Jung, Dr. Francine Brandenburg-Lustenberger and Dr. Alain Fontollet from the University of Zürich, Clinic for Fixed and Removable Prosthodontics and Dental Material Science, Switzerland.

Straumann does not give any guarantee that Straumann® Smart provides sufficient knowledge or instruction for the dental professional to become surgically active in the field of implantology. It is the dental professional's sole responsibility to ensure that he/she has the appropriate knowledge and instruction before placing dental implants.

Straumann® Smart does not replace a careful and thorough analysis of each individual patient by a dental professional. Further, it does not imply any guarantee or warranty with regard to completeness of the information provided to the patient. It does not replace the dental professional's duty to inform the patient about the treatment, the products and the risks involved and to receive the patient's informed consent. The dental professional is solely responsible for determining whether or not a treatment or product is suitable for a particular patient and circumstances. Knowledge of dental implantology and instruction in the handling of the relevant products is always necessary and the sole responsibility of the dental professional. The dental professional must always comply with the individual product's Instructions For Use as well as all laws and regulations.

STRAUMANN DISCLAIMS, TO THE EXTENT POSSIBLE BY LAW, ANY LIABILITY, EXPRESS OR IMPLIED, AND BEARS NO RESPONSIBILITY FOR ANY DIRECT, INDIRECT, PUNITIVE, CONSEQUENTIAL OR OTHER DAMAGES, ARISING OUT OF OR IN CONNECTION WITH ANY INFORMATION PROVIDED TO PATIENTS, ERRORS IN PROFESSIONAL JUDGMENT, IN PRODUCT CHOICES OR PRACTICE IN THE USE OR INSTALLATION OF STRAUMANN PRODUCTS.

All clinical content as well as clinical and radiographic images are provided by courtesy of Prof. Dr. Christoph Hämmerle, Prof. Dr. Ronald Jung, Dr. Francine Brandenburg-Lustenberger and Dr. Alain Fontollet from the University of Zürich, Clinic for Fixed and Removable Prosthodontics and Dental Material Science, Switzerland.

### **International Headquarters**

Institut Straumann AG

Peter Merian-Weg 12

CH-4002 Basel, Switzerland

Phone +41 (0)61 965 11 11

Fax +41 (0)61 965 11 01

[www.straumann.com](http://www.straumann.com)

© Institut Straumann AG, 2016. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.