

Stage 1 | Assessment and treatment planning

Step 2 Treatment planning

Assessment and treatment planning

Step 2 | Treatment planning

Overview

Assessment and treatment planning



Step 1 | Patient's expectations,
history and examination



Step 2 | Treatment planning



Step 3 | Consultation and consent



Step 4 | Fabrication of the surgical drill template

Surgical procedures



Step 1 | Implant surgery



Step 2 | Post-operative review and suture removal

7–10 days

Prosthetic procedures



Step 1 | Impression-taking



Step 2 | Fabrication of the final prosthesis

6–8 weeks



Step 3 | Insertion of the final prosthesis

Aftercare and maintenance



Step 1 | Review visit

2 weeks



Step 2 | Maintenance visit

3–6 months
(or as necessary)



In clinic with patient



Office / Lab work

Contents

Introduction	4
Learning objectives	5
1. Contacting the patient's physician	6
2. Study model analysis	6
2.1 Procedure	7
3. Radiographic planning	11
3.1 Outline important landmarks	11
3.2 Evaluate and assess	11
3.3 Measure the available vertical bone height and select the appropriate implant type, diameter and length	13
4. Risk assessment	16
5. Visit planning and cost proposal	18
5.1 Establish a diagnosis and prognosis of the remaining dentition to discuss with the patient.	18
5.2 Elaborate a treatment plan and calculate the amount of chairside time per visit and total number of visits required by the patient.	18
5.3 Calculate costs and write up a cost proposal	19
Appendix	20



Introduction

A thorough patient assessment is the basis of any dental treatment plan. You can finalize the diagnosis and prognosis of the oral situation, and undertake treatment planning in your office with the information gathered from the initial visit with the patient.

1. Contacting the patient's physician



Systematically carry out treatment planning by following points 1–5.

2. Study model analysis



3. Radiographic planning



4. Risk assessment



5. Visit planning and cost proposal



Click on the graphic to go directly to the chapter.



Assessment and treatment planning







Step 1 | History and examination

Learning
objectives



Learning objectives

Be able to use the information gathered from the patient's history and examination to:

-  Identify and select a patient who is suitable for straightforward implant treatment.
-  Complete the risk assessment considering the patient's medical or dental history.
-  Perform a thorough analysis of the study models and radiographs to finalize the diagnosis and prognosis of the dentition, determine which implant to use, and plan restorative options including any necessary adjunctive therapy before implant treatment.
-  Use the Straumann® X-ray Template and the formula to calculate the distortion of the radiographic image for accurate planning.
-  Be able to discuss with your dental technician the proposed case and any necessary restorative planning options.
-  Prepare the treatment plan, costs and alternative treatment options for discussion with the patient at the next visit.



Assessment and treatment planning

Step 2 | Treatment planning

1. Contacting the patient's physician

1. Contacting the patient's physician



The knowledge of former and current diseases, surgeries, and medications helps to identify patients who are at risk for complications of dental implant treatment. If there are any doubts, consult the physician or medical specialist treating the patient for further clarification¹.

Confirm the medical status with the patient's physician if necessary.

2. Study model analysis



A detailed and accurate space assessment is often difficult to do intraorally. It is recommended to obtain initial impressions to fabricate and mount study models, in order to properly perform diagnostic and planning steps.

The study models may be used for:

- Analysis of the occlusion
- Assessment of available space
- Assessment of other factors that might not be easily noted or recorded at the chair-side
- Creation of a diagnostic wax-up
- Creation of an X-ray template



[Video: Study model analysis](#)

Select the implant diameter, implant type, and position individually. Discuss with your [dental technician](#) on how to design the restoration so that the necessary oral hygiene measures can be carried out.

Use study models to properly analyze and plan the treatment.



Discuss with your dental technician about the design of the final restoration.



Assessment and treatment planning

Step 2 | Treatment planning

2. Study model analysis



A few important considerations to note about planning implants and their abutments:

- The implant abutments should always be axially loaded.
- The future implant should be planned with a restoratively-driven implant position in mind.
- The long axis of the implant should be aligned with the cusps of the opposing tooth.


To establish the ideal topographical situation, axial orientation, and choice of implants, we recommend the following to be fabricated by the [dental technician](#):

- A diagnostic wax-up in the edentulous area
- A resin-based surgical drill template to help determine the correct implant position during implant surgery

Important considerations when planning implants and final restorations

Ask your dental technician to create a diagnostic wax-up and surgical drill template.

2.1 Procedure



Apply the [Straumann® Diagnostic T](#) on the study model to obtain an initial measurement of the interproximal distance for the choice of the implant shoulder diameter and prosthetic reconstruction.

Use the Straumann® Diagnostic T to help choose the prosthetic platform of the future implant.

Analyze the implant position three-dimensionally on the study model.

- 2.1.1 Interproximal distance of bone
- 2.1.2 Bucco-lingual (or bucco-palatal) width of bone
- 2.1.3 Inter-occlusal distance

Assess the 3D space for the implant on the study model.



Assessment and treatment planning




Step 2 | Treatment planning

2. Study model analysis



2.1.1 Interproximal distance of bone

The interproximal bone availability is an important factor for choosing the implant type and diameter.



In a single tooth gap, a **minimal distance of 1.5 mm of bone** from the implant to the adjacent tooth **at bone level** (mesial and distal) is required. This will help to preserve bone and soft tissue to create a restoration with adequate emergence profile and soft tissue attachment to support oral hygiene measures and esthetics.

Reduced horizontal distance between tooth and implant can lead to bone resorption at adjacent teeth and reduced ability to clean around the restoration. If this minimum distance of 1.5 mm from the planned implant to the adjacent tooth at bone level is not possible, consider a smaller implant diameter or alternative treatment options (e.g., conventional prosthetics, orthodontics).

⚠ Caution: We do not recommend the use of a 3.3 mm endosteal diameter implant stand-alone in the posterior region due to the strength of masticatory forces in this area.

➔ For SP Implants, please refer to this chart of minimum bone widths.

➔ For BLT Implants, please refer to this chart of minimum bone widths.

Minimal distance of bone from implant to adjacent tooth at bone level: **1.5 mm.**

Chart for reference

Chart for reference



Assessment and treatment planning

Step 2 | Treatment planning

2. Study model analysis



2.1.2 Bucco-lingual (or bucco-palatal) width of bone

The facial and palatal bone wall must be **at least 1 mm** thick in order to ensure stable hard and soft tissue conditions.

For example, if the planned implant site has a bucco-lingual bone width of 6.5 mm, then you are able to use a 4.1 mm endosteal diameter Standard Plus or Bone Level Tapered Implant.


Minimum thickness of bucco-lingual or bucco-palatal bone wall in posterior areas: 1 mm.

→ For SP Implants, please refer to this chart of minimum bone widths.

Chart for reference

→ For BLT Implants, please refer to this chart of minimum bone widths.

Chart for reference



⚠ Caution: If the available orofacial bone is thinner than 1 mm or a layer of bone is missing on one or more sides around the planned implant site, a bone augmentation procedure is indicated. This technique should be employed only by dentists who have adequate experience in the use of augmentation procedures. If it only becomes apparent during surgery that this minimum distance is not present, it is recommended to close the flap and refer the patient to a specialist or colleague with adequate experience in GBR (Guided Bone Regeneration) procedures.

Refer cases which require GBR to more experienced colleagues.



Assessment and treatment planning


Step 2 | Treatment planning

2. Study model analysis



2.1.3 Inter-occlusal distance

The inter-occlusal distance is the distance from the future implant shoulder to the opposing tooth. The minimum distance is about **7 mm** for cement-retained restorations and about **4 mm** for screw-retained restorations.



Space requirements depend on the design of the restoration, including the selection of abutments and crowns. We recommend that you discuss this with your [dental technician](#) during treatment planning in order to have a restoratively-driven

orofacial implant position and axis. The presence of overeruption of the opposing dentition or the presence of cross-bites can reduce the space available for the final restoration.

⚠ Caution: If the available space is insufficient due to overeruption of the opposing teeth, depending on the extent of available space, minimal enameloplasty, orthodontic intervention, elective endodontics, crown lengthening, and crown in the opposing quadrant may be indicated. Such restorative cases should be handled only by clinicians who are experienced in treating advanced or complex cases.

Minimum inter-occlusal distance required:

≥ 7 mm for cement-retained restorations
≥ 4 mm for screw-retained restorations

For optimal results, discuss the final restoration with your dental technician.



3. Radiographic planning



The vertical bone availability determines the maximum allowable length of the implant that can be placed and can be assessed on a radiographic image.


3.1 Outline important landmarks

Outline important landmarks on the periapical or panoramic radiograph using a transparent foil or paper and a pencil.



[Video: Treatment planning with the Straumann® X-ray Template.](#)

Outline important landmarks on the radiographs.



3.2 Evaluate and assess

3.2.1 Adjacent anatomical structures when determining available bone height



Note the proximity of this implant to the mesial root of the molar

Inclination and condition of adjacent roots

Convergent adjacent roots may compromise available interproximal space. Respect this root anatomy to avoid injury or interference during implant placement. Injury of adjacent teeth during implant placement can result in root canal treatment and/or damage to the root that leads to further tooth loss.


Consider the anatomical shape of adjacent tooth roots.



Assessment and treatment planning

Step 2 | Treatment planning

3. Radiographic planning



Note the proximity of the sinus floor to the apex of the implant

Landmarks in the maxilla such as the sinus floor: A safety distance of **2 mm** is recommended from the floor of the sinus to the apex of the implant. Straumann® Smart users are advised to treat straightforward cases only. If this minimum distance of bone is not available, the patient should be referred to a clinician with adequate training and experience in sinus floor elevation (SFE) procedures.

Maintain a minimum safety distance of 2 mm from the implant to the maxillary sinus floor.




Note the anatomical landmarks in the mandible shown by the pencil markings

Landmarks in the mandible: Inferior alveolar nerve, mental foramen and submandibular fossa/concavity. A safety distance of **2 mm** of bone between apex of the implant and neurovascular structures should be respected.

Maintain a minimum safety distance of 2 mm from the implant to the mandibular canal or mental foramen.

⚠ Caution: The loss of bone height due to ridge flattening during surgery must always be considered.

3.2.2 Interproximal distance of bone at the crestal bone level (interradicular space)



Pay attention to the interproximal distance of bone at bone level.

3.2.3 Bone quality

- Less dense or Type-IV-bone may indicate a preference to use the [Bone Level Tapered Implant](#) for better primary stability.

Bone which appears less dense in the radiograph indicates soft bone quality.




Assessment and treatment planning

Step 2 | Treatment planning

3. Radiographic planning

3.3 Measure the available vertical bone height and select the appropriate implant type, diameter and length



This can be done with the help of the [Straumann® X-ray Reference Sphere](#) and [X-ray Templates](#) for Straumann® Standard Plus Implants and Straumann® Bone Level Tapered Implants.


Distortions can occur in X-rays, hence the implant dimensions shown on these individual templates are designed with the corresponding distortion factors (1:1 to 1.7:1).

Determining each magnification factor or scale is facilitated by showing the X-ray Reference Sphere on the template (next to the scale reference).



[Video: Treatment planning with the Straumann® X-ray Template](#)


Use the Straumann® X-ray Template corresponding to the type of implant chosen.



3.3.1 Compare the size of the X-ray Reference Sphere*

First, compare the size of the X-ray Reference Sphere* on the patient's X-ray with the size of the X-ray Reference Sphere on the template (marked with arrow). Superimpose the two pictures to find the correct scale. Then, determine the spatial relations around the implant.

*The X-ray Reference Sphere has a diameter of 5 mm. The image of the sphere on the X-ray provides the reference value for the magnification scale



(049.076V4) = Ø 6.0 mm
Straumann® Bone Level Tapered Implant



1.2 : 1

Example:

Scale 1.2 : 1 = Reference Sphere
Ø 6.0 mm

Use the Straumann® X-ray Template and compare with the Straumann® X-ray Reference Sphere in the radiograph to determine the correct magnification factor.

3.3.2 Bone availability (X-ray):

On the X-ray, measure the height of the available bone from the crest of the edentulous ridge to the anatomical landmarks that limit the placement of the implant. Make sure you do this measurement in line with the axis of the proposed implant.



Assessment and treatment planning

Step 2 | Treatment planning

3. Radiographic
planning



Calculation of bone availability

To calculate the effective vertical bone availability, use the following formula:

$$\frac{\text{X-ray Reference Sphere 5 mm} \times \text{bone availability (X-ray*)}}{\text{Reference Sphere diameter on the X-ray}} = \text{Effective bone availability}$$

Use this formula to calculate the available vertical bone.

* Taking into consideration all implant-related anatomic structures (e.g. mandibular canal, sinus maxillaris, etc.)


$$\frac{5 \text{ mm} \times 13 \text{ mm}}{6 \text{ mm}} = 10,8 \text{ mm}$$

For example, if the measured bone availability and Reference Sphere diameter on the X-ray are 13 mm and 6 mm (+ 20 % distortion), respectively:

The effective bone availability would be 10.8 mm, therefore you would choose an implant length of 10 mm.

⚠ Caution: The loss of bone height due to ridge flattening during surgery must always be considered.

3.3.3 Consider the additional length of the drill tip



Due to the construction and function of the drills, the drill tip is a maximum of 0.4 mm longer than the implant insertion depth. This additional length must be taken into consideration during the planning phase.

The tip of the drill is up to 0.4 mm longer than the implant insertion depth.




Assessment and treatment planning

Step 2 | Treatment planning

3. Radiographic planning




3.3.4 Consider the coronoapical (vertical) implant position



The [Standard Plus Implant](#) with its smooth neck collar of 1.8 mm should be submerged in the bone as far as the margin of the rough implant surface. It can be placed slightly deeper, but the preparation depth must be increased accordingly.

Consider the implant's vertical position in the bone.

The [Bone Level Tapered Implant](#) is best placed with the outer rim of the small 45° sloping edge (chamfer) at bone level.





Assessment and treatment planning

Step 2 | Treatment planning



4. Risk
Assessment



Straumann® Smart users are advised to consider these pre-requisites for a straightforward single tooth case:

- Single tooth gaps (including free-end situations) in the posterior non-esthetic zone
- Healthy patient (ASA²-1: a normal healthy patient; or ASA²-2, a patient with mild systemic disease) with undisturbed wound healing capability
- Good patient motivation and compliance, preferably a non-smoker
- No active periodontitis or occlusal parafunction
- Healed ridge (≥ 4 months post-extraction) with sufficient bone volume
- Conventional loading protocol (≥ 3 months after implant placement)
- Sufficient interproximal distance of bone for the regular platform (RN or RC) implants or wide neck platform (WN) implants
- Sufficient bucco-lingual or bucco-palatal bone width and height for placing implants without the need for bone augmentation
- Sufficient inter-occlusal distance for the prosthetic restoration
- Minimum vertical oral opening space of 30 mm to allow access with surgical instruments

Ensure that the pre-requisites of a straightforward implant case are fulfilled.

⚠ Caution: Any pathologic condition such as caries, periapical infections, temporomandibular joint disorders or oral mucosal lesions should be treated prior to implant placement.

In cases of residual probing depths (PD) ≥ 5 mm with concomitant bleeding on probing, full-mouth plaque scores $> 20\%$, and associated risk factors, pre-treatment and periodontal re-evaluation are mandatory before implant placement.





Assessment and treatment planning

Step 2 | Treatment planning


5. Visit planning and cost proposal



5.3 Calculate costs and write up a cost proposal.

- Analyze and calculate the approximate total costs for the surgical treatment in comparison to alternative treatment modalities.
 - Dental implants are a cost-effective alternative to conventional restorations and improve the quality of life of the patient³.
 - Especially for single-tooth replacements, an implant is to be regarded as a cost-effective treatment option in comparison to a traditional three-unit fixed dental prosthesis⁴.


A cost estimate is important for discussion with the patient.



Graph: After 7 years, the costs for a single-tooth restoration pay off and the long-term maintenance costs are lower than for a conventional bridge⁴.



Chart of minimum widths of bone for planning which SP Implant (RN/WN) to use











Implant type (endosteal diameter)		Shoulder diameter (mm)	Interproximal distance at bone level (mm)	Bucco-lingual or bucco-palatal width of bone (mm)	Recommended use for Straumann® Smart cases
SP Ø 3.3 mm RN		4.8	6.5	5.5	For narrow interdental spaces and narrow partially or fully edentulous bone ridges. Caution/Precaution: Small-diameter implants are not recommended for the posterior region. Specific indications for Titanium SLA® Ø 3.3 mm Standard Plus RN Implants: These implants are to be used only in cases for the following indications: • Partially dentate jaws with implant-borne, fixed constructions: Combine with implants Ø 4.1 mm and splint the suprastructure.
SP Ø 4.1 mm RN		4.8	7	6	For use in the maxilla and mandible, for restoration of partially or fully edentulous patients.
SP Ø 4.8 mm RN		4.8	8	7	For use in the maxilla and mandible, for restoration of partially or fully edentulous patients in wide interdental spaces and bony ridges.
SP Ø 4.8 mm WN		6.5	8	7	



Chart of minimum widths of bone for planning which BLT (NC/RC) Implant to use



For BLT Implants, the shoulder diameter is the same as the endosteal diameter at bone level.

Implant type (endosteal diameter)		Shoulder diameter (mm)	Interproximal distance at bone level (mm)	Bucco-lingual or bucco-palatal width of bone (mm)	Recommended use for Straumann® Smart cases
BLT Ø 3.3 mm NC		3.3	6.5	5.5	For narrow interdental spaces and narrow partially or fully edentulous bone ridges. Caution/Precaution: Small-diameter implants are not recommended for the posterior region.
BLT Ø 4.1 mm RC		4.1	7	6	For use in the maxilla and mandible, for restoration of partially or fully edentulous patients.
BLT Ø 4.8 mm RC		4.8	8	7	For use in the maxilla and mandible, for restoration of partially or fully edentulous patients in wide interdental spaces and bony ridges.



Assessment and treatment planning

Step 2 | Treatment planning

REFERENCES

- 1 Buser D, von Arx T, ten Bruggenkate C, Weingart D. Basic surgical principles with ITI implants. Clin Oral Impl Res 2000; 11 (Suppl.): 59–68. C Munksgaard 2000.
- 2 American Society of Anesthesiologists® Physical Status Classification System, <http://www.asahq.org/resources/clinical-information/asa-physical-status-classification-system>.
- 3 Vogel R, Smith-Palmer J, Valentine W. Evaluating the Health Economic Implications and Cost-Effectiveness of Dental Implants: A Literature Review. Int J Oral Maxillofac Implants 2013;28:343–356. doi: 10.11607/jomi.2921.
- 4 Priest, GF, Priest JE. The Economics Of Implants For Single Missing Teeth. Dental Economics 2004;94(5):130-138.



Assessment and treatment planning

Step 2 | Treatment planning

DISCLAIMER

Straumann® Smart is a blended training and education program focused on the education of general dentists who want to become surgically active in the field of dental implantology. The program is limited to information pertaining to straightforward implant cases and focuses on a reduced portfolio of products that are suitable for the treatment of such cases.

All clinical Straumann® Smart content – such as texts, medical record forms, pictures and videos – was created in collaboration with Prof. Dr. Christoph Hämmerle, Prof. Dr. Ronald Jung, Dr. Francine Brandenburg-Lustenberger and Dr. Alain Fontollet from the University of Zürich, Clinic for Fixed and Removable Prosthodontics and Dental Material Science, Switzerland.

Straumann does not give any guarantee that Straumann® Smart provides sufficient knowledge or instruction for the dental professional to become surgically active in the field of implantology. It is the dental professional's sole responsibility to ensure that he/she has the appropriate knowledge and instruction before placing dental implants.

Straumann® Smart does not replace a careful and thorough analysis of each individual patient by a dental professional. Further, it does not imply any guarantee or warranty with regard to completeness of the information provided to the patient. It does not replace the dental professional's duty to inform the patient about the treatment, the products and the risks involved and to receive the patient's informed consent. The dental professional is solely responsible for determining whether or not a treatment or product is suitable for a particular patient and circumstances. Knowledge of dental implantology and instruction in the handling of the relevant products is always necessary and the sole responsibility of the dental professional. The dental professional must always comply with the individual product's Instructions For Use as well as all laws and regulations.

STRAUMANN DISCLAIMS, TO THE EXTENT POSSIBLE BY LAW, ANY LIABILITY, EXPRESS OR IMPLIED, AND BEARS NO RESPONSIBILITY FOR ANY DIRECT, INDIRECT, PUNITIVE, CONSEQUENTIAL OR OTHER DAMAGES, ARISING OUT OF OR IN CONNECTION WITH ANY INFORMATION PROVIDED TO PATIENTS, ERRORS IN PROFESSIONAL JUDGMENT, IN PRODUCT CHOICES OR PRACTICE IN THE USE OR INSTALLATION OF STRAUMANN PRODUCTS.

All clinical content as well as clinical and radiographic images are provided by courtesy of Prof. Dr. Christoph Hämmerle, Prof. Dr. Ronald Jung, Dr. Francine Brandenburg-Lustenberger and Dr. Alain Fontollet from the University of Zürich, Clinic for Fixed and Removable Prosthodontics and Dental Material Science, Switzerland.

International Headquarters

Institut Straumann AG

Peter Merian-Weg 12

CH-4002 Basel, Switzerland

Phone +41 (0)61 965 11 11

Fax +41 (0)61 965 11 01

www.straumann.com

© Institut Straumann AG, 2016. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.