PT Válido para EUA / EN Valid for the USA / ES Válido para los Estados Unidos / IT Valido per gli Stati Uniti / DE Gültig für USA / FR Valable pour les États-Unis / DA Gyldigt for USA / HU Az USA-ra érvényes / RO Valabil pentru SUA / SV Giltig i USA / BG Валидно за САЩ / HR Vrijedi za SAD / CS Platné pro USA / NL Geldig voor de VS / ET Kehtib Ameerika Ühendriikides / FI Voimassa Yhdysvalloissa / EL Έγκυρο για τις ΗΠΑ / GA Baili do S.A.M / LV Derīgs ASV / LT Galioja Jungtinėms Amerikos Valstijoms / MT Validu għall-Istati Uniti / PL Obowiązuje w Stanach Zjednoczonych / SK Platí pre USA / SL Velja za ZDA

English

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The instruments that integrate the kit are sold separately. This device is intended for a specialized procedure, which must be performed by professionals qualified in Dental Implants and with knowledge of the Neodent Easy Guide technique. For optimum results, use the product knowing the appropriate techniques and always apply them under appropriate conditions, in an operating room.

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SANITATION

This product must be correctly cleaned after each use. Proceed as follows:

Manual cleaning and disinfection Cleaning

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2. Soak the disassembled instruments for at least 1 min in the cleaning solution (CIDEZYME®, 1.6 % v/v) so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments. Assist cleaning by careful brushing with a soft brush. Sway movable parts several times during cleaning.
If applicable, rinse all lumens of the instruments for, at least, five times using a single-use syringe (minimum volume of 10 mL).
3. Soak the disassembled instruments for 15 minutes in the cleaning solution (CIDEZYME®, 1.6 % v/v) with ultrasonic treatment so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments.
4. Remove the instruments from the cleaning solution and intensively post-rinse them for, at least, 3 times (for the minimum time of 1 minute) under running water.
If applicable, rinse all lumens of the instruments for, at least, five times at the beginning of the soaking time using a single-use syringe (minimum volume of 10 mL).

Disinfection

1. Soak the disassembled instruments for 12 minutes in the disinfectant solution (CIDEX® OPA – OPA Solution -, undiluted) so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments.
If applicable, rinse all lumens of the instruments for, at least, five times at the beginning of the soaking time using a single-use syringe (minimum volume of 10 mL).
2. Remove the instruments from the disinfectant solution and post-rinse them according to the instructions of the manufacturer of CIDEX® OPA – OPA Solution -:

Rinsing Instructions

- Following removal of the instruments from CIDEX® OPA – OPA Solution – Solution, thoroughly rinse the medical device by immersing it completely in a large volume of water. Use sterile water unless potable water is acceptable (maximum of 10 microorganisms/mL, maximum of 0.25 endotoxin/mL).
 - Keep the device totally immersed for, at least, 1 minute.
 - Manually flush all lumens with large volumes (not less than 100 mL) of rinse water.
 - Remove the device and discard the rinse water. Always use fresh volumes of water for each rinse. Do not reuse the water for rinsing or any other purpose.
 - Repeat the procedure for 2 additional times, concluding A TOTAL OF 3 RINSES, with large volumes of fresh water to remove CIDEX® OPA – OPA Solution – Solution residues. Residues may cause serious side effects.
3. Check and pack the instruments immediately after the removal.

Automated cleaning/disinfection (WD (Washer-Disinfector))

Use neodisher® MediZym.

1. Disassemble the instruments if possible (see specific disassemble instructions for each instrument, when applicable).
2. Transfer the disassembled instruments to the WD (pay attention that the instruments have no contact to each other).

3. Start the program.
4. Remove the instruments from the WD after the end of the program.
5. Check and pack the instruments immediately after the removal.

NOTE:

1. Pay attention to the following points during the selection of the WD:
 - approved efficiency of the WD (e.g. CE marking according to EN ISO 15883 or DGHM or FDA approval/clearance/registration);
 - possibility of an approved program for thermal disinfection (A0 value > 3000 or – in case of older devices – at least 5 minutes at 90 °C/194 °F; in case of dangerous chemical disinfection of disinfectant remnants on the instruments);
 - use appropriate program for instruments, as well as sufficient information on rinsing in the program;
 - post-rinsing only with sterile or low contaminated water (e.g. maximum of 10 microorganisms/mL, maximum of 0.25 endotoxin/mL);
 - use only filtered air (oil-free, low contamination with microorganisms and particles) for drying;
 - regular maintenance and check/calibration of the WD.
2. Please do not clean any instruments by use of using metal brushes or steel wool.
3. Check all instruments after cleaning and disinfection, on corrosion, damaged surfaces, and impurities. Do not use damaged instruments. Instruments that are still contaminated must be cleaned and disinfected again.
4. Packaging: insert the cleaned and disinfected instruments on the corresponding sterilization trays, in single-use sterilization packagings (single or double packaging) and/or sterilization containers, which fulfill the following requirements:
 - EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance);
 - suitable for steam sterilization ;
 - sufficient protection for instruments as well as for maintenance of sterilization packagings against mechanical damage;
5. After using the instruments, it is recommended to remove coarse impurities, performing the pre-treatment, before cleaning and disinfection (within a maximum deadline of 2 hours). The pre-treatment step must be performed for both cases of cleaning and disinfection (automated and manual).
 - a. Disassemble the instruments if possible;
 - b. Rinse the instruments for, at least, 1 minute under running water (temperature < 35 °C);
 - c. If applicable, rinse all lumens of the instruments five times per application using a single-use syringe (minimum volume of 10 mL). Sway movable parts several times during pre-treatment;
 - d. Remove manually all visible impurities by using a clean and soft brush (or a clean, soft, and lint-free cloth). In no case use metal brushes or steel wool;
 - e. Rinse again for, at least, 1 minute under running water.
6. If the cleaning/disinfection products mentioned cannot be found, make sure to use similar products to those indicated. Replacement is the owner's responsibility.

7. The drying of the parts is of utmost importance before storage and sterilization, because the accumulation of moisture on the products is harmful and may cause oxidation.

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This product is reusable and supplied non-sterile, being unitarily packaged. This product must be correctly cleaned and sterilized before each use. Sterilize the products on the previous day or on the day of the procedure. ATTENTION: This product cannot be autoclaved in its original packaging.

Please use for sterilization only the steam sterilization according to the following parameters:

	Fractionated Vacuum / Dynamic Air Removal ¹	Gravity Displacement
Sterilization Time	4 minutes	15 minutes
Sterilization Temperature	132° C / 270° F	132° C / 270° F
Drying Time	At least 20 minutes ²	At least 20 minutes ²

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PRECAUTIONS

- In case the user has any type of color blindness or presents any difficulties in differentiating colors, it is recommended to use the laser markings for identification.
- This product must be used sterile. Do not sterilize under dry heat (oven).
- The Kit Case cannot touch the walls of the autoclave, so its melting is avoided.
- Chemical sterilization is not recommended.
- Place the water measure for the autoclave as recommended by the manufacturer. Lack of water during the sterilization cycle may compromise the effectiveness of the autoclave and melt the tray, as well as damage the instruments.
- The contact between the tray and acrylic resin liquid results in damage to the structure of the Kit Case polymer.
- Be aware in cases of patients that present signs of allergy or hypersensitivity to chemical elements of the materials: polymer.
- For sterilization, the instruments and the Kit Case must be cleaned individually, according to their own instructions and accommodated in the Kit Case tray, which must be wrapped in surgical paper for autoclaving.
- Inappropriate planning can compromise the performance of the implant/prosthesis unit leading to failures of the system, as loss or fracturing of the implant, loosening or fracturing

of the prosthetic screws.

- Do not use the product if its packaging is damaged.
- Do not use the product with its validity expired.
- Make sure you have all the instruments needed to perform the procedure according to the surgical planning.
- Make sure to use the indicated drilling sequence according to the instructions for use of the planned implant.
- Before each procedure, make sure the pieces are properly seated.
- Before each procedure, check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, deleted markings, deformation or wear.
- Always use the Neodent product sequence. The use of prosthetic abutments and/or instruments from other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the professional's responsibility to use the Neodent products according to the instructions for use.
- In case there is an accident related to the device, inform the manufacturer using the communication channels available on the website and the local competent authorities.

ADVERSE EFFECTS

Adverse effects are not expected as long as the product is used in accordance with its instructions for use.

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The professional in charge is responsible for providing these guidelines to the patient. Instruct the patient as to the need for a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. Evaluate the need for follow-up with periodic clinical appointments and radiographic examinations.

STORAGE CONDITIONS

This product must be stored, in its original packaging, in a clean and dry location, at room temperature and protected from direct sunlight.

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Written on the label.

LIFE CYCLE

This product is recommended for up to 60 uses as long as Neodent recommended conditions of use are respected. Regardless of the amount of times the instrument has been used, the professional must always evaluate its conditions before/after each use.

Each instrument has its own recommendations of amount of cycles. Please consult the instructions for use from the instruments that integrate the

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GM Easy Guide Surgical Kit Case Narrow/Regular – Tray Layout



Table 1 – Maximum Load Configuration

1	103.583	NGS NARROW MUCOSA PUNCH
2	103.519	NGS NARROW BONE LEVELLING DRILL
3	103.545	NGS NARROW INITIAL DRILL
4	103.546	NGS NARROW TAPERED DRILL D3.5X8
5	103.547	NGS NARROW TAPERED DRILL D3.5X10
6	103.548	NGS NARROW TAPERED DRILL D3.5X11.5
7	103.549	NGS NARROW TAPERED DRILL D3.5X13
8	103.550	NGS NARROW TAPERED DRILL D3.5/3.75X8
9	103.551	NGS NARROW TAPERED DRILL D3.5/3.75X10
10	103.552	NGS NARROW TAPERED DRILL D3.5/375X11.5
11	103.553	NGS NARROW TAPERED DRILL D3.5/3.75X13
12	103.554	NGS NARROW TAPERED DRILL D3.75X8

13	103.555	NGS NARROW TAPERED DRILL D3.75X10
14	103.556	NGS NARROW TAPERED DRILL D3.75X11.5
15	103.557	NGS NARROW TAPERED DRILL D3.75X13
16	105.161	NGS GM NARROW DRIVER FOR CONTRA-ANGLE
17	105.162	NGS GM NARROW DRIVER FOR TORQUE WRENCH
18	---	---
19	104.060	NEO MANUAL SCREWDRIVER (MEDIUM)
20	105.160	NEO LONG SCREWDRIVER FOR CONTRA-ANGLE
21	---	---
22	125.170	NGS GM NARROW STABILIZER
23	125.170	NGS GM NARROW STABILIZER
24	---	---
25	125.170	NGS GM NARROW STABILIZER
26	---	---
27	---	---

28	103.395	GUIDED SURGERY DRILL 1.3
29	125.142	NGS GUIDE CLAMP
30	125.142	NGS GUIDE CLAMP
31	125.142	NGS GUIDE CLAMP
32	125.176	NGS PALATAL SETTER
33	103.558	DRILL FOR NGS PALATAL SETTER
34	104.050	TORQUE WRENCH
35	129.034	DEPTH PROBE

The instruments that integrate the kit are sold separately. This device is intended for a specialized procedure, which must be performed by professionals qualified in Dental Implants, with knowledge of the Neodent Easy Guide technique. For optimum results, use the product knowing the appropriate techniques and always apply them under appropriate conditions, including in an operating room.

DESCRIPTION

GM Easy Guide Surgical Kit Case Narrow/Regular are containers for storage of surgical instruments during their use, cleaning, and sterilization. It is manufactured in autoclavable polymer and enables the storage and safe fastening of each instrument. The Kit presents markings and color identification that guide the use of the instruments during the procedure. The kit may be equipped according to each procedure.

INDICATIONS FOR USE

Neodent Instrument Kits are intended to be used to enclose other medical devices that are to be sterilized by a health care provider. Neodent Instrument Kits are intended to allow sterilization of the enclosed medical devices. Neodent Instrument Kits require the use of FDA-cleared wrap to maintain the sterility of the enclosed devices.

The kits are to be enclosed in a sterilization wrap that is FDA-cleared for the indicated cycles, and moist heat (steam) sterilized using one of the following cycles:

Fractionated vacuum (pre-vacuum) – Exposure at 132 °C for 4 minutes, 20-minute dry time
Gravity displacement – Exposure at 132 °C for 15 minutes, 20-minute dry time.

Neodent Instrument Kits are intended for sterilization of non-porous loads.

The combined weight of the GM Easy Guide Surgical Kit Case Narrow/Regular and the associated instruments is 310,18 g. The weight of the empty Kit Case is 263,63 g.

Neodent Instrument Kit Cases should not be stacked during sterilization.

APPLICATIONS

The GM Easy Guide Surgical Kit Case Narrow/Regular is suitable for safe storage of surgical instruments, organizing the instruments in the sequence in which they will be used, as well as for sterilization of such instruments.

The GM Easy Guide Surgical Kit Case Narrow/Regular is recommended for the installation of GM implants with diameters 3.5 and 3.75, using the Neodent Easy Guide surgical technique.

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KIT CASE ASSEMBLY

To set up this kit, each reserved space is related to a number from the instrument table; refer to image of the tray layout. The instruments that go in the tray are sold separately and must be used according to their own Instructions for Use. The maximum load configuration is shown in Table 1. The combined weight of the GM Easy Guide Surgical Kit Case Narrow/Regular and the associated instruments is 310,18 g. The weight of the empty Kit Case is 263,63 g.

CLEANING

Clean this product according to the following instructions after each use:

1. Remove the instruments from the Kit Case and clean separately (see specific instrument instructions).
2. Disassemble the three parts of the Kit Case. Soak the lid, tray and base for at least 1 min in the cleaning solution (Enzol® enzymatic detergent, prepared per manufacturer directions) so that the components are sufficiently covered. Pay attention that there is no contact between the components. Assist cleaning by careful brushing with a soft brush, paying particular attention to the interface between the silicone block inserts and the tray or base. Do not use metal brushes or steel wool.
3. Soak the lid, tray and base for 15 min in the cleaning solution (Enzol® enzymatic detergent, prepared per manufacturer directions) with ultrasonic treatment so that the components are sufficiently covered. Pay attention that there is no contact between the components.
4. Remove the Kit Case components from the cleaning solution and rinse them at least 3 times intensively (at least 1 min) under running water (local potable water is sufficient).
5. Repeat the cleaning steps of components that remain visibly contaminated.
6. Check all components after cleaning for signs of corrosion or damage. Do not continue to use corroded or damaged components.
7. Dry the components thoroughly before storage and sterilization; the accumulation of moisture on the components and instruments may cause corrosion that can detrimentally affect instrument performance.
8. Insert the cleaned, disinfected and dried instruments in the corresponding locations of the Case Kit tray and base. Assemble the Kit Case by inserting the tray into the base and attaching the lid.
9. Protect the assembled Kit Case from moisture and dust if it will not be used for an extended time.

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This product is reusable and supplied non-sterile, being unitarily packaged. This product must be correctly cleaned and sterilized before each use. Sterilize the products on the previous day or on the day of the procedure.

ATTENTION: This product cannot be autoclaved in its original packaging.

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	Fractionated Vacuum / Dynamic Air Removal ¹	Gravity Displacement
Sterilization Time	4 minutes	15 minutes
Sterilization Temperature	132 °C	132 °C
Drying Time	20 minutes ²	20 minutes ²

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Note: After sterilization, pack the instruments at a dry and dust-free place.

PRECAUTIONS

- In case the user has any type of color blindness or presents any difficulties in differentiating colors, it is recommended to use the laser markings for identification.
- This product must be used sterile. Do not sterilize under dry heat (oven).
- The Kit Case cannot touch the walls of the autoclave, so its melting is avoided.
- Chemical sterilization is not recommended.
- Place the water measure for the autoclave as recommended by the manufacturer. Lack of water during the sterilization cycle may compromise the effectiveness of the autoclave and melt the tray, as well as damage the instruments.
- The contact between the tray and acrylic resin liquid results in damage to the structure of the Kit Case polymer.
- For sterilization, the instruments and the Kit Case must be cleaned individually, according to their own instructions and accommodated in the Kit Case tray, which must be wrapped in sterilizable wrap that is FDA-cleared for the indicated cycles.
- Inappropriate planning can compromise the performance of the implant/prosthesis unit leading to failures of the system, as loss or fracturing of the implant, loosening or fracturing of the prosthetic screws.
- Do not use the product if its packaging is damaged.
- Do not use the product with its validity expired.
- Make sure you have all the instruments needed to perform the procedure according to the surgical planning.
- Before each procedure, check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, deleted markings, deformation or wear.
- Always use the Neodent product sequence. The use of prosthetic abutments and/or instruments from other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the professional's responsibility to use the Neodent products according to the instructions for use, and to determine whether it suits the individual situation of each patient.

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Instruct the patient as to the need for a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. Evaluate the need for follow-up with periodic clinical appointments and radiographic examinations. The professional in charge is responsible for providing these guidelines.

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DISPOSAL OF MATERIAL

All products and consumables used during the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to observe the current legislation and adhere to it.

DATE OF EXPIRATION

Written on the label.

LIFE CYCLE

This product is recommended for up to 100 uses as long as Neodent recommended conditions of use are respected.

Regardless of the amount of times the instrument has been used, the professional must always evaluate its conditions before/after each use. Visually inspect the lid, tray, and case bottom to ensure there is no cracking, deformation, or other damage.

Visually inspect that all labeling printed on the lid and tray is clear and legible.

Verify that the lid, tray, and case bottom can be assembled and that the lid latches securely to the case bottom.

Do not use the Kit Case if any of the above or any other damage is observed, regardless of the amount of cycles of use.

Each instrument has its own recommendations of amount of cycles. Please consult the instructions for use from the instruments that integrate the Kit and replace them accordingly.

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

Symbols	Description
	Medical device
	Date of manufacture
	Manufacturer
	Use-by date
Qty	Quantity
	Consult instructions for use
PROIBIDO REPROCESSAR	Reprocessing forbidden
STERILE EO	Sterilized using ethylene oxide
STERILE R	Sterilized using irradiation
	Do not reuse
	Do not use if package is damaged
	Do not resterilize
	Fragile
	Upper limit of temperature
	Keep dry
	Keep away from sunlight
EC REP	Authorized representative in the European Community
EC IMP	Importer in the European Community
CH REP	Swiss Authorised Representative
	Single sterile barrier system
	Single sterile barrier system with protective packaging inside
Rx only	Notification required by FDA for United States market
	CE Mark
NON-STERILE 	Non-sterile

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ENZOL® is a registered trademark of Johnson & Johnson.

Federal (USA) law restricts this device to sale by or on the order of a licensed dentist or physician. Not all products are available in all countries. Please, contact the authorized distributor.

Symbols	Description
SIZE	Product size
REF	Catalog number
LOT	Batch code
SN	Serial number
MAT	Product manufacturing material
UDI	UDI code

1	103.583	EXTRATOR DE MUCOSA NGS ESTREITO
2	103.519	BROCA NIVELADORA DE OSSO NGS ESTREITA
3	103.545	BROCA INICIAL NGS ESTREITA
4	103.546	BROCA CONICA NGS ESTREITA D3.5X8
5	103.547	BROCA CONICA NGS ESTREITA D3.5X10
6	103.548	BROCA CONICA NGS ESTREITA D3.5X11.5
7	103.549	BROCA CONICA NGS ESTREITA D3.5X13
8	103.550	BROCA CONICA NGS ESTREITA D3.5/3.75X8
9	103.551	BROCA CONICA NGS ESTREITA D3.5/3.75X10
10	103.552	BROCA CONICA NGS ESTREITA D3.5/3.75X11.5
11	103.553	BROCA CONICA NGS ESTREITA D3.5/3.75X13
12	103.554	BROCA CONICA NGS ESTREITA D3.75X8

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

1	103.583	NGS NARROW MUCOSA PUNCH
2	103.519	NGS NARROW BONE LEVELLING DRILL
3	103.545	NGS NARROW INITIAL DRILL
4	103.546	NGS NARROW TAPERED DRILL D3.5X8
5	103.547	NGS NARROW TAPERED DRILL D3.5X10
6	103.548	NGS NARROW TAPERED DRILL D3.5X11.5
7	103.549	NGS NARROW TAPERED DRILL D3.5X13
8	103.550	NGS NARROW TAPERED DRILL D3.5/3.75X8
9	103.551	NGS NARROW TAPERED DRILL D3.5/3.75X10
10	103.552	NGS NARROW TAPERED DRILL D3.5/3.75X11.5
11	103.553	NGS NARROW TAPERED DRILL D3.5/3.75X13
12	103.554	NGS NARROW TAPERED DRILL D3.75X8
13	103.555	NGS NARROW TAPERED DRILL D3.75X10

14	103.556	NGS NARROW TAPERED DRILL D3.75X11.5
15	103.557	NGS NARROW TAPERED DRILL D3.75X13
16	105.161	NGS GM NARROW DRIVER FOR CONTRA-ANGLE
17	105.162	NGS GM NARROW DRIVER FOR TORQUE WRENCH
18	---	---
19	104.060	NEO MANUAL SCREWDRIVER (MEDIUM)
20	105.160	NEO LONG SCREWDRIVER FOR CONTRA-ANGLE
21	---	---
22	125.170	NGS GM NARROW STABILIZER
23	125.170	NGS GM NARROW STABILIZER
24	---	---
25	125.170	NGS GM NARROW STABILIZER
26	---	---
27	---	---
28	103.395	GUIDED SURGERY DRILL 1.3
29	125.142	NGS GUIDE CLAMP
30	125.142	NGS GUIDE CLAMP

31	125.142	NGS GUIDE CLAMP
32	125.176	NGS PALATAL SETTER
33	103.558	DRILL FOR NGS PALATAL SETTER
34	104.050	TORQUE WRENCH
35	129.034	DEPTH PROBE

The instruments that integrate the kit are sold separately.

This device is intended for a specialized procedure, which must be performed by professionals qualified in Dental Implants and with knowledge of the Neodent Easy Guide technique. For optimum results, use the product knowing the appropriate techniques and always apply them under appropriate conditions, in an operating room.

DESCRIPTION

Kit Case for storage of surgical instruments for use in implant installation procedures using the Neodent Easy Guide technique. Made of autoclavable polymer, the case has silicone supports for safe conditioning and fixation of each instrument during sanitation, sterilization and surgical procedures. It has markings and color identification that organize and guide its use. The Kit Case can be equipped according to each procedure.

INDICATIONS FOR USE

Product indicated for conditioning and safe fastening surgical and/or prosthetic instruments during their use and sterilization. The use of this product rationalizes the storage and organization of instruments during and after the dental procedure.

APPLICATIONS

The kit is used for the safe storage of surgical instruments, as well as for support during sterilization. The GM Easy Guide Surgical Kit Case Narrow/Regular is recommended for the installation of GM implants with diameters 3.5 and 3.75, using the Neodent Easy Guide Technique. Its target population is presented dental surgeons who perform dental implant installation procedures.

WARNING

As the clinical outcome of dental treatment is influenced by multiple variables, the residual risks described below can occur even if the product is used according to the instructions for use. These are: allergy, discomfort, pain, edema, treatment failure, local infection, inflammation, local irritation, pyrogenic reaction, loss of product function, prosthetic and/or aesthetic and/or biomechanical issues.

CONTRAINDICATIONS

This product presents no contraindications as long as it is used according to its instructions for use. Limitations of use for each instrument must be observed in its corresponding instruction for use.

HANDLING

Remove the lid to access the instruments by gently pressing on the latch system between the case bottom and the lid. To close the Kit Case, place the lid on the case bottom and press it down. During the surgical procedure, the tray must be placed on a sterile surface and handled with a sterile glove.

KIT CASE ASSEMBLY

To set up this kit case, each reserved space is related to a number from the instrument table; refer to image of the tray layout. The instruments that go in the cassette are sold separately and must be used according to their own Instructions for Use.

SANITATION

This product must be correctly cleaned after each use.

Proceed as follows:

Manual cleaning and disinfection

Cleaning

1. Disassemble the instruments if possible (see specific disassemble instructions for each instrument, when applicable).

2. Soak the disassembled instruments for at least 1 min in the cleaning solution (CIDEZYME®, 1.6 % v/v) so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments. Assist cleaning by careful brushing with a soft brush. Sway movable parts several times during cleaning.

If applicable, rinse all lumens of the instruments for, at least, five times using a single-use syringe (minimum volume of 10 mL).

3. Soak the disassembled instruments for 15 minutes in the cleaning solution (CIDEZYME®, 1.6 % v/v) with ultrasonic treatment so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments.

4. Remove the instruments from the cleaning solution and intensively post-rinse them for, at least, 3 times (for the minimum time of 1 minute) under running water.

If applicable, rinse all lumens of the instruments for, at least, five times at the beginning of the soaking time using a single-use syringe (minimum volume of 10 mL).

Disinfection

1. Soak the disassembled instruments for 12 minutes in the disinfectant solution (CIDEX® OPA – OPA Solution -, undiluted) so that the instruments are sufficiently covered. Pay attention that there is no contact between the instruments. If applicable, rinse all lumens of the instruments for, at least, five times at the beginning of the soaking time using a single-use syringe (minimum volume of 10 mL).

2. Remove the instruments from the disinfectant solution and post-rinse them according to the instructions of the manufacturer of CIDEX® OPA – OPA Solution -:

Rinsing Instructions

• Following removal of the instruments from CIDEX® OPA – OPA Solution – Solution, thoroughly rinse the medical device by immersing it completely in a large volume of water. Use sterile water unless potable water is acceptable (maximum of 10 microorganisms/mL, maximum of 0.25 endotoxin/mL).

• Keep the device totally immersed for, at least, 1 minute.

• Manually flush all lumens with large volumes (not less than 100 mL) of rinse water.

• Remove the device and discard the rinse water. Always use fresh volumes of water for each rinse. Do not reuse the water for rinsing or any other purpose.

• Repeat the procedure for 2 additional times, concluding A TOTAL OF 3 RINSES, with large volumes of fresh water to remove CIDEX® OPA – OPA Solution – Solution residues. Residues may cause serious side effects.

3. Check and pack the instruments immediately after the removal.

Automated cleaning/disinfection

(WD (Washer-Disinfector))

Use neodisher® MediZym.

1. Disassemble the instruments if possible (see specific disassemble instructions for each instrument, when applicable).

2. Transfer the disassembled instruments to the WD (pay attention that the instruments have no contact to each other).

3. Start the program.

4. Remove the instruments from the WD after the end of the program.

5. Check and pack the instruments immediately after the removal.

NOTE:

1. Pay attention to the following points during the selection of the WD:

- approved efficiency of the WD (e.g. CE marking according to EN ISO 15883 or DGHM or FDA approval/clearance/registration);

- possibility of an approved program for thermal disinfection (A0 value > 3000 or – in case of older devices – at least 5 minutes at 90 °C/194 °F; in case of dangerous chemical disinfection of disinfectant remnants on the instruments);

- use appropriate program for instruments, as well as sufficient information on rinsing in the program;
- post-rinsing only with sterile or low contaminated water (e.g. maximum of 10 microorganisms/mL, maximum of 0.25 endotoxin/mL);

- use only filtered air (oil-free, low contamination with microorganisms and particles) for drying;
- regular maintenance and check/calibration of the WD.

2. Please do not clean any instruments by use of using metal brushes or steel wool.

3. Check all instruments after cleaning and disinfection, on corrosion, damaged surfaces, and impurities. Do not use damaged instruments. Instruments that are still contaminated must be cleaned and disinfected again.

4. Packaging: insert the cleaned and disinfected instruments on the corresponding sterilization trays, in single-use sterilization packagings (single or double packaging) and/or sterilization containers, which fulfill the following requirements:

- EN ISO/ANSI AAMI ISO 11607 (for USA: FDA clearance);

- suitable for steam sterilization ;

- sufficient protection for instruments as well as for maintenance of sterilization packagings against mechanical damage;

5. After using the instruments, it is recommended to remove coarse impurities, performing the pre-treatment, before cleaning and disinfection (within a maximum deadline of 2 hours).

The pre-treatment step must be performed for both cases of cleaning and disinfection (automated and manual).

a. Disassemble the instruments if possible;

b. Rinse the instruments for, at least, 1 minute under running water (temperature < 35 °C);

c. If applicable, rinse all lumens of the instruments five times per application using a single-use syringe (minimum volume of 10 mL). Sway movable parts several times during pre-treatment;

d. Remove manually all visible impurities by using a clean and soft brush (or a clean, soft, and lint-free cloth). In no case use metal brushes or steel wool;

e. Rinse again for, at least, 1 minute under running water.

- If the cleaning/disinfection products mentioned cannot be found, make sure to use similar products to those indicated. Replacement is the owner's responsibility.
- The drying of the parts is of utmost importance before storage and sterilization, because the accumulation of moisture on the products is harmful and may cause oxidation.

NOTE: During cleaning, try to avoid contact between cutting tools and other tools so the cutting power is not harmed.

PRESENTATION AND STERILIZATION

For traceability, each kit case presents on its tray the laser engraving of the UDI code (Unique Device Identification).

This product is reusable and supplied non-sterile, being unitarily packaged. This product must be correctly cleaned and sterilized before each use. Sterilize the products on the previous day or on the day of the procedure. ATTENTION: This product cannot be autoclaved in its original packaging. Please use for sterilization only the steam sterilization according to the following parameters:

	Fractionated Vacuum / Dynamic Air Removal ¹	Gravity Displacement
Sterilization Time	4 minutes	15 minutes
Sterilization Temperature	132 °C / 270 °F	132 °C / 270 °F
Drying Time	At least 20 minutes ²	At least 20 minutes ²

¹ At least three vacuum steps.

² The effectiveness required in drying time depends directly on the parameters of sole responsibility of the user (density and load configuration, sterilizing conditions, which must be determined by the user). Nevertheless, a drying time shorter than 20 minutes cannot be applied.

NOTE: After sterilization, pack the instruments at a dry and dust-free place.

PRECAUTIONS

- In case the user has any type of color blindness or presents any difficulties in differentiating colors, it is recommended to use the laser markings for identification.
- This product must be used sterile. Do not sterilize under dry heat (oven).
- The Kit Case cannot touch the walls of the autoclave, so its melting is avoided.
- Chemical sterilization is not recommended.
- Place the water measure for the autoclave as recommended by the manufacturer. Lack of water during the sterilization cycle may compromise the effectiveness of the autoclave and melt the tray, as well as damage the instruments.
- The contact between the tray and acrylic resin liquid results in damage to the structure of the Kit Case polymer.
- Be aware in cases of patients that present signs of allergy or hypersensitivity to chemical elements of the materials: polymer.
- For sterilization, the instruments and the Kit Case must be cleaned individually, according to their own instructions and accommodated in the Kit Case tray, which must be wrapped in surgical paper for autoclaving.
- Inappropriate planning can compromise the performance of the implant/prosthesis unit leading to failures of the system, as loss or

fracturing of the implant, loosening or fracturing of the prosthetic screws.

- Do not use the product if its packaging is damaged.
- Do not use the product with its validity expired.
- Make sure you have all the instruments needed to perform the procedure according to the surgical planning.
- Make sure to use the indicated drilling sequence according to the instructions for use of the planned implant.
- Before each procedure, make sure the pieces are properly seated.
- Before each procedure, check the conditions of the Neodent surgical instruments, always respecting their useful life. Replace the instruments if there is damage, deleted markings, deformation or wear.
- Always use the Neodent product sequence. The use of prosthetic abutments and/or instruments from other manufacturers does not ensure the perfect function of the Neodent Implant System and exempts any product warranty.
- It is the professional's responsibility to use the Neodent products according to the instructions for use.

ADVERSE EFFECTS

Even if the product is used according to the instructions for use, the residual risks described under "Warning" can occur.

POST-OPERATIVE PRECAUTIONS AND MAINTENANCE

The professional in charge is responsible for providing these guidelines to the patient.

Instruct the patient as to the need for a professional medical monitoring after the surgery and to obey the guidelines regarding the precautions, hygiene and prescription of drugs. Evaluate the need for follow-up with periodic clinical appointments and radiographic examinations.

STORAGE CONDITIONS

Until its first use, the product must be stored in its original packaging, in a clean and dry place, at room temperature and protected from direct sunlight. For the storage of reusable items after the first use, follow the guidelines in the "Sanitation" and "Presentation and Sterilization" sections of these instructions for use.

DISPOSAL OF MATERIAL

All products and consumables used during the surgery for the installation of dental implants may endanger the health of anyone who handles them. Before discarding them into the environment, it is recommended to observe the current legislation and adhere to it. In case there is no current legislation, pack them in sharps collector and dispose them in medical waste.

DATE OF EXPIRATION

Written on the label.

LIFE CYCLE

This case is recommended for up to 60 uses, as long as the conditions of use recommended by Neodent are observed. Regardless of how many times the instruments have been used, the professional should always assess the condition of the instruments before and after each use. Each instrument has its own recommended number of uses. Please refer to the instructions

for use of the instruments that make up the case and replace them according to the guidelines.

SYMBOLS

The table describes the symbols that may be printed on the product labeling. Please, refer to the physical labeling to see which symbols are applicable to the product.



Symbols	Description
	Product size
	Catalog number
	Batch code
	Serial number
	Product manufacturing material
	UDI code
	Medical device
	Date of manufacture
	Manufacturer
	Use-by date
	Quantity
	Consult instructions for use
	Reprocessing forbidden
	Sterilized using ethylene oxide
	Sterilized using irradiation
	Do not reuse
	Do not use if package is damaged
	Do not resterilize
	Fragile
	Upper limit of temperature
	Keep dry
	Keep away from sunlight
	Authorized representative in the European Community
	Importer in the European Community
	Swiss Authorised Representative
	Single sterile barrier system
	Single sterile barrier system with protective packaging inside
	Notification required by FDA for United States market
	CE Mark

Symbols	Description
NON-STERILE 	Non-sterile

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Not all products are available in all countries. Please contact your authorized distributor.

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© 2021 – JJGC Indústria e Comércio de Materiais Dentários S.A. All rights reserved.
Neodisher® MediZym is a registered trademark of Chemische Fabrik Dr. Weigert GmbH & Co. KG. CIDEZYME® and CIDEX® OPA are a registered trademark of Johnson & Johnson.
© 2023 – JJGC Indústria e Comércio de Materiais Dentários S.A. All rights reserved.
© 2023 – JJGC Indústria e Comércio de Materiais Dentários S.A. Id-drittijiet kollha miżmuma. Neodisher® MediZym hija trademark reġistrata ta' Chemische Fabrik Dr. Weigert GmbH & Co. KG. CIDEZYME® u CIDEX® OPA huma trademarks reġistrati ta' Johnson & Johnson.
© 2023 – JJGC Indústria e Comércio de Materiais Dentários S.A. Vse pravice pridržane.
Neodisher® MediZym je registrirana blagovna znamka družbe Chemische Fabrik Dr. Weigert GmbH & Co. KG.