Sterile Decellularized Dermis

DESCRIPTION

Straumann AlloGraft products are manufactured by LifeNet Health, 1864 Concert Drive, Virginia Beach, Virginia 23453, and distributed in the US by Straumann USA, LLC and in Canada by Straumann Canada Limited. Decellularized dermis was processed from donated human tissue, resulting from the generous gift of an individual or his/her family. The dermis was processed using a proprietary technology, which safely renders the dermal matrix acellular without compromising the biomechanical properties for its intended surgical applications. Decellularized dermis achieves a sterility assurance level (SAL) of 10⁻⁶ via gamma irradiation and is preserved with glycerol.

INDICATIONS FOR USE

Decellularized dermis serves as a scaffold which is suitable for the reinforcement of damaged or inadequate integumental tissue at the surgical site.

CONTRAINDICATIONS

The contraindications include, but are not limited to:

- Use in any patient who has a known or suspected allergy to any of the antibiotics and/or processing reagents listed in this package insert.

WARNINGS AND PRECAUTIONS

The same medical/surgical conditions or complications that apply to any surgical procedure may occur during or following implantation. The surgeon is responsible for informing the patient of the risks associated with their treatment and the possibility of complications or adverse reactions. As with any allograft bio-implant, the potential for transmission of infectious agents exists.

This bio-implant may contain residuals of antibiotics (Gentamicin, Lincomycin, Polymyxin B Sulfate, and/or Vancomycin), N-Lauroyl Sarcosinate (detergent), Benzoylase® or Denarase® (recombinant endonuclease), and/or glycerol. Caution should be exercised if the patient has a known sensitivity to any of these antibiotics and/or reagents.

POTENTIAL ADVERSE EVENTS

Potential adverse events or outcomes include, but are not limited to, infection, allograft tissue rejection, allergic reaction to residual processing reagents, loss of bio-implant structural integrity, reoperation and/or death.

Promptly report any adverse event(s) or outcome(s) potentially attributable to the decellularized dermis (See COMPLAINTS AND RETURNS section).

DONOR SCREENING AND TESTING

All donors have been screened and tissues recovered, processed, stored, tested and distributed in accordance with current U.S. federal regulations as promulgated in 21 CFR 1270 and 1271, current Standards for Tissue Banking set forth by the American Association of Tissue Banks (AATB) and international laws and regulations as required.

This allograft bio-implant was deemed suitable for implantation by LifeNet Health. A physician medical director evaluated the following donor variables to determine donor suitability: infectious disease test results, current donor medical history, behavioral risk assessment interview, physical assessment, relevant medical records, including previous medical history, laboratory test results, and autopsy or coroner reports (if performed).

All donors are tested for relevant infectious diseases. Testing is performed by laboratories that are registered with the U.S. Food and Drug Administration (FDA) and certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) and 42 CFR 493. Test methods that are FDA-licensed, approved, or cleared for donor screening are used as available. The following test criteria were met for the donor of this allograft bio-implant:

Required Infectious Disease Testing	
Test	Acceptance Criteria
HBcAb: Hepatitis B Total Core Antibody	Negative/Non-Reactive
HBsAg: Hepatitis B Surface Antigen	Negative/Non-Reactive
HCV NAT: Hepatitis C Virus Nucleic Acid Test	Negative/Non-Reactive
HCvAb: Hepatitis C Antibody	Negative/Non-Reactive
HBV NAT: Hepatitis B Virus Nucleic Acid Test*	Negative/Non-Reactive
HIV-1 NAT: Human Immunodeficiency Virus Type 1 Nucleic Acid Test	Negative/Non-Reactive
HIV 1/2 Ab: Human Immunodeficiency Virus Types 1/2 Antibody	Negative/Non-Reactive
RPR/STS or Equivalent: Syphilis	Confirmatory Negative/Non-Reactive
HTLV I/II Ab: Human T-Lymphotropic Virus Types I/II Antibody**	Negative/Non-Reactive

*Not required for donors recovered after December 16, 2016. Performed as required by international laws and regulations.

**Not required for donors recovered prior to March 31, 2010. Performed as required by international laws and regulations.

STORAGE REQUIREMENTS

The distributor, intermediary, and/or end-user clinician or facility is responsible for storing decellularized dermis under appropriate conditions prior to further distribution or implantation. Dermis must be stored as listed in the table below.

Storage Temperature	Special Conditions
Refer to the label	Do not freeze or refrigerate
	Store in its original cardboard sleeve
	Minimize excessive exposure to light and protect from excessive heat

The packaging may contain a temperature sensitive dot that will turn from white to pink or red if the upper temperature limit has been exceeded. Do not use the decellularized dermis if the temperature dot appears to be a color other than white.

DOCUMENTATION REQUIREMENTS

Recipient records must be maintained by the practitioner for the purpose of tracking tissue post-transplantation. Please make note of the S/N number as indicated on the package in the patient's chart. Labels are provided with each product for your convenience. Straumann® also provides an online Patient Tracking System. To access, go to <https://www.straumann.com/us/en/dental-professionals/allograft.html> and complete the Patient Tracking Form as needed. This online form does not replace the need for the practitioner to track patient recipient records of allograft tissue use.

INSTRUCTIONS FOR USE

It is important to read and understand the following instructions prior to clinical use. Improper preparation technique may adversely affect handling properties and/or performance.

PLEASE READ THESE INSTRUCTIONS AS THEY ARE DIFFERENT THAN THOSE FOR TRADITIONAL CRYOPRESERVED SKIN ALLOGRAFTS.

GENERAL INSTRUCTIONS:

- Use on a single occasion for a single patient only.
- Once the packaging is opened, the dermis must be used for the current procedure or discarded.
- Any unused dermis must be discarded.
- Inspect the dermis, inner and outer packaging, and labels carefully:
 - Do not use past the expiration date as indicated on the label.
 - Do not use if the dermis is damaged or the packaging integrity is compromised.
 - Do not use if there are discrepancies in label information.
 - When the temperature dot is present, do not use the dermis if the dot appears to be a color other than white.
- Use aseptic technique at all times.
- Do not re-sterilize.
- Keep the dermis stored according to recommended storage instructions until preparing it for implantation.

PREPARATIONS FOR USE:

ORIENTATION: Decellularized dermis has two physically distinct sides; a reticular side and a papillary side. In general, when applied, the papillary side will face up while the reticular side is placed against the surgical wound or the most vascularized tissue. The dermis is packaged with the papillary side visible through the clear side of the packaging.

- Non-Sterile Team Member:** Open the cardboard sleeve and retrieve the pouch from within.
- Aseptically open the outer peel pack and present inner pouch to the **Sterile Team Member**.
- Sterile Team Member:** To maintain orientation of the dermis, the papillary side should be marked with a sterile marker immediately after opening the inner pouch. **The dermis is packaged with the papillary side visible through the clear side of the packaging.** Dermis: Open the inner peel pouch and remove the dermis along with its slip sheet. Remove the slip sheet prior to application. Biowashers: Open the inner peel pouch and remove the dermis while maintaining orientation of the individual grafts. Biowashers do NOT have a slip sheet.
- NOTE: Rinsing is not required prior to application; however, it may improve handling.** If a rinse is preferred by the physician, continue to the rinse instructions below.

If not used immediately, keep dermis moist until implantation.

RINSE INSTRUCTIONS (Optional)

- Non-Sterile Team Member:** Prepare a sterile rinse basin with enough sterile isotonic solution (e.g., sterile saline) to completely cover the dermis. **CAUTION: Ensure the rinse solution does not exceed 42°C as this may damage the dermis.**
- Sterile Team Member:** After opening the packaging per the instructions above, remove the dermis from the slip sheet and immerse the dermis in sterile isotonic solution for a minimum of 1 minute. Ensure the dermis is completely submerged in solution during the rinse.
- Keep the dermis completely submerged in sterile isotonic solution until needed.

The maximum sterile isotonic solution exposure time for decellularized dermis is 4 hours.

TRACEABILITY

Regulations pertaining to human tissue products in the United States and Canada only require Straumann to track distributed tissue to the customer facility. It is the responsibility of the implanting facility to maintain a record of the specific patient who has received a human tissue implant. This Tissue Tracking System is provided as a courtesy to our customer. It is not a substitute for your own patient records. Straumann accepts no liability for maintaining, verifying or otherwise collecting the information identified above.

COMPLAINTS AND RETURNS

Please contact Customer Service for information regarding Straumann's Return Policy.

WARRANTY STATEMENT

Due to the inherent variability of allograft tissue, biological and biomechanical properties cannot be guaranteed by LifeNet Health or SWAL.

Straumann Canada Limited CTO Establishment #100190
Source Establishment: LifeNet Health CTO #100038

Manufactured by:
LifeNet Health
1864 Concert Drive
Virginia Beach, VA 23453 U.S.A.

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Straumann Canada Limited CTO Établissement n°100190