## STERILE FREEZE-DRIED ALLOGRAFT BIO-IMPLANT

### DESCRIPTION

Straumann® AlloGraft products are manufactured by LifeNet Health, 1864 Concert Drive, Virginia Beach, Virginia, 23453, and distributed in the US by Straumann USA LLC and in Canada by Straumann Canada Limited. This allograft bio-implant was processed from donated human tissue, resulting from the generous gift of an individual or his/her family. This bio-implant was processed using LifeNet Health's patented and proprietary Allowash XG® technology. The tissue is first disinfected through the Allowash® process, and then terminally sterilized via gamma radiation.

### INDICATIONS FOR USE

Straumann® AlloGraft is indicated for filling and/or augmenting intraoral/maxillo-facial osseous defects, such as intrabony periodontal osseous and furcation defects, augmentation of bony defects of the alveolar ridge, filling tooth extraction sites, and sinus elevation grafting.

### CONTRAINDICATIONS

This bio-implant may contain residuals of antibiotics (Gentamicin and/or Vancomycin), alcohol, surfactants, and/or glycerol. Caution should be exercised if the patient has a known sensitivity to any of these antibiotics.

### WARNINGS AND PRECAUTIONS

- Federal law (USA) restricts this tissue for use by a licensed clinician only.
- Tissue may transmit infectious agents.
- As individual patient responses to receipt of allograft will vary, LifeNet Health makes no claims concerning the biologic or biomechanical properties of allograft tissue.
- For single patient use only.
- For use on a single occasion only.
- Once the packaging has been opened, the tissue must either be transplanted or discarded.
- Do not use if package integrity has been compromised.
- Do not re-sterilize.

### POTENTIAL ADVERSE EVENTS

Any adverse events or outcomes must be reported promptly. Please call our Regulatory Affairs department and have the S/N number available.

### TISSUE RETURNS

Please contact Customer Service for information regarding Straumann's Return Policy.

### DONOR SCREENING AND TESTING

Donor-eligibility determination has been performed by LifeNet Health, 1864 Concert Drive, Virginia Beach, Virginia, 23453; LifeNet Health is a full service, not-for-profit tissue bank, an accredited member of AATB, and an FDA registered and ISO 13485:2016 certified company. All tissues have been screened, recovered, manufactured and distributed according to current FDA and American Association of Tissue Banks (AATB) standards and regulations. This human tissue has been determined to be suitable for transplantation by LifeNet Health. A physician medical director has evaluated the following to determine suitability: infectious disease test results, current donor medical history and behavioral risk assessment interview, physical assessment, relevant medical records, including previous medical history, laboratory test results, and autopsy or coroner reports (if performed).

All donors are tested for relevant infectious diseases. Testing is performed by laboratories that are registered with the U.S. Food and Drug Administration (FDA) and certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) and 42 CFR 493. Test methods that are FDA-licensed, approved, or cleared for donor screening are used as available. The following test criteria were met for the donor of this allograft bio-implant:

Required Infectious Disease Testing	
Test	Acceptance Criteria
HBcAb: Hepatitis B Total Core Antibody	Negative/Non-Reactive
HBsAg: Hepatitis B Surface Antigen	Negative/Non-Reactive
HCV NAT: Hepatitis C Virus Nucleic Acid Test	Negative/Non-Reactive
HCVAb: Hepatitis C Antibody	Negative/Non-Reactive
HBV NAT: Hepatitis B Virus Nucleic Acid Test*	Negative/Non-Reactive
HIV-1 NAT: Human Immunodeficiency Virus Type 1 Nucleic Acid Test	Negative/Non-Reactive
HIV 1/2 Ab: Human Immunodeficiency Virus Types 1/2 Antibody	Negative/Non-Reactive
RPR/STS or Equivalent: Syphilis	Confirmatory Negative/Non-Reactive
HTLV I/II Ab: Human T-Lymphotropic Virus Types I/II Antibody**	Negative/Non-Reactive

\*Not required for donors recovered after December 16, 2016. Performed as required by international laws and regulations.

\*\*Not required for donors recovered prior to March 31, 2010. Performed as required by international laws and regulations.

### STORAGE REQUIREMENTS

- Freeze-Dried - store at ambient temperature.
- Do not freeze.
- The tissue dispensing service and/or the end-user clinician or facility is responsible for storing this product under appropriate conditions prior to transplantation.

### DOCUMENTATION REQUIREMENTS

Recipient records must be maintained by the practitioner for the purpose of tracking tissue post-transplantation. Please make note of the S/N number as indicated on the package in the patient's chart. Labels are provided with each product for your convenience. Straumann® also provides an online Patient Tracking System. To access, go to <https://www.straumann.com/us/en/dental-professionals/allograft.html> and complete the Patient Tracking Form as needed. This online form does not replace the need for the practitioner to track patient recipient records of allograft tissue use.

### STRAUMANN'S COMMITMENT TO QUALITY

Straumann® is committed to providing high quality products that are produced under a rigorous quality assurance program. As such, Straumann® has partnered with LifeNet Health, a known leader in the allograft tissue market, to provide this product to you. If you have any questions or comments, we would like to hear from you. Please contact our Customer Service department to speak with one of our representatives.

### PREPARATIONS FOR USE:

#### Preparation Notes:

Use aseptic technique at all times.

#### Non-Sterile Team Member:

1. Visually inspect packaging to ensure that it is intact and that its integrity has not been breached. If the packaging is damaged, the enclosed graft may be contaminated and should not be used.
2. Peel open outer tray foil lid stock and present contents to the Sterile Team Member.
3. If tissue is packaged in a jar, firmly grasp the jar and remove from outer tray. If tissue is packaged in a plastic tray, firmly grasp the "Lift Here" tab and remove from outer tray.
4. Visually inspect the container to ensure that it is intact and that its integrity has not been breached. If the packaging is damaged, the enclosed graft may be contaminated and should not be used.
5. If rehydrating, hydrating media may include antibiotic solution, sterile saline, I.V. fluids, plasma, bone marrow, or other specific blood components. Rehydrate until required consistency and handling are achieved as per clinician preference.
6. The graft is ready for implantation. Do not allow graft to dry. Keep graft completely submerged in sterile solution until surgically implanted.

If the surgery is postponed, discard the allograft.

For patent information, please visit: [www.lifenethhealth.org/patents](http://www.lifenethhealth.org/patents)

Straumann Canada Limited CTO Establishment #100190  
Source Establishment: LifeNet Health CTO #100038

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Straumann Canada Limited CTO Établissement n°100190