

# Technical Sheet

## GLO Take-Home Whitening Kit GLO Science, Inc.



### PRODUCT INDICATION

|              |                         |
|--------------|-------------------------|
| Product type | Dental esthetic device  |
| Indication   | Teeth whitening at home |

### PRODUCT PROPERTIES

|                              |  |
|------------------------------|--|
| Origin                       | USA  |
| Storage temperature          | Room temperature   |
| Product specifications       | <ul style="list-style-type: none"><li>• Take-Home teeth whitening kit</li><li>• Contains patented G.L.O (Guided Light Optics) warming heat &amp; light technology in a universally-sized mouthpiece with 6 % hydrogen peroxide whitening gel vials</li><li>• One patient treatment 4 x 8min = 32 min with 1 vial</li><li>• Whitening mouthpiece is comprised of a flexible circuit board containing light- and heat-emitting elements, encased within a silicone resin overmold</li><li>• Mouthpiece warms up to 52°C to activate the whitening gel and accelerate the formation of whitening oxygen ions</li><li>• GLO Control battery pack provides power to whitening mouthpiece via Lithium Ion polymer battery</li><li>• GLO Control contains an automatic 8-minute off-switch/timer; users administer up to four 8-minute applications at once</li><li>• Includes the same professional teeth whitening GLO control device, which is used in the office for GLO chairside whitening treatments and is available to patients exclusively through dental professionals</li></ul> |
| Component parts              | <p>Each kit includes :</p> <ul style="list-style-type: none"><li>• Professional GLO Control Device and Mouthpiece</li><li>• 10 Professional-strength Whitening gel GLO Vials 6%</li><li>• Mouthpiece case</li><li>• USB cable &amp; power adaptor</li><li>• Travel case</li><li>• User manual</li></ul>  |
| Sterilization / Disinfection | <ul style="list-style-type: none"><li>• Rinse home kit mouthpiece with water to clean, not autoclavable</li><li>• Wipe GLO control clean</li></ul>   |
| Shelf life                   | 2 years  |
| Regulatory status            | Cosmetic<br>CE   |
| Entry on the market          | 2014   |

## MANUFACTURING

### Packaging

Cardboard boxes

## AVAILABLE PRODUCT LINE

|  | Description / Packaging  | Art.         |
|--|--|--------------|
|  | 6 boxes GLO Take-Home Whitening Kits 6% HP (HP-Hydrogen Peroxide)<br>Weight (kg): 3.06<br>Dimensions (cm): 41.3(L) x 25.4(W) x 25.46(H)      | 311106EU     |
|  | 1 box GLO Take-Home Whitening Kit 6% HP (HP-Hydrogen Peroxide)<br>Weight (kg): 0.51<br>Dimensions (cm): 20.32(L) x 20.32(W) x 7.62(H)        | 458113EU     |
|  | 6 boxes GLO Vials Refill (x10), 6% HP (HP-Hydrogen Peroxide) + Lip Care<br>Weight (kg): 2.7<br>Dimensions (cm): 23.1(L) x 20.3(W) x 17.8(H)  | 2581052106EU |
|  | 1 box GLO Vials Refill (x10), 6% HP (HP-Hydrogen Peroxide) + Lip Care<br>Weight (kg): 0.45<br>Dimensions (cm): 15.24(L) x 17.78(W) x 1.91(H) | 258105210EU  |

## HOTLINE

### UK

• (+44) 01 293 651230

### GERMANY

• (+49) 761 450 13 50

### ITALY

(+39) 800 810 500

## International Headquarters

Institut Straumann AG  
Peter Merian-Weg 12  
CH-4002 Basel, Switzerland  
Phone +41 (0)61 965 11 11  
Fax +41 (0)61 965 11 01  
www.straumann.com

© Institut Straumann AG, 2019. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

