








## Explanation of coding and types

Color coding		
	blue	Endosteal implant diameter 2.9 mm
	yellow	Endosteal implant diameter 3.3 mm
	red	Endosteal implant diameter 4.1 mm
	green	Endosteal implant diameter 4.8 mm

Color coding – Straumann® SNOW		
	yellow	Endosteal implant diameter 3.6 mm
	red	Endosteal implant diameter 4.0 mm
	green	Endosteal implant diameter 5.0 mm

Implant types		
TL	S	Standard Implant
	SP	Standard Plus Implant
	MI	Mini Implant
	CIM	Ceramic Implant Monotype
	CI	Ceramic Implant
BL	BL	Bone Level Implant
	BLT	Bone Level Tapered Implant
	SNOW BL	SNOW BL Ceramic Implant

Connection types	
NNC: Narrow Neck CrossFit® Ø 3.5 mm	Ø 3.5 mm 
RN: Regular Neck Ø 4.8 mm	Ø 4.8 mm 
WN: Wide Neck Ø 6.5 mm	Ø 6.5 mm 
MI: Mini Implant System with Optiloc® Retentive System	Ø 3.4 mm 
SC: Small CrossFit® Ø 2.9 mm	Ø 2.9 mm 
NC: Narrow CrossFit® Ø 3.3 mm	Ø 3.3 mm 
RC: Regular CrossFit® Ø 4.1 and Ø 4.8 mm	Ø 4.1 mm    Ø 4.8 mm 
ND: Narrow Diameter Ø 3.5 mm	Ø 3.5 mm 
RD: Regular Diameter Ø 4.8 mm	Ø 4.8 mm    Ø 4.8 mm 

Surgery