

axiom[®]
tissue**LEVEL**



**COMPONENTS
OVERVIEW.**

ANTHOGYR 
A Straumann Group Brand

AXIOM® TISSUE LEVEL SURGICAL COMPONENTS

HEALFIT® SH ANATOMICAL AND SCANNABLE HEALING SCREWS				
Manual ≤10 N.cm STERILE Manual hexagonal surgical wrench OPCS100	Shape	Coronary height	PLATFORM R/4.8 mm	PLATFORM N/4.0 mm
	A	2.5 mm	-	TSHSA-N200
		4.5 mm	-	TSHSA-N400
	B	2.5 mm	TSHSB-R200	TSHSB-N200
		4.5 mm	TSHSB-R400	TSHSB-N400
	C	2.5 mm	TSHSC-R200	TSHSC-N200
		4.5 mm	TSHSC-R400	TSHSC-N400
	D	2.5 mm	TSHSD-R200	-
		4.5 mm	TSHSD-R400	-
	E	2.5 mm	TSHSE-R200	-
		4.5 mm	TSHSE-R400	-

Healfit® SH

Anatomical and scannable healing solution with biomimetic design

- 3-in-1 solution: Heal – Shape – Take the impression.
- 5 anatomical shapes based on the study of hundreds of tooth profiles.
- Open software solution available in major planning softwares, with a direct access to CAD libraries on Anthogyr website.





DISCOVER MORE



Axiom® TL countersinks for optimized bone crest preparation		
	PLATFORM R/4.8 mm	PLATFORM N/4.0 mm
Ø 2.4 mm	TOPFF24R	TOPFF24N
Ø 3.0 mm	TOPFF30R	TOPFF30N
Ø 3.6 mm	TOPFF36R	TOPFF36N
Ø 4.2 mm	TOPFF42R	-

HEALING SCREWS				
 STERILE  Manual hexagonal surgical wrench  OPCS100	Conical	mm	PLATFORM R/4.8 mm	PLATFORM N/4.0 mm
		2		
			THS-R210	THS-N210
		4		
		THS-R410	THS-N410	
	Straight	2		
		THS-R200	THS-N200	

IMPRESSION TAKING								TEMPORIZATION			
SINGLE RESTORATION								SINGLE RESTORATIONS			
TRANSFER	CONVENTIONAL				DIGITAL			Laboratory scan adapter	Indexed temporary abutment	PLATFORM R/4.8 mm	PLATFORM N/4.0 mm
	PLATFORM R 4.8 mm		PLATFORM N 4.0 mm		PLATFORM R 4.8 mm		PLATFORM N 4.0 mm				
	Pop-in transfer indexed				Digital transfer for Intra Oral Scanner impression						
Screw included	Short	Long	Short	Long					 M1.6 titanium screw included		
					PEEK	Metal	PEEK				
	TT200-R	TT200L-R	TT200-N	TT200L-N							
					156-02-DT	156-02-S-MDT	156-01-DT	156-01-S-MDT	156-0X-SAA	TC100-R	TC100-N
Long and short screws included	Pick-up transfer indexed										
											
	TT300-R			TT-300-N							
ANALOG	Implant analog				Analog for printed model			Wrench for scan adapter			
											
	TA100-R (x1) TA100-R-4 (x4)	TA100-N (x1) TA100-N-4 (x4)		Wrench for analog for printed model							
					PA-TOOL-01			SATOOL-01			
PLURAL RESTORATIONS								PLURAL RESTORATIONS			
TRANSFER	CONVENTIONAL				DIGITAL			Laboratory scan adapter	Non indexed temporary abutment	PLATFORM R/4.8 mm	PLATFORM N/4.0 mm
	PLATFORM R 4.8 mm		PLATFORM N 4.0 mm		PLATFORM R 4.8 mm		PLATFORM N 4.0 mm				
	Pop-in transfer non indexed				Digital transfer for Intra Oral Scanner impression						
Screw included									 M1.6 titanium screw included		
				PEEK	Metal	PEEK	Metal				
	ILT200-R-4 (x4)	ILT200-N-4 (x4)									
				156-02-DT-IL	156-02-S-MDT	156-01-DT-IL	156-01-S-MDT	156-0X-SAA	TC100-R-P	TC100-N-P	
Long and short screws included	Pick-up transfer non indexed										
											
	ILT100-R (x1) ILT100-R-4 (x4)	ILT100-N (x1) ILT100-N-4 (x4)									
ANALOG	Implant analog				Analog for printed model			Wrench for scan adapter			
											
	TA100-R (x1) TA100-R-4 (x4)	TA100-N (x1) TA100-N-4 (x4)		Wrench for analog for printed model							
					PA-TOOL-01			SATOOL-01			

AXIOM® TISSUE LEVEL PROSTHETIC COMPONENTS

ABUTMENT					
SINGLE & PLURAL CEMENTED RETAINED RESTORATIONS					
Esthetic abutments	Angulation	PLATFORM R/4.8 mm		PLATFORM N/4.0 mm	
		Ø 5.0 mm	Ø 6.0 mm	Ø 4.0 mm	Ø 5.0 mm
	0°				
	7°				
	15°				
	23°				

Axiom® prosthetic kit

INMODOPP3



TI-BASE AND PRE-MILLED FOR LABSIDE CAD/CAM																	
SINGLE CEMENTED OR SCREW RETAINED RESTORATIONS																	
Ti-base  M1.6 black screw included  T3160  Hexagonal wrench	CEREC® <table><tr><td>Ø 5.0 mm Size L Block</td><td>Ø 4.0 mm Size S Block</td></tr><tr><td>PLATFORM R/4.8 mm</td><td>PLATFORM N/4.0 mm</td></tr><tr><td> TBASEC-R-L</td><td> TBASEC-N-S</td></tr></table>	Ø 5.0 mm Size L Block	Ø 4.0 mm Size S Block	PLATFORM R/4.8 mm	PLATFORM N/4.0 mm	 TBASEC-R-L	 TBASEC-N-S										
	Ø 5.0 mm Size L Block	Ø 4.0 mm Size S Block															
	PLATFORM R/4.8 mm	PLATFORM N/4.0 mm															
	 TBASEC-R-L	 TBASEC-N-S															
	CUSTOMIZED TITANIUM ABUTMENTS																
PRE-MILLED  Definitive titanium M1.6 screw T3161 included   Hexagonal wrench	<table><tr><th>PreFace abutments</th><th>TI-Forms-Abutments forM-Series</th></tr><tr><td>MEDENTIKA® holder</td><td>Amann Girschbach® holder</td></tr><tr><td>PLATFORM R/4.8 mm</td><td>PLATFORM N/4.0 mm</td></tr><tr><td>PLATFORM R/4.8 mm</td><td>PLATFORM N/4.0 mm</td></tr><tr><td></td><td></td></tr><tr><td></td><td></td></tr><tr><td>Ø 11.5 mm:</td><td>TPM-M-TI-R11 TPM-M-TI-N11 TPM-C-TI-R11 TPM-C-TI-N11</td></tr><tr><td>Ø 16.0 mm:</td><td>TPM-M-TI-R16 TPM-M-TI-N16 - -</td></tr></table>	PreFace abutments	TI-Forms-Abutments forM-Series	MEDENTIKA® holder	Amann Girschbach® holder	PLATFORM R/4.8 mm	PLATFORM N/4.0 mm	PLATFORM R/4.8 mm	PLATFORM N/4.0 mm					Ø 11.5 mm:	TPM-M-TI-R11 TPM-M-TI-N11 TPM-C-TI-R11 TPM-C-TI-N11	Ø 16.0 mm:	TPM-M-TI-R16 TPM-M-TI-N16 - -
	PreFace abutments	TI-Forms-Abutments forM-Series															
	MEDENTIKA® holder	Amann Girschbach® holder															
	PLATFORM R/4.8 mm	PLATFORM N/4.0 mm															
	PLATFORM R/4.8 mm	PLATFORM N/4.0 mm															
																	
																	
Ø 11.5 mm:	TPM-M-TI-R11 TPM-M-TI-N11 TPM-C-TI-R11 TPM-C-TI-N11																
Ø 16.0 mm:	TPM-M-TI-R16 TPM-M-TI-N16 - -																

Torq Control®

TS501

Universal torque wrench




AXIOM® TISSUE LEVEL PROSTHETIC COMPONENTS

X-Base®

Labside innovative range for greater flexibility

- Patented laser surface treatment for increased adhesion and time saving.
- Large choice of coronary heights and gingival profiles.
- Intuitive & accurate index.
- Angulated Access up to 25°.

INSPIRED BY SPACE EXPLORATION PATENTED



X-BASE® FOR CAD/CAM LABSIDE					
SINGLE SCREW-RETAINED RESTORATIONS					
X-Base	Ch*	PLATFORM R/4.8 mm		PLATFORM N/4.0 mm	
		TFLEX-R4-S	TFLEX-N4-S	TS160	TS163 short
	4				
	6				
	Cast coping				TS162 long
With Angulated Access					
X-Base	Ch*	TFLEX-R4-SAA / TFLEX-R4-SAAU		TFLEX-N4-SAA / TFLEX-N4-SAAU	
		4			TS163P-2 short
	5				TS162P-2 long
	6			TFLEXS-AA	TFLEXSL-AA
	Cast coping AA 25°				
PLURAL SCREW-RETAINED RESTORATIONS					
X-Base®	Ch*	TFLEX-R4-P		TFLEX-N4-P	
		4			TS163P-2 short
	5				TS162P-2 long
	6			TFLEXS-AA	TFLEXSL-AA
	Cast coping				
With Angulated Access / FOR PLURAL					
X-Base® AA*	Ch*	TFLEX-R4-PAA		TFLEX-N4-PAA	
		4			TFLEXS-PAA
	5				TFLEXSL-PAA-2 short
	6				
	Cast coping AA 25°				

AXIOM® TISSUE LEVEL PROSTHETIC COMPONENTS

AXIN® SCREW RETAINED TOOTH WITH AA*			AXIN® INSTRUMENTS & SCREWS			
TEMPORIZATION			AXIN® wrench	AXIN-TOOL-AML		
AXIN® Base	PLATFORM R/4.8 mm	PLATFORM N/4.0 mm		 AXINLABKITAML		
	 AXIN156-02-B	 AXIN156-01-B				
	 AXIN156-0X-S					
AXIN® temporary coping	AXIN-C-48-10	AXIN-C-40-10	AXIN® laboratory kit	 AXINLABKITAML		
	 AXIN-C-48-15	AXIN-C-40-15				
	 AXIN-C-48-25	AXIN-C-40-25				
AXIN® channel rod Ø 2.0 mm For temporary or permanent	AXIN-C-48-10	AXIN-C-40-10	AXIN® clinical screw	 AXIN156-0X-S		
	AXIN-C-48-15	AXIN-C-40-15				
	AXIN-C-48-25	AXIN-C-40-25				
AXIN® protection cap	AXIN-PCC-48	AXIN-PCC-40	AXIN® laboratory screw	 AXIN156-0X-SL		
LABORATORY COMPONENTS			MULTI-UNIT LABORATORY COMPONENTS			
AXIN® laboratory bases	PLATFORM R/4.8 mm	PLATFORM N/4.0 mm	Protective analog	Ø 4.8 mm / R	Ø 4.0 mm / N	
	 AXIN156-02-L	 AXIN156-01-L		 MUA200 (x1) MUA200-4 (x4)	 MUNA200	
	 AXIN156-0X-SL			 CR25-4 (x4)		
TITANIUM BASES FOR MULTI-UNIT			Channel rod Ø 2.5 mm for MUA			
Straight						
PLATFORM R/4.8 mm	PLATFORM N/4.0 mm					
Ch: 5 mm	MUFLEX-5	Ch: 4 mm	MUNFLEX-4			
MUFLEX		MUNFLEX				
X-Base with Angulated Access			Temporary Abutments	Ø 4.8 mm / R	Ø 4.0 mm / N	
Ch: 5 mm	MUFLEX-5-AA	Ch: 4 mm		MUNFLEX-4-AA	MUC100	MUNC100
MUFLEX-AA25		MUNFLEX-AA25			MUC200	MUNC200
			Hexagonal wrench MU titanium screw included	MUC200	MUNC200	
				Peek	MUC200	MUNC200
				with Angulated Access	MUC200	MUNC200
			Ball wrench	MUC200	MUNC200	
				Angulation 0°	MUC200	MUNC200
				10°	MUC200	MUNC200
			M1.4 Titanium screw with AA included MUA141	15°	MUC200	MUNC200
				25°	MUC200	MUNC200
					MUC200	MUNC200

AXIOM® TISSUE LEVEL PROSTHETIC COMPONENTS

INSTRUMENTS			
Hexagonal wrench	SHORT	LONG	XL
	INCHECV (14 mm)	INCHELV (21.5 mm)	INCHEXLV (29 mm)
	INMHECV (20 mm)	INMHELV (27 mm)	INMHEXLV (35 mm)
	INBW100S (23 mm) INBW100S-3 (X3)	INBW100L (29 mm) INBW100L-3 (X3)	INBW100XL (37 mm) INBW100XL-3 (X3)
	INBM100S (20.5 mm) INBM100S-3 (X3)	INBM100L (27 mm) INBM100L-3 (X3)	INBM100XL (35 mm) INBM100XL-3 (X3)
Hexagonal mandrel	MUM100 (16 mm)	MUM100L (25.5 mm)	MUM100XL (35 mm)
Ball wrench			
Ball mandrel			
Mandrel for straight MU			
Wrench for straight MU abutments			
Short prehensile tool for MU			
SCREWS			
Black M1.6 prosthetic screw (single)			TS160
Titanium M1.6 prosthetic screw (single)			TS161
Short LAB M1.6 screw (single)			TS163
Long LAB M1.6 screw (single)			TS162
Black M1.6 prosthetic screw (plural)			TS160P
Titanium M1.6 prosthetic screw (plural)			TS161P
Black M1.6 prosthetic screw with AA for Simeda® (plural)			TSAA160P
Titanium M1.6 prosthetic screw with AA for Simeda® (plural)			TSAA161P
Short LAB M1.6 screw (plural)			TS163P-2 (x2)
Long LAB M1.6 screw (plural)			TS162P-2 (x2)
Short laboratory M1.6 screw with AA for Simeda® (plural)			TSAA163P-2 (x2)
AA screw gripper			AATool / AATool-4 (x4)
M1.4 MU black screw for ti cocr prosthesis			MUA140
M1.4 MU titanium screw for zirconium prosthesis			MUA141
M1.4 MU black screw for ti cocr prosthesis			MU140
M1.4 MU titanium screw for zirconium prosthesis			MU140Z
M1.4 MU titanium lab screw blue			MU141
MU short lab screw			MUT101
MU long lab screw			MUT102

AXIOM® TISSUE LEVEL PROSTHETIC COMPONENTS

ABUTMENT					
Definitive abutment	REMOVABLE OVERDENTURE				
	NOVALOC®		LOCATOR®		
	GINGIVAL HEIGHT	PLATFORM R/4.8 mm	PLATFORM N/4.0 mm	GINGIVAL HEIGHT	PLATFORM R/4.8 mm
	0.0 mm	TOPNOVA000-R	TOPNOVA000-N	0.5 mm	O2830
	1.0 mm	TOPNOVA010-R	TOPNOVA010-N	1.5 mm	O2831
Impression components	2.0 mm	TOPNOVA020-R	TOPNOVA020-N	2.5 mm	O2832
	3.0 mm	TOPNOVA030-R	TOPNOVA030-N	3.5 mm	O2833
	4.0 mm	TOPNOVA040-R	TOPNOVA040-N	4.5 mm	O2834
Secondary parts	Impression coping	2010.722-NOV (x4)	Transfer	O8505 (x4)	
	Model analog	2010.721-NOV (x4)	Analog	O8530 (x4)	
	Block out spacer	2010.723-NOV (x4)			
	Processing collar	2010.724-NOV (x10)			
	Processing insert	2010.725-NOV (x4)			
Instruments	Matrix housing	2010.701-NOV (x4)	«Pink» light retention insert	O8527 (x4)	Insert, Mid Arch
	Matrix housing - PEEK	2010.702-NOV (x4)	«Blue» extra light retention insert	O8529 (x4)	Denture attachment housing
	Matrix housing - extended	2010.703-NOV (x4)	«Grey» 0-retention insert	O8558 (x4)	Processing Package Anterior/Posterior
	Extra-light retention insert, red	2010.710-NOV (x4)	«Red» extra light retention insert	O8548 (x4)	Processing Package Mid Arch
	Light retention insert, white	2010.711-NOV (x4)	«Orange» light retention insert	O8915 (x4)	Seating and Removal Tool 046.828
Instruments	Medium retention insert, yellow	2010.712-NOV (x4)	«Green» medium retention insert	O8547 (x4)	Systems Core Tool
	Strong retention insert, green	2010.713-NOV (x4)			
	Extra-strong retention insert, blue	2010.714-NOV (x4)			
	Ultra-strong retention insert, black	2010.715-NOV (x4)			
	Processing package Titanium	2010.601-NOV (x2)	Locator® torque wrench insert latch	O8913	Seating Tip Assembly
Instruments	Processing package PEEK	2010.611-NOV (x2)			
	Equipment box with 3 instruments	2010.101-NOV	Locator® torque wrench insert driver	O8260	Removal Tip Assembly 10156-01
	Processing insert removal instrument	2010.731-NOV			Tip Wrench 10155-01
	Retention insert instrument	2010.741-NOV			Optional Removal Tip 10095
	Matrix housing extraction instrument	2010.751-NOV	Locator® 3 in 1 core tool	O8393	Systems Core Tool Collet 8348