

Sina XT-T: Discover the latest Simea[®] zirconia

Sina XT-T is a new hybrid zirconia with gradient of translucency for a smooth transition from incisal to cervical.

Its mechanical strength has been validated to be used for all types of Simea[®] restorations in zirconia performed at Anthogyr Center.

It is a new choice in the Anthogyr Web Order, in addition to the current zirconias Sina Z, T, ML. Prosthesis in Sina XT-T will be at the same price as the equivalent prosthesis in Sina ML

BENEFITS



+ Flexibility and time saving :

- Suitable for homothetical, cut-back and anatomical designs, according to indications
- Less time of veneering, according to the design

+ Know-how and quality :

- Implant bridge without Ti-base, with Angulated Access solution on Multi-Unit abutments
- High mechanical strength from 1000 MPa in the incisal area to 1200 MPa in the cervical area
- Incisal part milled in the most translucent part of the disk
- More than 10 years of expertise in zirconia
- Serenity[®] warranty : 5 years

+ Wide choice of indications :

- Abutment, partial to full bridge on Anthogyr, Straumann[®], Neodent[®] and MPS platforms
- Single and plural tooth supported restorations
- Choice between 16 shades according to VITA shade guide and 1 Bleach shade with a gradient of translucency from 44% to 51% in the incisal part

ANTHOGYR SA

8, um Miercherbiereg P.B 184,
L-7526 Mersch – Luxembourg

Follow us



*Multi-Platform Solutions