

simeda[®]



Customised
prosthetic
solutions




Anthogyr
A Straumann Group Brand

Anthogyr CAD-CAM solutions

Anthogyr produces customised prostheses adapted to all indications and workflows.

- A single partner
- A machining centre of excellence



CAD-CAM made accessible to all

A **Simeda**® solution for any laboratory organisation:

- Preparation of prostheses from your STL files
- Design of prostheses from digitised or non-digitised models
- Digitisation of your models and design of your prostheses

simeda®


Why choose Anthogyr CAD-CAM solutions?

Single **practitioner** - prosthetist dual partner

- High precision
- Innovation
- Expertise in CAD-CAM machining
- Premium materials
- Complete prosthetic range



Simeda® an added value service offering

- Technical support: personalised assistance with expert dental technicians 
- Training
- Order online through Web Order
- **Serenity**®* warranty programme

*See the general warranty conditions for the Anthogyr Serenity® Programme.

Simeda®

High precision and innovative CAD-CAM prosthetic solutions for all types of single or multiple restorations.


INNOVATION

inLink®

A unique lock system fully integrated into the prosthesis, specifically for screw-retained multiple restorations on Axiom® Multi Level®.

Angulated Access

Integrated inLink® lock system



Angulated Access on Multi-Unit Axiom®BL & TL abutments and Connect+* MPS platforms

Angulated Access that threads screws through the screw channel.
Angulations:
10°, 15°, 20°, 25°.

*Full list of compatible platforms available at www.anthogyr.com

Customised multiple-unit prostheses

- Cemented and screw-retained
- Bridges, suprastructures, bars
- Passive solutions
- Easy to maintain

Premium materials

Simeda®, specialising in ceramics, also offers a wide range of materials.

- Sina Z®, Sina T®, Sina ML zirconias
- Titanium and Cobalt-chromium
- VITA SUPRINITY® PC, VITA ENAMIC® and IPS e.max® CAD

Customised single-unit prostheses

INNOVATION

AxiN®

The biological safety of a customised **screw-retained zirconia** for all Axiom® Multi Level® implant sectors.


Abutments and screw-retained teeth

Improved design of the crown for greater durability and better cosmetic support.



DigitALL All inclusive

A fully digital workflow from an optical impression taken with an intra-oral scanner.



High precision M.P.S (Multi-Platform Solutions) prosthetic solutions.

More than 250 implant platforms produced by the Anthogyr CAD-CAM machining centre including some of the most popular ones on the market:

- Straumann®
 - Nobel Biocare®
 - Zimmer Biomet®
 - Dentsply Sirona®
- Full list available at www.anthogyr.com

Solutions SCAN-CAD

- High performance MEDIT® Identica T500 & T300 scanners.
- Exocad® design software.



Solutions Lab**side**



Labside library

A library of CAD software to create your customised prostheses in the laboratory on the **Axiom® Multi Level®** platform.

- Axiom® Multi Level® scan-adapters
- for laboratory scanners
- Axiom® Multi Level® digital transfers
- for intraoral scanners
- Axiom® Multi Level® analogues
- for impression models.

 dental wings

3shape 

exocad

C 184_EN - 2019-12

Photo credits: Anthogyr, Getty images, Julien Cury
All rights reserved, non-binding pictures.



ANTHOGYR SA

8, um Miercherbiërg P.B 184,
L-7526 Mersch - Luxembourg

www.anthogyr.com

 **Anthogyr**
A Straumann Group Brand