

The background of the entire image is a vibrant, deep blue space filled with numerous white stars of varying sizes. Swirling throughout the scene are ethereal nebulae in shades of red, pink, and purple, some appearing as wispy clouds and others as more dense, filamentary structures. In the upper center, a large, spherical planet with a blue and white swirling pattern is partially visible. In the lower right, a ringed planet, similar to Saturn, is seen from a distance. The central focus of the image is the Axiom X3 spacecraft, a large, dark grey, multi-segmented structure with a complex, angular design. It is oriented diagonally, with its front end pointing towards the upper right. A bright, white light source, possibly a star or a distant planet, is positioned behind the spacecraft, creating a strong backlighting effect and casting a long, bright glow across the scene. The overall composition is dynamic and futuristic, emphasizing the theme of space exploration.

axiomx3[®]

BEYOND FRONTIERS

Anthogyr
A Straumann Group Brand



**DISCOVER
A NEW KIND
OF IMPLANT
SOLUTION...**

axiom^{x3}

Explore new territories, take up new implant challenges to better meet the needs of all your patients.

A holistic response adapted to clinical diversity, while respecting and preserving the most precious environment: bone.

Always bring you more advantages: simplify your actions, offer agile and adapted protocols according to indications and in particular immediate treatments.

At your side, whatever your daily challenge, an effective clinical response, a new experience, beyond your frontiers.

Anthogyr
A Straumann Group Brand

01

Universal bone anchorage

Enjoy a high performance design for **universal anchorage** in all bone densities

- Designed for less constraints in hard bone with minimal instrument sequence
- Flexible drilling protocol for additional stability in soft bone
- Adapted protocols to suit immediate treatment

ENTER NEW WORLDS

02

Bone preservation

More space for bone, decreasing insertion efforts in hard bone without sacrificing stability and efficient osteotomy

- Axiom X3® allows to save 52% bone compared to a classical protocol*
- Minimally invasive and flexible drilling protocols, less drills, NO tap:

BONE IS GOLD

*In soft bone average of 52% of bone preserved compared to Axiom® REG protocol - internal test.

03

Full control guidance

More than a powerful insertion: go where you want to go!

- Make the choice of optimal guidance and insertion
- Designed for dynamic engagement and optimized stability, suitable in all situations, including immediate treatments



axiomx3®

B E Y O N D F R O N T I E R S

Slim threads

Easier insertion

Alternative thread length

*Bone preservation
Lower insertion efforts
Lateral guidance*

3 bi-directional cutting flutes

*Dynamic engagement all over the implant height
Insertion axis correction*

Thin apex

Efficient engagement in slim implant beds

BONE CARE DESIGN
PATENTED

Single conical connection

*Common to
Axiom® Bone Level - Ø 2.7 mm*

Reduced neck diameter

Less cortical stress

Tapered profile

*Optimal primary stability
in soft bone*

Rounded apex

Sinus membrane protection




04

Frictionless integration

Enjoy the Axiom® ecosystem, designed for efficient implantology

- A unique conical connection to simplify prosthetic steps
- A single surgery kit, common to all Axiom® implants
- Integration into most digital platforms and flows
- Access to Simeda® customized prosthesis



Axiom® confidence, almost 15 years of clinical experiences and success

05

Answer patients needs

Meet the requirements of your patients, offer **modern therapeutic options**, whatever their clinical situation

- Chair time saving oriented, minimally invasive protocols, short to **immediate treatments**
- Designed to answer the diversity of indications



Anthogyr
A Straumann Group Brand

C220_GB - 2022-02 Photos credits: Anthogyr. All rights reserved, actual products may vary

C220_GB - 2022-02

MORE THAN AN IMPLANT
A NEW UNIVERSE



EXPLORE
NEW WORLDS
NOW

STAY TUNED!



anthogyr.com

Anthogyr
A Straumann Group Brand