



Connected Implantology motor



Onthogyr Performance

at the service of patients

Connected Implantology motor

With **Xpert Unit**[®], enjoy the convenience of using a simple and effective implantology motor. Discover the 1st connected motor for true innovation in terms of services and patient follow-up.

SIMPLICITY

The features are mastered from the first use.

Performance

Effective manoeuvres during surgery thanks to the high torque, up to 80 N.cm.



Wireless pedal Freedom of movement



Anthogyr Xpert Unit[®] - Connected Implantology motor

Total freedom of movement



PERFORMANCE

Pedal command during surgery:

- Navigation through different programmes.
- Reverse operation.
- Irrigation management.
- Torque increase and decrease.



3 FREEDOM

Wireless pedal: One less constraint in your working environment!

Appealing and powerful

SIMPLICITY

- Large touchscreen: intuitive navigation.
- Direct use from the first screen.
- Implantology sequence, simple and ready to use.
- Surgery programme.

ERGONOMICS

- Compact unit: space optimisation.
- Streamlined design: easy maintenance.
- Easy placement of irrigation lines.

FLEXIBILITY

 All parameters are easy to change: torque, speed, irrigation flow.

COMFORT AND SAFETY

- Large screen.
- Displaying key information during surgery: excellent legibility.

An application to make your life easier

INNOVATION

TRACEABILITY

- Synchronisation of surgery data. • Organisation of data by patient to
- facilitate information management.

EASY ACCESS TO INFORMATION

- User manual, video tutorials, FAQs, technical support

The extra service



PEACE OF MIND

Up to 1000 hours of operation: **Extended operation**

DURABILITY

Rechargeable battery: practical and durable.



Mont-Blanc XS

A concentration of power in your hands!

ULTRA LIGHT

Easy to hold Micromotor and Mont Blanc XS contra-angle set.

EFFICIENCY

N cm

- Speed: 15 to 2000 rpm.
- Torque: 5 to 80 N.cm: the pleasure of a high tightening torque.

WORKING COMFORT

- Easy to hold.
- Light: more visibility in the mouth.
- Light, flexible cord: easier handling.

SAFETY AND DURABILITY

- Thermodisinfectable and autoclavable micromotor: easier maintenance.
- Removable, thermodisinfectable and autoclavable contra-angle: adapted to the strictest maintenance protocols.

COMPUS



Training center in implantology

- 1 operating theatre, 1 dental prosthesis laboratory
- 800 dental professionals trained each year
- Conferences, live surgeries, workshops
- At the foot of the Mont-Blanc, between Chamonix and Geneva

C 190_GB 2022-09

Photo credits : Anthogyr, Godefroy de Maupeou, Getty images . All rights reserved, actual products may vary.

ANTHOGYR 2237, avenue André Lasquin 74 700 Sallanches - France www.anthogyr.co.uk



f 🛎 in 💟