

- Guided surgery

## The FIRST DRILL solution

- + High precision instruments
- + Compatibility with various kinds of implant planning software
- + Compact and progressive kit



Let yourself be guided with total peace of mind.



# CAMPUS



### Training center in implantology

- 1 operating theatre, 1 dental prosthesis laboratory
- 800 dental professionals trained each year
- Conferences, live surgeries, workshops
- At the foot of the Mont-Blanc, between Chamonix and Geneva

C196\_GB - 2022-11

Photo credits : Anthogyr, Getty images, Godefroy de Maupeou.  
All rights reserved, actual products may vary.

STRAUMANN LTD  
3 Pegasus Place  
Crawley  
RH10 9AY  
United Kingdom  
[www.anthogyr.com](http://www.anthogyr.com)



**Anthogyr**  
A Straumann Group Brand

## FIRST DRILL Guided surgery



**Anthogyr**  
A Straumann Group Brand



## PRECISION



Precise guidance thanks to drilling perfectly adjusted to Anthogyr dedicated metallic sockets.

## FREEDOM

• Compatible with different kinds of implant planning software.



- Suitable for all kinds of indications.
- Designed for **Axiom® Bone Level** and **Tissue Level** implants.

*\*Future compatibility.*

**+** Differentiate your practice by integrating digital tools.  
Result: a reassuring explanation, calm patients!

## ERGONOMICS

- Direct drill guiding with no intermediate tools (no spoons, no sliders).
- Compact size surgical guide thanks to the narrow sockets.
- Placement facilitated by a compact and functional kit, consistent with the **Axiom® Multi Level®** solution.



Pointer drill with bit

Ø2.0 mm drill with depth markers

Ø2.0 mm anti-rotation guiding socket



Drill pin with bit



Stabilisation pin for surgical guides (mucosal support)

Pin socket

With the FIRST DRILL guided surgery solution,

- Start your surgery with total peace of mind.
- Ensure the accurate placement of the implant axis and of the implant emergence point as planned, thanks to the guided passage of the initial drills.
- Complete the surgical protocol and implant placement quickly and confidently.

## INTEGRATION IN THE DIGITAL COLLABORATIVE WORKFLOW



**FIRST DRILL** Guided surgery complements Anthogyr's digital solutions and assists you in the implementation of a digital workflow.