

integral



**Guided
surgery.**

ANTHOGYR

A Straumann Group Brand




Unleash the power of digital.

Integral is a fully guided surgery solution developed to address all your **Axiom**® guided cases with the flexibility you need.



Integral is an open solution that allows you to integrate major workflows. Enjoy the predictability, safety and minimally invasive procedures of guided surgery. Choose **Integral** and its distinctive design developed to increase your surgical efficiency with **Axiom**® solutions.



integral

axiom®

ANTHOGRYR 
A Straumann Group Brand

integral

ERGONOMICS AND PERFORMANCE.

KEYLESS SYSTEM


Offering optimal comfort and peace of mind during surgery

- Self-guiding drills and instruments, no wrench, no spoon required
- Design allowing single-handed procedures



Unleash the power of digital.

ACCURACY.




SEAMLESS ROUTE

Intuitive color coding and smart marking for easy identification of the sequence

- Clear organization of instruments
- 1 single kit during surgery


V-SHAPE SLEEVES WITH IMPLANT INDEXATION

- Optimized design for interdental spaces
- Unique placement of the sleeves in the surgical guide.



ACCURATE COMPONENTS ADJUSTMENT FOR HIGH PRECISION GUIDANCE

- Achievement of desired accuracy.
- Markers for appropriate final implant placement in depth and axis.
- Designed to provide comfortable handling.




FLEXIBILITY.

Enjoy the flexibility of an exhaustive range of **4 KITS** to cover main clinical indications with **Axiom®** implants.

- Axiom® Bone Level and Tissue Level implants
- X3, REG and PX profiles
- Ø3.4 mm to Ø4.6 mm
- 6.5 mm to 14 mm lengths

axiom®

♥ OPTIMIZED KIT CONFIGURATION covering 85% of indications




1

Kit Ø3.6 - Ø4.2

Main clinical situations

- Axiom® BL Ø3.4 / 4.0 mm
- Axiom® TL N Ø3.4 / 4.0 mm




2

Kit Ø3.6

Narrow interdental spaces

- Axiom® BL Ø3.4 mm




3

Kit Ø4.2

Intermediate interdental spaces

- Axiom® BL Ø3.4 / 4.0 mm
- Axiom® TL N Ø3.4 / 4.0 mm



4

Kit Ø5.0

Large interdental spaces

- Axiom® BL Ø3.4 / 4.0 / 4.6 mm
- Axiom® TL N Ø3.4 / 4.0 mm
- Axiom® TL R Ø3.4 / 4.0 / 4.6 mm

FREEDOM.

OPEN COMPATIBILITY

- Seamless integration into dental offices and optimized communication with dental labs.
- Freedom for Labside design and printing of surgical guides.

Compatible with major planning software*

coDiagnostiX®



smop
powered by swissmeda

3shape
Implant Studio

STAY TUNED!



See exhaustive compatibilities

*All trademarks are the property of their respective owners.

SMILE IN A BOX® SERVICES

Focus on your passion, we take care of the rest.

- Comprehensive and flexible digital workflow.
- From planning to temporary restoration, through surgical guide.
- Including design and manufacturing.



FOLLOW US.



www.anthogyr.com

ANTHOGYR 
A Straumann Group Brand