

WALK ON A TRUSTED PATH

INTEGRAL

Guided surgery

 **Anthogyr**
A Straumann Group Brand



INTEGRAL

Guided surgery

WALK ON A TRUSTED PATH


INTEGRAL is a fully guided surgery solution developed to address all your guided cases with the flexibility you need.

The **INTEGRAL** range gives you the freedom to match your clinical needs as closely as possible. Smart kits and instruments are designed to enhance the accuracy and ergonomics of your surgical

procedures. **INTEGRAL** is an open solution that allows you to integrate major workflows.

Enjoy the predictability, safety and minimally invasive procedures of guided surgery. Choose **INTEGRAL** and its distinctive design developed to increase your surgical efficiency.






ERGONOMICS AND PERFORMANCE

Keyless system offering optimal comfort and peace of mind during surgery



- Self-guiding drills and instruments, no wrench, no spoon required
- Design allowing single-handed procedures



Follow the path!

Intuitive color coding and smart marking for easy identification of the path

- Clear organization of instruments
- 1 single kit during surgery



Optimized kit configuration covering 85% of indications


FLEXIBILITY

Enjoy the **flexibility** of an exhaustive range of 4 kits, to cover main clinical indications with Axiom® Multi Level® implants.



- Axiom® Bone Level and Tissue Level implants
- X3, REG and PX profiles
- Ø3.4 mm to Ø4.6 mm
- 6.5 mm to 14 mm lengths

Ø 3.6
Ø 4.2




- **Main clinical situations** ♥
- Axiom® BL Ø3.4 / 4.0 mm
- Axiom® TL N Ø3.4 / 4.0 mm

Ø3.6



- **Narrow interdental spaces**
- Axiom® BL Ø3.4 mm

Ø4.2



- **Intermediate interdental spaces**
- Axiom® BL Ø3.4 / 4.0 mm
- Axiom® TL N Ø3.4 / 4.0 mm

Ø5.0



- **Large interdental spaces**
- Axiom® BL Ø3.4 / 4.0 / 4.6 mm
- Axiom® TL N Ø3.4 / 4.0 mm
- Axiom® TL R Ø3.4 / 4.0 / 4.6 mm

FREEDOM

Open compatibility with major planning software

- Seamless integration into dental offices and optimized communication with dental labs
- Freedom for Labside design and printing of surgical guides

Now compatible with* :



New compatibilities upcoming:
STAY TUNED!

*All trademarks are the property of their respective owners.


ACCURACY

V-shape sleeves with implant indexation

- Optimized design for interdental spaces
- Unique placement of the sleeves in the surgical guide

Accurate components adjustment for high precision guidance

- Achievement of desired accuracy
- Markers for appropriate final implant placement in depth and axis
- Designed to provide comfortable handling



Flat sides

Unique setting with constant distance of 9 mm above implant body for a secured protocol



INTEGRAL

Guided surgery

READY FOR

axiomx3[®]

Photos credits: Anthogyr, Unsplash. All rights reserved, actual products may vary.

C222X_EN - 2022-06

STAY TUNED!



Follow us!



www.anthogyr.com

Anthogyr

A Straumann Group Brand