

LABSIDE compatibility list

Our libraries are compatible with CAD software : Exocad®, 3Shape® .

They can be downloaded from www.anthogyr.com section DIGITAL / LABSIDE.

Before any installation, it is important to identify :

- . The CAD software and its version,
- . The component references for Axiom® BL and Axiom® TL
- . The scanning accessories

Before any installation, it is important to uninstall previous versions of libraries that are no longer useful.

02/11/2022

C208 - IndC

Implant Platform		Milling options		Workflow for conventional impression		Workflow for digital impression with IOS (Intra-Oral Scanner)		Exocad® All versions			3Shape® Scanners and CAD software				
Components references	Components type	Indexed	Rotational	Scan-Adapters for lab scanners	Analog for plaster model	Digital transfers for IOS	Analog for printed model	Lab scanner Imetric®	Lab scanner Medit®	Others configurations					
								version 5 and upper	version 2	all opened Exocad® versions	versions 2017 and 2018	version 2019	version 2020	version 2021	version 2022
Name of the associated library files															
OPPM-M-TI-11 OPPM-M-TI-16	PreFace Axiom BL	✓	-	152-27-SAA 152-27-SAA	OPIA100 OPIA100	152-27-DT 152-27-DT	152-27-PA 152-27-PA	Anthogyr_Labside_Imetric_V2.02.02	Anthogyr_Labside_exoCAD_V2.02.02	Anthogyr_Labside_exoCAD_V2.02.02	NA	NA	Anthogyr_Labside_2020_3shape_V2.03.02	Anthogyr_Labside_2021_3shape_V2.03.02	Anthogyr_Labside_2022_3shape_V2.03.02
TPM-M-TI-N11 TPM-M-TI-N16 TPM-M-TI-R11 TPM-M-TI-R16	PreFace Axiom TL	✓	-	156-0X-SAA 156-0X-SAA 156-0X-SAA 156-0X-SAA	TA100-N TA100-N TA100-R TA100-R	156-01-DT 156-01-DT 156-02-DT 156-02-DT	156-01-PA 156-01-PA 156-02-PA 156-02-PA								
OPFLEX413/423/433 OPFLEX513/523/533	FlexiBase® Axiom BL	✓	-	152-27-SAA 152-27-SAA	OPIA100 OPIA100	152-27-DT 152-27-DT	152-27-PA 152-27-PA								
MUNFLEX MUFLEX	FlexiBase® Multi-Unit	-	✓	151-04-SAO 151-03-SAO	MUNA100 MUA100	151-04-DT-MUN 151-03-DT-MU	151-04-PA 151-03-PA								
TFLEX-N TFLEX-R	FlexiBase® Axiom TL	✓	-	156-0X-SAA 156-0X-SAA	TA100-N TA100-R	156-01-DT 156-02-DT	156-01-PA 156-02-PA								
TFLEX-N-P TFLEX-R-P		-	✓	156-0X-SAO 156-0X-SAO	TA100-N TA100-R	156-01-DT-IL 156-02-DT-IL	156-01-PA 156-02-PA								
OPFLEX403 OPFLEX503	FlexiBase® Axiom BL	✓	-	152-27-SAA 152-27-SAA	OPIA100 OPIA100	152-27-DT 152-27-DT	152-27-PA 152-27-PA	Anthogyr_Labside-OPFLEX403_503-SA-DT-PA-Exocad-Imetric	Anthogyr_Labside-OPFLEX403_503-SA-DT-PA-Exocad-Medit	Anthogyr_Labside-OPFLEX403_503-SA-DT-PA-Exocad	Anthogyr_Labside-OPFLEX403_503-SA-DT-PA-3Shape2018	Anthogyr_Labside-OPFLEX403_503-SA-DT-PA-3Shape2019	Anthogyr_Labside-OPFLEX403_503-SA-DT-PA-3Shape2020	NA	NA