

axiom®



Axiom®
Components
Overview



Anthogyr
A Straumann Group Brand

AXIOM® SURGICAL COMPONENTS

		axiomX3®						axiom®PX													
IMPLANT RANGE																					
STERILE																					
Cover screw included																					
OPIM100																					
OPCS100																					
mm																					
length (mm)																					
6.5		-	XT40060	XT46060	XT52060	XT58060	XT64060	-	-	PX46060	PX52060										
8		XT34080	XT40080	XT46080	XT52080	XT58080	XT64080	PX34080	PX40080	PX46080	PX52080										
10		XT34100	XT40100	XT46100	XT52100	XT58100	XT64100	PX34100	PX40100	PX46100	PX52100										
12		XT34120	XT40120	XT46120	XT52120	XT58120	XT64120	PX34120	PX40120	PX46120	PX52120										
14		XT34140	XT40140	XT46140	-	-	-	PX34140	PX40140	PX46140	-										
16		XT34160	XT40160	-	-	-	-	PX34160	PX40160	-	-										
18		XT34180	XT40180	-	-	-	-	PX34180	PX40180	-	-										
Bur / Drill		Round bur (INFB20)				Ø 1.5 Pointer drill (OPPO15150)		Ø 2.0 Lindermann bur (OPR20)													
BONE DENSITY		D1	D2/3	D4	D1	D2/3	D4	D1	D2/3	D4	D1	D2/3	D4	D1	D2/3	D4	D2/3/4	D2/3/4	D2/3/4	D2/3/4	
Initial drill Gauge	Ø 2.0																				
Ø: mm																					
Stepped drill	Ø 2.0 2.4																				
	Ø 2.4 3.0																				
	Ø 3.0 3.6																				
	Ø 3.6 4.2																				
Ø 4.2 4.8																					
Ø 4.8 5.4																					
Ø 5.4 6.0																					
OPTIONAL GAUGES																					
Calibrating film		Axiom® X3 XTLFC_NOT										Calibrating film									
												Axiom® PX PXFC-NOT									

AXIOM® PROSTHETIC COMPONENTS

Mont-Blanc
from the Aiguille du Midi
Chamonix, France

A modern, angular building with large glass windows and a balcony, situated on a steep, rocky mountain slope. The building is surrounded by snow-covered mountain peaks and glaciers under a clear blue sky.