


Straumann Group News Fall Edition 2020

#Care. Innovate. BounceBack.



**DISCOVER INNOVATIONS,
CHANGES, AND TRENDS
IN ESTHETIC DENTISTRY.**

#Excellence in Immediacy

#Edentulous Patients

#Prosthetic Efficiency

#Esthetic Dentistry

#Enablement&Education

#Straumann Group Digital Solutions

#Orthodontics

#Straumann Group Brands

Some products and services may not be available in all countries. Please contact your country sales representative for details.

straumanngroup

#Care. Innovate. BounceBack.

Everybody is aware of the impact Covid-19 has on our private lives, businesses and societies. More than ever, our key focus lies in supporting people now – those who are impacted and those who care for them.

But what's next?

Let's adapt to the new reality – with a mindset of innovation, by doing things differently and even better than before. We want to secure the future of families and the dental economy – practices, labs and health

care organizations – by sympathetically looking after our customers and bouncing back smartly.

Be invited to learn more about our innovative approaches on how to deal with the current situation by employing means that are safe, reasonable and decent.

This is our contribution to the future of dentistry, and your invitation to bounce back with us.



TABLE OF CONTENT



EXCELLENCE IN IMMEDIACY | 8

- Straumann® ImmediaXy™
- BLX Implant System
- VeloDrill™/Modular Cassette



EDENTULOUS PATIENTS | 16

- Edentulous Solutions Overview
- Mini Implants
- Novaloc®
- Pro Arch
- BL Prosthetics
- Straumann® Zygomatic Implant System



PROSTHETIC EFFICIENCY | 30

- Angled Solutions
- CARES® Implant Prosthetics
- Straumann® Prosthetic Selection Guide
- Ceramic Healing Abutments
- nIce®



ESTHETIC DENTISTRY | 42

- Straumann® Emdogain®
- Straumann® Emdogain® FL
- Straumann® XenoGraft/XenoFlex
- Jason® membrane
- PURE Ceramic Implant System
- Straumann Group Implant Health



ENABLEMENT & EDUCATION | 54

- Smart
- Courses
- ITI
- Confidence. Supported by Scientific Evidence



**STRAUMANN GROUP
DIGITAL SOLUTIONS | 58**

- Solutions for Dentists
- Solutions For Labs
- Smile in a Box™
- Virtuo Vivo™
- P series



ORTHODONTICS | 72

- ClearCorrect
- DenToGo™



STRAUMANN GROUP | 82

- Neodent®
- Medentika®
- Anthogyr®
- Career
- Innovations
- WIN



GLOBAL

The global leader in tooth replacement with the vision to become the leading global provider in dentistry.



FOUNDATION

Iconic Swiss engineering united with world-class clinicians and technicians for enhanced dental care.



straumanngroup

United for global
excellence in dentistry.

More than creating smiles – restoring
confidence. We want to be the partner
of choice in esthetic dentistry.



INNOVATION

Pioneering powerhouse
turning ideas into meaningful
and disruptive
innovation.



QUALITY

Committed to delivering
uncompromising quality
through rigorous
standards.



CULTURE

Purpose driven, passionate
high performance organization
with a goal of improving the
dental world.





**NO. 1 IN GLOBAL
IMPLANTOLOGY**

The gold standard in
evidence-based implant and
esthetic dentistry.



**INNOVATIVE
BY NATURE**

Industry benchmark for
disruptive innovation backed
by scientific evidence.



Straumann® – A world-leading brand for confidence in esthetic dentistry.



EDUCATION & SCIENTIFIC NETWORK

Largest global scientific network and total commitment to education in dentistry.



LEGENDARY SWISS ENGINEERING

Swiss premium quality, precision and pioneering innovative technologies.



INTEGRATED DIGITAL PERFORMANCE

Efficient and validated workflows, high-tech dental equipment and premium materials.



TRUSTED SERVICE

Global leader through highest service standards and enabling services focused on customer and patient care.

EXCELLENCE IN



IMMEDIACY

Straumann® ImmediaXy™ Taking immediacy to the next level

With ImmediaXy™, Straumann takes immediate treatments to the next level. The Straumann® ImmediaXy™ ecosystem of solutions covers the whole workflow, from data acquisition to restoration, and was designed to enable clinicians to:

- Predictably and confidently perform immediate treatments, even in challenging indications
- Gain time and cost efficiency
- Accelerate revenue generation

This results in improved quality of life for patients and increased growth and efficiency for the dental practice.



END-TO-END

Straumann® ImmediaXy™ covers and seamlessly connects all stages of the workflow: Data Acquisition, Planning & Design, Surgery and Restoration.



CONFIDENCE

Straumann® ImmediaXy™ capitalizes on Straumann's scientific approach, proven expertise and evidence-based solutions.^{1,2}



EFFICIENCY

Straumann® ImmediaXy™ puts operational simplicity and workflow optimization at the core of its solutions, whether digital, surgical or restorative.

Taking immediacy to the next level.

Straumann® ImmediaXy™ builds on innovation at every stage:

- Seamless and flexible digital workflows: chairside, centralized and everything in-between
- Implants and instruments perfected for immediacy
- Highly effective biomaterials to meet any challenge
- Efficient prosthetic options for all indications



COMPREHENSIVE PARTNERSHIP

As a true partner, Straumann offers service, training and practice growth initiatives in combination with excellent products and solutions.



Data
Acquisition

Planning &
Design

Surgery

Restoration

“Push boundaries predictably and confidently in immediacy.”

Straumann® BLX Implant System

Confidence beyond Immediacy.



BLX combines progressive functional design with our high-performance Roxolid® material and the clinically proven SLActive® surface, designed to provide you with confidence for all clinical situations. Innovations like the VeloDrill™ system, Straumann® Dynamic Bone Management and our Esthetic Ease Concept are aimed at significantly improving surgical and restorative workflows.

Contact your local Straumann representative or visit: www.confidence-in-you.com



DYNAMIC BONE MANAGEMENT

Intelligent implant concept supporting immediate protocols in all bone types.



ESTHETIC EASE CONCEPT

One connection and undercontoured prosthetics to deliver esthetics with ease.

**Straumann® TLX
COMING FALL 2021!**



**Iconic Tissue Level
meets Immediacy.**



**REAL
CONFIDENCE**

Backed by long-term scientific evidence for the Roxolid® material and SLActive® surface technology.



Discover BLX:
www.straumann.com/blx

Straumann® VeloDrill™ System

High performance. Low temperature.

The Straumann® VeloDrill™ System shortens oral surgery chair time, minimizes heat generation and delivers high drilling stability.

Suitable for all surgical protocols, the system lowers setup cost by maximizing synergy between instrumentation for freehand and guided surgery.

For flexibility and convenience, the VeloDrill™ System is stored in the Straumann® Modular Cassette.



FLEXIBLE

Shared instrumentation for freehand and guided surgery.



FAST

Shorter chair-time with a simplified drilling protocol.



LOW TEMPERATURE DRILLING

Avoids overheating surrounding structure.





Straumann® Modular Cassette

A cassette that grows with you.

For dental professionals who are starting to work with Straumann® implant systems, the Straumann® Modular Cassette is a compact yet customizable cassette for all workflows, designed to blend into your clinic's daily workflow, allowing surgeons and assistants to work together efficiently.

Start small and build up your cassette system over time, depending on your needs.



CUSTOMIZE

Assemble your kit according to your needs.



COLLABORATE

Designed to blend into your daily clinical workflow.



EFFICIENT

Streamline your workflow.



reddot design award
winner 2019

EDENTULOUS PATIENTS



EDENTULOUS PATIENTS

Straumann® Edentulous Solutions
Because we are all unique.

Out of millions of people facing dental problems, there are no two cases the same. Having the right solution for each individual patient is a challenge every clinician faces.

With Straumann® Edentulous Solutions you have a whole portfolio designed to meet expectations of different patient groups and treat more patients.

Straumann® Edentulous Solutions

Patient-centric approach.

Out of millions of people facing dental problems, there are no two cases the same. Having the right solution for each individual patient is a challenge every clinician faces.

With Straumann® Edentulous Solutions you have a whole portfolio designed to meet expectations of different patient groups and treat more patients.

STRAUMANN® REMOVABLE EDENTULOUS SOLUTIONS



Straumann® Mini Implant System



Straumann® Novaloc® Retentive System






STRAUMANN® FIXED EDENTULOUS SOLUTIONS



**Straumann® Pro Arch with
BLT and BLX systems**



Straumann® TL and BL implant systems



**Straumann® Zygomatic Implants
combined with BLX**





Straumann® Mini Implants Mini Implant. Maximum trust.

Smiling, laughing and eating naturally is a quality of life every edentulous patient would love to maintain.

Now you can provide edentulous patients presenting reduced horizontal bone availability with a less invasive*, immediate removable fixation of their overdenture.

You and your patients can rely on a trusted brand, premium quality and long-term scientific evidence.

*If GBR can be avoided



SLA® – LONG-TERM SCIENTIFIC EVIDENCE

- High and consistent survival rates between 95.1% and 98.8% documented by different studies after 5- and 10-year follow-up.³⁻⁹
- Very low prevalence of peri-implantitis (1.8%) over the 10-year follow-up period.⁴



ROXOLID® – PROVEN QUALITY

- Higher mechanical strength compared to titanium.¹
- The successful use of Roxolid® has been documented in numerous clinical trials with up to 5-year follow-ups.²



OPTILOC® – DURABILITY AND EFFICIENCY

- Reduced maintenance, excellent resistance to wear, space-saving design.



Contact your local Straumann representative or visit www.straumann.com/mini-implants.



Straumann® Novaloc®

A reliable connection that endures.

When durability and endurance count.

There are moments in life when you want a connection to be extremely reliable. Until now, hybrid denture attachment systems may have reached their limit in challenging implant situations. Now there's Novaloc®.

Straumann® Novaloc® for dental professionals, who want to offer edentulous patients a long-lasting, yet cost efficient removable restoration to give them back quality of life.



HIGH DURABILITY

Combination of ADLC abutment coating and PEEK retention inserts designed to resist wear for long-term performance of the prosthetic fixation and therefore minimized maintenance.



HIGH FLEXIBILITY

Situations with high implant divergences can be restored thanks to the Straumann® Novaloc® matrix design and to the Straumann® Novaloc® angulated abutments.



Discover more: www.straumann.com/novaloc

The Straumann® Novaloc® Retentive System is engineered for low maintenance and for high patient comfort. It offers excellent wear resistance for a durable high attachment performance. The angled abutments allow additional compensation of divergence between implants, providing high prosthetic flexibility.



PATIENT COMFORT

The Straumann® Novaloc® Retentive System is designed for convenient insertion and removal of the prosthesis and for optimal adjustment of the denture retention according to patient's needs.



QUALITY OF LIFE

An implant supported removable denture is a cost-efficient implantable solution that can bring back quality of life to edentulous patients by improving comfort, function, esthetics, speech and self-esteem.





Straumann® Pro Arch Tailored to fit. Designed to last.

Straumann® Pro Arch offers different treatment options for implant-borne fixed full-arch restorations. Help your patients regain lasting quality of life.

Proven and predictable

Results from retrospective study Straumann® Pro Arch 99.7% of 440 arches were loaded at the day of surgery. The overall implant survival rate was 98.1% after a median follow-up time of 1,5 years.¹⁰

TREATMENT OPTIONS

Variety of treatment options for specific indications and different patient needs.^{11,12}

IMMEDIATE LOADING

Straumann® BLT and BLX implants designed for reliable primary stability.

REDUCED INVASIVENESS WITH ROXOLID®

Preserves bone, if bone grafting can be avoided, and allows more treatment options with narrow and short implants.^{13,14}



Now with Smile in a Box™



SLACTIVE® FOR COMPROMISED PATIENTS

Peace of mind even in challenging situations.¹⁵⁻²⁰

RESTORATIVE FLEXIBILITY

Prosthetic portfolio addresses patient's esthetic expectations within their financial resources.

PRACTICE GROWTH

Practice development and patient communication tools support your practice growth.



Discover more: www.straumann.com/proarch

Straumann® BL Prosthetics

Continuous improvements and portfolio completions.

Listening to your feedback has inspired us to take the Screw-retained Abutments (SRA) portfolio to the next level, and to add new options to our BL prosthetics offering.

Users especially value SRAs as part of immediate loading protocols in Straumann® Pro Arch. However, their versatility combined with a broad range of accessories also makes them an attractive solution for bridges and crowns.



Enhanced user experience and peace of mind:

- Straight Screw-retained abutments: New one-piece design
- Angled Screw-retained abutments: New edge-free founded design

New options to optimally match your workflow:

- New SRA protective caps and temporary copings
- Additional BL/BLT impression post for implant level impression, non-engaging, to facilitate the construction of implant-level superstructures



Straumann® Zygomatic Implant System Designed by experts with the patient in mind.

The Straumann® Zygomatic Implant System provides an immediate, predictable fixed solution for edentulous and patients with terminal dentition.^{21,22,24}

Our unique implant designs and straight-forward prosthetic concept provide versatile treatment options with immediate protocols and predictable fixed solutions, even in the most demanding clinical situations. The Straumann® Zygomatic Implant System helps you to tailor the solution to your patients anatomy types and bone deficits.²⁷

The Straumann® Zygomatic Implant System completes the Straumann edentulous portfolio with a premium solution so that you have all the tools to treat full arch indications from one source.



**FROM EXPERTS
FOR EXPERTS**

Bringing experts together to provide a treatment you can trust in.

**IMMEDIATE
LOADING^{21,22,25}**

Provides primary stability and allows immediate loading.

**ANATOMICALLY
BASED DESIGNS**

Two designs with a unique combination of features to respect patients' anatomical structures.²⁶

**STRAIGHT-
FORWARD, VERSATILE
PROSTHETICS²³**

Versatile prosthetic portfolio with a single implant connection reduces complexity; compatible with Straumann® Edentulous Solutions.

PROSTHETIC EFF



EFFICIENCY



Straumann® Prosthetic Solutions
Performance. Partnership. Perfection.

For labs and dental professionals alike, Straumann stands for partnership in prosthetic efficiency. The Straumann® Prosthetic Solutions with integrated digital workflows provide a complete portfolio of proven systems designed to cover all indications, serve patients' needs and support professional preferences.

Straumann® Angled Solutions

Flexibility and precision for a great performance.

What makes a great performance?

Expand your options to leverage the quality of Straumann® Screw-retained restorations. With flexibility and precision from a different angle.





Flexibility sets the stage for a great performance. The Angled Solutions combine angled screw-channels and a best-in-class self-retaining screwdriver. They offer treatment options for a wide variety of indications in the anterior and posterior zone and give the option for either a conventional or digital design.

Precision and attention to detail make the difference. We have refined the angled screw-channel of the restoration to help achieve the best esthetic and functional outcome. No matter the implant, the screw-channel can be tilted up to 25° in all directions. Yet its diameter has been designed to be as small as possible to leave more space for a restoration. The proven Straumann® conical screw connection and insertion torque of 35 Ncm provide strong retentive power.



EASY TO USE

Straight-forward workflow, easy handling.



ESTHETIC

Up to 25° of screw-channel angulation. Small screw-channel diameter.



FLEXIBLE

All indications, platforms, digital and conventional.



ORIGINAL

Leverage the precise fit and the quality of the original Straumann® connection.



Discover Straumann® Variobase® for Crown AS: www.straumann.com/variobase-as

Straumann® CARES® Implant Prosthetics

Ultimate restorative flexibility.

With the CARES® workflow, you have access to a comprehensive portfolio of customized implant prosthetics for multiple brands, through a common digital workflow.

Benefit from decades of research, development and experience that Straumann® and its partner companies, such as Createch, a pioneer in the field of screw retained bars and bridges, have built into their innovative products and solutions.





PERFORMANCE

Reached through best-in-class digital equipment, technologies, materials, services and know-how.



FLEXIBILITY

A comprehensive portfolio of restorative options and a variety of workflows.



SIMPLICITY

The one-stop shop for all your CAD/CAM restorative needs for multiple implant brands.



Straumann® Prosthetic Selection Guide

Make a difference by working in the most efficient way.

To help you master the complexity you have to deal with, Straumann introduces the Prosthetic Selection Guide, a unique digital tool that gives you access to all the components you need for your case.



Download from
Google Play



Download from
App Store



Download as
desktop version

Straumann® Ceramic Healing Abutments

The first step to harmonious soft-tissue healing.

The Straumann® Ceramic Healing Abutments ensure favorable conditions for soft-tissue attachment, thereby supporting a healthy peri-implant environment.

Their well-proven zirconia material also helps surgeons and prosthodontists looking for less plaque attachment and enhanced soft tissue healing from the day of surgery.





**FAVORABLE
SOFT-TISSUE
ATTACHMENT**

With blood circulation similar to that around natural tooth.^{1-2*}



**DESIGNED
FOR HEALTHY
PERI-IMPLANT
ENVIRONMENT**

Less plaque attachment.^{2-5*}



**ESTHETICS FROM
THE DAY OF
SURGERY**

Ceramic for healing phase and final restoration.



EASE OF USE

Aspiration-security and color-coding.

*references on www.straumann.com/ceramic-healing-abutments

Straumann® n!ce®
Fully crystallized glass-ceramic.
n!ce® to meet you.






New high-tech, glass-ceramic solution from Straumann® that delivers fast, esthetic single tooth restorations



Discover more: www.straumann.com/nice

ESTHETIC DENT



ISTRY

Straumann® Esthetic Dentistry
Let your patients enjoy life with a
beautiful smile.

Looking for state-of-the-art products that meet the increasing demand for esthetic solutions? We have the portfolio that will support you in achieving and maintaining your patients' beautiful smiles.

Straumann® Emdogain® – 25 years Mastering periodontal regeneration and oral wound healing.



Our biological solution that induces regeneration of periodontal tissues and significantly enhances the oral wound healing is celebrating 25 years of clinical success.

25 years passed since the introduction of Straumann® Emdogain®. Over time, our biological material became a gold standard in periodontology due to its extensive documentation and solid clinical evidence.



STABILITY OVER LONG-TERM

Documented in over 1000 publications including 10 years follow-up studies.^{1,2}



MORE PATIENT COMFORT

Enhances the wound healing and leads to less pain and less post-operative swelling.³⁻⁸



BETTER ESTHETIC OUTCOME

Provides enhanced esthetic results: keratinized tissue and improved recession coverage.^{9,10}



Straumann® Emdogain® FL Cultivating periodontal regeneration.



For more information on the clinical success visit:
www.straumann.com/emdogain-fl

Differentiate your practice by offering a flapless regenerative periodontal treatment with Straumann® Emdogain® FL.

For the first time, true periodontal regeneration can be achieved in a minimally invasive way: without surgical flap.^{11,12} Based on decades of clinical experience in regenerative flap surgery, periodontal pockets up to 9 mm deep can now be regenerated and eliminated as part of flapless periodontal procedures with Straumann® Emdogain® FL.



MORE EFFECTIVE

Eliminate more pockets.¹³



MORE EFFICIENT

Achieve results comparable to periodontal flap surgery.¹⁴



INCREASE PATIENT COMFORT

Render your procedures milder and more tolerable.¹⁵



Straumann® XenoGraft and XenoFlex Mastering everyday cases.

Following the well-established treatment protocol using demineralized bovine bone granules for the efficient treatment of bone defects.

Straumann® XenoGraft and XenoFlex:

- Are easy to handle
- Have long-term volume stability
- Already successfully applied in over 500,000 cases worldwide



VERSATILITY



VOLUME STABILITY



EASY HANDLING

Jason® membrane Mastering performance.

The heart beats 100,000 times a day,
meaning 4 million times a year

For every beat, the heart requires a similar effort as needed to squeeze a tennis ball with one hand. The pericardium surrounds the heart and needs to constantly allow and resist the acting forces.

Jason® membrane exploits the power of the heart: thanks to its native collagen structure of porcine pericardium, this membrane has set a new standard in bone augmentation procedures. Its excellent biocompatibility¹⁶ and the extended barrier function¹⁷ ensure a reliable performance in all kinds of augmentations, including complex cases and thin biotypes.



TRUST IN NATURE

Over 10 years of clinical success
and scientific publications.



TENSILE STRENGTH

Tear resistance thanks to the
native pericardium structure.



EASY HANDLING

Very thin yet not sticky to
instruments.





Straumann® PURE Ceramic Implant System

The natural and strong solution.
Makes patients love their smiles.





Nothing is more winning than a light-hearted and happy smile. With the Straumann® PURE Ceramic Implant System even very demanding patients can smile with confidence without compromising on esthetics, reliability and the most natural choice of material.

With the Straumann® PURE Ceramic Implant System you can grant the best esthetic, natural and strong treatment to your patients. They will benefit from all the highly esthetic advantages of a natural ceramic implant, ivory-colored like a natural tooth root and not shining through, even in cases of thin gingival biotypes. Yet you can rely on our high-performance zirconia ceramic material, which is even stronger than the gold standard, grade-4 titanium implants.



Discover natural PURE white:
www.straumann.com/pure



OUTSTANDING ESTHETICS

Favorable soft tissue attachment, high-end esthetic restorations.^{18,19}



PROVEN QUALITY

High performance zirconia ceramic, 100% proof tested.²⁰



EXCEPTIONAL SURFACE

ZLA® surface with revolutionary osseointegration features.²¹⁻²³



FLEXIBILITY

Flexible treatment protocols, digital pre-operative implant planning.



INNOVATIVE

A new system that helps you expand your patient pool.



Straumann Group Implant Health Fit for life.



HEALTHY PERI-IMPLANT CONDITIONS²⁵⁻⁴⁸

Ensure predictable, long-term treatment success, benefit from clinically proven solutions.



BROADENED TREATMENT OPTIONS

Offer a broad range of dental treatment solutions and services to prevent, diagnose and treat peri-implant complications.

PREVENTION. DIAGNOSIS. TREATMENT.

Peri-implant complications are inflammatory conditions affecting the soft tissue and bone around dental implants. Recent data²⁴ indicates that one in five patients suffer from peri-implantitis 10 years after implant placement.

Peri-implantitis and its precursor peri-implant mucositis create suffering, discomfort and high costs for patients. If not treated early enough or in the right manner, peri-implantitis can lead to implant loss.

Straumann Group Implant Health is a holistic service offering designed to successfully prevent, diagnose and treat peri-implant complications.



TIME AND COST SAVING

Reduce complexity with a one-stop shop solution and treatment guidance.
Facilitate greater patient acceptance for recalls and home care.



ENABLEMENT &



EDUCATION

Straumann® Enablement and Education
Your success is our success.

Technologically, dentistry is evolving at a rapid pace. At the same time, “soft factors” and marketing activities are becoming more important in gaining new patients. We provide support that goes beyond products and technical skills. Whether you are experienced or just about to start your career as a dental professional: Straumann offers a variety of dedicated programs, as well as practice management and marketing tools. Furthermore, an extensive offering in training and continuous education will ensure that you and your practice stay at the forefront of dentistry.

Straumann® Smart
Immerse to succeed.
Discover the world of implant
dentistry with us.



**CLASSROOM
TRAINING**



**ONLINE
TRAINING**



MENTORING



**SIMPLE PRODUCT
PORTFOLIO**



**LAB
COMMUNICATION**



**PRACTICE
AND PATIENT
MARKETING**

Straumann Courses.

The Straumann Course website provides an overview of local and international courses from Straumann, botiss and ITI. Courses range from all levels of complexity with hands-on sessions and live surgeries.



Check local and international
courses here:
www.straumann.com/courses



For more information visit
www.straumann.com/smart





Your confidence.
Supported by
scientific evidence.



THE ITI – A STRONG COMMITMENT TO EVIDENCE-BASED PRACTICE

For over 40 years, Straumann has benefited from a close and fruitful relationship with the International Team for Implantology (ITI), the largest international organization in implant dentistry. The ITI is a driving force in education and the dissemination of knowledge in this field, focusing on the development of well-documented, evidence-based treatment guidelines and the consolidation of long-term data.

Almost 20,000 members from all over the world participate in ITI education and research activities and exchange their expertise to continuously improve their skills and treatment provided to the patients.

Straumann products are highly reliable and clinically safe, and this can only be achieved through a strong dedication to Science.

The dental community highly respects our scientific outlook and commitment, which is reinforced by the vast body of evidence that our products and solutions offer. Year after year, hundreds of peer-reviewed publications by independent researchers validate their excellent and long-term performance. This is extremely relevant for dentists and patients who choose products that are strongly supported by evidence and provide predictable outcomes.

Since the establishment of Straumann, we have closely collaborated with the leading clinicians, research institutes and universities worldwide to provide the clinical data and substantiation for our products that support their exceptional standards. Our products are continuously documented in different treatment scenarios. We always ascertain we have evidence before launching a product in market.

DIGITAL SOLUTIONS

A black and white photograph of a woman with her arm raised in a crowd, with a blue horizontal bar at the top containing the text 'DIGITAL SOLUTIONS'. The woman is wearing a dark, sleeveless top and has her right arm raised high. The background is dark with some blurred lights, suggesting a concert or event setting.

ONS



Straumann Group Digital Solutions Orchestrating digital efficiency.

Welcome to Digital Solutions, the integrated digital performance offering and brands of the Straumann Group. Efficient and validated workflows, high-tech dental equipment and premium materials.

SOLUTIONS FOR DENTISTS.

The Straumann Group Digital Solutions offering for dentists entails cost-effective, efficient, and easy-to-use hardware for intraoral scanning, printing and in-house milling. Dentists, lab technicians, as well as patients – they all benefit from a digital workflow that delivers time savings and allows for more time in communication and patient care.

Scanning — Planning — Designing — Milling —



3D-Printing

P series



P wash and P cure

Consumables

P Pro Resins



Discover more: www.straumann.com/cares-dentists

SOLUTIONS FOR LABS. PLAYING TOGETHER SEAMLESSLY.

Our efficient, flexible and collaborative workflow solution will make daily life in your lab easier and more predictable. Straumann Group Digital Solutions for dental labs means validated, digital workflows from scan to manufacture. It's a comprehensive package of first class hardware and software dedicated to support you in meeting all the upcoming challenges of digital dentistry.

Scanning — Planning — Designing — Milling



T500



7 series



coDiagnostiX™
CARES® Designing Software

D series and M series





3D-Printing



P series

P wash and P cure

Consumables

P Pro Resins



Discover more: www.straumann.com/cares-labs

CUSTOMIZED PROSTHETIC SOLUTIONS

Our Centralized Milling Centers utilize high-speed cutting (HSC) milling technology to manufacture high-quality restorations in a validated workflow. Our state-of-the-art milling machines combine a mixture of size, speed and agility to deliver predictable and consistent prosthetic quality.





FLEXIBILITY

Extend the indication range for your milling equipment.



TIME

Benefit from a faster and controlled turnaround time.



PEACE OF MIND

Straumann® Guarantee and support.




PROFITABILITY

Maximize your return on investment.



PRECISION

Perfectly matching components.



Smile in a Box™

Focus on your passion – we take care of the rest.

Smile in a Box™ is a flexible treatment planning and manufacturing service to gain efficiency and practice growth. You select the services you want to focus on - we deliver everything else you need for your treatment in one box.



**IMPROVE PATIENT EXPERIENCE
AND TREATMENT ACCEPTANCE**

by faster time to teeth and effective
communication.



GAIN EFFICIENCY IN YOUR PRACTICE

with immediate protocols and flexible workflow outsourcing.



MEET PATIENT EXPECTATIONS

with a predictable guided surgery workflow. Benefit from all-in-one delivery.



GROW YOUR BUSINESS IN A FLEXIBLE WAY

that suits your situation. Establish digital workflows without investment.



Get access: www.smileinabox.com

Virtuo Vivo™

Capture each note.

Accurate, full arch, powder-free scanning in real-life colors is just the beginning of your digital journey. With Virtuo Vivo™ you have access to a full digital ecosystem.

With Virtuo Vivo™ you can quickly and easily acquire digital impression data.







P series

Strike a chord with precision and productivity.

Our P series from Rapid Shape deliver a new dimension in fast and professional 3D printing to allow you to produce top-quality dental products in-house and on-demand. The P series provides new functions to improve quality, flexibility, productivity and efficiency.



P Pro Resins deliver consistent, high-quality products you can rely on, with Straumann® service and after-sales support. Flexible enough to meet the demand for ease-of use in clinical practice and high productivity in dental laboratories. Suited for a broad variety of planned and released indications as well as applications: models, splints, gingiva masks, temporary crowns and bridges, drill guides, casts, denture bases, impression trays.



ORTHODONTIC





S

Straumann® Orthodontics
The next generation of integrated
orthodontics.

We are witnessing a paradigm shift created by changing patient behaviours and groundbreaking innovations in diagnostics, procedures and products. Be it artificial intelligence, digital equipment, software or clear aligners, the line between dentists and orthodontists is becoming increasingly blurred. We understand: both patient and dentist expect integrated, patient-friendly and effective solutions that are accessible and affordable.



**Clear.
Simple.
Friendly.**



clearcorrect

A Straumann Group Brand

In 2006, one doctor went the extra mile for his patients: he founded his own company to make aligners for them. One by one, other dentists signed up to try this new clear aligner solution, and ClearCorrect grew from there.

Today, ClearCorrect has served tens of thousands of dentists all over the world and has delivered literally millions of aligners to their patients. Through it all, we've made it our priority to listen to our dentists, and to put their needs first. That's what's kept us in business for over 12 years, and that's how we know we'll always be here.

To correct malocclusion, a growing number of patients and doctors prefer clear aligners to traditional wires and brackets, due to advantages such as:

- Nearly invisible appliances
- No dietary restrictions
- No obstructions to flossing or brushing
- Controlled single-tooth movement
- Easier remote monitoring

With a low barrier to entry and fast return on investment, clear aligners are an excellent way to expand your practice by offering treatment that patients actively seek.

ClearCorrect aligners:

Precision engineering. Over a decade of experience.

STRONG, COMFORTABLE POLYURETHANE

Supports trimline to achieve
challenging movements

AUTOMATED THERMOFORMING & TRIMMING

Ensures consistent, highly
detailed aligners

STAIN & CRACK RESISTANT

Tested thoroughly

LASER MARKED

To ensure traceability

Four good reasons to choose ClearCorrect



Attractive and accessible

With a strong price/performance ratio and smart training options, ClearCorrect makes aligners accessible to more dentists and patients.




Founded by doctors

Dentist-friendly by nature, ClearCorrect has been listening to dentists and collaborating with them to create solutions they love since 2006.

Proven to optimize retention.

ClearCorrect aligners feature a smooth, unscaloped trim line that extends past the gingival margin. Research shows that aligners trimmed with this technique are significantly more retentive than scalloped aligners, which can reduce the need for attachments. Fewer attachments means the aligners fit the teeth tighter, and that means better esthetics, which means happier patients.

Study results | retention yielded from trim line



Source: "The effect of gingival-margin design on the retention of thermoformed aligners" by Daniel P. Cowley, James Mah, and Brendan O'Toole in the Journal of Clinical Orthodontics: JCO 11/2012; 46(11):697-702.



Quality: born in the USA

Engineered in the USA with care and pride at our Center of Excellence in Texas, ClearCorrect aligners are made from material optimized for stress retention and clarity while resisting stains and cracks.



Top-notch support

Dealing with us is easy and uncomplicated. We offer real, human service from fully committed and highly-trained support specialists who go the extra mile to help dentists succeed.

DenToGoTM Powered By
DENTAL
MONITORING


BECAUSE YOU CAN'T ALWAYS
BE THERE



Straumann® Virtual Clinic > DenToGo™

Making clear aligner therapy more efficient.

Straumann is excited to announce DenToGo™, one of the world's first mobile orthodontics monitoring solutions, which allows for more efficient patient scheduling, increased patient compliance, and personalized messaging. DenToGo™ helps you keep a close eye on your patients' progress, with dynamic wear schedules, and helps improve patient compliance, while allowing dentists to schedule appointments more efficiently.



Efficient Scheduling

DenToGo™ lets you keep a closer eye on your patients' ClearCorrect treatment, so you can schedule appointments more efficiently. Your patients scan their teeth at home, and the DenToGo™ app will help you assess whether they're ready for the next step.

Increased Compliance

The DenToGo™ app helps patients take a more active role in wearing their aligners, monitoring their progress, and working towards their goals.

Personalized Messages

With personalized messages, you can engage even more frequently with your patients without additional chair time, and without sacrificing the unique, human approach your patients want.

STRAUMANN G



GROUP





GLOBAL BRAND

Available in more than 70 countries with more than 1,6 millions implant sold annually.



PROGRESSIVE TREATMENT CONCEPTS

Modern and reliable solutions.



FOUNDED BY A DENTIST TO DENTISTS

A legacy of more than 25 years focused on ease of use.



THERAPY ACCESS AND AFFORDABILITY

Accessibility to proven and affordable solutions.





NEODENT® – Creating new smiles everyday. A leader in immediate treatment.

Greatness is an achievement
NEODENT® GRAND MORSE®

The Neodent® Grand Morse® implant system is the achievement of 25 years of experience in implant dentistry, and shared experiences with many clinicians worldwide. Continuing with a unique purpose to always deliver high quality treatment options that change patients' lives, the Grand Morse® implant system is the Neodent® evolution. Rooted in our philosophy of respecting mechanical and biological principles, this makes it THE implant of choice in dental implant therapy.

Helix GM[®] Unbeatable Versatility.



Screw-retained restoration		Screw or Cement-retained restoration		
Titanium Temporary Abutment	Pro-Peek Abutment	Titanium Base	Titanium Base C	Titanium Base AS
				
Temporary Single or Multiple-unit	Temporary Single-unit	Single-unit		



Grand Esthetic

Deliver immediate natural esthetics.



Grand Stability

Optimized implant design for high initial stability in different bone types and techniques.



Grand Simplicity

Ease of use at its best.



Grand Reliability

Stable and strong foundation designed for long-term success.

Cement-retained restoration

Screw-retained restoration

Overdenture

Screw-retained prosthesis

Titanium Base for Bridge

Titanium Block

CoCr Abutment

Anatomic Abutment

Universal Abutment

Abutment

Angled Mini Conical Abutment

Novaloc

Mini Conical Abutment

Micro Abutment



(AG or Medentika Holder)

(straight and angled)

(straight and angled)

(straight and angled)

Multiple-unit

Single or Multiple-unit

Single-unit

Multiple-unit

Single or Multiple-unit



A SMILE FOR EVERYONE NEODENT® NEOARCH®

Increasing expectations for shortened treatment duration represent a significant challenge for dental professionals especially in patients with anatomical deficiencies. The Neodent® Implant System offers an optimized solution for immediate fixed treatment protocols in edentulous patients even with severe atrophic maxilla. Neodent® NeoArch® allows to significantly improve patient satisfaction and quality of life by immediately restoring function and esthetics.⁽¹⁾



Immediate function resulting in shorter treatment times.



Immediate natural-looking esthetics with versatile restorative options.



Immediate peace of mind thanks to a stable foundation.

BONE RESORPTION



Helix GM®



Helix GM® Long



Zygoma GM™



DIFFERENTIATE YOUR PRACTICE

NEODENT® DIGITAL SOLUTIONS



NEODENT® TITANIUM BASES ANGLED SOLUTION

- For single-unit restoration with angled screw channel;
- Several gingival heights (from 0.8 to 2.5 mm);
- Several abutment diameters (4.0, 4.5 and 5.5 mm);
- Allows maximum angulation of 25° depending on the gingival height and cementable area.



NEODENT® TITANIUM BASES FOR BRIDGES

- For multi-unit restoration;
- Several gingival heights (from 0.8 to 4.5 mm);
- Several abutment diameters (3.5, 4.5 and 5.5 mm);
- Allows maximum divergence of 10° for ø3.5 and 16° for ø4.5 and ø5.5.



NEODENT® TITANIUM BASES FOR CROWN

- For single-unit restoration;
- Several gingival heights (from 0.8 to 4.5 mm);
- Several abutment diameters (3.5, 4.5, 5.5 and 6.5 mm);
- Allows crown customization and maximum angulation of 30°.



NEODENT® TITANIUM BLOCK

- For Single or multi-units restorations;
- For Medentika holder in two different diameters: 11.5 and 15.8 mm;
- For Amann Girrbach holder in one diameter: 12 mm;
- Allows maximum angulation of 30°.



NEODENT® TITANIUM BASES C

- For single-unit restoration;
- For CEREC workflow;
- Several gingival heights (from 0.8 to 5.5 mm);
- Allows crown customization and maximum angulation of 20°.



ONE STEP HYBRID COPING

- For abutment level multi-unit or edentulous restoration;
- For Screw Retained Bar and bridges
- Provides a passive fit;
- Allows maximum divergence of 28°.



NEODENT® UNIVERSAL ABUTMENT

- For single-unit restorations;
- Straight and angled abutments (17° and 30°);
- Several gingival heights (from 0.8 to 5.5 mm for straight abutment and from 1.5 to 3.5 mm for angled abutments);
- Abutment diameters (3.3 and 4.5 mm).



TITANIUM COPING

- Coping to be used between the customized crown and Neodent® abutments;
- For screw-retained single-unit prosthesis;
- Available for GM Abutment and GM Micro Abutment;
- Allows digital workflow.



Neodent® libraries are available for the main CAD softwares:



CARES Visual



3Shape



Dental Wings



Exocad



LEADING MANUFACTURER

With a variety of connections
- compatible with the major
implant systems - one of the
leading supplier of compatible
prosthetics



MADE IN GERMANY

A long tradition of precision
craftsmanship and
industrial production.



PASSION FOR PRECISION

Dedication and ambition to
perfect products and enable
access to supreme technology.



MEDENTiKA[®] – Prosthetics
perfected – implants
engineered.

IPS IMPLANT SYSTEM

MPS MULTI PLATFORM SYSTEM



MICROCONE



- High-precision, implant-abutment connection
- Promotes preservation of the crestal bone* thanks to the innovative (continuous) micro-macro thread*



With Microcone, Medentika® has developed a unique implant system demonstrating excellent esthetic long-term results and optimal functionality.

QUATTROCONE



- Optimized primary stability
- Bone friendly
- Particularly suitable for:
 - Immediate placement
 - Immediate loading
 - QuattroFix concept



Quattrocone is the most innovative implant system from Medentika®. It has been developed by implantologists for implantologists. Quattrocone30 is designed to enhance the range especially for the QuattroFix concept.



* Based on the principle of crestal bone preservation by a micro thread.

** Camlog products are supplied by Camlog Biotechnologies AG, who have no legal or commercial relation with Medentika®.

PROCONE



- Proven implant-abutment connection
- Promotes preservation of the crestal bone*
- High primary stability



With Procone, Medentika® offers a proven implant specifically developed to meet the surgical requirements to enable crestal bone* preservation and provide high primary stability even in demanding situations. The implant connection is compatible with the Camlog Tube-in-Tube™ implant-abutment connection**.

MINICONE




- A slim and permanent implant solution
- Suitable in situations where bone availability is greatly reduced
- ADLC (amorphous diamond-like carbon) coating sets new standards in durability



The Medentika® Minicone implant has been developed for intraosseous insertion in the upper or lower jaw. The Minicone is a one-piece, root-shaped implant with a sand-blasted and acid-etched, rough, residue-free surface, including an integrated prosthetic Optiloc® abutment.



Quattrocone



QUATTROCONE30

Quattrocone30 is specially developed and patented for the QuattroFix treatment concept and all indications with angulated implant placement.

QUATTROCONE

Quattrocone is a completely new implant concept. It is not only suitable for use in the QuattroFix treatment concept combined with immediate function/loading, but also in soft bone with direct implant placement in extraction sockets and in the aesthetically critical region.

QUATTROFIX

Fixed restoration for atrophic ridges allows for a comprehensive treatment plan for edentulous patients, of full-arch immediate restoration, using just two straight and two 30° angulated Quattrocone Implants restored with Multi-unit Abutments.

Prosthetic Components

NEW



MULTI-UNIT SCREW PATRIX MEDENTILOC AND NOVALOC

The new multi-unit screw patrix gives you another option for restoration. It enables you to switch easily from a screw-retained restoration to a removable overdenture construction.



NOVALOC®

Novaloc provides maximum hardness with optimum sliding Characteristics

The surface quality of the ADLC coating (amorphous diamond-like carbon) sets new standards.

The Novaloc matrix system with its well established technology is a prefabricated connector for retaining removable restorations on Novaloc abutments.



MULTI-UNIT ABUTMENT

The Multi-unit abutment supports a wide range of prosthetic restorations. It is ideal for creating patient -specific, individual hybrid restorations or as the basis for a tailored restoration in the esthetic zone.



IPS IMPLANT SYSTEM

MPS MULTI PLATFORM SYSTEM

Titanium base ASC Flex **NEW** Platform switch and GH 2.5 mm



OPTIMIZED EMERGENCE PROFILE

The optimized, even slimmer emergence profile with a gingival height of 2.5mm supports and protects the soft tissue. Available for conical implant-abutment connections



PLATFORM SWITCHING DESIGN

The platform switch design for the Camlog®* and Procone implant-abutment connection is used to support the hard and soft tissues.

With conical implant systems, the platform switch effect is produced by the design of the implant-abutment connection.

Titanium base ASC Flex rotating for bridges

NEW



ROTATING

Without indexing in the
implant-abutment connection

The Titanium base ASC Flex are available for multi-unit bridge restorations in esthetically demanding situations.


The screw channel can be located in the non-esthetic area. The Titanium base ASC Flex rotating can also be used for a bridge restoration with a straight screw channel.





GENUINE EAGERNESS TO PROGRESS

We have a passion for disruptive advancements and dynamic concepts in dental implantology.



HIGH-TECH EXPERTISE

More than 70 years of French engineering and manufacturing know-how in medical devices.



PATIENTS ESTHETICS FIRST

Implantology concepts with prosthetic and integrated digital solutions in mind.



PERFORMANCE SOLUTIONS

A smart and user-friendly product line for all treatment preferences – simple and demanding procedures.



Anthogyr - Performance in dental implantology driven by ingenuity.



Anthogyr designs, manufactures and delivers a complete range of implants and instruments to support dental health professionals in treating patients across the world. Both innovative and reactive, we assist implantologists and labs in their dental challenges.


axiom[®]
MULTI LEVEL[®]

Axiom[®] Multi level[®] is an implant concept that offers a smart range of surgical solutions, complementary philosophies with Bone and Tissue level implants, fully driven by prosthetics outcome for patients.

Axiom[®] Bone Level

Designed for biological integration

With more than 30 years of experience in implantology, Anthogyr launched the Axiom[®] implant system 10 years ago to improve access to implantology by offering innovative and accessible solutions, a greater comfort for practitioners and performance in their everyday practice.



+ Available in
2 body profiles
REG and **PX**

+ Sealed internal conical connection

- Morse taper type connection.
- Alveolar bone is preserved.

+ Platform switching

- Development of a gingival sleeve.
- Soft tissue stabilisation.

+ Subcrestal positioning

- Adapted response to all bone situations.



A SIMPLE RANGE TO COVER ALL INDICATIONS


2 profiles - REG PX
4 diameter, 7 lengths



REG

PX

SINGLE
DIAMETER
PROSTHETIC
CONNECTION




+ Prosthetic flexibility

- Prosthetic components, exhaustive and rational.
- Labside CAD CAM, a currently expanding range.



Axiom[®] Tissue Level

Transgingival implant designed for biological safety



OPTIMAL TISSUE MANAGEMENT

+ Maximised tissue volume



+ Design favorable to soft tissues

+ Safe prosthetic approach at gingival level

+ Epithelium and connective tissue attachment preserved

PROSTHETIC COMFORT

+ Double connection

- **inLink®** innovative connection for multiple screw-retained restorations direct on **Axiom® TL**.



+ Possible evolution from single to plural

A COMPREHENSIVE RANGE TO COVER ALL INDICATIONS

- 2 profiles: REG and PX
- 4 diameters, 6 lengths
- 2 platforms, 3 neck heights



1 SINGLE COMPACT SURGICAL KIT

REG and PX profiles
Both Axiom® BL-TL implants

+ Simplicity

+ Freedom and flexibility

- Possibility to adapt the implant choice during surgery

+ Efficient protocols

+ High performance cutting drills

Intuitive color codes guide



Simeda[®] CAD/CAM prosthesis

High precision and innovative customized solutions for Axiom[®] implants and MPS connections*

AxiN[®] INNOVATION


Simeda[®] CAD/CAM screw-retained tooth for all sectors

AxiN[®] base

Specific integrated retainer

Narrow Angulated Access

Machined AxiN[®] tooth Simeda[®]
Sina Z, T and ML zirconia



Simeda®, specialist in ceramics, also offers premium materials for long term confidence

SIMEDA® CAD/CAM PROSTHESIS OFFER

- + Customized single and multiple unit prostheses
- + Angulated access on Multi-Unit Axiom® and MPS platforms
- + Freedom of design with Simeda® Additive



WHAT MAKES AxIN® UNIQUE?

+ Biological safety

- No glue, no sealing cement.

+ Respect of the natural tooth anatomy

- Channels not visible from the outside.
- Angulation up to 25°.
- Narrow diameter Ø2.0 mm.

+ Temporary angulated components to prepare soft tissues

+ Easy handling

- Removable titanium base.
- Screw integrated in the restoration.

+ Lower maintenance

- Channels outside the fragile areas.
- Lower risk of chipping.

+ Peace of mind

- Full CAD/CAM Simeda® prosthesis.

CAMPUS

Anthogyr Education Center



Develop your surgical and prosthetics skills with Anthogyr Campus located in the Mont Blanc Valley close to Chamonix.

- + A 2-day-Anthogyr journey to discover our roots, DNA, values, people and behind the scene, our production site.
- + Be surprised by the outstanding benefits about the **Axiom**® concept and our **Simeda**® CAD CAM innovations.
- + Let's be practical! Participate to Campus hands-on sessions. In the heart of your activity, benefit from a surgical day with a Live Surgery.



**PARTICIPATE IN AN AMAZING
INTERNATIONAL SEMINAR!**

900

Participants
per year

33

Nationalities

50

International
implantology
instructors

Col du Plan and La Mer de Glace (Chamonix France)

The Straumann Group. We are #ChangeMakers

The Straumann Group is a powerful combination of brands and partners driving innovation in dentistry. United to build trust, engage, collaborate and create opportunities with a customer-focused mindset to serve patients around the world. Our goal has always been to create more than just smiles, but lasting confidence for dentists, labs, patients, employees, partners, dental networks and the global community.

Game changing companies like ours are built on #ChangeMakers with a high performance mindset and values that allow individuals and teams to flourish.



STRAUMANN GROUP FACTS

> 7,000

Employees
worldwide

#1

Implant
manufacturer
worldwide

79

Nationalities work
at Straumann Group

100+

Countries are
served

>1 million

Doctors trained
since our
foundation

1 second

Every second a
Straumann Group
product, service or
component is used

We are #ChangeMakers. A team of players and learners.

At Straumann Group we offer the perfect playing field for #ChangeMakers. We empower our employees to perform with us, make an impact and to question the status quo.

What does changemaking mean for us? It starts small, in our daily business, changing little things and it expands to exploring and entering new grounds that lead to industry-changing solutions.

We learn what is needed, review what we have done, identify what we can improve.

Take action. Initiate new things. Embrace change. This is what we have in common. We call it the **Player-Learner-Mindset** – we cultivate change rapidly to nurture a culture of high performance. We never stop learning and exploring. We create win-win situations to make things better.

We are always looking for talented professionals to join our team. Are you ready to make an impact? Find your opportunities at www.straumann.com/jobs



ASPIRATION

Agile & Independent.
Be energized and embrace change to shape our future.



AMBITION

Innovative & Scientific.
Initiating and harnessing scientific and technological innovation.



ATMOSPHERE

Colorful & Informal.
A truly international company with an informal atmosphere.



Turning ideas into innovation.

To create the dentistry of tomorrow, we strive to be your partner of choice in turning your ideas into innovative products and solutions. Our Straumann Group Idea Portal is open to clients, researchers, clinicians, companies and anyone interested in collaborating with us. We look forward to hearing your ideas within the 8 innovation clusters indicated below. It only takes a few minutes to share your idea with us in confidentiality and potentially start a long-lasting collaboration.




MATERIALS

- Materials for 3D printing
- Biomaterials



SURFACE & COATINGS

- Soft tissue attachment



DESIGN CONCEPTS

- Immediacy
- Soft tissue management



Innovation allows us to dive into new promising workflows, technologies and solutions. Please join us on our journey and work with us to create the dentistry of tomorrow



<https://innovation.straumann.com>



TREATMENT CONCEPTS

- Peri-implantitis
- Immediate treatment workflows



DIGITAL PLANNING & WORKFLOWS

- Patient and treatment simulations
- Artificial intelligence



MANUFACTURING TECHNOLOGIES

- Immediate and cost-effective dentures




BUSINESS INNOVATION

- Customer analytics
- Machine learning

We are a global network of dental professionals:

Engaging women in implantology
Empowering women to follow their ambitions
and advance into leadership roles



Why?

Women account for an increasing percentage of dental school graduates, and only a minority become active in implant dentistry. If this trend continues, patient access to implant therapy and maintenance is at risk. There is also a shortage of female role models: surgeons, speakers, key opinion leaders and practice owners. We aim to engage women at all stages of their career in implantology and inspire them to advance as leaders and role models.

By building a strong global community, supported by the biggest players in Implantology, we will make a change.

How?

1. **Raising awareness** of the situation and featuring female talents.
2. **Networking opportunities** at WIN events, dental congresses and online.
3. **State-of-the-art implant dentistry** via the WIN classroom, summaries of scientific evidence and skill-building workshops.
4. **Mentorship** with WIN Circles, online mentoring, and leadership tips.

Your benefits as a WIN member



Annual WIN events*
Such as WIN Summits and
WIN get togethers.



Scientific Highlights
Updates on the latest findings
in implant dentistry.



Professional development
Skill-building workshops
Visibility and lecturing
opportunities for role models



Free access to the
WIN Classroom
High quality clinical practice
recommendations from the
ITI Academy.



WIN online community
Where you can network and
interact with peers as well as
benefit from mentoring in a
dedicated online group.



20% reduction on first year
ITI membership
Digital membership (not
cumulative with other
discounts).

*Might be free of charge or with additional fee



Local WIN
Additional offerings from your local WIN Community such as
WIN Talks, Implantology Day and more.

Patients worldwide need more women to be engaged in implant dentistry, and female dentists need the strength and encouragement that WIN provides.
Charlotte Stilwell, ITI President Elect

supported by:



Scan to
join WIN
today.



win@straumann.com
www.straumann-group.com/win

REFERENCES

Excellence in Immediacy

1 Chiapasco M, Casentini P, Zaniboni M, Corsi E, Anello T. Titanium-zirconium alloy narrow-diameter implants (Straumann Roxolid[®]) for the rehabilitation of horizontally deficient edentulous ridges: prospective study on 18 consecutive patients. *Clin Oral Implants Res.* 2012 Oct;23(10):1136-41. 2 Quirynen M, Al-Nawas B, Meijer HJA, Razavi A, Reichert TE, Schimmel M, Storelli S, Romeo E. Small-diameter titanium Grade IV & titanium-13zirconium implants in edentulous mandibles: 3-year results from a double-blind, randomized controlled trial. *Clin Oral Implants Res.* 2015 Jul;26(7):831-40.

Edentulous Patients

11 Norm ASTM F67 (states min. tensile strength of annealed titanium). 2 <http://www.straumann.com/science-roxolid.html>. 3 Fischer K, Stenberg T.: Prospective 10-year cohort study based on a randomized controlled trial (RCT) on implant-supported full-arch maxillary prostheses. Part 1: sandblasted and acid-etched implants and mucosal tissue. *Clin Implant Dent Relat Res.* 2012 Dec;14(6):808-15. 4 van Velzen FJ, Ofec R, Schulten EA, Ten Bruggenkate CM. 10-year survival rate and the incidence of peri-implant disease of 374 titanium dental implants with an SLA surface: a prospective cohort study in 177 fully and partially edentulous patients. *Clin Oral Implants Res.* 2015 Oct;26(10):1121-8 5 Cochran DL, Jackson JM, Bernard JP, ten Bruggenkate CM, Buser D, Taylor TD, Weingart D, Schoolfield JD, Jones AA, Oates TW Jr. A 5-year prospective multicenter study of early loaded titanium implants with a sandblasted and acid-etched surface. *Int J Oral Maxillofac Implants.* 2011 Nov-Dec;26(6):1324-32. 6 Cochran D, Oates T, Morton D, Jones A, Buser D, Peters F. Clinical field trial examining an implant with a sand-blasted, acid-etched surface. *J Periodontol.* 2007 Jun;78(6):974-82. 7 Bornstein MM, Schmid B, Belser UC, Lussi A, Buser D. Early loading of non-submerged titanium implants with a sandblasted and acid-etched surface. 5-year results of a prospective study in partially edentulous patients. *Clin Oral Implants Res.* 2005 Dec;16(6):631-8. 8 Rocuzzo M, Aglietta M, Bunino M, Bonino L. Early loading of sandblasted and acid-etched implants: a randomized-controlled double-blind split-mouth study. Five-year results. *Clin Oral Implants Res.* 2008 Feb;19(2):148-52. 9 Derks J, Schaller D, Håkansson J, Wennström JL, Tomasi C, Berglundh T. Effectiveness of Implant Therapy Analyzed in a Swedish Population: Prevalence of Peri-implantitis. *J Dent Res.* 2016 Jan;95(1):43-9. doi 10 Calvo-Guirado JL, Lopez Torres JA, Dard M, Javed F, Perez-Albacete Martinez C, Mate Sanchez de Val JE. Evaluation of extrashort 4-mm implants in mandibular edentulous patients with reduced bone height in comparison with standard implants: a 12 month results. *Clin Oral Impl Res* 27, 2016, 867-874. 11 Piano S, Romeo E, Sbricoli L, Pisoni G, Cea N, Lops D. Simplified procedure for the immediate loading of a complete fixed prosthesis supported by four implants in the maxillary jaw: a 2-year prospective study. *Clin. Oral Impl. Res.* 00, 2015, 1-7. 12 Müller F et al. Small-diameter titanium grade IV and titanium-zirconium implants in edentulous mandibles: five-year results from a double-blind, randomized controlled trial. (2015). *BMC Oral Health.* 2015 Oct 12;15(1):123. 13 Al-Nawas B, Domagała P, Fragola G, Freiburger P, ^ Ortiz-Vigón A, Rousseau P, Tondela J. A prospective non-interventional study to evaluate survival and success of reduced diameter implants made from titanium-zirconium alloy. *J Oral Impl*, submitted. PMID:24666383. 14 Nicolau P, Guerra F, Reis R, Krafft T, Benz K, Jackowski J. 10-year results from a randomized controlled multicenter study with immediately and early loaded SLActive implants in posterior jaws. Presented at 25th Annual Scientific Meeting of the European Association of Osseointegration – 29 Sep – 1 Oct 2016. 15 Paris Heberer S, Kilic S, Hossamo J, Raguse J-D, Nelson K. Rehabilitation of irradiated patients with modified and conventional sandblasted, acid-etched implants: preliminary results of a split-mouth study. *Clin. Oral Impl. Res.* 22, 2011; 546–551. 16 Nelson K, Stricker A, Raguse J-D and Nahles S (2016). Rehabilitation of irradiated patients with chemically modified and conventional SLA implants: a clinical clarification. *J Oral Rehabil.* 43: 871–872. doi:10.1111/joor.12434. 17 Machuca G, Cabrera JJ. “A prospective, case-control clinical study of titanium-zirconium alloy implants with hydrophilic surface in patients with Type 2 diabetes mellitus” Manuscript submitted. 18 Buser D et al. (2017). Development of Implant Stability Quotient values of implants placed with simultaneous sinus floor elevation - results of a prospective study with 109 implants. *Clin Oral Implants Res.* 2017 Jan;28(1):109-115. doi: 10.1111/clr.12768. Epub 2016 Jan 16. A Proline-Rich Peptide Mimic Effects of EMD in Rat Oral Mucosal Incisional Wound Healing. 19 Villa O, Wohlfahrt JC, Mdlá I, Petzold C, Reseland JE, Snead ML, Lyngstadaas SP. *J Periodontol.* 2015 Dec;86(12):1386-95. 20 Eckert SE, Hueler G, Sandler N, Elkattah R, McNeil DC. Immediately Loaded Fixed Full-Arch Implant-Retained Prosthesis: Clinical Analysis When Using a Moderate Insertion Torque. *Int J Oral Maxillofac Implants.* 2019 Jan 31. 21 Tuminelli FJ, Walter LR, Neugarten J, Bedrossian E. Immediate loading of zygomatic implants: A systematic review of implant survival, prosthesis survival and potential complications. *Eur J Oral Implantol.* 2017;10 Suppl 1:79-87. 22 Fatigue Testing Report – ZAGA Implant System. 23 Boyes-Varley JG, Howes DG, Lowrie JF, Blackbeard GA. Surgical Modifications to the Branemark Zygomaticus Protocol in the treatment of the severely resorbed Maxilla: A Clinical Report. *Int. J Oral Maxillofac. Implants* 2003;18:232–237. 24 Wang F, Monje A, Lin GH, Wu Y, Monje F, Wang HL, Davó R. Reliability of four zygomatic implantsupported prostheses for the rehabilitation of the atrophic maxilla: a systematic review. *Int J Oral Maxillofac Implants.* 2015 Mar-Apr;30(2):293-8. 25 Jung RE, Al-Nawas B, Araujo M, Avila-Ortiz G, Barter S, Brodala N, Chappuis V, Chen B, De Souza A, Almeida RF, Fickl S, Finelle G, Ganeles J, Gholami H, Hammerle C, Jensen S, Jokstad A, Katsuyama H, Kleinheinz J, Kunavisarut J, Mardas N, Monje A, Pappaspyridakos P, Payer M, Schiegnitz E, Smeets R, Stefanini M, Ten Bruggenkate C, Vazouras K, Weber HP, Weingart D, Windisch P. Group 1 ITI Consensus Report: The influence of implant length and design and medications on clinical and patient-reported outcomes - *Clin Oral Implants Res.* 2018 Oct;29 Suppl 16:69-77. 26 Aparicio C, López-Piriz R, Albrektsson T. ORIS Criteria of Success for the Zygoma-Related Rehabilitation: The (Revisited) Zygoma Success Code. *Int J Oral Maxillofac Implants* 2020;35:366–378. 27 Aparicio C, Balshi T J, Chow J, David L, Davo R, McLellan G, Nicolopoulos C, Peñarocha M, Pikos M A, Polido W, Zarrinkelk H, Malevez C. The Zygomatic Osteotomy Revisited. ZAGA Consensual Guide for Zygomatic Surgery Decision Taking. *Int J Oral Maxillofac Implants* 2020 (accepted).

Prosthetic Efficiency

1 Kajiwara N, Masaki C, Mukaibo T, Kondo Y, Nakamoto T, Hosokawa R (2015). Soft tissue biological response to zirconia and metal implant abutments compared with natural tooth: microcirculation monitoring as a novel bioindicator. *Implant Dent* 24(1):37-41 2 Degidi M, Artese L, Scarano A, Perrotti V, Gehrke P, Piattelli A. Inflammatory infiltrate, microvessel density, nitric oxide synthase expression, vascular endothelial growth factor expression, and proliferative activity in peri-implant soft tissues around titanium and zirconium

oxide healing caps. *J Periodontol.* 2006 Jan;77(1):73-80

3 Liñares A, Grize L, Muñoz F, Pippenger BE, Dard M, Domken O, Blanco-Carrión J. Histological assessment of hard and soft tissues surrounding a novel ceramic implant: a pilot study in the minipig. *J Clin Periodontol.* 2016 Jun;43(6):538-46.

4 S. Roehling, M. Astasov-Frauenhoffer, I. Hauser-Gerspach, O. Braissant, H. Engelhardt, T. Waltimo, M. Gahlert. In vitro biofilm formation on commercially available machined and micro-roughened titanium and zirconia implant surfaces. *Clin. Oral Impl. Res.* 26 (Suppl. 12), 2015

5 Scarano A, Piattelli M, Caputi S, Favero GA, Piattelli A. Bacterial adhesion on commercially pure titanium and zirconium oxide disks: an in vivo human study. *J Periodontol.* 2004 Feb;75(2):292-6.

Esthetic Dentistry

1 Sculean A, Kiss A, Miliauskaitė A, Schwarz F, Arweiler NB, Hannig M. Ten-year results following treatment of intra-bony defects with enamel matrix proteins and guided tissue regeneration. *J Clin Periodontol.* 2008 Sep;35(9):817-24.

2 McGuire MK, Scheyer ET, Nunn M. Evaluation of human recession defects treated with coronally advanced flaps and either enamel matrix derivative or connective tissue: comparison of clinical parameters at 10 years. *J Periodontol.* 2012 Nov;83(11):1353-62.

3 Jepsen S, Heinz B, Jepsen K, Arjomand M, Hoffmann T, Richter S, Reich E, Sculean A, Gonzales JR, Bödeker RH, Meyle J. A randomized clinical trial comparing enamel matrix derivative and membrane treatment of buccal Class II furcation involvement in mandibular molars. Part I: Study design and results for primary outcomes. *J Periodontol.* 2004 Aug;75(8):1150-60.

4 Tonetti MS, Fourmousis I, Suvan J, Cortellini P, Bragger U, Lang NP. Healing, postoperative morbidity and patient perception of outcomes following regenerative therapy of deep intrabony defects. *J Clin Periodontol.* 2004;31:1092-1098.

5 Ozcelik O, Haytac MC, Seydaoglu G. Immediate post-operative effects of different periodontal treatment modalities on oral health-related quality of life: a randomized clinical trial. *J Clin Periodontol.* 2007 Sep;34(9):788-96.

6 Wennström JL, Lindhe J. Some effects of enamel matrix proteins on wound healing in the dento-gingival region. *J Clin Periodontol.* 2002 Jan;29(1):9-14.

7 Lee JH, Kim DH, Jeong SN. Comparative assessment of anterior maxillary alveolar ridge preservation with and without adjunctive use of enamel matrix derivative: A randomized clinical trial. *Clin Oral Implants Res.* 2020;31(1):1-9. doi:10.1111/clr.13530.

8 Shin SH, Cueva MA, Kerns DG, Hallmon WW, Rivera-Hidalgo F, Nunn ME. A comparative study of root coverage using acellular dermal matrix with and without enamel matrix derivative. *J Periodontol.* 2007;78(3):411-421. doi:10.1902/jop.2007.060170.

9 Pilloni A, Paolantonio M, Camargo PM. Root coverage with a coronally positioned flap used in combination with enamel matrix derivative: 18-month clinical evaluation. *J Periodontol.* 2006 Dec;77(12):2031-9.

10 Górski B, Górska R, Wysiokińska-Miszczuk J, Kaczyński T. Tunnel technique with enamel matrix derivative in addition to subepithelial connective tissue graft compared with connective tissue graft alone for the treatment of multiple gingival recessions: a randomized clinical trial [published online ahead of print, 2020 May 7]. *Clin Oral Investig.* 2020;10.1007/s00784-020-03312-6. doi:10.1007/s00784-020-03312-6.

11 Shujaa Addin A, Akizuki T, Matsuura T, Hoshi S, Ikawa T, Maruyama K, Ono W, Fukuba S, Izumi Y. Histological healing after nonsurgical periodontal treatment with enamel matrix derivatives in canine experimental periodontitis. *Odontology.* 2018 Jul;106(3):289-296.

12 Mellonig JT, Valderrama P, Gregory HJ, Cochran DL. Clinical and histologic evaluation of non-surgical periodontal therapy with enamel matrix derivative: a report of four cases. *J Periodontol.* 2009 Sep;80(9):1534-40.

13 Graziani F, Gennai S, Petrini M, Bettini L, Tonetti M. Enamel matrix derivative stabilizes blood clot and improves clinical healing in deep pockets after flapless periodontal therapy: A Randomized Clinical Trial. *J Clin Periodontol.* 2019 Feb;46(2):231-240.

14 Aimetti M, Ferrarotti F, Mariani GM, Romano F. A novel flapless approach versus minimally invasive surgery in periodontal regeneration with enamel matrix derivative proteins: a 24-month randomized controlled clinical trial. *Clin Oral Investig.* 2017 Jan;21(1):327-337.

15 Wennström JL, Lindhe J. Some effects of enamel matrix proteins on wound healing in the dento-gingival region. *J Clin Periodontol.* 2002 Jan;29(1):9-14.

16 Rothamel et al. 2012, Biocompatibility and Biodegradation of a Native, Porcine Pericardium Membrane. Results from in vitro/in vivo Examination. *J Int J Oral Maxillofac Implants.* 27(1):146-54.

17 Barbeck et al. 2015, Porcine dermis and pericardium-based, non-cross-linked materials induce multinucleated giant cells after their in vivo implantation: A physiological reaction? *J Oral Implantol.* 41(6):e267-81.

18 Kajiwara N, Masaki C, Mukaibo T, Kondo Y, Nakamoto T, Hosokawa R. Soft tissue biological response to zirconia and metal implant abutments compared with natural tooth: microcirculation monitoring as a novel bioindicator. *Implant Dent.* 2015;24(1):37-41.

19 Glauser S, Sailer I, Wohlwend A, Studer S, Schibli M, Schärer P. Experimental zirconia abutments for implant-supported single-tooth restorations in esthetically demanding regions: 4-year results of a prospective clinical study. *Int J Prosthodont.* 2004 May-Jun;17(3):285-90.

20 Institut Straumann AG, Straumann® PURE Ceramic Implants: A 100 % proof test ensures reliable implant strength. 03/14 490.028/en/B/001. www.straumann.com.

21 Bormann KH, Gellrich NC, Kniha H, Dard M, Wieland M, Gahlert M. Biomechanical evaluation of a microstructured zirconia implant by a removal torque comparison with a standard Ti-SLA implant. *Clin Oral Implants Res.* 2012 Oct;23(10):1210-6. doi: 10.1111/j.1600-0501.2011.02291.x. Epub 2011 Nov 14.

22 Gahlert M, Roehling S, Sprecher CM, Kniha H, Milz S, Bormann K. In vivo performance of zirconia and titanium implants: a histomorphometric study in mini pig maxillae. *Clin Oral Implants Res.* 2012 Mar;23(3):281-6. doi: 10.1111/j.1600-0501.2011.02157.x.

23 Gahlert M1, Röhling S, Wieland M, Eichhorn S, Küchenhoff H, Kniha H. A comparison study of the osseointegration of zirconia and titanium dental implants. A biomechanical evaluation in the maxilla of pigs. *Clin Implant Dent Relat Res.* 2010 Dec;12(4):297-305. Doi: 10.1111/j.1708-8208.2009.00168.x.

24 Lee C T et al., *Journal of Dentistry.* 2017;62:1-12

25 Gahlert M et al. *Clin Oral Implants Res.* 2012 Mar;23(3):281-6.

26 Bormann KH et al. *Clin. Oral Impl. Res.* 23, 2012, 1210-1216

27 Gahlert M et al. *Clin Implant Dent Relat Res.* 2010 Dec;12(4):297- 305.

28 Gahlert M et al. *Clin Oral Implants Res.*, Volume 27 (12) – Dec 1, 2016.

29 Janner SFM et al. *Clin Oral Implants Res.* 2018 Mar;29(3):277- 289.

30 Bormann KH et al. *BMC Oral Health.* 2018 Nov 1;18(1):181.

31 Roehling S et al. *J of Periodontol* 2017 Mar;88(3):298-307.

32 Linares A et al. *J Clin Periodontol* 2016; 43: 538-546

33 Buser D et al. *Clin Implant Dent Relat Res.* 2012 Dec;14(6):839- 51.

34 Rocuzzo M et al. *Clin Oral Implants Res.* 2014 Oct;25(10):1105-12

35 van Velzen FJ et al *Clin Oral Implants Res.* 2015 Oct;26(10):1121-8.

36 Fischer K et al. *Clin Implant Dent Relat Res.* 2012 Dec;14(6):808-15

37 J. Derks et al. *J Dent Res.* 2015 Mar;94(3 Suppl):445-515.

38 Bormann KH et al. Publication in preparation

39 Alomrani AN et al. *Clin Oral Implants Res.* 2003 Jun;14(3):251-62

40 Gigandet M et al. *Clin Implant Dent Relat Res.* 2014 Apr;16(2):303-11.

41 Joda T et al. *Int J Oral Maxillofac Implants.* 2015 Nov-Dec;30(6):1272-9.

42 Kajiwara N et al. *Implant Dent* 2015;24(1):37-41

43 Heberer S et al. *Clin. Oral Impl. Res.* 22, 2011; 546-551.


44 C. Nack et al. *Journal of Oral Rehabilitation* 2015 42; 57- 64.

45 Nelson, K et al. *J Oral Rehabil.* 43: 871-872.

46 Tonetti MS et al. *J Clin Periodontol.* 2002 Apr;29(4):317-25.

47 Pilloni A et al. *J Periodontol.* 2006 Dec;77(12):2031-9

48 McGuire MK et al. *J Periodontol.* 2012 Nov;83(11):1353-62.



International Headquarters

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone +41 (0)61 965 11 11
Fax +41 (0)61 965 11 01
www.straumann.com

© Institut Straumann AG, 2020. All rights reserved.
Straumann® and/or other trademarks and logos from
Straumann® mentioned herein are the trademarks or
registered trademarks of Straumann Holding AG and/or
its affiliates.

Straumann News
[490.602/en/A/00](https://www.straumann.com/490.602/en/A/00) 09/20