

Neodent DirectFit Screw

FIT AND FINISH, DIGITALLY

Fit the prosthesis directly on the Multi-Unit Abutments with the Neodent DirectFit Screw, delivering unparalleled accuracy and reliability — the precision essential for digital workflows.



ACCURACY AND RELIABILITY THROUGH ENGINEERED EVIDENCE:

- 0.8 mm thickness designed to promote resistance and reliability over time;
- Precise workflow from planning to milling, available for exocad®, 3Shape®, and CARES® Visual.

FAST-TRACK DIGITAL RESTORATIONS WITH IMPROVED WORKFLOW EFFICIENCY:

- Reduce the number of clinical and lab steps along the way;
- Reduce chairtime and turnaround time when compared to existing workflows with titanium copings linked to the restoration.

SCALE UP YOUR BUSINESS WITHOUT THE COMPLEXITY OF TRADITIONAL METHODS:

- Expand the potential of your financial resources: predictable restorations, with less manual adjustments and no coping debonding failures;
- Explore different raw materials combined with a wide range of Neodent® implant lines.



RESTORATION MATERIALS:



PMMA



Zirconia



SRBB



Discover how Neodent
DirectFit Screw can be used.
Scan the QR code!

PLATFORMS:





AVAILABLE DIGITAL WORKFLOWS:

Chairside or In Lab workflow

GM Mini Conical Abutment			
0.8 mm	1.5 mm	2.5 mm	
115.243	115.244	115.245	
3.5 mm	4.5 mm	5.5 mm	
115.246	115.247	115.248	

GM Exact Mini Conical Abutment 17°/ 30°/ 45°/ 45° slim/ 52°/ 60°

	17°	30°	45°
1.5 mm	115.275	115.278	115.281
2.5 mm	115.276	115.279	115.282
3.5 mm	115.277	115.280	
	45° slim	52°	60°
1.5 mm	115.302	115.300	115.285
2.5 mm	115.303	115.301	115.286





With hybrid workflow (physical model)

GM Mini Conical Abutment			
0.8 mm	1.5 mm	2.5 mm	
115.243	115.244	115.245	
3.5 mm	4.5 mm	5.5 mm	
115.246	115.247	115.248	

GM Exact Mini Conical Abutment 17°/ 30°/ 45°/ 45° slim/ 52°/ 60°

	17°	30°	45°
1.5 mm	115.275	115.278	115.281
2.5 mm	115.276	115.279	115.282
3.5 mm	115.277	115.280	
	45° slim	52°	60°
1.5 mm	115.302	115.300	115.285
2.5 mm	115.303	115.301	115.286



With IOS + Straumann Exact

GM Mini Conical Abutment			
0.8 mm	1.5 mm	2.5 mm	
115.243	115.244	115.245	
3.5 mm	4.5 mm	5.5 mm	
115.246	115.247	115.248	

GM Exact Mini Conical Abutment 17°/ 30°/ 45°/ 45° slim/ 52°/ 60°

	17°	30°	45°
1.5 mm	115.275	115.278	115.281
2.5 mm	115.276	115.279	115.282
3.5 mm	115.277	115.280	
	45° slim	52°	60°
1.5 mm	115.302	115.300	115.285
2.5 mm	115.303	115.301	115.286

