

etkoniDent



Conversion Tables

		etkon® iDent Series		AA				
		MEDENTIKA® Series		E				
		Compatibility		Nobel Biocare®/ NobelReplace® Tapered				
		Implant Connection Diameter (mm)		NP 3.5	RP 4.3	WP 5.0	D 6.0	
		Gingiva height - GH (mm)		0.3	0.3	0.3	0.3	
		Recommended Torque		35 Ncm	35 Ncm	35 Ncm	35 Ncm	
Stock	Ti base	Ti base 3.5	etkon iDent	1AA010.00	1AA010.01	1AA010.02	1AA010.03	
			Medentika	E 1000	E 1010	E 1020	E 1030	
		Ti base 5.5	etkon iDent	1AA012.00	1AA012.01	1AA012.02	1AA012.03	
			Medentika	E 1100	E 1110	E 1120	E 1130	
	Analog	Analog	etkon iDent	1AA021.00	1AA021.01	1AA021.02	1AA021.03	
			Medentika	E 50	E 51	E 52	E 53	
		Repo Analog	etkon iDent	AA0	AA1	AA2	AA3	
			Medentika	E 80	E 81	E 82	E 83	
		Tool for Repo Analog	etkon iDent	AA1.00	AA1.01	AA1.02	AA1.03	
			Medentika	E 80 T	E 81 T	E 82 T	E 83 T	
	Scanbody	Scanbody	etkon iDent	AA0.00	AA0.01	AA0.02	AA0.03	
			Medentika	E 1400	E 1410	E 1420	E 1430	
	Screw	Abutment Screw	etkon iDent	1AA015.00	1AA015.01			
			Medentika	E 60	E 61			
		Lab Screw	etkon iDent	1AA016.00	1AA016.01			
			Medentika	E 60 L	E 61 L			
	Abutment Blank Titanium	Ø 11.5 mm	etkon iDent	1AA017.00	1AA017.01	1AA017.02		
			Medentika	E 9000-R	E 9010-R	E 9020-R		
		Ø 15.8 mm	etkon iDent		1AA018.00	1AA018.01		
			Medentika		E 9010-16R	E 9020-16R		
		Abutment Blank CoCr	Ø 11.5 mm	etkon iDent	1AA019.00			
				Medentika	E 9800-R			
	Ø 15.8 mm		etkon iDent		1AA020.00	1AA020.01		
			Medentika		E 9810-16R	E 9820-16R		
Digital	Customized abutment	Customized Abutment Ti	etkon iDent	2AA023.00	2AA023.01	2AA023.02		
			Medentika	E 9000	E 9010	E 9020		
		Customized Abutment CoCr	etkon iDent	2AA024.00	2AA024.01	2AA024.02		
			Medentika	E 9000	E 9010	E 9020		
Stock	Screwdriver	Handpiece Long	etkon iDent	1ZA027.00				
			Medentika	M 03-5				
		Handpiece Short	etkon iDent	1ZA027.01				
			Medentika	M 09-5				
		Ratchet Long	etkon iDent	1ZA030.00				
			Medentika	M 10-5				
		Ratchet Short	etkon iDent	1ZA030.01				
			Medentika	M 11-5				

Note: Digital Ti-bases and Screws article numbers:
Medentika: add -ET at the end of stock art.nr.
etkon® ident: first digit of art. nr. is 2 instead of 1

		etkon® iDent Series	AB					
		MEDENTIKA® Series	F					
		Compatibility	Nobel Biocare®/ NobelActive®, NobelReplace® Conical					
		Implant Connection Diameter (mm)	D 3.0	NP 3.5	NP 3.5	RP 4.3/5.0	RP 4.3/5.0	
		Gingiva height - GH (mm)	0.65	0.65	1.15	0.65	1.15	
		Recommended Torque	15 Ncm	35 Ncm	35 Ncm	35 Ncm	35 Ncm	
Stock	Ti base	Ti base 3.5	etkon iDent	1AB010.02	1AB010.00	1AB010.03	1AB010.01	1AB010.04
			Medentika	F 1020	F 1000	F 1200	F 1010	F 1210
		Ti base 5.5	etkon iDent	1AB012.02	1AB012.00	1AB012.03	1AB012.01	1AB012.04
			Medentika	F 1120	F 1100	F 1300	F 1110	F 1310
	Analog	Analog	etkon iDent	1AB021.02	1AB021.00		1AB021.01	
			Medentika	F 52	F 50		F 51	
		Repo Analog	etkon iDent	AB2	AB0		AB1	
			Medentika	F 82	F 80		F 81	
		Tool for Repo Analog	etkon iDent	AB1.02	AB1.00		AB1.01	
			Medentika	F 82 T	F 80 T		F 81 T	
	Scanbody	Scanbody	etkon iDent	AB0.02	AB0.00		AB0.01	
			Medentika	F 1420	F 1400		F 1410	
	Screw	Abutment Screw	etkon iDent	1AB015.02	1AB015.00		1AB015.01	
			Medentika	F 66	F 60		F 61	
		Lab Screw	etkon iDent	1AB016.02	1AB016.00		1AB016.01	
			Medentika	F 66 L	F 60 L		F 61 L	
	Abutment Blank Titanium	Ø 11.5 mm	etkon iDent	1AB017.02	1AB017.00		1AB017.01	
			Medentika	F 9020-R	F 9000-R		F 9010-R	
		Ø 15.8 mm	etkon iDent		1AB018.00		1AB018.01	
			Medentika		F 9000-16R		F 9010-16R	
	Abutment Blank CoCr	Ø 11.5 mm	etkon iDent	1AB019.00				
			Medentika	F 9820-R				
		Ø 15.8 mm	etkon iDent		1AB020.00		1AB020.01	
			Medentika		F 9800-16R		F 9810-16R	
Digital	Customized abutment	Customized Abutment Ti	etkon iDent	2AB023.02	2AB023.00		2AB023.01	
			Medentika	F 9020	F 9000		F 9010	
		Customized Abutment CoCr	etkon iDent	2AB024.02	2AB024.00		2AB024.01	
			Medentika	F 9820	F 9800		F 9810	
Stock	Screwdriver	Handpiece Long	etkon iDent	1ZA027.00				
			Medentika	M 03-5				
		Handpiece Short	etkon iDent	1ZA027.01				
			Medentika	M 09-5				
		Ratchet Long	etkon iDent	1ZA030.00				
			Medentika	M 10-5				
		Ratchet Short	etkon iDent	1ZA030.01				
			Medentika	M 11-5				

Note: Digital Ti-bases and Screws article numbers:
Medentika: add -ET at the end of stock art.nr.
etkon® ident: first digit of art. nr. is 2 instead of 1

		etkon® iDent Series		AE		
		MEDENTIKA® Series		K		
		Compatibility		Nobel Biocare® / Brånemark®		
		Implant Connection Diameter (mm)		NP 3.5	RP 4.1	WP 5.1
		Gingiva height - GH (mm)		0.6	0.6	0.6
		Recommended Torque		35 Ncm	35 Ncm	35 Ncm
Stock	Ti base	Ti base 3.5	etkon iDent	1AE010.00	1AE010.01	1AE010.02
			Medentika	K 1000	K 1010	K 1020
		Ti base 5.5	etkon iDent	1AE012.00	1AE012.01	1AE012.02
			Medentika	K 1100	K 1110	K 1120
	Analog	Analog	etkon iDent	1AE021.00	1AE021.01	1AE021.02
			Medentika	K 50	K 51	K 52
		Repo Analog	etkon iDent	AE0	AE1	AE2
			Medentika	K 80	K 81	K 82
		Tool for Repo Analog	etkon iDent	AE1.00	AE1.01	AE1.02
			Medentika	K 80 T	K 81 T	K 82 T
	Scanbody	Scanbody	etkon iDent	AE0.00	AE0.01	AE0.02
			Medentika	K 1400	K 1410	K 1420
	Screw	Abutment Screw	etkon iDent	1AE015.00	1AE015.01	1AE015.02
			Medentika	K 60	K 61	K 62
		Lab Screw	etkon iDent	1AE016.00	1AE016.01	1AE016.02
			Medentika	K 60 L	K 61 L	K 62 L
	Abutment Blank Titanium	Ø 11.5 mm	etkon iDent	1AE017.00	1AE017.01	
			Medentika	K 9000-R	K 9010-R	
		Ø 15.8 mm	etkon iDent			1AE018.00
			Medentika			K 9020-16R
	Abutment Blank CoCr	Ø 11.5 mm	etkon iDent	1AE019.00	1AE019.01	
			Medentika	K 9800-R	K 9810-R	
		Ø 15.8 mm	etkon iDent			1AE020.00
			Medentika			K 9820-16R
Digital	Customized abutment	Customized Abutment Ti	etkon iDent	2AE023.00	2AE023.01	2AE023.02
			Medentika	K 9000	K 9010	K 9020
		Customized Abutment CoCr	etkon iDent	2AE024.00	2AE024.01	2AE024.02
			Medentika	K 9800	K 9810	K 9820
Stock	Screwdriver	Handpiece Long	etkon iDent	1ZA027.00		
			Medentika	M 03-5		
		Handpiece Short	etkon iDent	1ZA027.01		
			Medentika	M 09-5		
		Ratchet Long	etkon iDent	1ZA030.00		
			Medentika	M 10-5		
		Ratchet Short	etkon iDent	1ZA030.01		
			Medentika	M 11-5		

Note: Digital Ti-bases and Screws article numbers:
Medentika: add -ET at the end of stock art.nr.
etkon® ident: first digit of art. nr. is 2 instead of 1

		etkon® iDent Series		AC			
		MEDENTIKA® Series		H			
		Compatibility		BIOMET 3i® / Certain®			
		Implant Connection Diameter (mm)		D 3.4	D 4.1	D 5.0	
		Gingiva height - GH (mm)		0.3	0.3	0.3	
		Recommended Torque		20 Ncm	20 Ncm	20 Ncm	
Stock	Ti base	Ti base 3.5	etkon iDent	1AC010.00	1AC010.01	1AC010.02	
			Medentika	H 1000	H 1010	H 1020	
		Ti base 5.5	etkon iDent	1AC012.00	1AC012.01	1AC012.02	
			Medentika	H 1100	H 1110	H 1120	
	Analog	Analog	etkon iDent	1AC021.00	1AC021.01	1AC021.02	
			Medentika	H 50	H 51	H 52	
		Repo Analog	etkon iDent	AC0	AC1	AC2	
			Medentika	H 80	H 81	H 82	
		Tool for Repo Analog	etkon iDent	AC1.00	AC1.01	AC1.02	
			Medentika	H 80 T	H 81 T	H 82 T	
	Scanbody	Scanbody	etkon iDent	AC0.00	AC0.01		
			Medentika	H 1400	H 1410		
	Screw	Abutment Screw	etkon iDent	1AC015.00			
			Medentika	H 60			
		Lab Screw	etkon iDent	1AC016.00			
			Medentika	H 60 L			
	Abutment Blank Titanium	Ø 11.5 mm	etkon iDent	1AC017.00	1AC017.01		
			Medentika	H 9000-R	H 9010-R		
		Ø 15.8 mm	etkon iDent		1AC018.00	1AC018.01	
			Medentika		H 9010-16R	H 9020-16R	
	Abutment Blank CoCr	Ø 11.5 mm	etkon iDent	1AC019.00			
			Medentika	H 9800-R			
		Ø 15.8 mm	etkon iDent		1AC020.00	1AC020.01	
			Medentika		H 9810-16R	H 9820-16R	
	Digital	Customized abutment	Customized Abutment Ti	etkon iDent	2AC023.00	2AC023.01	2AC023.02
				Medentika	H 9000	H 9010	H 9020
			Customized Abutment CoCr	etkon iDent	2AC024.00	2AC024.01	2AC024.02
				Medentika	H 9800	H 9810	H 9820
Stock	Screwdriver	Handpiece Long	etkon iDent	1ZA025.01			
			Medentika	M 03-2			
		Handpiece Short	etkon iDent	1ZA025.03			
			Medentika	M 09-2			
		Ratchet Long	etkon iDent	1ZA028.01			
			Medentika	M 10-2			
		Ratchet Short	etkon iDent	1ZA028.03			
			Medentika	M 11-2			

Note: Digital Ti-bases and Screws article numbers:
Medentika: add -ET at the end of stock art.nr.
etkon® ident: first digit of art. nr. is 2 instead of 1

		etkon® iDent Series		AD			
		MEDENTIKA® Series		I			
		Compatibility		BIOMET 3i®/ External Hex			
		Implant Connection Diameter (mm)		D 3.4	D 3.4	D 3.4	
		Gingiva height - GH (mm)		0.6	0.6	0.6	
		Recommended Torque		35 Ncm	35 Ncm	35 Ncm	
Stock	Ti base	Ti base 3.5	etkon iDent	1AD010.00	1AD010.01	1AD010.02	
			Medentika	I 1000	I 1010	I 1020	
		Ti base 5.5	etkon iDent	1AD012.00	1AD012.01	1AD012.02	
			Medentika	I 1100	I 1110	I 1120	
	Analog	Analog	etkon iDent	1AD021.00	1AD021.01	1AD021.02	
			Medentika	I 50	I 51	I 52	
		Repo Analog	etkon iDent	AD0	AD1	AD2	
			Medentika	I 80	I 81	I 82	
		Tool for Repo Analog	etkon iDent	AD1.00			
			Medentika	I 80 T			
	Scanbody	Scanbody	etkon iDent	AD0.00	AD0.01		
			Medentika	I 1400	I 1410		
	Screw	Abutment Screw	etkon iDent	1AD015.00			
			Medentika	I 61			
		Lab Screw	etkon iDent	1AD016.00			
			Medentika	I 61 L			
	Abutment Blank Titanium	Ø 11.5 mm	etkon iDent	1AD017.00	1AD017.01		
			Medentika	I 9000-R	I 9010-R		
		Ø 15.8 mm	etkon iDent			1AD018.00	
			Medentika			I 9020-16R	
	Abutment Blank CoCr	Ø 11.5 mm	etkon iDent	1AD019.00	1AD019.01		
			Medentika	I 9800-R	I 9810-R		
		Ø 15.8 mm	etkon iDent			1AD020.00	
			Medentika			I 9820-16R	
	Digital	Customized abutment	Customized Abutment Ti	etkon iDent	2AD023.00	2AD023.01	2AD023.02
				Medentika	I 9000	I 9010	I 9020
			Customized Abutment CoCr	etkon iDent	2AD024.00	2AD024.01	2AD024.02
				Medentika	I 9800	I 9810	I 9820
Stock	Screwdriver	Handpiece Long	etkon iDent	1ZA025.00			
			Medentika	M 03-1			
		Handpiece Short	etkon iDent	1ZA025.02			
			Medentika	M 09-1			
		Ratchet Long	etkon iDent	1ZA028.00			
			Medentika	M 10-1			
		Ratchet Short	etkon iDent	1ZA028.02			
			Medentika	M 11-1			

Note: Digital Ti-bases and Screws article numbers:
Medentika: add -ET at the end of stock art.nr.
etkon® ident: first digit of art. nr. is 2 instead of 1

		etkon® iDent Series		AF			
		MEDENTIKA® Series		R			
		Compatibility		Zimmer® Dental / Tapered Screw-Vent®			
		Implant Connection Diameter (mm)		D 3.5	D 4.5	D 5.7	
		Gingiva height - GH (mm)		0.3	0.3	0.3	
		Recommended Torque		30 Ncm	30 Ncm	30 Ncm	
Stock	Ti base	Ti base 3.5	etkon iDent	1AF010.00	1AF010.01	1AF010.02	
			Medentika	R 1000	R 1010	R 1020	
		Ti base 5.5	etkon iDent	1AF012.00	1AF012.01	1AF012.02	
			Medentika	R 1100	R 1110	R 1120	
	Analog	Analog	etkon iDent	1AF021.00	1AF021.01	1AF021.02	
			Medentika	R 50	R 51	R 52	
		Repo Analog	etkon iDent	AF0	AF1	AF2	
			Medentika	R 80	R 81	R 82	
		Tool for Repo Analog	etkon iDent	AF1.00	AF1.01	AF1.02	
			Medentika	R 80 T	R 81 T	R 82 T	
	Scanbody	Scanbody	etkon iDent	AF0.00	AF0.01	AF0.02	
			Medentika	R 1400	R 1410	R 1420	
	Screw	Abutment Screw	etkon iDent	1AF015.00			
			Medentika	R 60			
		Lab Screw	etkon iDent	1AF016.00			
			Medentika	R 60 L			
	Abutment Blank Titanium	Ø 11.5 mm	etkon iDent	1AF017.00	1AF017.01	1AF018.02	
			Medentika	R 9000-R	R 9010-R		
		Ø 15.8 mm	etkon iDent	1AF018.00	1AF018.01	1AF018.02	
			Medentika	R 9000-16R	R 9010-16R	R 9020-16R	
	Abutment Blank CoCr	Ø 11.5 mm	etkon iDent				
			Medentika				
		Ø 15.8 mm	etkon iDent	1AF020.00	1AF020.01	1AF020.02	
			Medentika	R 9800-16R	R 9810-16R	R 9820-16R	
	Digital	Customized abutment	Customized Abutment Ti	etkon iDent	2AF023.00	2AF023.01	2AF023.02
				Medentika	R 9000	R 9010	R 9020
			Customized Abutment CoCr	etkon iDent	2AF024.00	2AF024.01	2AF024.02
				Medentika	R 9800	R 9810	R 9820
Stock	Screwdriver	Handpiece Long	etkon iDent	1ZA025.00			
			Medentika	M 03-1			
		Handpiece Short	etkon iDent	1ZA025.02			
			Medentika	M 09-1			
		Ratchet Long	etkon iDent	1ZA028.00			
			Medentika	M 10-1			
		Ratchet Short	etkon iDent	1ZA028.02			
			Medentika	M 11-1			

Note: Digital Ti-bases and Screws article numbers:
Medentika: add -ET at the end of stock art.nr.
etkon® ident: first digit of art. nr. is 2 instead of 1

		etkon® iDent Series		AG				
		MEDENTIKA® Series		S				
		Compatibility		DENTSPLY® Implants/ ASTRA TECH OsseoSpeed® TX				
		Implant Connection Diameter (mm)		D 3.0	D 3.5/4.0	D 3.5/4.0	D 4.5/5.0	D 4.5/5.0
		Gingiva height - GH (mm)		0.6	0.6	1.1	0.6	1.1
		Recommended Torque		15 Ncm	25 Ncm	25 Ncm	25 Ncm	25 Ncm
Stock	Ti base	Ti base 3.5	etkon iDent	1AG010.02	1AG010.00	1AG010.03	1AG010.01	1AG010.04
			Medentika	S 1030	S 1000	S 1200	S 1020	S 1220
		Ti base 5.5	etkon iDent	1AG012.02	1AG012.00	1AG012.03	1AG012.01	1AG012.04
			Medentika	S 1130	S 1100	S 1300	S 1120	S 1320
	Analog	Analog	etkon iDent	1AG021.02	1AG021.00		1AG021.01	
			Medentika	S 53	S 50		S 52	
		Repo Ana-log	etkon iDent	AG2	AG0		AG1	
			Medentika	S 83	S 80		S 82	
		Tool for Repo Ana-log	etkon iDent	AG1.02	AG1.00		AG1.01	
			Medentika	S 83 T	S 80 T		S 82 T	
	Scanbody	Scanbody	etkon iDent	AG0.02	AG0.00		AG0.01	
			Medentika	S 1430	S 1400		S 1420	
	Screw	Abutment Screw	etkon iDent	1AG015.02	1AG015.00		1AG015.01	
			Medentika	S 67	S 60		S 61	
		Lab Screw	etkon iDent	1AG016.02	1AG016.00		1AG016.01	
			Medentika	S 67 L	S 60 L		S 61 L	
	Abutment Blank Titanium	Ø 11.5 mm	etkon iDent	1AG017.02	1AG017.00		1AG017.01	
			Medentika	S 9030-R	S 9000-R		S 9020-R	
		Ø 15.8 mm	etkon iDent		1AG018.00		1AG018.01	
			Medentika		S 9000-16R		S 9020-16R	
	Abutment Blank CoCr	Ø 11.5 mm	etkon iDent	1AG019.00				
			Medentika	S 9830-R				
		Ø 15.8 mm	etkon iDent		1AG020.00		1AG020.01	
			Medentika		S 9800-16R		S 9820-16R	
Digital	Customized abutment	Customized Abutment Ti	etkon iDent	2AG023.02	2AG023.00		2AG023.01	
			Medentika	S 9030	S 9000		S 9020	
		Customized Abutment CoCr	etkon iDent	2AG024.02	2AG024.00		2AG024.01	
			Medentika	S 9830	S 9800		S 9820	
Stock	Screwdriver	Handpiece Long	etkon iDent	1ZA025.00				
			Medentika	M 03-1				
		Handpiece Short	etkon iDent	1ZA025.02				
			Medentika	M 09-1				
		Ratchet Long	etkon iDent	1ZA028.00				
			Medentika	M 10-1				
Ratchet Short	etkon iDent	1ZA028.02						
Medentika	M 11-1							

Note: Digital Ti-bases and Screws article numbers:
Medentika: add -ET at the end of stock art.nr.
etkon® ident: first digit of art. nr. is 2 instead of 1

		etkon® iDent Series		AH			
		MEDENTIKA® Series		T			
		Compatibility		DENTSPLY® Implants/ XIVE® S			
		Implant Connection Diameter (mm)		D 3.4	D 3.8	D 4.5	D 5.5
		Gingiva height - GH (mm)		0.3	0.3	0.6	0.6
		Recommended Torque		25 Ncm	25 Ncm	25 Ncm	25 Ncm
Stock	Ti base	Ti base 3.5	etkon iDent	1AH010.00	1AH010.01	1AH010.02	1AH010.03
			Medentika	T 1000	T 1005	T 1010	T 1020
		Ti base 5.5	etkon iDent	1AH012.00	1AH012.01	1AH012.02	1AH012.03
			Medentika	T 1100	T 1105	T 1110	T 1120
	Analog	Analog	etkon iDent	1AH021.00	1AH021.03	1AH021.01	1AH021.02
			Medentika	T 50	T 55	T 51	T 52
		Repo Analog	etkon iDent	AH0	AH3	AH1	AH2
			Medentika	T 80	T 85	T 81	T 82
		Tool for Repo Analog	etkon iDent	AH1.00	AH1.03	AH1.01	AH1.02
			Medentika	T 80 T	T 85 T	T 81 T	T 82 T
	Scanbody	Scanbody	etkon iDent	AH0.00	AH0.01	AH0.02	
			Medentika	T 1400	T 1405	T 1410	
	Screw	Abutment Screw	etkon iDent	1AH015.00			
			Medentika	T 60			
		Lab Screw	etkon iDent	1AH016.00			
			Medentika	T 60 L			
	Abutment Blank Titanium	Ø 11.5 mm	etkon iDent	1AH017.00	1AH017.01	1AH017.02	
			Medentika	T 9000-R	T 9005-R	T 9010-R	
		Ø 15.8 mm	etkon iDent		1AH018.00	1AH018.01	1AH018.02
			Medentika		T 9005-16R	T 9010-16R	T 9020-16R
	Abutment Blank CoCr	Ø 11.5 mm	etkon iDent	1AH019.00			
			Medentika	T 9800-R			
		Ø 15.8 mm	etkon iDent		1AH020.00	1AH020.01	1AH020.02
			Medentika		T 9805-16R	T 9810-16R	T 9820-16R
Digital	Customized abutment	Customized Abutment Ti	etkon iDent	2AH023.00	2AH023.01	2AH023.02	2AH023.03
			Medentika	T 9000	T 9005	T 9010	T 9020
		Customized Abutment CoCr	etkon iDent	2AH024.00	2AH024.01	2AH024.02	2AH024.03
			Medentika	T 9800	T 9805	T 9810	T 9820
Stock	Screwdriver	Handpiece Long	etkon iDent	1ZA025.01			
			Medentika	M 03-2			
		Handpiece Short	etkon iDent	1ZA025.03			
			Medentika	M 09-2			
		Ratchet Long	etkon iDent	1ZA028.01			
			Medentika	M 10-2			
		Ratchet Short	etkon iDent	1ZA028.03			
			Medentika	M 11-2			

Note: Digital Ti-bases and Screws article numbers:
Medentika: add -ET at the end of stock art.nr.
etkon® ident: first digit of art. nr. is 2 instead of 1

Straumann North American Headquarters

Straumann USA, LLC

60 Minuteman Road

Andover, MA 01810

Phone 800/448 8168 (US)

Fax 978/747 2490

www.straumann.us

ifu.straumann.com