

STRAUMANN® **NICE®** ZIRCONIA

Where innovation
meets esthetics.





ESTHETICS

High color accuracy with built-in shade matching and multi-layered feature designed to deliver natural- looking prosthetic outcomes.



VERSATILITY

Providing for your patient cases a wide range of indications from single units to multi-unit bridges and full arches.*



EFFICIENCY

Just one step to finish the restoration with full contour design and an efficient digital manufacturing workflow.

*Product offering may vary from country to country. Not all products and workflows are available in all countries. Indication and material availability dependent upon registrations. Not all shades and/or disc heights offered in all material and/or translucency types. Product status and availability subject to change.

WHERE INNOVATION MEETS ESTHETICS

Zirconia can be applied to a wide range of indications. It lends itself well to an efficient digital manufacturing workflow.

Why Zirconia? Zirconia is the strongest ceramic, prosthetic-milling material available for dental restorations. It combines durability and high strength with esthetics. It provides a high-level of natural looking color accuracy and shade-matching.

Why n!ce®? With the introduction of n!ce Zirconia, Straumann® now offers a material portfolio designed to include a wide range of treatment options. Straumann n!ce Zirconia discs feature user-friendly packaging and labeling to aid in facilitating a digital CAD/CAM manufacturing workflow.

Why Straumann? Straumann stands for premium Swiss quality, precision and cutting-edge innovation that delivers confidence in esthetic dentistry.



FLEXIBILITY

Synergies between digital workflows, Straumann® offers outsourcing options on demand.*



HIGH QUALITY

Premium Swiss quality designed for predictable prosthetic outcomes you can trust.



CONVENIENCE

User-friendly packaging and labeling help to facilitate digital manufacturing workflow.

*Product offering may vary from country to country. Not all products and workflows are available in all countries. Indication and material availability dependent upon registrations. Not all shades and/or disc heights offered in all material and/or translucency types. Product status and availability subject to change.

ESTHETICS



The Straumann® nIce® Zirconia range delivers premium esthetics with innovative color technology and accurate shade matching.

The multi-layered design with integrated color layers ensures a seamless color gradient for natural looking results, without the need of an additional veneering or staining step.

The nIce® Zirconia portfolio includes a wide range of translucencies: LT (Low Translucency), HT (High Translucency), and XT (Extra High Translucency) as well as monochromatic and multi-layered VITA shades for esthetic restorations.

Shade evaluation was performed by dental professionals to aid in developing accurate color-matching and the natural-looking prosthetic outcomes you expect from Straumann®.

- Natural ceramic material
- Heights and shades available to assist with achieving the optimal esthetic result
- Accurate color and shade-matching technology





“You probably know Straumann® as the market leader in dental implants. But did you know that we’re changing the face of dental prosthetics with our n!ce® Zirconia solutions and next generation digital workflow?”

Florian J. Kirsch, Head Digital Business Unit

HIGH QUALITY

Straumann delivers premium Swiss quality, precision and cutting-edge innovation that delivers confidence and predictable prosthetic outcomes you can trust.

Zirconia is the strongest ceramic, prosthetic-milling material available for dental restorations. It combines durability and high strength with natural-looking esthetics..

The n!ce® Zirconia range:

- Esthetic prosthetic outcomes
- Validated process and precision engineered
- Swiss-quality production

CONVENIENCE



Straumann® n!ce® Zirconia discs feature user-friendly packaging and labeling help to facilitate digital manufacturing workflow.

Our new packaging is user-friendly to aid in more clearly identifying material. Our disc labelling concept assists with providing critical information required to mill and sinter, such as shade and enlargement factor, remain clearly visible on the disc itself even if the disc is partly milled. This ensures that the disc can be re-positioned and re-used later to produce subsequent restorations.

Our n!ce Zirconia discs can be used as part of your existing zirconia process without the need to adapt your established workflow. New processing steps or hardware should not be required thanks to the standard disc geometry of the n!ce Zirconia materials.

Furthermore, we use minimal packaging for storage optimization. Our packaging uses recycled, and environmentally friendly FSC* cartons.

*FSC: The Forest Stewardship Council is an international non-profit organization with a mission to create and sustain responsible forestry.




Just one step to finish the restoration with full contour design to limit post processing and facilitate an efficient digital manufacturing workflow.

EFFICIENCY

With our nlce® material range we are harmonizing and simplifying our milling materials portfolio for centralized and in-house milling for enhanced efficiency.

With Straumann® digital solutions you have the potential to increase productivity and profitability from planning to production and finishing. Improve your efficiency and shorten turnaround times with a streamlined and reliable workflow.



Monochromatic discs available in various shades.



Multi-layered discs for esthetic restorations.



FLEXIBILITY

Straumann® offers a range of outsourcing options with synergies between digital workflows.* Manufacture in-house or gain access to milled restorations from our centralized service. The choice is yours.

Straumann in-house options	SCAN	PLAN & DESIGN	MILL	FINISH
	SCAN & SHAPE		CENTRALIZED MILLING	FINISH

From our centralized milling services we can manufacture the restoration according to your design. Or we can deliver the design with our Scan & Shape service, a cost-effective option to in-house milling without compromising on quality or esthetics. Straumann delivers premium Swiss quality restorations.

Or benefit from synergy with your inLab solution. The complete portfolio of n!ce® Zirconia discs, shades and translucency levels are available for use with in-house production process*. Enjoy the flexibility to select the workflow that works for you.

*Product offering may vary from country to country. Not all products and workflows are available in all countries. Indication and material availability dependent upon registrations. Not all shades and/or disc heights offered in all material and/or translucency types. Product status and availability subject to change.





VERSATILITY

Premium treatment options

With the introduction of nIce® Zirconia, Straumann® now offers a material portfolio designed to include a wide range of treatment options. With our portfolio, nIce Zirconia provides for your patient cases indications from a single crown to multi-unit bridges and full arches.

Shade matching

nIce Zirconia is available in up to 16 shades corresponding to the VITA tooth shade guide and bleach shades to facilitate post processing to enable you to deliver the final restoration in the desired shade.*

- Disc heights range from 12 mm to 30 mm
- Restoration design options from framework to full contour monolithic restorations
- Processing options: glazing, staining, layering/veneering

An offering range for esthetic restorations

	APPLICATION	
Low Translucency	Framework	Up to full arch
High Translucency	Full contour	Up to full arch
Extra High Translucency	Full contour	Up to 3-unit bridges

*nIce® ZirconiaLT is available in white and 5 mono shades (L1, L2, M1, M2, and I1). nIce® ZirconiaHT is available in all 16 VITA shades (A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3 and D4) and 3 Bleach shades (BL1 for mono and multi, BL2 for multi, and BL3 for multi). nIce® ZirconiaXT is available in 6 VITA shades (A1, A2, A3, A3.5, A4, and B1) and 2 Bleach shades (BL1 for multi and BL2 for multi).

THE NEXT GENERATION IN DIGITAL INNOVATION. POWERED BY STRAUMANN®.

Innovation drives productivity. And productivity delivers growth. The next generation of end-to-end digital manufacturing solutions is changing the face of dentistry.

With Straumann® digital solutions you have the potential to increase productivity and profitability from planning to production and finishing.

SCAN

Discover our broad spectrum of intraoral and lab scanners. Choose the digital solution that best fits your needs from plug & play simplicity to empower innovation that redefines what is possible. The perfect starting point in your digital workflow, scanning technology empowers you to obtain and transmit accurate, digital impressions seamlessly.




Virtuo Vivo™ and 3Shape TRIOS®

PLAN & DESIGN

Connectivity, piece of mind, and great results are a reality with coDiagnostiX® and CARES® Visual.


- CoDiagnostiX: State-of-the-art dental treatment planning software for predictable results.
- CARES Visual: User-friendly software with a complete and comprehensive dental interface, giving you access to a variety of specific services and design solutions that can enhance your daily practice experience and deliver esthetic solution for your patient.



coDiagnostiX and CARES Visual

MILL & SINTER

With multiple workflows options available, you have the flexibility to choose which option suites your needs. Whether that is implementing an in-house milling and sintering workflow or outsourcing to Straumann's centralized milling facility during peaks in work or complex cases.



Roland Mills



n!ce® Zirconia discs

RESTORE

Patient receives a final zirconia restoration created from n!ce Zirconia.



n!ce Zirconia

STRAUMANN® N!CE® ZIRCONIA – WHERE INNOVATION MEETS ESTHETICS.



- A premium range of treatment options
- Esthetic prosthetic outcomes
- User-friendly packaging and labelling solution
- Digital manufacturing workflow



Straumann North American Headquarters

Straumann USA, LLC
60 Minuteman Road
Andover, MA 01810
Phone 800/448 8168 (US) · 800/363 4024 (CA)
Fax 978/747 2490
www.straumann.us · www.straumann.ca

ifu.straumann.com
ifu.dentalwings.com

© Straumann USA, LLC 2023. All rights reserved.
Straumann® and/or other trademarks and logos from Straumann® that are mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates. All rights reserved.

