



EFFORTLESS **CHAIRSIDE** EXCELLENCE

**DENTAL 3D PRINTING HAS NEVER
BEEN THIS SMART AND SIMPLE**


**MEET STRAUMANN SIGNATURE MIDAS
BY SPRINTRAY**

The world's first digital restoration press and the simplest and smartest method of 3D printing dental restorations.

Midas is an all-new solution for fabricating restorations chairside. Patented Digital Press Stereolithography combines the tank, resin, and platform into a single capsule, allowing it to print with materials up to 10x more viscous than conventional 3D printers.

Up to Three Capsules

Midas accommodates up to three capsules at once. Print a range of dental restoratives in 10 minutes total.



INTEGRATED SEAMLESSLY


Straumann Signature Midas is seamlessly integrated with **Straumann SIRIOS™** and open to any intraoral scanner, go from scan to smile in one efficient, intuitive workflow via **Straumann AXS™** – while you focus on what matters most: **your patients.**



INTEGRATED WITH AXS

Designed to make your workflow flow.

Upload your scans to AXS and send your Smart Design to your printer for seamless manufacturing.



End-to-end integration via Straumann AXS™



VALIDATED MATERIALS

Features SprintRay Ceramic Crown, cleared for use in 3D printing permanent single tooth restorations.

THE FUTURE OF DENTAL RESTORATIONS IS PRINTED



Scan to Learn More

Discover how Midas is set to transform chairside restorations. Scan the QR code or visit: straumann.us/chairside



Straumann North American Headquarters

Straumann USA, LLC
60 Minuteman Road
Andover, MA 01810
Phone 800/448 8168 (US)
800/363 4024 (CA)

www.straumann.us
www.straumann.ca

2710 Media Center
Drive, Suite 100
Los Angeles, CA 90065

www.sprintray.com
info@sprintray.com

© Straumann USA, LLC 2025. All rights reserved. Straumann® and/or other trademarks and logos from Straumann® that are mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates. All rights reserved.