

Straumann® nice® Zirconia
Where innovation
meets esthetics.

COMING
SOON..!



A wide selection of zirconia for all your patient cases.



n!ce® Zirconia ^{LT} Low Translucency ^{1,2}					
	<ul style="list-style-type: none"> White and mono shades in a low translucency disc For final restorations For framework and up to full arch 				
	↑ Translucency	Shades	Disc Diameter	Disc Heights	↓ Flexural Strength ³
	Low	White, L1, L2, M1, M2, I1	98.5 mm	12, 14, 16, 18, 20 mm	≥ 1,200 MPa

n!ce Zirconia ^{HT} High Translucency ^{1,2}					
	<ul style="list-style-type: none"> White, mono, multilayered, and bleach shades in a high translucency disc For final restorations For full contour and up to full arch 				
	↑ Translucency	Shades	Disc Diameter	Disc Heights	↓ Flexural Strength ³
	High	White, 16 VITA Classical shades and 2 Bleach shades ⁴	98.5 mm	12, 14, 16, 18, 20, 25, 30 mm	≥ 1,000 MPa

n!ce Zirconia ^{XT} Extra Translucency ^{1,2}					
	<ul style="list-style-type: none"> White, multilayered and bleach shades in an extra high translucency disc For final restorations For full contour and up to 3-unit restorations 				
	↑ Translucency	Shades	Disc Diameter	Disc Heights	↓ Flexural Strength ³
	Extra High	White, 16 VITA Classical shades and 2 Bleach shades ⁴	98.5 mm	12, 14, 16, 18, 20, 25mm	≥ 500 MPa

¹ Product offering may vary from country to country. Not all products and workflows are available in all countries. Indication and material availability dependent upon registrations. Not all shades and/or disc heights offered in all material and/or translucency type.

² Not intended for use with a Variobase

³ Converted three-point bending

⁴ n!ce Zirconia is available in all 16 VITA shades (A1, A2, A3, A3.5, A4, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3 and D4) and 2 Bleach shades (BL1 and BL2)

Straumann North American Headquarters

Straumann USA, LLC
 60 Minuteman Road
 Andover, MA 01810
 Phone 800/448 8168 (US) • 800/363 4024 (CA)
 Fax 978/747 2490
 www.straumann.us • www.straumann.ca

ifu.straumann.com