

More than in-house milling.

Adding value to your lab.

With Straumann® Pre-milled Abutment Blanks you can now produce original Straumann® one-piece customized titanium abutments with your in-house milling equipment. Boost your efficiency by using your in-house milling machine and, additionally, you benefit from the original Straumann implant-abutment connection.

SIMPLY EFFICIENT

- Extend the indication range for your milling equipment
- Maximize your return on investment
- Benefit from a faster and controlled turn-around time



JUSTIFIED TRUST WITH THE ORIGINAL IMPLANT-ABUTMENT CONNECTION

- Perfectly matching components
- Excellent service and support
- Straumann Guarantee*

*For applicable terms and conditions of the Straumann Guarantee (including the guarantee territory), please refer to the Straumann Guarantee brochure (152.360/en)

Your way to your in-house milled original Straumann® one-piece titanium abutment

Milling prerequisites

• Straumann® Pre-milled Abutment Blanks are compatible with Medentika abutment blank holders which are available for a steadily growing number of milling machines:

- Dental Concept DC5®
- Imes-icore®
- Datron D5®
- VHF
- Wissner Gamma 202®
- Röders RXD®
- MB Maschinen Cobra Mill®

Available
June, 2015

Coming soon

Please contact your milling machine manufacturer for further information on compatibility and integration.

Software design

- Download the Straumann® Pre-milled Abutment Blanks Implant Kit for Dental Wings, 3shape and Exocad from the Straumann Website.
- Design your one-piece titanium abutment according to your requirements.

Ordering information

The Pre-milled Abutment Blanks are available in two blank sizes for the following implant platforms:

		Art.-Nr.	Article	Quantity	Price
Bone Level		022.0028	NC Pre-milled Abutment Blank for Medentika Holder Ø 11.5		
		022.0029	NC Pre-milled Abutment Blank for Medentika Holder Ø 15.8		
		022.0030	RC Pre-milled Abutment Blank for Medentika Holder Ø 11.5		
		022.0031	RC Pre-milled Abutment Blank for Medentika Holder Ø 15.8		
Tissue Level		022.0034	RN Pre-milled Abutment Blank for Medentika Holder Ø 11.5		
		022.0035	RN Pre-milled Abutment Blank for Medentika Holder Ø 15.8		
		022.0036	WN Pre-milled Abutment Blank for Medentika Holder Ø 11.5		
		022.0037	WN Pre-milled Abutment Blank for Medentika Holder Ø 15.8		