

STRAUMANN® CARES® DIGITAL SOLUTIONS

CARES® Implant Prosthetics

Ultimate restorative flexibility.



Customized restorations of exceptional quality from one trusted partner.



PERFORMANCE

Reached through best-in-class digital equipment, technologies, materials, services and knowhow.




FLEXIBILITY

A comprehensive portfolio of restorative options and a variety of workflows.



SIMPLICITY

The one-stop-shop for all your CAD/CAM restorative needs for multiple implant brands.



With the CARES® workflow, you have access to a comprehensive portfolio of customized implant prosthetics for multiple brands, through a common digital workflow.

CARES® Implant Prosthetics combines carefully selected dental equipment with the latest digital technology, premium materials, customized services into a best in class portfolio of tailor-made prosthetic solutions, that is continuously expanding.

CARES® Implant Prosthetics enables you to increase your competitiveness through access to a complete portfolio of technologies, tools and materials. You can leverage the full extent of your individual digital set-up and workflow preferences, enabling you to maximize your productivity and efficiency with the peace of mind that you are delivering restorations of exceptional quality.

When deciding for CARES® Implant Prosthetics, you have the freedom of choice in every situation. You can produce your restorations in-house or use our centralized milling services enabling you to maximize your productivity and efficiency.

A portfolio to outperform your expectations.

CARES® Implant Prosthetics deliver high performance through a unique combination of best-in-class digital equipment, technologies, materials, and services.

With CARES® Implant Prosthetics you benefit from decades of research, development and experience that Straumann® and its partner companies, such as Createch, a pioneer in the field of screw retained bars and bridges, have built into their innovative products and solutions.




CARES® IMPLANT PROSTHETICS WORKFLOW STRICTLY CONTROLLED MANUFACTURING PROCESS

Order → Production → Quality Control → Delivery



From scan to design.



State-of-the-art milling machines utilizing high-speed cutting (HSC) milling technology.



Thorough finishing and quality control of all parts.



A seamless workflow until delivery to your office door.

Experience the combination of excellence in design engineering together with unrivalled precision manufacturing from our CARES® centralized production facilities. It all starts with the careful selection of raw materials, each component is manufactured under strictly controlled conditions, including a meticulous quality control process to guarantee the highest standards. We are determined in our commitment to providing customers and patients with reliability, quality and peace of mind.

Lean on a comprehensive CAD/CAM portfolio.

No matter if you restore a straightforward single unit or a complex full arch rehabilitation, CARES® Implant Prosthetics is always the right choice when it comes to best-in-class tailor made restorative solutions for virtually any indication.

FOR IN HOUSE MILLING



CARES® X-Stream™

Maximize efficiency

Minimize turnaround time and costs by having your abutment and the corresponding, excellent fitting coping or crowns manufactured and delivered together. One scan, one design, one delivery from one source.



For your crown and bridge prosthetics, you have the flexibility to select from a set of cement- and screw- retained solutions. For your edentulous cases you can opt from a broad range of removable or fixed solutions.

CENTRALIZED MILLING SERVICES



Angled screw-channel solutions

Uncompromised esthetics and functionality

Ensures that for your screw-retained restorations, the position of the screw exit hole is always in the perfect occlusal position, even if implants are placed in an unfavorable position.*







*not available for all indications in all markets

Flexible workflows for unlimited freedom.


Access a comprehensive portfolio of customized implant prosthetics through your preferred workflow. Whatever your level of equipment, CARES® will assist you in enhancing your individual lab set-up.

Be it a scan, a design or a complete production, you can let Straumann® do the whole prosthetic work for you or only transfer single steps to us and do the remaining work yourself. If you prefer producing your restorations on your own, we offer all the components you need to create best-in-class CAD/CAM implant prosthetics. The choice is yours.

	Lab set-up	Scan	Design	Production
				
WITH SCANNER	Straumann® CARES® Scanner	Straumann® CARES® Lab Scanner	Straumann® CARES® Visual Software	<ul style="list-style-type: none"> • CARES® centralized production <ul style="list-style-type: none"> – Straumann® – Createch • Inlab production <ul style="list-style-type: none"> – CARES® inlab milling equipment²
	3 rd party scanner with connectivity (DWOS, 3shape, exocad) ¹	3 rd party scanner system <ul style="list-style-type: none"> • Dental Wings • 3shape • exocad 	3 rd party CAD/CAM software plus <ul style="list-style-type: none"> • Straumann® CARES® App for DWOS • Straumann® DME for 3shape • Straumann® libraries for exocad 	<ul style="list-style-type: none"> • CARES® centralized production <ul style="list-style-type: none"> – Straumann® – Createch • Inlab production <ul style="list-style-type: none"> – CARES® inlab milling equipment² – 3rd party inlab milling equipment
WITHOUT SCANNER	No equipment required (file upload, model or wax-up)	Straumann® CARES® Scan & Shape Service	Straumann® CARES® Scan & Shape Service	<ul style="list-style-type: none"> • CARES® centralized production <ul style="list-style-type: none"> – Straumann® – Createch

¹ Dental Wings, 3-shape, exocad with respective Straumann libraries / Createch libraries (for details see technical brochure 490.318/en)

² Please contact your respective Straumann partner on how to connect Straumann inlab equipment to your individual lab setup



Everything you need from one
trusted partner to make life easier.

Today's dental lab reality is that you have to deliver prosthetic solutions for a variety of implant systems produced from different companies.

With CARES® Implant Prosthetics, you benefit from a seamless ordering process. It is your one-stop-shop for all your CAD/CAM restorative needs and customer services across multiple implant systems and brands, all this from a single source you can rely on.

- One order
- One shipment
- One invoice

CARES® Implant Prosthetics covers all Straumann Group brands: Straumann®, Neodent® and MEDENTiKA®. In addition you get access to compatible prosthetics for the major implant brands³.



³ For more details see technical brochure 490.318/en (CARES® abutments, Pre-Milled Abutment Blanks, Screw-Retained Bars and Bridges)

Quality and service for your peace of mind.

When you choose CARES® Implant Prosthetics you can have absolute trust that you will receive a product you can rely on and manufactured to the highest standards. In fact, we are so confident of the quality of our CARES® Implant Prosthetics that we offer our customers unprecedented guarantees on our components, and in some cases even a Lifetime Plus guarantee*.

In addition, you benefit from our excellence in service and support through our highly skilled personnel:

- Hotline support from our call centers around the world
- Online technical support for your digital cases
- Training and education through our support teams



*For more details please refer to the respective guarantee brochures:
Straumann®, Neodent®, MEDENTIKA®, Createch



Contact your local Straumann representative now or visit
www.straumann.com.



490.317/en/B/00 12/18

International Headquarters

Institut Straumann AG

Peter Merian-Weg 12

CH-4002 Basel, Switzerland

Phone +41 (0)61 965 11 11

Fax +41 (0)61 965 11 01

www.straumann.com

© Institut Straumann AG, 2018. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

