



















CARES® X-Stream™ restorative options

MATERIAL OVERVIEW

- 1 Bridges and bars are available with up to 10 implant abutment connections.
- 2 Molar bridge excluded
- 3 Not available for Bars
- 4 Only available for Fixed bars
- 5 Not available for NNC and WN
- 6 Not available for NNC
- 7 Restorations only available without screw channel hole
- 8 Restorations only available with screw channel hole

		Tooth replacement options Full-contour or frameworks for single tooth, bridges and bars													
		Zirconia							Metal		Polymer				
		nIce®	IPS e.max® CAD	zerion® GI	zerion LT® 5	zerion® ML	3M™ ESPE™ Lava™ Plus Zirconia ⁴	3M™ Lava™ Esthetic	zerion® UTML ^{2,3}	coron®	ticon®	polycon® ae ³	VITA® CAD-Temp	JUVORA™ PEEK	
															
PURE Ceramic	Implant Monotype		Single	Single		Single	Single	Single	Single		Single	Single	Single	Single	
CARES® Abutments	CARES® Abutment, Zirconium dioxide ^{5,6}		Single	Single	Full	Full	Full	Full	≤ 3 U	≤ 3 U			Full	Full	Full
	CARES® Abutment, Titanium ^{6,7}		Single	Single	Full	Full	Full	Full	≤ 3 U	≤ 3 U	Full	Full	Full	Full	Full
	CARES® Abutment, CoCr ^{6,7}		Single	Single	Full	Full	Full	Full	≤ 3 U	≤ 3 U	Full	Full	Full	Full	Full
Variobase®	Variobase® for Crown		Single	Single	Single	Single	Single	Single	Single	Single	Single		Single	Single	Single
	SC Variobase® for Crown				Single	Single	Single	Single	Single	Single	Single		Single	Single	Single
	Variobase® for Crown AS					Single		Single			Single		Single		
	Variobase® for Bridge/Bar Cylindrical ⁸				Full	Full	Full	Full	≤ 3 U	≤ 3 U	Full	Full	Full	Full	Full
	Variobase® for Bridge/Bar Cylindrical Coping for Screw-retained Abutments ⁸				Full	Full	Full	Full	≤ 3 U	≤ 3 U	Full	Full	Full	Full	Full
Stock Abutments	BL Cementable Abutments		Single	Single	≤ 3 U	≤ 3 U	≤ 3 U	≤ 3 U	≤ 3 U	≤ 3 U	≤ 3 U	≤ 3 U	≤ 3 U		
	TL synOcta® Cementable Abutments		Single	Single	≤ 3 U	≤ 3 U	≤ 3 U	≤ 3 U	≤ 3 U	≤ 3 U	≤ 3 U	≤ 3 U	≤ 3 U		

Single = Single tooth | Full = Full arch | U = Units

Name		Material	Shades
n!ce®		Lithium aluminosilicate ceramic reinforced with lithium disilicate	Available in 16 shades
IPS e.max® CAD		Lithium- disilicate	Available in 45 shades
zerion® GI		Gingiva shaded zirconia	Available in 2 shades
zerion® LT		High-strength translucent zirconia	Available in 9 shades
zerion® ML ⁴		High translucent multi-layered zirconium dioxide	Available in 5 shades
3M™ ESPE™ Lava™ Plus Zirconia ⁴		High translucent zirconium dioxide	Available in 16 shades
3M™ Lava™ Esthetic		Multi-layered ultra-high translucent zirconia	Available in 5 shades
zerion® UTML ^{2,3}		Multi-layered ultra-high translucent zirconia	Available in 7 shades
coron®		Cobalt Chromium	n/a
ticon®		Titanium Grade 4	n/a
polycon® ae ⁴		PMMA	Available in 5 shades
VITA® CAD-Temp		Monochromatic or multi-layer cross-linked acrylate polymer with microfiller	Available in 4 monochromatic and 3 multi-layered shades
JUVORA™ PEEK		High performance polymer	Available in 2 shades

Note:

Application and material availability might differ from country to country. Please contact your local sales representative for a detailed overview of the available applications and prosthetic lines.

International Headquarters

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone +41 (0)61 965 11 11
Fax +41 (0)61 965 11 01
www.straumann.com

3M™, ESPE™, Lava™ are trademarks of 3M or 3M Deutschland GmbH.

Used under license in Canada. IPS e.max® is a registered trademark of Ivoclar Vivadent AG, Liechtenstein.

VITA® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & C o. KG, Germany.

© Institut Straumann AG, 2020. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.