

Technical Sheet

Straumann® Emdogain® FL
Institut Straumann AG



PRODUCT INDICATION

Product type	Adjunct to scaling and root planing and dental hygiene procedures
Indication	<p>Emdogain® FL is intended for subgingival and topical application in conjunction with scaling and root planing procedures to provide periodontal pocket regeneration in a minimally invasive way: without surgical flap.</p> <ul style="list-style-type: none">• Emdogain® FL has been shown to be effective in residual pockets with probing depths from 5 mm to 9 mm with no furcation involvement in patients with adequate plaque control.• Emdogain® has also been shown to enhance the early healing of periodontal soft tissue wounds resulting from the instrumentation of periodontal pockets.

PRODUCT PROPERTIES

Origin	Sweden
Storage temperature	2°C–8°C
Product specifications	<ul style="list-style-type: none">• Straumann® Emdogain® consists of one unit of Straumann® Emdogain® (0.15 ml or 0.3 ml) and one unit of Straumann® PrefGel® (0.6 ml) and 2 application syringes• Straumann® Emdogain® contains 30 mg/ml Enamel Matrix Derivative• Straumann® PrefGel® contains Ethylenediaminetetraacetate (EDTA) 24 %
Accessories	N/A
Sterilization / Disinfection	Sterile & Single use
Shelf life	2 years from the date of production
Regulatory status	CE-marked
Entry on the market	March 2019

MANUFACTURING

Production process

Manufactured from porcine unerupted tooth buds according to the highest pharmaceutical standards under current good manufacturing practice (cGMP). Complies with all relevant ISO and FDA standard and guidelines.

Packaging

Cardboard packaging

AVAILABLE PRODUCT LINE

Product	Type/Size/Amount	Art.-No.
Straumann® Emdogain® FL 0.15 ml	Contains 0.15 ml Emdogain®, 0.6 ml PrefGel® and 2 application cannulas	075.130
Straumann® Emdogain® FL 0.3 ml	0.3 ml Emdogain®, 0.6 ml PrefGel® and 2 application cannulas	075.131

HOTLINE

UK: (+44) 01 293 651230

Italy: (+39) 800 810 500

Germany: (+49) 761 450 13 50

International Headquarters

Institut Straumann AG

Peter Merian-Weg 12

CH-4002 Basel, Switzerland

Phone +41 (0)61 965 11 11

Fax +41 (0)61 965 11 01

www.straumann.com

© Institut Straumann AG, 2018. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

