

STRAUMANN® N!CE® ZIRCONIA

Where innovation meets esthetics.



WHERE **INNOVATION** MEETS ESTHETICS.

Zirconia can be applied to a wide range of indications. It lends itself well to an efficient digital manufacturing workflow.

Why Zirconia? Zirconia is the strongest ceramic, prosthetic-milling material available for dental restorations. It combines durability and high strength with unsurpassed, natural esthetics. It provides a high-level of natural looking color accuracy and shade-matching.

Why n!ce®? With the introduction of n!ce® Zirconia, Straumann® takes Zirconia to the next level of innovation with a unique portfolio designed to provide an unparalleled range of treatment options that deliver natural, long-lasting results. Our n!ce® discs feature an ingenious, simple and user-friendly packaging and labelling solution facilitating a seamless full-digital manufacturing workflow.

Why Straumann®? Straumann® stands for premium Swiss quality, precision and cutting-edge innovation that delivers confidence in esthetic dentistry backed by the largest global scientific network.

Simply the natural choice.





ESTHETICS

High color accuracy with built-in shade matching and multi-layered feature designed to deliver natural-looking prosthetic outcomes.



VERSATILITY

All you need for your patient cases for a wide range of indications from single units to multi-unit bridges and full arches.



EFFICIENCY

Just one step to finish the restoration with full contour design and an efficient digital manufacturing workflow.



FLEXIBILITY

Synergies between all digital workflows, Straumann® offer simple outsourcing options on demand.



HIGH QUALITY

Premium Swiss quality designed for predictable prosthetic outcomes you can trust.



CONVENIENCE

An ingenious, simple and user-friendly packaging and labelling solution facilitating a seamless full-digital manufacturing workflow.

ESTHETICS

The Straumann® n!ce® Zirconia range delivers premium, natural esthetics with innovative color technology and accurate shade matching.



The multi-layered design with integrated color layers ensures a seamless color gradient for natural looking results, without the need of an additional veneering or staining step.

The n!ce® Zirconia portfolio includes a wide range of translucencies: LT (Low Translucency), HT (High Translucency), XT (eXtra high Translucency) and a new material generation GT (Gradient Translucency), as well as monochromatic and multi-layered VITA shades for esthetic restorations.

Shade evaluation was performed by dental professionals to guarantee accurate color-matching and the natural-looking prosthetic outcomes you expect from Straumann®.

- Natural ceramic material
- All heights and shades available to achieve the optimal esthetic result
- Accurate color and shade-matching technology
- True-to-life multi-layered gradient



“You probably know Straumann® as the market leader in dental implants. But did you know that we’re changing the face of dental prosthetics with our n!ce® Zirconia solutions and next generation digital workflow?”

Florian J. Kirsch, Head Connected Customer Solutions.

HIGH QUALITY

Straumann® delivers premium Swiss quality, precision and cutting-edge innovation that delivers confidence and predictable prosthetic outcomes you can trust.

Backed by the largest global scientific network, Straumann® has taken Zirconia to the next level of innovation. Zirconia is the strongest ceramic, prosthetic-milling material available for dental restorations. The n!ce® Zirconia portfolio combines reliability, durability and high strength with an unparalleled range of treatment options that deliver natural, reliable, long-lasting results.

The n!ce® Zirconia range:

- Delivering high-quality results you can rely on
- Validated process and precision engineered
- Swiss-quality production, that delivers peace-of-mind

CONVENIENCE

Straumann® n!ce® Zirconia discs feature an ingenious, simple and user-friendly packaging and labelling solution that facilitates a seamless full-digital, CAD/CAM manufacturing workflow. Streamlined for simplicity, convenience and efficiency.



Our new convenient packaging is userfriendly to ensure clear and simple material identification. Our ingenious disc labelling concept ensures that all critical information required to mill and sinter, such as shade and enlargement factor, remain clearly visible on the disc itself even if the disc is partly milled. This ensures that the disc can easily be re-positioned and re-used later to produce subsequent restorations.

Our n!ce® Zirconia discs can be used as part of your existing zirconia process without the need to adapt your established workflow. No need to learn new processing steps or purchase new hardware. Thanks to the standard disc geometry, our n!ce® Zirconia materials are optimized to work with your existing Straumann® milling machines and furnaces as well as third-party hardware.

Furthermore, we use minimal packaging for storage optimization. All our packaging uses recycled, and environmentally friendly FSC* carton instead of traditional foam in line with our sustainability goals.



Our n!ce® discs feature an ingenious, simple and customer friendly packaging & labelling solution.



*FSC: The Forest Stewardship Council is an international non-profit organization with a mission to create and sustain responsible forestry.

EFFICIENCY

Just one step to finish the restoration with full contour design to limit post processing and facilitate an efficient digital manufacturing workflow.

With our nlce® material range we are harmonizing and simplifying our milling materials portfolio for centralized and in-house milling for enhanced efficiency.

With Straumann® digital solutions you can increase productivity and profitability from planning to production and finishing. Improve your efficiency and shorten turnaround times with a streamlined and reliable workflow.

Save time with improved efficiency and high predictability.



Monochromatic discs available in various shades.



Multi-layered discs for esthetic restorations.



FLEXIBILITY

Straumann® offers a range of outsourcing options with synergies between all digital workflows. Manufacture in-house or gain access to milled restorations from our centralized service. The choice is yours.

Straumann® in-house options	SCAN	PLAN & DESIGN	MILL	FINISH
Straumann® services	SCAN & SHAPE		CENTRALIZED MILLING	FINISH

From our centralized milling services we can manufacture the restoration according to your design. Or we can deliver the design with our Scan & Shape service, a cost-effective option to in-house milling without compromising on quality or esthetics. Straumann® delivers affordable premium Swiss quality restorations.

Or benefit from full synergy with your inLab solution. The complete portfolio of 400 n!ce® Zirconia discs, shades and translucency levels are available and fully compatible with your established in-house production process. Enjoy the flexibility to select the workflow that works for you so that you can maximize resource management.





VERSATILITY

Premium treatment options

With the introduction of n!ce® Zirconia, Straumann® takes Zirconia to the next level of innovation with a unique portfolio designed to provide an unparalleled range of treatment options that deliver natural, long-lasting results. With our broad portfolio, n!ce® Zirconia provides all you need for your patient cases for all indications from a single crown to multi-unit bridges and full arches.

Excellent shade matching

n!ce® Zirconia is available in 16 shades corresponding to the VITA tooth shade guide and two bleach shades to facilitate post processing to enable you to deliver the final restoration in the desired shade.

- Disc heights range from 12 mm to 30 mm
- Restoration design options from framework to full contour monolithic restorations
- Processing options: glazing, staining, layering/veneering

A complete offering for esthetic restorations

	Application	
Low Translucency	Framework	Up to full arch
High Translucency	Full contour	Up to full arch
eXtra high Translucency	Full contour	Up to 3-unit bridges
Gradient Translucency	Full contour	Up to full arch

THE **NEXT GENERATION** IN DIGITAL INNOVATION. POWERED BY STRAUMANN.®

Innovation drives productivity. And productivity delivers growth. The next generation of end-to-end digital manufacturing solutions is changing the face of dentistry.

With Straumann® digital solutions you can increase productivity and profitability from planning to production and finishing. Improve your efficiency and shorten turnaround times with a streamlined and reliable workflow. We help you do more, faster and better. So that you can focus on what matters. Your patients.



SCAN

Discover our broad spectrum of intraoral and lab scanners. Choose the digital solution that best fits your needs from plug&play simplicity to empowering innovation that redefines the boundaries of possibility. The perfect starting point in your digital workflow. Our innovative scanning technology empowers you to scan and transmit accurate, digital impressions seamlessly.




Virtuo Vivo™, 3Shape TRIOS®, Medit scanners



PLAN & DESIGN

Connectivity, piece of mind and great results with coDiagnostiX® and CARES® Visual. CoDiagnostiX® – Our state-of-the-art dental treatment planning software for predictable results.

CARES® Visual – The user-friendly software with a complete and comprehensive dental interface, giving you access to a variety of specific services and design solutions that will enhance your daily practice experience and deliver the best esthetic solution for your patient.



coDiagnostiX® and CARES® Visual



nIce® Zirconia discs



MILL

The M series/Motion 2 is a universal, cost-efficient production system that can handle the majority of cases. With multiple workflow options, you can enjoy the full flexibility of using one system regardless of what challenges arise. And for work peaks or complex cases, our centralized milling facility operates as an extension to your lab.



M series/Motion 2



SINTER


The Therm 3 furnace maximizes process reliability to deliver homogeneous, distortion-free restorations. With constant temperature control and distribution, and operational simplicity ensure consistent performance.



Therm 3

STRAUMANN® NICE® ZIRCONIA – WHERE INNOVATION MEETS ESTHETICS.

- A premium range of treatment options
- Natural, long-lasting results
- User-friendly packaging and labelling solution
- Seamless, digital manufacturing workflow



490.770/en/B/00 11/23

International Headquarters

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone +41 (0)61 965 11 11
Fax +41 (0)61 965 11 01
www.straumann.com

© Institut Straumann AG, 2023. All rights reserved.
Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

 **straumann**