


A collection of various dental prostheses, including metal and ceramic partial dentures, complete dentures, and individual dental implants, displayed on a green and black background.

SCALE UP WITH STRAUMANN UN!Q™

An outsourced, open-platform design and production service for patient-specific prosthetics.

In a time of digital transformation, when resources are scarce and competition is fierce, Straumann® UN!Q empowers you with a premium service to outsource the planning, design and manufacturing of your custom implant prosthetics.

EXCEL WITH CUSTOM PROSTHETICS



Scalable: grow your business without further investment in CAD resources or CAM equipment.



Reliable: deliver consistent, predictable results by leveraging the Straumann® 70-year legacy of Swiss engineering.



Comprehensive: extensive portfolio and wide compatibility with over 50 implant platforms and design options.



Flexible: pick and choose the on-demand workflow that is tailored to your needs.

70

Years of
Straumann Excellence
and Innovation

2 MILLION+

Smiles Restored

125

Implant Platforms Covered



RESTORATIVE FLEXIBILITY

With original systems, unequalled precision and consistent high quality, the Straumann® prosthetic portfolio addresses patients' esthetic expectations with the flexibility to fit their budget.

STRAUMANN UN!Q™ **FIXED** FULL ARCH PORTFOLIO HIGHLIGHTS



STRAUMANN® COREBAR™ TO SUPPORT ZIRCONIA RESTORATIONS

- Provides larger surface area for bonding vs. titanium-base coping and enhanced distribution of forces
- Lighter than restorations produced with solid zirconia
- Available through design service
- Efficient workflow – Download the core file for the bar from the AXS platform and begin milling your zirconia sleeve while the bar is manufactured



THIMBLE BRIDGE

- Designed for easy maintenance and repair
- Angled screw channels for enhanced esthetics
- Individual crowns allow for a combination of esthetic and high-strength materials



BASIC AND ADVANCED HYBRIDS

- Acrylic hybrids are an ideal solution for ongoing tissue management with easy tissue relining
- Design versatility for retention
- Circumferential ledge provides additional support for teeth to prevent shearing, increasing durability and stability

STRAUMANN UN!Q™ REMOVABLE PORTFOLIO HIGHLIGHTS



NOVALOC® LOCATOR® BAR

- Provides more resistance against mastication forces compared to abutment-supported attachments
- Aligns attachments into a common insertion axis
- Customizable height and width
- 0°, 4°, 6° and 8° wall taper



2 IN 1 TELESCOPIC

- Provides a combination of the passive fit of a fixed restoration and the cleanability of a removable restoration
- Circumferential seal, created between the gingiva and the secondary structure, promotes hygiene and gingival health



FULL RANGE OF REMOVABLE PROSTHETICS


- The Straumann UN!Q™ portfolio includes multiple options for removable prosthetic attachment systems including:
 - Dolder® U Shape & Egg-Shape Bars
 - Preci Vertex and Horix Bar
 - Ackermann® Bar
 - Round® Bar
 - MP-Clip® Bar

"Working with your service feels effortless, the flow is very fast and simple. I love the design approval process and the design quality is great."


YDL Dental, Texas

DESIGN & MANUFACTURING SERVICES

Outsourcing your design work is only a few clicks away. Take advantage of our premium on-demand service via Straumann AXS™, the online case ordering and management portal.



www.straumann.com/axs



Key benefits – Straumann AXS™ Portal

- Simple and easy to use
- Track & trace your cases all in one place
- Seamless core file workflow
- Quick access to digital training & education opportunities
- Services easily accessible via our e-shop

TURNAROUND TIMES

DESIGN TURNAROUND TIMES

**SINGLE UNITS
(ABUTMENTS)
& 3D MODELS**



24 hours

Cut-off time:



2 pm CET

**MULTI-UNITS
(BRIDGES AND BARS)**



48 hours

Cut-off time:



2 pm CET

PRODUCTION TURNAROUND TIMES


CUSTOM ABUTMENT TURNAROUND (Titanium and Cobalt Chrome*)

Case received or design approved prior to 2 pm		Case Shipped	
Monday	→	Tuesday	
Tuesday	→	Wednesday	
Wednesday	→	Thursday	
Thursday	→	Friday	
Friday	→	Monday	

BAR TURNAROUND

Case received or design approved prior to 2 pm		Case Shipped	
Monday	→	Wednesday	
Tuesday	→	Thursday	
Wednesday	→	Friday	
Thursday	→	Monday	
Friday	→	Tuesday	

*Turnaround times for Ti and CoCr abutments are shown above.



International Headquarters

Institut Straumann AG

Peter Merian-Weg 12

CH-4002 Basel, Switzerland

Phone +41 (0)61 965 11 11

Fax +41 (0)61 965 11 01

www.straumann.com

© Institut Straumann AG, 2025. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

