



STRAUMANN SIRIOS™ X3

Capture Excellence.




Product overview

STRAUMANN SIRIOS™ X3 **CAPTURE** **EXCELLENCE**

Straumann SIRIOS™ X3 is your gateway to the ultimate dental experience: unlocking tailored treatment journeys and several possibilities for every unique clinical need.

Experience next-level scanning, and step into the Straumann® integrated digital ecosystem, where connectivity is effortless and collaboration is easy, redefining what's possible in your practice. Powered by world-class education and support, Straumann SIRIOS™ X3 is built to transform your digital dentistry experience.





UNLOCK TAILORED TREATMENT JOURNEY EXPERIENCE

As the gateway to the Straumann® integrated digital ecosystem, Straumann SIRIOS™ X3 unlocks multiple treatment possibilities.



NEXT-LEVEL SCANNING EXPERIENCE

Straumann SIRIOS™ X3 transforms the intraoral scanning experience – making it faster and smarter.



ELEVATED STRAUMANN EXPERIENCE

From Straumann®, you can expect nothing less than premium service and customer experience that follow you since the beginning of your journey with us.

UNLOCK TAILORED TREATMENT JOURNEY EXPERIENCE

As the gateway to the Straumann® integrated digital ecosystem, Straumann SIRIOS™ X3 unlocks multiple treatment possibilities.

- Access all workflows and applications directly from your Straumann SIRIOS™ X3 software.
- Gain instant access to the Straumann® integrated digital ecosystem, thanks to automatic upload and synchronization with Straumann AXS™.



START YOUR TAILORED TREATMENT JOURNEY DIRECTLY FROM YOUR STRAUMANN SIRIOS™ X3

Select your preferred workflow and enjoy seamless integration with intuitive, AI-powered guidance – all from your Straumann SIRIOS™ X3.

IMPLANT-BORNE AND TOOTH-BORNE WORKFLOWS

Straumann SIRIOS™ X3 features automatic scanbody recognition to streamline your workflow and provide real-time validation of scan quality.

Dedicated and validated full-arch workflow is designed to achieve high accuracy for full-arch implant restorative treatment, all combined in one single scanning case. Experience superior workflow efficiency when using Straumann EXACT™, our validated solution for full-arch restoration.

SURGICAL PLANNING WORKFLOW

Plan implant cases and design surgical guides easily. Export IO scans to coDignostiX® installed on your laptop.

CHAIRSIDE 3D PRINTING WORKFLOW

Produce your temporary or permanent single-tooth restorations effortlessly thanks to seamless integration with Straumann Signature Midas by SprintRay.

ORTHODONTIC WORKFLOW

Treat aligner, revision and retainer cases from the Straumann SIRIOS™ software. Enjoy direct access to the scan comparison app along with seamless patient case connectivity to Dr. Portal from ClearCorrect®.



APPLICATIONS

Gain direct access to design capabilities and patient communication tools:

Design



Model Builder



Chairside Solutions

Patient Communication



Scan Compare


STREAMLINE YOUR TAILORED TREATMENT JOURNEY AND ENHANCE COLLABORATION WITH STRAUMANN AXS™


Straumann SIRIOS™ X3 connects seamlessly to the Straumann® integrated digital ecosystem via Straumann AXS™ – your central hub for end-to-end digital workflows.

Straumann AXS™ is a cloud-based platform that enables secure, open, and modular workflow integration across treatments.

With automatic scan uploads and synchronized patient data, you eliminate duplicate entry and gain instant access to the next step – making your workflow simpler, faster, and more connected.

This level of integration enables greater efficiency, consistency, and precision – so you can deliver optimal outcomes for every patient.





SEND TO LAB WORKFLOW

Streamline your practice and communication easily sharing cases with your lab in only one click.

SMILE IN A BOX® WORKFLOW

Outsource guided surgery cases easily with the streamlined Smile in Box® case submission, directly from Straumann AXS™.

STRAUMANN UN!Q™ WORKFLOW

Outsource design and production of patient-specific restorations easily on Straumann AXS™.

GET PEACE-OF-MIND THANKS TO A SECURE PLATFORM


Store, process and share patient information safely and in full compliance with regulations in the cloud in one single place.

FLEXIBILITY TO SELECT YOUR PREFERRED WORKFLOW WITH AN OPEN-SYSTEM SOLUTION


Enjoy the freedom to import and export files and integrate with different hardwares and softwares. Work with your preferred workflows and tailor your process to your specific needs.




Smile in a Box®



Straumann
SIRIOS™ X3



Straumann
EXACT™



Straumann
UN!Q™

STRAUMANN AXS™

■ SURGICAL ■ RESTORATIVE

NEXT-LEVEL SCANNING EXPERIENCE

Straumann SIRIOS™ X3 transforms the intraoral scanning experience – making it faster and smarter. It builds on trusted performance and takes every feature a step further, to help you capture excellence every time.

PEN-GRIP PREMIUM DESIGN

This sleek and compact wireless scanner combines cutting-edge technology with a sophisticated, pen-like design. It provides a natural grip and maximum comfort – for both clinician and patient.

LIGHTWEIGHT AT JUST 240 g

Our lightweight solution makes scanning effortless and helps to reduce fatigue coming from an extended use.

ONE OPTIMIZED TIP SIZE

One optimized tip size that suits your clinical needs. It ensures premium experience for your patients while providing an excellent field of view, even in challenging areas.

LED RING

The built-in LED ring offers you real-time visual guidance throughout the scanning process, enhancing usability and helping you focus on what matters most – providing the best patient experience.





INTELLIGENT ACTION BUTTON

Gain instant access to streamlined workflows and your most-used functions or shortcuts, streamlining every step of your scanning experience.

FLEXIBLE SCANNING STRATEGY

Enjoy flexibility in your scanning process. Adapt easily to each clinical situation, without the need for movement restriction or mandatory scanning strategy.

HIGH SPEED SCANNING

Without compromising great outcomes, Straumann SIRIOS™ X3 delivers seamless performance: whether you're scanning a single tooth or a full arch. With fast-responsive technology and real-time rendering, every detail is captured smoothly and accurately, ensuring a faster, more efficient clinical workflow.

FAST WIRELESS CHARGING

Straumann SIRIOS™ X3 is designed to deliver uninterrupted superior performance. It reaches 50% charge in just 60 minutes, providing over 2 hours of uninterrupted scanning on a full charge. With a 6-month standby time, it's always ready when you are.

CALIBRATION FREE

Engineered for lasting precision, it's ready to scan from day one. No need for costly calibration tools or time-consuming steps.



THE ELEVATED STRAUMANN® EXPERIENCE

From Straumann®, you can expect nothing less than premium service and customer experience that follow you since the beginning of your journey with us.

SUPERIOR TRAINING AND EDUCATION


Our commitment to education is a cornerstone of effective patient care and practice success. We offer a wide range of educational programs, including classroom courses, eLearning, webinars – all integrated into one digital platform – worldwide.

SUPPORT YOU CAN RELY ON

Have access to a customer support program that offers quick issue resolution, availability in multiple geographic locations, and local language assistance. Benefit from a comprehensive customer care package that includes extended warranty, full accidental coverage, and all the valuable Straumann® customer care services you already know and trust.

STRAUMANN® IS THE FLAGSHIP BRAND OF THE WORLD'S LARGEST IMPLANT MANUFACTURER GROUP – THE STRAUMANN GROUP.

Renowned as one of the most iconic brands in dentistry, we have a long-standing legacy of pioneering significant advancements in the field. We are at the forefront of expanding the scope of dental care by relentlessly and consistently providing excellence and transformational innovation to a global community of forward-thinking dental professionals.



STRAUMANN SIRIOS™ X3 CAPTURE **EXCELLENCE.**



491.046/en/A/00 08/25

International Headquarters

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone +41 (0)61 965 11 11
Fax +41 (0)61 965 11 01
www.straumann.com

© Institut Straumann AG, 2025. All rights reserved.
Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

