

DIGITAL MATERIALS OVERVIEW

MATERIALS ¹

	MILLING				SERVICE / SOLUTION			APPLICATION ¹											
	CARES® C Series	CARES® M series	CARES® D series	Centralized milling	CARES® X-Stream™	CARES® Scan and Shape	Connectivity	Tooth restoration			Others	Direct connection to implants							
								Inlays ² / Onlays ² / Veneers	Copings / Crowns	Bridges / Bars (maximum units)	Splints	CARES® Customized Abutments and Pre-Milled Blanks	etkon® iDent Customized Abutments and Pre-Milled Blanks	CARES® Screw-Retained Bridges and Bars					
Glass-ceramics	Translucency																		
nIce®	LT, HT																		
IPS e.max® CAD	MO, LT, HT																		
VITA® Suprinity	LT, HT																		
VITA® Mark II	HT																		
VITA® TriLuxe / TriLuxe Forte	HT																		
Hybrids	Translucency																		
3M™ ESPE™ Lava™ Ultimate	LT, HT																		
VITA® Enamic	LT, HT																		
Ceramics	Translucency																		
zerion® UTML	Ultra HT																		
zolid® SHT ML	Super HT																		
zolid® SHT / SHT preshade	Super HT																		
zerion® ML	HT																		
3M™ ESPE™ Lava™ Plus Zirconia	HT																		
zolid® HT / HT preshade	HT																		
zerion® LT	LT																		
zi	LT																		
zerion® GI	LT																		
Zirconium Dioxide for CARES®	LT																		
Metal																			
coron® Cobalt-Chromium																			
Cobalt-Chromium for etkon® iDent																			
Sintron®																			
TAN for CARES®																			
ticon® Titanium Grade 4																			
Titanium for CARES®																			
Titanium for etkon® iDent																			
Polymer																			
polycon® ae																			
VITA® CAD-Temp																			
Temp / Temp multilayer																			
PEEK																			
Splintec																			
polycon® cast																			
PMMA (burn out polymer)																			
Wax																			

¹ application and material availability might differ from country to country.

² veneerable inlay/onlay for 3M™ ESPE™ Lava™ Plus Zirconia, zerion®, ticon®, coron®, polycon® ae and polycon® cast materials.

³ material listed and mill-able

⁴ CARES® Abutments available, no Pre-Milled Abutments Blanks available

MO: Medium Opacity LT: Low Translucency HT: High Translucency GI: Gingiva shaded

VITA® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & Co. KG, Germany.

3M™, ESPE™ and Lava™ are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada.

JUVORA™ is a trademark of JUVORA Ltd.

IPS e.max® is a registered trademark of Ivoclar Vivadent AG, Liechtenstein.