










# CARES® X-Stream™ restorative options

## MATERIAL OVERVIEW

- 1 Bridges and bars are available with up to 10 implant abutment connections.
- 2 Molar bridge excluded
- 3 Not available for Bars
- 4 Only available for Fixed bars
- 5 Not available for NNC and WN
- 6 Not available for NNC
- 7 Restorations only available without screw channel hole
- 8 Restorations only available with screw channel hole

		Tooth replacement options								
		Full-contour or frameworks for single tooth, bridges and bars								
		Ceramic					Metal		Polymer	
		nIce®	IPS e.max® CAD	zerion® UTML <sup>2,3</sup>	zerion® ML	3M™ ESPE™ Lava™ Plus Zirconia <sup>4</sup>	zerion® <sup>4</sup>	coron®	ticon®	polycon® ae <sup>3</sup>
PURE Ceramic	Implant Monotype									
CARES® Abutments	CARES® Abutment, Zirconium dioxide <sup>5,6</sup>			≤ 3 units	Full arch	Full arch	Full arch			Full arch
	CARES® Abutment, Titanium <sup>6,7</sup>			≤ 3 units	Full arch	Full arch	Full arch	Full arch	Full arch	Full arch
	CARES® Abutment, CoCr <sup>6,7</sup>			≤ 3 units	Full arch	Full arch	Full arch	Full arch	Full arch	Full arch
Variobase®	Variobase® for Crown			Single tooth	Single tooth	Single tooth	Single tooth	Single tooth		Single tooth
	Variobase® for Crown AS <sup>8</sup>			Single tooth	Single tooth	Single tooth	Single tooth	Single tooth		Single tooth
	Variobase® for Bridge/Bar			≤ 3 units	Full arch	Full arch	Full arch	Full arch	Full arch	Full arch
	Variobase® for Bridge/Bar Coping for Screw-retained Abutments <sup>8</sup>			≤ 3 units	Full arch	Full arch	Full arch	Full arch	Full arch	Full arch
	Variobase® for Bridge/Bar Cylindrical <sup>8</sup>			≤ 3 units	Full arch	Full arch	Full arch	Full arch	Full arch	Full arch
	Variobase® for Bridge/Bar Cylindrical Coping for Screw-retained Abutments <sup>8</sup>			≤ 3 units	Full arch	Full arch	Full arch	Full arch	Full arch	Full arch
Stock Abutments	BL Cementable Abutments			≤ 3 units	≤ 3 units	≤ 3 units	≤ 3 units	≤ 3 units	≤ 3 units	≤ 3 units
	TL synOcta® Cementable Abutments			≤ 3 units	≤ 3 units	≤ 3 units	≤ 3 units	≤ 3 units	≤ 3 units	≤ 3 units

Name		Material	Shades
nIce®		Lithium aluminosilicate ceramic reinforced with lithium disilicate	12 shades
IPS e.max® CAD		Lithium- disilicate	45 shades
zerion® UTML <sup>3,4</sup>		Ultra translucent multi-layered zirconium dioxide	4 VITA shades and 1 Bleach shade
zerion® ML <sup>5</sup>		High translucent multi-layered zirconium dioxide	3 shades
3M™ ESPE™ Lava™ Plus Zirconia <sup>5</sup>		High translucent zirconium dioxide	16 VITA shades and 2 Bleach shades
zerion® <sup>5</sup>		Low translucent zirconium dioxide	9 shades
coron®		Cobalt Chromium	n/a
ticon®		Titanium Grade 4	n/a
polycon® ae <sup>4</sup>		PMMA	5 shades

**Note:**

Application and material availability might differ from country to country. Please contact your local sales representative for a detailed overview of the available applications and prosthetic lines.

**International Headquarters**

Institut Straumann AG  
 Peter Merian-Weg 12  
 CH-4002 Basel, Switzerland  
 Phone +41 (0)61 965 11 11  
 Fax +41 (0)61 965 11 01  
 www.straumann.com

3M™, ESPE™, Lava™ are trademarks of 3M or 3M Deutschland GmbH.  
 Used under license in Canada. IPS e.max® is a registered trademark of Ivoclar Vivadent AG, Liechtenstein.  
 VITA® is a registered trademark of VITA Zahnfabrik H. Rauter GmbH & C o. KG, Germany.

© Institut Straumann AG, 2018. All rights reserved.  
 Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.