

# Technical Sheet

## CURODONT™ REPAIR Credentis



### PRODUCT INDICATION

<b>Product type</b>	Biomimetic mineralization of initial carious lesions and similar defects
<b>Indication</b>	<ul style="list-style-type: none"><li>• CURODONT™ REPAIR is recommended for Guided Enamel Regeneration of initial carious lesions (proximal, buccal and occlusal) and similar defects</li><li>• Guided Enamel Regeneration (GER) for the treatment of initial carious lesions</li></ul>

### PRODUCT PROPERTIES

<b>Origin</b>	Switzerland
<b>Storage temperature</b>	Up to 25°C
<b>Product Specifications</b>	1 box contains 10 Applicators
<b>Component Parts</b>	Self-Assembling Peptide P11-4, Synthetic Product One CURODONT™ REPAIR applicator unit contains both CUROLOX® Technology for Guided Enamel Regeneration and water for activation of the device
<b>Sterilization / Disinfection</b>	n/a
<b>Shelf life</b>	4 years
<b>Regulatory status</b>	Medical Device, Class IIa
<b>Entry on the market</b>	2012 1 <sup>st</sup> gen, 2015 2 <sup>nd</sup> gen

## MANUFACTURING

<b>Product</b>	Applicators are coated with the active substance (FMOC-peptide synthesis)
----------------	---

<b>Packaging</b>	Applicator in one bag, 10 bags in one box
------------------	---

## AVAILABLE PRODUCT LINE

Product	Type/Size/Amount	Art.-No.
---------	------------------	----------

CURODONT™ REPAIR – Distributed by Straumann	1 applicator in one bag, 10 bags in Box	P540.S
--	---	--------

## HOTLINE

**UK:** (+44) 01 296 651230  
**Italy:** (+39) 800 810 500  
**Germany:** (+49) 761 450 13 50

### International Headquarters

Institut Straumann AG  
Peter Merian-Weg 12  
CH-4002 Basel, Switzerland  
Phone +41 (0)61 965 11 11  
Fax +41 (0)61 965 11 01  
[www.straumann.com](http://www.straumann.com)