

Technical Sheet

CARISOLV®
RLS Global AB



PRODUCT INDICATION

Product type	Caries Removal Gel
Indication	Softens demineralized carious dentine prior to mechanical removal during chemo-mechanical minimal invasive treatment of patients diagnosed with caries

PRODUCT PROPERTIES

Origin	Sweden
Storage temperature	Keep the product refrigerated (2–8°C) in its original packaging during storage
Product specifications	CARISOLV® consists of two solutions provided in a double syringe with a closable cap. The solutions are mixed to form a gel at the time of use. Each package consists of 2 syringes packaged in aluminium pouches, 2 plungers, 10 static mixers and 10 intra-oral tips. The mixing of the solutions forms a clear gel, with a slight chlorine odour
Ingredients	A clear, viscous solution consisting of water, carboxymethyl cellulose, sodium hypochlorite, sodium chloride, and amino acids (glutamic acid, leucine, lysine) with alkaline pH
Sterilization / Disinfection	N/A
Shelf life	15 months
Regulatory status	EU: Medical Device Class I
Entry on the market	1998

PACKAGING

Packaging	2 pieces of Carisolv units are packaged together with 2 plungers, 10 static mixers, 10 intra-oral tips and an IFU in a shelf box
------------------	--

AVAILABLE PRODUCT LINE

Product	Type/Size/Amount	Art.-No.
Carisolv®	2 pcs twin syringes containing 2.5 ml gel, 10 pcs mixing elements and 10 pcs applicator tips For up to 8 treatments depending on size of the cavity	10491

CUSTOMER CONTACT

RLS Global AB
Customer contact

Tel.
+46 31-780 68 20
E-mail
customer@rlsglobal.se

UK: (+44) 01 296 651230
Italy: (+39) 800 810 500
Germany: (+49) 761 450 13 50

International Headquarters

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone +41 (0)61 965 11 11
Fax +41 (0)61 965 11 01
www.straumann.com