

Technical Sheet

PerioPatch®

Izun Pharmaceuticals Corp.



PRODUCT INDICATION

Product type	Ultra-thin topical oral Patch containing Botanical Ingredients
Indication	<ul style="list-style-type: none">• Gingival inflammation• Oral wounds

PRODUCT PROPERTIES

Origin	Botanicals
Storage temp.	15–30°C (room temperature)
Product specifications	Adhesive matrix containing Echinacea purpurea, Sambucus nigra, Centella asiatica
Component Parts	Ethylcellulose, polyacrylic acid, Sambucus nigra, castor oil, acacia gum, methyl hydroxypropyl cellulose, glycerol, strawberry flavor, Centella asiatica, titanium dioxide, Echinacea purpurea, polysorbate 80, CI 77491
Sterilization / Disinfection	Non-sterile product
Shelf life	3 years
Regulatory status	<ul style="list-style-type: none">• CE Europe• FDA 510(k) cleared USA
Entry on the market	2016 (USA)

MANUFACTURING

Production process

Thin film extrusion technology to produce bilayer film on PET liner

Packaging

- Patch on PET liner in individually heat-sealed pouch
- 6 patches on each card; one card per envelope - represents individual patient unit of use
- 5 envelopes per display box

AVAILABLE PRODUCT LINE

Product	Type/Size/Amount	Art.-No.
PerioPatch® Envelope	6 individually pouched patches	5115018
PerioPatch® Display Box	5 envelopes (total of 30 patches)	5115056

HOTLINE

UK: (+44) 01 296 651230

Italy: (+39) 800 810 500

Germany: (+49) 761 450 13 50

International Headquarters

Institut Straumann AG

Peter Merian-Weg 12

CH-4002 Basel, Switzerland

Phone +41 (0)61 965 11 11

Fax +41 (0)61 965 11 01

www.straumann.com

© Institut Straumann AG, 2018. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

