

EXCELLENCE IN IMMEDIACY

Straumann® BLX
Confidence beyond Immediacy.



The next generation solution for confidence beyond immediacy.

Straumann® BLX has been perfected for immediacy and is also an excellent solution for all other indications to suit the dentist's preferred treatment protocol – ranging from immediate to conventional placement and loading.

Straumann® BLX has been developed to improve the predictability and utilization of immediacy as an option for every-day applications. It instils confidence in and for dentists through the combination of legacy, experience and innovation. Straumann® BLX is the culmination of Swiss premium quality, precision and pioneering innovation combined with enabling education in dentistry backed by one of the largest global scientific networks, the ITI.

Straumann® BLX System lets you discover and capitalize on new business opportunities by using chair-time more efficiently to treat more patients with confidence.

“For immediate protocols I want products that provide me with peace of mind in demanding clinical situations. Straumann® BLX Implants with Roxolid® and SLActive® give me that confidence.


BLX is exciting because it extends the treatment options we can offer with Straumann® products. It's a new era of implant treatment.”

Dr. Eirik Salvesen
Stavanger, Norway



IMMEDIATE PROTOCOLS:

- Immediate procedures are clinically well documented
- High survival rates, comparable to conventional protocols¹



FOR THE PATIENT THIS MEANS:

- Fewer dentist visits
- Shorter treatment cycle

FOR THE CLINICIAN THIS MEANS:

- Develop and capitalize on new business opportunities by addressing increasing patient demand
- Shorter treatment times

*Example scenario

The efficient workflow protocol that leads to shorter treatment cycles for the patient.

Straumann® BLX

Confidence beyond Immediacy.

The high-performance material Roxolid® combined with the clinically proven surface, SLActive® provides confidence beyond immediacy even in the most demanding clinical situations.

TorcFit™ Connection

Hybrid internal conical connection with high flexibility and strength.

| Esthetic Ease Concept

One prosthetic line for enhanced simplicity and versatility to deliver natural esthetics.



Straumann® CARES®

Digital Solutions

- Fully digitally integrated workflows and seamless connectivity



| Roxolid®

Reducing invasiveness with smaller implants

- Preserves vital structures and vascularization^{2, 3}
- Reduces invasive grafting procedures^{2, 3}
- Increases treatment options in challenging anatomical situations and narrow interdental spaces^{3, 4}

Simplified drill protocol

Less steps and flexible sequence with a new drill design for minimized heat generation.

| Ø 3.75 mm for all indications

Slim core and fully tapered implant design enabling bone-preserving protocols with a Ø 3.75 mm implant applicable for all indications.

| Dynamic Bone Management™

Straumann® BLX Implant features unique properties that enable enhanced control over insertion torque to achieve optimal primary stability and ensure confidence in immediate protocols.

| SLActive®

Designed to maximize treatment success and predictability

- Reduce healing time to 3–4 weeks^{5, 6}
- Deliver high long-term predictability in immediate loading⁷
- Provide outstanding success rates in compromised patients^{8, 9, 10, 11}
- Promote bone regeneration, even at compromised sites^{12, 13, 14, 15}

Controlled bone management
in all surgical situations.



The intelligent implant design allows for Dynamic Bone Management™ and makes immediate protocols achievable and predictable in all bone types.



Dynamic Bone Management.™

- The unique “Cut, Collect, Condense”-properties redistribute native bone along the implant leading to high bone-to-implant contact
- Bi-directional cutting elements enable control over insertion torque to achieve optimal primary stability in all bone types
- Reduced implant neck designed for stress-free crestal bone
- Enhanced drill design with minimized heat generation means fewer drilling steps





Esthetic Ease Concept.




Simplified but versatile portfolio, featuring one connection and under-contoured prosthetics to deliver esthetics with ease.

SIMPLICITY MEETS VERSATILITY

- One connection for simplified restorative workflows
- Slim and under-contoured abutment designs for predictable esthetic outcomes
- Consistent emergence profile for easy soft-tissue conditioning



ONE CONNECTION





Real Confidence.

To respond to market needs for shorter treatment times and simplified workflows, Straumann® BLX delivers outstanding versatility, reliability and predictable outcomes for immediate function.



INNOVATION BACKED BY SCIENTIFIC EVIDENCE DELIVERS RESULTS YOU CAN TRUST

- Proven technologies, Roxolid® material and SLActive® surface
- Scientific backing and educational network
- Swiss precision engineering and quality
- Fully digitally integrated workflows
- Complete range of biomaterials

Perfected for immediacy and also an excellent solution for all other indications. The Straumann® BLX intelligent implant design makes immediate protocols achievable and predictable in all bone types.

Immediacy from the
confidence brand.



Contact your local Straumann representative now or visit:
www.confidence-in-you.com

Straumann® BLX for confidence beyond Immediacy.

- The winning combination of innovative design and surface technology with high-performance material
- Combines strength, healing power and optimal primary stability
- Flexible and simplified workflows
- Predictable outcomes with results you can trust

REFERENCES

1 Moraschini V, et al. Immediate versus conventional loaded single implants in the posterior mandible: a meta-analysis of randomized controlled trials. *Int J Oral Maxillofac Surg* 2016 Jan;45(1):85-92. Literature Review, full list available. 2 Ioannidis A, Gallucci GO, Jung RE, Borzangy S, Hämmerle CH, Benic GI. Titanium-zirconium narrow-diameter versus titanium regular-diameter implants for anterior and premolar single crowns: 3-year results of a randomized controlled clinical study. *J Clin Periodontol*. 2015 Nov;42(11):1060-70. doi:10.1111/jcpe.12468. Epub 2015 Nov 14. 3 Al-Nawas B, Domagala P, Fragola G, Freiburger P, Ortiz-Vigón A, Rousseau P, Tondela J. A Prospective Noninterventive Study to Evaluate Survival and Success of Reduced Diameter Implants Made From Titanium-Zirconium Alloy. *J Oral Implantol*. 2015 Aug;41(4):e118-25. doi: 10.1563/AAID-JOI-D-13-00149. Epub 2014 Mar 25. 4 Altuna P, Lucas-Taulé E, Gargallo-Albiol J, Figueras-Álvarez O, Hernández-Alfaro F, Nart J. Clinical evidence on titanium-zirconium dental implants: a systematic review and meta-analysis. *Int J Oral Maxillofac Surg*. 2016 Jul;45(7):842-50. doi: 10.1016/j.ijom.2016.01.004. Epub 2016 Feb 3. 5 Lang NP, Salvi GE, Huynh-Ba G, Ivanovski S, Donos N, Bosshardt DD. Early osseointegration to hydrophilic and hydrophobic implant surfaces in humans. *Clin Oral Implants Res*. 2011 Apr;22(4):349-56. doi:10.1111/j.1600-0501.2011.02172.x 6 Oates TW, Valderrama P, Bischof M, Nedir R, Jones A, Simpson J, Toutenburg H, Cochran DL. Enhanced implant stability with a chemically modified SLA surface: a randomized pilot study. *Int J Oral Maxillofac Implants*. 2007 Sep-Oct;22(5):755-60 7 Nicolau P, Guerra F, Reis R, Krafft T, Benz K, Jackowski J. 10-year outcomes with immediate and early-loaded implants with a chemically modified SLA surface, Quintessence Int. (article in press). 8 Nelson, K., Stricker, A., Raguse, J.-D. Nahles, S. (2016), Rehabilitation of irradiated patients with chemically modified and conventional SLA implants: a clinical clarification. *J Oral Rehabil*, 43: 871–872. doi:10.1111/joor.12434. 9 Heberer S, Kilic S, Hossamo J, Raguse J-D, Nelson K. Rehabilitation of irradiated patients with modified and conventional sandblasted, acid-etched implants: preliminary results of a split-mouth study. *Clin. Oral Impl. Res.* 22, 2011; 546–551. 10 Cabrera-Domínguez J, Castellanos-Cosano L, Torres-Lagares D, Machuca-Portillo G. A Prospective Case-Control Clinical Study of Titanium-Zirconium Alloy Implants with a Hydrophilic Surface in Patients with Type 2 Diabetes Mellitus. *Int J Oral Maxillofac Implants*. 2017 Sep-Oct;32(5):1135-1144. doi: 10.11607/jomi.5577; 11 Cabrera-Domínguez J. A prospective, two-year clinical trial of titanium-zirconium alloy implants (Roxolid® Straumann®) with hydrophilic surface (SLActive®) in patients with Type 2 Diabetes Mellitus. presented during 26th Annual Scientific Meeting of the European Association of Osseointegration – 5-7 Oct 2017, Madrid, Spain 12 Kuchler U, Chappuis V, Bornstein mm, Siewczyk M, Gruber R, Maestre L, Buser D. Development of Implant Stability Quotient values of implants placed with simultaneous sinus floor elevation – results of a prospective study with 109 implants. *Clin. Oral Impl. Res.* 28, 2017, 109–115 13 Schwarz, F., et al., Bone regeneration in dehiscence-type defects at chemically modified (SLActive®) and conventional SLA® titanium implants: a pilot study in dogs. *J Clin.Periodontol.* 34.1 (2007): 78–86 14 Schwarz, F., et al., Bone regeneration in dehiscence-type defects at non-submerged and submerged chemically modified (SLActive®) and conventional SLA titanium implants: an immunohistochemical study in dogs. *J Clin.Periodontol.* 35.1 (2008): 64–75. 15 El Chaar E, Zhang L, Zhou Y, et al. (Superior) osseointegration of superhydrophilic implants placed in defect grafted bones. *Int J Oral Maxillofac Implants* 2018. (article in press).



International Headquarters

Institut Straumann AG
Peter Merian-Weg 12
CH-4002 Basel, Switzerland
Phone +41 (0)61 965 11 11
Fax +41 (0)61 965 11 01
www.straumann.com